

BODILINK® HEAD SUPPORT PREMIUM ORDER FORM



HCPCS CODING: E0955 - CUSHIONED HEADREST
 HCPCS CODING: E1028 - REMOVABLE HARDWARE
 HCPCS CODING: E1028 - FLIP-DOWN HARDWARE

Since U.S. Medicare coding is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process.



PT Max 8 Hardware
Shown on Acta-Back®

Company Name/ACCT #: _____

P.O. Number: _____

Requested By: _____

Phone: _____ Fax: _____

Email: _____

Ship To: _____

Client Reference: _____

For best results, do not fill out in your browser. Interactive form should be completed using Adobe Reader after saving to your local drive. Then email or print and fax to Customer Support at orders.comfort@permobil.com or 1.406.522.8563.

DO NOT SEND PROTECTED HEALTH INFORMATION. IT IS NOT NEEDED TO MAKE THE PRODUCT YOU ARE REQUESTING.

1 BODILINK HEAD SUPPORT PAD

- You may make one selection per section.
- Choose pad firmness and size in 1.1.
- Choose a cover material in 1.2.

1.1

PREMIUM PAD SIZE

HCPCS CODING: E0955 MSRP \$249.00

Length (L) refers to the actual size dimension of the support from bottom to top edge.

		Pad Firmness	
L	6W x 3L	Soft	BL-HSP1-6W3L
		Extra Soft	BL-HSP2-6W3L
L	8W x 4L	Soft	BL-HSP1-8W4L
		Extra Soft	BL-HSP2-8W4L
L	10W x 5L <i>with adjustable wings</i>	Soft	BL-HSP1-10W5L
		Extra Soft	BL-HSP2-10W5L
L	12W x 5L <i>with adjustable wings</i>	Soft	BL-HSP1-12W5L
		Extra Soft	BL-HSP2-12W5L
L	14W x 5L <i>with adjustable wings</i>	Soft	BL-HSP1-14W5L
		Extra Soft	BL-HSP2-14W5L
L	18W x 5L <i>with adjustable wings</i>	Soft	BL-HSP1-18W5L
		Extra Soft	BL-HSP2-18W5L

1.2

PREMIUM PAD COVER MATERIAL

	FOR SKIN PROTECTION & SHEAR REDUCTION. NOT FLUID RESISTANT.
	FOR FLUID PROTECTION & AN EASILY CLEANED SURFACE.
	FOR PATIENT COMFORT.

	GlideWear® (GLIDEWEAR) (add \$26.00)
	Comfort-Tek® (COMFORT-TEK)
	Stretch-Air® (STRETCH-AIR)



Premium Head Support Pads angles can be formed by adjusting the plates in the pad.

2 BODILINK™ HEAD SUPPORT HARDWARE

- Make one selection per column in subsection A or B.
- If A is chosen, choose a stem length option.
- Choose which mounting block option you require in 2.2.

2.1 PT HARDWARE

RANGE OF MOTION

PT MODEL	TOTAL # OF LINKS	TOTAL LINK LENGTH	MAX ANTERIOR REACH	MAX POSTERIOR REACH	MAX LATERAL REACH
MID 6	2	6"	5.5"	1"	N/A
LAT 6	3	6"	5.5"	1"	6"
MID 8	2	8"	7.5"	1"	N/A
LAT 8	3	8"	7.5"	1"	8"
MAX 8	4	8"	7.5"	3"	8"

A. PT HARDWARE, REMOVABLE

HCPCS E1028



SHAPE	MSRP	PART NUMBER
MID 6	\$235.00	(BL-HS-PT1RMMID6)
LAT 6	\$273.00	(BL-HS-PT1RMLAT6)
MID 8*	\$235.00	(BL-HS-PT1RMMID8)
LAT 8*	\$273.00	(BL-HS-PT1RMLAT8)
MAX 8**	\$344.00	(BL-HS-PT1RMMAX8)

*Recommended for Acta-Relief® back supports. See options for mounting hardware selections.
 **Provides the greatest amount of articulation and range of motion.

PT HARDWARE, REMOVABLE, STEM LENGTH OPTION



STEM LENGTH	MSRP
STANDARD 7" STEM	\$0.00
***3" STEM	\$0.00

PART NUMBER
(HS-STEM-7)
(HS-STEM-3)

***3" Stem is recommended for use with power mounting where there may be special constraints with power chair functions.

B. PT HARDWARE, FLIP-DOWN

HCPCS E1028

Flip-Down Hardware has 7" Stem



SHAPE	MSRP	PART NUMBER
MID 6	\$287.00	(BL-HS-PT1FDMID6)
LAT 6	\$326.00	(BL-HS-PT1FDLAT6)
MID 8*	\$287.00	(BL-HS-PT1FDMID8)
LAT 8*	\$326.00	(BL-HS-PT1FDLAT8)
MAX 8**	\$375.00	(BL-HS-PT1FDMAX8)

*Recommended for Acta-Relief® back supports. See options for mounting hardware selections.
 **Provides the greatest amount of articulation and range of motion.

2.2 HEAD SUPPORT MOUNTING OPTION



MOUNTING BLOCK-POSTERIOR A or B
 A: For Acta-Back®, Acta-Embrace®, Visco-Back, and Acta-Relief with Power Mount Only



MOUNTING BLOCK-ANTERIOR A
 For ROHO® AGILITY and other OEM back supports



MOUNTING BLOCK-QUANTUM A
 For Quantum Tru-Comfort back supports

B: For Radius Back® or OEM wood/t-nut back supports



MOUNTING BLOCK-PERMOBIL A
 For Permobil Corpus 3G Back Supports (23"- 28"H)



MOUNTING BLOCK-PERMOBIL B
 For Permobil Corpus 3G Back Supports (20"H)







BRACKET WITH MOUNTING BLOCK-ACTA-RELIEF
 18W Shown Above
 For Non-Power Mount Only.

MOUNTING BLOCK OR BRACKET WITH MOUNTING BLOCK	MSRP	PART NUMBER
POSTERIOR A MOUNTING BLOCK	\$0.00	(HS-MB-POSTA)
POSTERIOR B MOUNTING BLOCK	\$0.00	(HS-MB-POSTB)
ANTERIOR A MOUNTING BLOCK	\$0.00	(HS-MB-ANTA)
QUANTUM A MOUNTING BLOCK	\$90.00	(HS-MB-PW1A)
PERMOBIL A MOUNTING BLOCK	\$90.00	(HS-MB-PW2A)
PERMOBIL B MOUNTING BLOCK	\$90.00	(HS-MB-PW2B)
*ACTA-RELIEF 14W BRACKET WITH MOUNTING BLOCK	\$72.00	(HS-MB-AR14W)
*ACTA-RELIEF 16W BRACKET WITH MOUNTING BLOCK	\$72.00	(HS-MB-AR16W)
*ACTA-RELIEF 18W BRACKET WITH MOUNTING BLOCK	\$72.00	(HS-MB-AR18W)
*ACTA-RELIEF 20W BRACKET WITH MOUNTING BLOCK	\$72.00	(HS-MB-AR20W)


*For Acta-Relief, includes mounting bracket and Posterior A block.
 If you already have an Acta-Relief® head support bracket, order Posterior A Block

3 HEAD SUPPORT ADDITIONAL MOUNTING HARDWARE


3.1 EXTRA LINKS & WING LOCKS

	LINK LENGTH	MSRP	EXTRA PART					
	EXTRA 2" LINK	\$60.00	(BL-HS-PT1LINK2)					
	LINK LENGTH	MSRP	EXTRA PART					
	EXTRA 3" LINK	\$60.00	(BL-HS-PT1LINK3)					
	LINK LENGTH	MSRP	EXTRA PART					
	EXTRA 4" LINK	\$60.00	(BL-HS-PT1LINK4)					
	TOOL-LESS WING LOCK	MSRP	EXTRA PART		EXTRA PART	EXTRA PART	EXTRA PART	
	1 LOCK	\$23.00/each	(BL-HS-52WINGA1)		(BL-HS-52WINGA1)	(BL-HS-52WINGA1)	(BL-HS-52WINGA1)	


3.2 EXTRA MOUNTING HARDWARE




MOUNTING BLOCK-POSTERIOR A or B
A: For Acta-Back®, Acta-Embrace, Visco-Back, and Acta-Relief with Power Mount Only
B: For Radius Back® or OEM wood/t-nut back supports




MOUNTING BLOCK-ANTERIOR A
 For ROHO® AGILITY and other OEM back supports




MOUNTING BLOCK-QUANTUM A
 For Quantum Tru-Comfort back supports



MOUNTING BLOCK-PERMOBIL B
 For Permobil Corpus 3G Back Supports (20"H)



BRACKET-ACTA-RELIEF
 18W Shown Above



MOUNTING BLOCK-PERMOBIL A
 For Permobil Corpus 3G Back Supports (23"- 28"H)

MOUNTING BLOCK OR BRACKET	MSRP	PART NUMBER
POSTERIOR A MOUNTING BLOCK	\$55.00	(BL-HS-MBPOSTA)
POSTERIOR B MOUNTING BLOCK	\$55.00	(BL-HS-MBPOSTB)
ANTERIOR A MOUNTING BLOCK	\$55.00	(BL-HS-MBANTA)
QUANTUM A MOUNTING BLOCK	\$95.00	(BL-HS-MBPW1A)
PERMOBIL A MOUNTING BLOCK	\$95.00	(BL-HS-MBPW2A)
PERMOBIL B MOUNTING BLOCK	\$95.00	(BL-HS-MBPW2B)
*ACTA-RELIEF 14W BRACKET	\$160.00	(BL-HS-AR14W)
*ACTA-RELIEF 16W BRACKET	\$160.00	(BL-HS-AR16W)
*ACTA-RELIEF 18W BRACKET	\$160.00	(BL-HS-AR18W)
*ACTA-RELIEF 20W BRACKET	\$160.00	(BL-HS-AR20W)

*Bracket only, does not include mounting block.

ADDITIONAL	MSRP	EXTRA PART
1/4" SPACER PLATE	\$39.00/each	(BL-HS-SPCRA)

The following are trademarks and registered trademarks of Permobil Seating and Positioning: ACTA-BACK®, ACTA-EMBRACE®, AGILITY®, BodiLink®, Comfort-Tek®, RADIUS BACK®, and Stretch-Air®. Glidewear® is a registered trademark of Tamarack Habilitation Technologies, Inc.