

Power wheelchairs


K450 MX



permobil

High or low, explore more

The power seat-to-floor function allows easy access to the ground for story time, play, or transfers



Specifications

Maximum user weight	125 lb
Maximum speed	5 mph
Range	17 mi*
Base width	25"
Base length	35"
Minimum turning radius	28"
Ground clearance	2.5"
Recommended battery type	Group 34
Weight including batteries	337 lb
Drive electronics	90A R-net
Seat-to-floor height	3" – 26" (elevator)
Optional seat elevator	8"
Power tilt	0° – 45°
Manual recline	85° – 105°
Manual legrest elevation	90° – 125°
Seat depth	10" – 18" (by 1" incr.)



Seat width 10", 12", 14", 16"

Backrest height 14" – 22.5"

**Actual driving range will vary based on driving and battery conditions.*

Specifications may vary based upon configuration. Please refer to your owners manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions.

Color options

