

F3 Corpus

AADL #W688

permobil



Dealer Information

* **Contact:** _____
Dealer Code: _____
* **Dealer Name:** _____
Address: _____
* **City:** _____
* **Province/Postal Code:** _____
* **Phone#:** _____
* **Fax#:** _____
PO#: _____
Email Address: _____


Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* **First Name:** _____
* **Last Name:** _____
Diagnosis: _____
Funding Source: _____
Client Age: _____

Client Measurements

* **Weight:** _____
* **Height:** _____
A) Top of Shoulders: _____
B) Chest Depth: _____
C) Chest Width: _____
D) Seat Depth: _____
E) Top of Head: _____
F) Elbow to Hand: _____
G) Seat Pan to Elbow: _____
H) Hip Width: _____
* **I) Knee to Foot:** _____



Three column pricing on order/quote:

Please send order/quote to fax#: (800) 231-3256
Email to: sales@permobil.com

Prices effective **October 1, 2021** and payable in Canadian dollars.

Permobil Ltd.
12 - 75 Mary St.
Aurora, Ontario
L4G 1G3
Tel: (800) 736-0925
Fax: (800) 231-3256
www.permobil.ca

✓

Base Options

I107433-99-0	Set Seat Height at 18.5" <i>If ordering Power Adjustable Seat Height, 1" of elevator travel is used to set seat at 18.5".</i>	No Charge
I107434-99-0	Set Seat Height at 19.5" <i>If ordering Power Adjustable Seat Height, 2" of elevator travel is used to set seat at 19.5".</i>	No Charge

MUST CHOOSE BATTERIES:

Part Number	Description	AADL - CA
BATTERY OPTIONS:		
I107031-99-0	Batteries, Grp 24 (72Ah) Sealed Gel, x2, Installed	965.00
I107032-99-0	Batteries, Grp 34 (60Ah) Sealed Gel, x2, Installed	953.00

OPTIONAL ITEMS:

Part Number	Description	AADL - CA
I107459-99-0	Essentials Box with Reflectors (Client Cost) <i>Please Note: The total backrest angle will be limited to 150° to avoid contact between backrest components and the Essentials Box.</i>	552.00
I108417-99-0	Essentials Box with Lights (Client Cost) <i>Please Note: The total backrest angle will be limited to 150° to avoid contact between backrest components and the Essentials Box. Must also select Light Kit.</i>	919.00
I108418-99-0	LED Light Kit for F3 (Client Cost)	1,195.00
✓ I11143	Enhanced Steering Perform Unit (ESP) R-net	1,223.00
I10715	Compact Laptop Ventilator Tray - Corpus Seat <i>Mounts to rear of Corpus backrest hinge. Accommodates Pulmonetic Systems / CareFusion LTV Series and similar sized Laptop Ventilators. The tray is 330mm Tall, adjustable Width between 200mm - 330mm, and adjustable Depth between 50mm - 165mm. Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option. Additional drive restrictions may apply.</i>	1,732.00
I10787	Battery Tray for Laptop Ventilator Tray <i>Mounts to rear of Laptop Ventilator Tray. Measures 290mm W x 170mm H x 42mm D. Accommodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	277.00
I10788	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 175mm Tall. Width infinitely adjustable 90mm - 170mm. Depth infinitely adjustable 95mm - 125mm. Custom Option.</i>	391.00
ID51113	Programming Key for R-net - Permobil Version	540.00

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	AADL - CA
PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i>		
I107779-99-0	Color Joystick Module, (CJSM) <i>This R-net color joystick features: A Large High-Resolution 3.5" VGA color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	1,375.00
I11186	Permobil Joystick Module <i>Designed exclusively for Permobil, this R-net color joystick features: a 2.8" color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	No Charge
I107829-99-0	3.5mm Y-Splitter Cable for PJSM/CJSM <i>For use with the PJSM or CJSM Mode/Profile jack. Allows 2 mechanical switches to control 4 separate functions. Requires PC Programming to assign the short hold and long hold commands to each switch. Switches sold separately.</i>	43.00
I105700-99-0	R-net Remote Colour JS w/ Lights, Heavy Duty Gimbal <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Colour LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	2,103.00

Joystick Options

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	AADL - CA
I10045	Cone Shaped Joystick Knob - Standard	No Charge
I100960-99-0	Chin Cup For Joystick Knob	123.00
I100961-99-0	"T" Handle For Joystick Knob	181.00
I100962-99-0	"Large Ball" For Joystick Knob	123.00
I100963-99-0	"Softball" For Joystick Knob	165.00
I100964-99-0	"Mushroom" For Joystick Knob	165.00
I100965-99-0	"Stick" For Joystick Knob	165.00
I105442-99-0	1.75" Black Foam Ball For Joystick Knob	123.00
I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	158.00
I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	158.00
I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	158.00

Part Number	Description	AADL - CA
I10664	Joystick Mount Right - Fixed	No Charge
I10665	Joystick Mount Left - Fixed	No Charge
I108037-99-0	Slimline Retractable Joystick Mount - Right	466.00
I108038-99-0	Slimline Retractable Joystick Mount - Left	466.00
I11232	Adjustable Height Panel Bracket <i>Provides 1.77" of vertical adjustment below the armrest and 25° of angle adjustment. Compatible with all Permobil joystick mounts except for the M1 Fixed Joystick Mount.</i>	361.00
I105210-99-0	Mount Standard Joystick on Corpus Back Hinge <i>Standard R-net Joystick is mounted to backrest hinge. Can be chosen instead of PN I105408-99-0. Mount on Left: _____ Mount on Right: _____</i>	323.00
I108268-99-0	Bodypoint Tri-Lock Midline Mount Assembly <i>Includes L-shaped joystick arm and height adjustable mounting clamp. Joystick mount sold separately (select from the mounts listed below). Please indicate if you want midline mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	711.00
I108280-99-0	Bodypoint Tri-Lock Mount for R-net JS	291.00
I108279-99-0	Bodypoint Tri-Lock Mount for Compact Joystick	291.00
I108310-99-0	Bodypoint Tri-Lock Mount for Omni	332.00

Part Number	Description	AADL CA
I10842	PERMOFix Mount for Mini/Compact Joystick <i>Includes a shim that allows the PERMOFix joint to be mounted on the joystick arm. Provides an adjustable mount for a Mini Joystick or Compact Joystick.</i>	291.00

Part Number	Description	AADL - CA
✓ I10411	Expandable Controller - R-net <i>Must be ordered with all expandable/remote joysticks.</i>	No Charge
✓ I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller.</i>	No Charge
I105408-99-0	Attendant Control R-net	1,346.00
I10582	Output Module (IOM) for R-net Control System (Client Cost) <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	1,801.00

Seating System

CHOOSE POSTERIOR TILT/RECLINE PACKAGE:

Part Number	Description	AADL - CA
	PLEASE NOTE: <i>A half speed limitation is imposed when the Total Backrest Angle (TBA) > 150°. TBA = Tilt Angle + Backrest Recline Angle. Example 1: 45° Tilt + 100° Recline = 145° TBA (Full Speed). Example 2: 35° Tilt + 120° Recline = 155° TBA (Half Speed).</i>	

Seating System

I11131	Corpus Seat w/ No Tilt & 180° Power Recline	AADL #W993 for Power Tilt	9,437.00
I11132	Corpus Seat w/ 50° Tilt & Manual Recline <i>Manual Recline angle limited to 120°.</i>		8,461.00
I11133	Corpus Seat w/ 50° Tilt & 180° Power Recline		13,941.00
I11211	Corpus Seat w/ Manual Recline		2,495.00

ACTIVE REACH PACKAGES:

Part Number	Description	AADL - CA
	PLEASE NOTE: <i>Active Reach should be recommended by a licensed clinical professional and/or ATP. While operating Active Reach, positioning equipment such as lap belts and chest supports should be considered to meet the user's specific safety needs. Note: A low speed limitation is imposed from 5° to 10° of Active Reach, an extra low speed limitation is imposed from 10° to 20° of Active Reach, and a drive lockout is imposed above 20° of Active Reach.</i>	
I11141	5° Active Reach (Anterior Tilt) <i>When ordered with a Power Adjustable Seat Height, tilt range will be -5° to +50°. When ordered without a Power Adjustable Seat Height, tilt range will be -5° to +45°. If this option is selected with a fixed seat, the minimum seat to floor height at 0° of tilt will be approximately 19.5". Power Legrest recommended for full functionality. Minimum legrest angle will be ~ 97° when manual legrest is selected.</i>	No Charge
I107090-99-0	10° Active Reach (Anterior Tilt) <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height and Power Tilt. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Power Legrest recommended for full functionality. Minimum legrest angle will be ~ 103° when manual legrest is selected.</i>	No Charge
I107091-99-0	20° Active Reach (Anterior Tilt) (Client Cost) <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height, Power Tilt, Power Recline and Power Legrest Elevation. Must select Corpus VS Power Elevating Legrest, Corpus VS Footplates and Corpus VS Knee Support from VS Legrest Options section. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Not compatible with retractable positioning belts.</i>	3,047.00
I11100	30° Active Reach w/ Transfer (Anterior Tilt) (Client Cost) <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height, Power Tilt, Power Recline and Power Legrest Elevation. Must select Corpus VS Power Elevating Legrest, Corpus VS Footplates and Corpus VS Knee Support from VS Legrest Options section. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Not compatible with retractable positioning belts.</i>	3,426.00

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	AADL - CA
I10418	Single Seat Function Control Switch for R-net <i>Does not include thru drive control operation of power seat functions.</i>	No Charge
I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options.</i>	1,711.00
I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options.</i>	2,851.00

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	AADL - CA
I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge
I102799-99-0	Toggle Switches on Seat Function Control Kit	305.00
I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge

Seating System Options

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	AADL - CA
I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge
I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge
I106256-99-0	Mount ICS Switchbox on Backrest (Client Cost) <i>Mount on Left: _____ Mount on Right: _____</i>	171.00

OPTIONAL ITEMS:

Part Number	Description	AADL - CA
I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,346.00
I10750	Memory Seat Program, Corpus w/ R-net	125.00
I106209-99-0	Legrest Movement Linked to Recline Electronically (Client Cost) <i>Power Legrest will raise when Power Recline moves backward and will lower when back moves forward. If Power Legrest is operated, Power Recline will NOT automatically reposition. Feature works through ICS Switchbox and Driver Control.</i>	50.00
I10898	Enable Independent Repositioning Mode <i>IRM allows for sequential movement of power seat functions into a posterior position (35° tilt, 115° recline, 135° legrest) and returns to upright position (5° tilt, 100° recline, 110° legrest). This can be accessed using a standard drive control or a single switch. On the Permobil Joystick Module, IRM is accessible through the mode/profile jack (long switch hold). Please see the Alt Drive Accessories section for switch options.</i>	No Charge

Backrest Options

CHOOSE SLIDING BACKREST KIT:

Part Number	Description	AADL - CA
I105067-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No Charge

CHOOSE BACKREST TYPE:

Part Number	Description Price	AADL - CA
	BACKREST HEIGHT NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Backs include lateral wedges and lumbar pads.</i>	

Part Number	Description	AADL - CA
I106545-99-0	Solid Corpus Back - 14"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	No Charge
I106546-99-0	Solid Corpus Back - 16"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	No Charge
I106547-99-0	Solid Corpus Back - 18"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	No Charge
I106548-99-0	Solid Corpus Back - 20"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	No Charge

PLEASE NOTE:

Taller backrests listed below include: lateral wedges, foam lumbar pads, UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided. * Foam lumbar pads not included with Solid Backs.*

I106532-99-0	Solid Corpus Back - 14"W x 23" - 28"T (C3GERGBCA)	No Charge
I11107	Solid Corpus Back - 16"W x 23" - 28"T (C3GERGBCA)	No Charge
I11108	Solid Corpus Back - 18"W x 23" - 28"T (C3GERGBCA)	No Charge
I11109	Solid Corpus Back - 20"W x 23" - 28"T (C3GERGBCA)	No Charge

MUST ALSO SELECT BACK HEIGHT:

Set Corpus Backrest at 23" _____
Set Corpus Backrest at 24" _____

Backrest Options

- Set Corpus Backrest at 25" _____
- Set Corpus Backrest at 26" _____
- Set Corpus Backrest at 27" _____
- Set Corpus Backrest at 28" _____

CHOOSE BACKREST CUSHION:

Part Number	Description	AADL - CA
I11335	Leatherette Cushion - Corpus Back <i>Includes Dual Density Backrest Foam.</i>	No Charge
I11336	Stretch-Air Cushion - Corpus Back <i>Includes Dual Density Backrest Foam with Stretch-Air breathable fabric cover.</i>	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description (Seating Benefits)	AADL - CA
	ROHO AGILITY NOTE: <i>The ROHO AGILITY Back Systems are available in widths of 14", 16", 18" and 20" and shell heights of 16", 18" and 20". Direct Mount system includes backrest shell, cushion, air inserts, UT rail for headrest and Direct Mount adapter for Corpus.</i>	
I11258	ROHO AGILITY Max Contour Back System - 16" Tall <i>Backrest can be installed at a height (seat pan to top of AGILITY) between 22" and 25". Specify Shell Width (14", 16", 18", 20"W): _____</i>	1,155.00
I11259	ROHO AGILITY Max Contour Direct Mount - 18" Tall <i>Backrest can be installed at a height (seat pan to top of AGILITY) between 22" and 26". Specify Shell Width (14", 16", 18", 20"W): _____</i>	1,155.00
I11260	ROHO AGILITY Max Contour Direct Mount - 20" Tall <i>Backrest can be installed at a height (seat pan to top of AGILITY) between 24" and 27". Specify Shell Width (14", 16", 18", 20"W): _____</i>	1,155.00
	MUST ALSO SELECT BACK HEIGHT:	
	Set AGILITY Backrest at 22" _____	
	Set AGILITY Backrest at 23" _____	
	Set AGILITY Backrest at 24" _____	
	Set AGILITY Backrest at 25" _____	
	Set AGILITY Backrest at 26" _____	
	Set AGILITY Backrest at 27" _____	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description (Seating Benefits)	AADL - CA
I107488-99-0	Permobil 3" x 5" Lateral Supports - Pair <i>Must also select I107490-99-0. Only compatible with Corpus backrest shells. Backrest Lateral Wedges will be omitted.</i>	473.00
I107489-99-0	Permobil 5" x 7" Lateral Supports - Pair <i>Must also select I107490-99-0. Only compatible with Corpus backrest shells. Backrest Lateral Wedges will be omitted.</i>	473.00
I107490-99-0	Adj/Removable Lateral Link Hardware - Pair	648.00

OPTIONAL ITEMS:

Part Number	Description	AADL - CA
I105782-99-0	Tubular Backrest Frame - Corpus (1" Canes) <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Canes are attached with 2" long UniTrack linkage disks that can rotate to provide space for backrest hardware. The tops of the canes can be adjusted up to a height of 27" above the seat pan. Depending on the aftermarket backrest, several inches of seat depth may be lost when mounting this system.</i>	491.00
I106072-99-0	Direct Backrest Frame - Corpus <i>Used to mount ROHO Agility, Jay 3, Matrx Elite, Varilite Icon and Ride Java backrests ranging in width from 14" to 20". Slotted brackets and 1" UniTrack linkage disks provide a close, adjustable mount. Backrest pitch can be tuned by rotating the linkage disks. Some seat depth may be lost due to the position and orientation of the backrest.</i>	491.00

Seat Sizing

CHOOSE SEAT WIDTH:

Part Number	Description	AADL - CA
I107023-99-0	Set Seat Width at 17": _____	No Charge
I107024-99-0	Set Seat Width at 19": _____	No Charge
I107025-99-0	Set Seat Width at 21": _____	No Charge
I107026-99-0	Set Seat Width at 23": _____	No Charge

CHOOSE SEAT DEPTH:

Part Number	Description	AADL - CA
	PLEASE NOTE: <i>There are three seat depth kits, each providing three seat depth settings: (1) 14" - 16", (2) 17" - 19" and (3) 20" - 22". Please select the required seat depth and the appropriate kit will be provided.</i>	
I107013-99-0	Seat Depth 14"	No Charge
I107014-99-0	Seat Depth 15"	No Charge
I107015-99-0	Seat Depth 16"	No Charge
I107016-99-0	Seat Depth 17"	No Charge
I107017-99-0	Seat Depth 18"	No Charge
I107018-99-0	Seat Depth 19"	No Charge
I107019-99-0	Seat Depth 20"	No Charge
I107020-99-0	Seat Depth 21"	No Charge
I107021-99-0	Seat Depth 22"	No Charge

CHOOSE SEAT CUSHION:

Part Number	Description (Seating Benefits)	AADL - CA
	NOTE: <i>A Corpus II Ergonomic Seat Cushion will be provided for select seat sizes.</i>	
I11337	Corpus Ergo Cushion - Leatherette Cover	623.00
I11338	Corpus Ergo Cushion - Stretch-Air Cover	623.00
I10965	Corpus Ergo Air Cushion (C3GRHOSCA) <i>Features ROHO Dry Floatation technology. The removable ROHO cell pad provides adjustable skin protection while the contoured foam and thigh support wedges help maintain proper positioning and stability. Includes removable stretch fabric cover. Available in widths of 17", 19" and 21" and depths of 17", 19" and 21", as well as 17" x 23". Not Available in 21" x 17"</i>	661.00
I11113	Incontinence Cover for Ergo Cushion <i>Made from a stretch Dartex fabric and suitable for all Ergo cushions including the Corpus Ergo Air. The cover is cut to the contours of the Ergo cushion and is designed to be used as an outer layer of protection. Cover size will be matched to cushion size. Sent as a part. Part Number will be provided on the quote / order confirmation based on the Seat Width and Seat Depth selected.</i>	198.00
I11166	ROHO QUADTRO SELECT HIGH PROFILE	769.00
I11168	ROHO HIGH PROFILE Single Compartment	630.00
I11170	ROHO HIGH PROFILE Sensor Ready <i>Single compartment cushion designed to work with Smart Check Device I107505-99-0. Not available in all seat widths and depths.</i>	630.00
I107505-99-0	ROHO Smart Check Device <i>Sensing device that allows users to setup and maintain proper inflation pressure in ROHO Sensor Ready cushions. Only available when ordered with a ROHO Sensor Ready cushion.</i>	295.00
I11173	Incontinence Cover for ROHO Cushion <i>Cover size will be matched to ROHO cushion size. Sent as a part.</i>	144.00
I105118-99-0	Omit Seat Cushion	No Charge

Armrest Options

CHOOSE ARMBAR KIT AND MOUNTING POSITION:

Part Number	Description	AADL - CA
	PLEASE NOTE: <i>Narrow Armbar Kit can be adjusted from 15"W to 19"W between the armrest pads without any additional parts. Wide Armbar Kit can be adjusted from 19"W to 23"W between the armrest pads without any additional parts. Please select both Armbar Kit and Width Setting below.</i>	
I106543-99-0	Narrow Armbar Kit, Inf Adjustable 15"W - 19"W	No Charge
I106544-99-0	Wide Armrest Kit, Inf Adjustable 19"W - 23"W <i>Required when selecting a 20" wide backrest. Set the Width Between the Armrests to: _____ Must be wider than backrest width to ensure proper fit. Please round to nearest 1".</i>	No Charge

CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number	Description	AADL - CA
✓ I105069-99-0	Height Adjustable Armrests SET ARMREST HEIGHT AT: _____ <i>Corpus Armrests are infinitely adjustable from 8" - 15", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.</i>	No Charge

CHOOSE ARMREST SWIVEL FIXED OR ENABLED:

Part Number	Description	AADL - CA
I105071-99-0	Left Armrest Bracket - Fixed <i>Medial swivel cannot be enabled when selecting this option. Must be selected when ordering arm trough adapters.</i>	No Charge
I105070-99-0	Right Armrest Bracket - Fixed <i>Medial swivel cannot be enabled when selecting this option. Must be selected when ordering arm trough adapters.</i>	No Charge
I105073-99-0	Left Armrest Bracket - Medial Swivel Enabled	No Charge
I105072-99-0	Right Armrest Bracket - Medial Swivel Enabled	No Charge

CHOOSE ARMREST PAD RIGHT:

Part Number	Description	AADL - CA
I105100-99-0	4"x10" Corpus Armrest Pad Right - Leatherette	No Charge
I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge
I105105-99-0	4"x16" Corpus Armrest Pad Right - Leatherette	No Charge
I105999-99-0	4"x18" Corpus HD Armrest Pad Right - Leatherette	No Charge

CHOOSE ARMREST PAD LEFT:

Part Number	Description	AADL - CA
I105098-99-0	4"x10" Corpus Armrest Pad Left - Leatherette	No Charge
I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge
I105107-99-0	4"x16" Corpus Armrest Pad Left - Leatherette	No Charge
I105998-99-0	4"x18" Corpus HD Armrest Pad Left - Leatherette	No Charge

OPTIONAL ITEMS:

Part Number	Description (Seating Benefits)	AADL - CA
I105463-99-0	Corpus Arm Trough Adapter 1" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option.</i>	357.00
I105462-99-0	Corpus Arm Trough Adapter 7/8" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option.</i> Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____	357.00
I107339-99-0	Armrest Pouch - Right	260.00
I107338-99-0	Armrest Pouch - Left	260.00

Armrest Options

I107840-99-0	UT Armrest Side Support - Right Outer <i>Includes 2”T x 6.5”W tapered pad.</i>	372.00
I107863-99-0	UT Armrest Side Support - Left Outer <i>Includes 2”T x 6.5”W tapered pad.</i>	372.00
I107864-99-0	UT Armrest Side Support - Right Inner <i>Includes 2”T x 6.5”W tapered pad. May need to increase width between the armrests to ensure armrests will flip-up.</i>	372.00
I107865-99-0	UT Armrest Side Support - Left Inner <i>Includes 2”T x 6.5”W tapered pad. May need to increase width between the armrests to ensure armrests will flip-up.</i>	372.00
I107839-99-0	UT Elbow Support - Right <i>Must also select PN I10881.</i>	200.00
I107838-99-0	UT Elbow Support - Left <i>Must also select PN I10882.</i>	200.00
I10881	Adjustable Elbow Support Hardware - Right	291.00
I10882	Adjustable Elbow Support Hardware - Left	291.00

Legrest Options

CORPUS LEGREST OPTIONS:

Part Number	Description	AADL - CA
	CHOOSE LEGREST TYPE:	
I107764-99-0	Corpus Manual Elevating Legrest <i>No charge when ordering power tilt and/or power recline.</i>	No Charge
I107762-99-0	Corpus Power Elevating Legrest - 170° (Client Cost) <i>Includes mechanical articulation.</i>	3,655.00
	CHOOSE FOOTPLATE STYLE:	
I105437-99-0	2pc Wide Corpus Footplates - 9”D x 7.5”W	No Charge
I105436-99-0	2pc Narrow Corpus Footplates - 9”D x 6”W	390.00
I106760-99-0	1pc Narrow Corpus Footpl - 9.75”D x 12.5”W <i>Footplate is 12.5”W in front, 8”W in back.</i>	390.00
	CALF SUPPORT OPTIONS:	
I105252-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5”W x 4.5”T Pad.</i>	No Charge

CORPUS LEGREST OPTIONS:

Part Number	Description	AADL - CA
I105258-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5”W x 6.5”T Pad.</i>	No Charge
I105250-99-0	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5”W x 6”T Pad.</i>	No Charge
I105256-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7”W x 7”T Curved Pad.</i>	404.00
I105254-99-0	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8”W x 8”T Pad.</i>	404.00

CORPUS VS LEGREST OPTIONS (FOR SELECT ACTIVE REACH PACKAGES): **(Client Cost)**

Part Number	Description	AADL - CA
	CORPUS VS LEGREST AND KNEE SUPPORT: <i>NOTE: The Corpus VS Legrest is only available with 20°, 30° and 45° of Active Reach.</i>	
I107759-99-0	Corpus VS Power Elevating Legrest - 180° <i>Includes power articulation.</i>	3,655.00
I107038-99-0	Knee Support - Corpus VS <i>Select when ordering Corpus VS Legrest. Must also select I11129.</i>	215.00
I11129	Adj Removable Knee Support Hardware (1 Each)	314.00

Legrest Options

	CHOOSE FOOTPLATE STYLE:	
I106969-99-0	1pc Corpus VS Footplate - 9"D x 16.5"W	No Charge
I106971-99-0	2pc Wide Corpus VS Footplates - 9"D x 7.5"W	390.00
I106967-99-0	2pc Narrow Corpus VS Footplates - 9"D x 6"W	390.00
	CALF SUPPORT OPTIONS:	
I107211-99-0	VS Calf Support Kit (J) w Mesh Cover & Adj Hardware 5"W x 6.5"T Pad. Compatible with footplate heights of 10" - 18".	No Charge
I106190-99-0	VS Calf Support Kit (A) w Mesh Cover & Adj Hardware 5.5"W x 6"T Pad. Compatible with footplate heights of 10" - 18".	No Charge
I106188-99-0	VS Calf Support Kit (H) w Mesh Cover & Adj Hardware 7"W x 7"T Curved Pad. Compatible with footplate heights of 11" - 18".	433.00
I106186-99-0	VS Calf Support Kit (D) w Mesh Cover & Adj Hardware 8"W x 8"T Pad. Compatible with footplate heights of 12" - 18".	433.00

FOOTPLATE HEIGHT POSITION:

Part Number	Description	AADL - CA
	SET FOOTPLATE HEIGHT AT: _____ <i>This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90°. DO NOT INCLUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments, please indicate desired footplate height. NOTE: Corpus footplate height range is 6"-17" and Corpus VS height range is 9"-17". Height may need to be reduced for shipping. Calf Pads are not available at heights less than 9".</i>	

OPTIONAL ITEMS:

Part Number	Description (Seating Benefits)	AADL - CA
I106530-99-0	Stealth Gel Pads - Narrow Footplates (Pair)	411.00
I106225-99-0	Stealth Gel Pads - Wide Footplates (Pair)	423.00

Accessories

I104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge
--------------	----------------------------------------	-----------

HEADREST OPTIONS:

Part Number	Description (Seating Benefits)	AADL - CA
I107470-99-0	Permobil 7"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	452.00
I107469-99-0	Permobil 10"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	452.00
I107468-99-0	Permobil 14"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	452.00
I107480-99-0	Adj/Removable Headrest Link Hardware - Standard	314.00
I107481-99-0	Adj/Removable Headrest Link Hardware - Extended <i>Includes an additional 3" link to extend headrest pad range.</i>	339.00
I11305	Swing Away Egg Switch Assembly <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	614.00
I11306	Swing Away Microlight Assembly <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	705.00
I11307	Swing Away Facial Pad Assembly 4"L x 2"T <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	500.00
I11308	Swing Away Facial Pad Assembly 6"L x 2"T <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	500.00
I11276	Provide Headrest Receiver for AGILITY Backrest <i>Select for AGILITY backs and third party backs ordered off chair.</i>	No Charge
I106331-99-0	Permobil 10"W Headrest, Leatherette <i>Not available with ROHO AGILITY Back Systems. Must also select PN I10552.</i>	640.00
I10552	Adj/Removable Headrest Hardware	137.00

Accessories

I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	258.00
--------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------

OPTIONAL ITEMS:

Part Number	Description (Seating Benefits)	AADL - CA
I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	384.00
I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	384.00
I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left: _____ Mount Amp Suppt on Right: _____</i>	137.00
I106104-99-0	UT Hip Suppt (B) w/ Mesh Cover - High (Pair) <i>4.5"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	313.00
I106103-99-0	UT Hip Suppt (A) w/ Mesh Cover - High (Pair) <i>6"W x 5.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	313.00
I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	508.00
I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	508.00
I102977-99-0	UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	508.00
I10550	Adj Removable - Thigh Support Hardware (Pair) (Seating Benefits)	137.00
I106762-99-0	Abductor Kit Corpus - Leatherette (Seating Benefits) <i>Must also select PN I10555.</i>	380.00
I10555	Adj Removable - Abductor Hardware (1 Each) (Seating Benefits)	137.00
I105260-99-0	Transfer Handles, Corpus Seat (Seating Benefits)	528.00
I105292-99-0	Push Handles, Corpus Seat	No Charge
I106210-99-0	RAM Self-Leveling Cup Holder (Client Cost)	91.00
I106806-99-0	RAM X-Grip Phone Holder (Client Cost)	91.00
I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount (Client Cost) <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount. Note: Customization may be required if ordered along with an Omni.</i>	504.00
I107222-99-0	USB Charger (Client Cost) <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Please select armrest for mounting, keeping in mind that charger may need to be relocated. Mount on left armrest: _____ Mount on right armrest: _____</i>	215.00
I105805-99-0	Medical Necessities Bag Hooks (Client Cost) <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall backrest. Custom required if ordered with a ROHO AGILITY backrest.</i>	99.00
I102059-99-0	Medical Necessities Bag (Client Cost) <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall backrest. Custom required if ordered with a ROHO AGILITY backrest.</i>	280.00
I105840-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit (Seating Benefits) <i>18.25"L x 4"W Pad Size. Item will be shipped as a part.</i>	193.00
I105839-99-0	Bodypoint Monoflex Large Chest Support Belt Kit (Seating Benefits) <i>21.5"L x 4.5"W Pad Size. Item will be shipped as a part.</i>	193.00
I105850-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit (Seating Benefits) <i>24.5"L x 5"W Pad Size. Item will be shipped as a part.</i>	193.00
I105262-99-0	Joystick Cutout Tray - 16" W (Seating Benefits) <i>Must also select PN I10556.</i>	469.00
I105263-99-0	Joystick Cutout Tray - 18" W (Seating Benefits) <i>Must also select PN I10556.</i>	469.00

Accessories

I10843	Small Angle Adjustable Tray (Seating Benefits) <i>Must also select PN I10556. Mount on Left: _____ Mount on Right: _____</i>	469.00
I10556	Adj Removable - Tray Hardware (1 Each)	137.00
I108230-99-0	Permobil Adjustable Tray - 13"W (Seating Benefits)	345.00
I108228-99-0	Permobil Adjustable Tray - 15"W (Seating Benefits)	345.00
I108227-99-0	Permobil Adjustable Tray - 17"W (Seating Benefits)	345.00
I108229-99-0	Permobil Adjustable Tray - 19"W (Seating Benefits)	345.00
I108232-99-0	Permobil Tray Hardware Right Side Adjustable (Seating Benefits)	250.00
I108233-99-0	Permobil Tray Hardware Left Side Adjustable (Seating Benefits)	250.00
I105224-99-0	Oxygen holder (Seating Benefits)	341.00
I105223-99-0	Crutch holder (Seating Benefits)	354.00

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	AADL - CA
I10751	Mini-Joystick for Omni - (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	3,122.00

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	AADL - CA
I103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,523.00
I103894-99-0	Compact Joystick LITE for Omni - (9-pin connector) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,743.00

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	AADL - CA
	NOTE: <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options.</i>	
I104690-99-0	Alt. JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	296.00
I107135-99-0	Alt JS/Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod. Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	599.00
I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	959.00
I10752	Haims Harness for Chin Control Mount, Large <i>For neck width up to 19cm.</i>	685.00
I10753	Haims Harness for Chin Control Mount, Small <i>For neck width up to 14cm.</i>	685.00

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	AADL - CA
I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode & On/Off. Switch operation force = 130 g Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	284.00

Alternative Drive Systems

I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 130 g</i>	517.00
--------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------

SIP & PUFF SYSTEMS:

Part Number	Description	AADL - CA
I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Includes breath tube kit.</i>	1,831.00
I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Includes breath tube kit.</i>	2,437.00
I106807-99-0	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Includes breath tube kit.</i>	1,011.00
I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>	296.00

PROPORTIONAL FOOT CONTROL:

Part Number	Description	AADL - CA
I10767	Proportional Foot Control <i>Not available with Active Reach Packages. Mount on Left: _____ Mount on Right: _____</i>	3,215.00
I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	296.00

Total Control Head Array

CONFIGURE HEAD ARRAY:

Part Number	Description	AADL - CA
I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module separately.</i>	5,650.00
I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module separately.</i>	No Charge
	CHOOSE OCCIPITAL PAD SIZE (REQUIRED)	
I105920-99-0	Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge
I105921-99-0	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge
	CHOOSE OCCIPITAL PAD MOUNT (REQUIRED):	
I10847	Fixed Mount for Occipital Pad - Head Array	No Charge
I105936-99-0	Adj. Swivel Mount for Head Array Occipital Pad <i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i>	506.00
	CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):	
I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links <i>Please select the length of the vertical stem below: Standard length (10"): _____ Short length (6"): _____</i>	314.00
I105941-99-0	Stealth TWB480 Backrest Mnt Hdwr with Flipback	519.00
I105925-99-0	Stealth HMO475 MultiAxis Backrest Mounting Hdwr	266.00
	SELECT LEFT ARM PACKAGE (OPTIONAL)	
I10976	Add Left Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	428.00
I10977	Add Left Lateral Arm w Egg Switch	285.00

Total Control Head Array

	Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____	
I10978	Add Left Lateral Arm w Pad (No Sensor or Switch)	174.00
	Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____	
	SELECT RIGHT ARM PACKAGE (OPTIONAL):	
I10979	Add Right Lateral Arm w Proximity Sensor	428.00
	Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____	
I10980	Add Right Lateral Arm w Egg Switch	285.00
	Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____	
I10981	Add Right Lateral Arm w Pad (No Sensor or Switch)	174.00
	Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____	

Alt Drive Systems Setup

CHOOSE OMNI (REQUIRED):

Part Number	Description	AADL - CA
I108500-99-0	Omni2 Control Module for R-net - Corpus <i>Two-section design features 3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology for control of up to 4 Bluetooth devices. Must also select mounting hardware below. Please specify which type of alternative drive input device you will be using. Input Device: _____</i>	1,995.00
I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	290.00

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	AADL - CA
I10080	Mount Seat Control Switchbox on Backrest	No Charge

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	AADL - CA
I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	318.00
I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	559.00
I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	806.00
I105907-99-0	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	416.00
I105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	277.00
I105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	169.00
I106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	146.00

Alt Drive Accessories

OPTIONAL ITEMS:

Part Number	Description	AADL - CA
I101769-99-0	Microlight Switch Qty: _____	183.00
I101767-99-0	Buddy Button Qty: _____	138.00
I106293-99-0	Egg Switch Black Qty: _____	152.00
I101771-99-0	Mini Cup Switch Qty: _____	135.00
I107500-99-0	Proximity Switch Qty: _____	640.00
I107495-99-0	Mini Joystick 5 Switch OC for Omni - (9-pin)	1,087.00
I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.	570.00
I10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ Mount on Left: _____ Mount on Right: _____ w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.	855.00

OPTIONAL ITEMS:

Part Number	Description	AADL - CA
I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	39.00
I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	43.00
I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	47.00
I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	52.00
I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.	43.00

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, LTD. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgment:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____
 Printed Name: _____
 Title: _____
 Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.