

| Customer and Order Information           |                  |                    |           |
|------------------------------------------|------------------|--------------------|-----------|
| Purchase Order *                         |                  | Tag                | v.0910    |
| Dealer/Provider Name *                   |                  |                    |           |
| Purchaser Name and Contact Information * |                  |                    |           |
| Shipping Address *                       |                  |                    |           |
| User Information (if available)          | Weight _____ lbs | Height _____ ft/in | *required |
| Order Notes/Instructions                 |                  |                    |           |

| Base Chair & Frame Options                                        |                                                                                                                                                                                                                                                                                            |                                                     |                                                                              |
|-------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------|------------------------------------------------------------------------------|
| <input checked="" type="checkbox"/> 1001000                       | Bentley 20° Manual Tilt Wheelchair                                                                                                                                                                                                                                                         |                                                     | W322 \$ 2,685                                                                |
| <b>Chair Weight Capacity</b>                                      |                                                                                                                                                                                                                                                                                            |                                                     |                                                                              |
| <input type="checkbox"/> 1002000                                  | 250 lbs Capacity Available for 14" to 22" Seat Width.                                                                                                                                                                                                                                      |                                                     |                                                                              |
| <input type="checkbox"/> 1002050                                  | 300 lbs Capacity Available for 14" to 22" Seat Width.                                                                                                                                                                                                                                      |                                                     | Prior Approval \$ 200                                                        |
| <input type="checkbox"/> 1005006                                  | 450 lbs Capacity Available for 20" to 32" Seat Width. Not compatible with Transit Tie Down System.                                                                                                                                                                                         |                                                     | Prior Approval \$ 500                                                        |
| <b>User Weight Range</b> Required for setup of tilt uplift force. |                                                                                                                                                                                                                                                                                            |                                                     |                                                                              |
| <input type="checkbox"/> 1002430                                  | Up to 75 lbs                                                                                                                                                                                                                                                                               | <input type="checkbox"/> 1002431 75 lbs to 100 lbs  | <input type="checkbox"/> 1002432 100 lbs to 115 lbs                          |
| <input type="checkbox"/> 1002433                                  | 115 lbs to 150 lbs                                                                                                                                                                                                                                                                         | <input type="checkbox"/> 1002434 150 lbs to 200 lbs | <input type="checkbox"/> 1002435 200 lbs to 250 lbs                          |
| <input type="checkbox"/> 1002436                                  | 250 lbs to 300 lbs                                                                                                                                                                                                                                                                         | <input type="checkbox"/> 1002437 300 lbs to 375 lbs | <input type="checkbox"/> 1002438 375 lbs to 450 lbs                          |
| <input type="checkbox"/> 1002300                                  | <b>Transit Tie Down System (TTDS)</b><br>TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.<br>Not compatible with Seat Widths greater than 22". Tested to max. Weight Capacity of 300 lbs. |                                                     |                                                                              |
| <b>PDG Premium Colours</b>                                        |                                                                                                                                                                                                                                                                                            |                                                     |                                                                              |
| <input type="checkbox"/> 1005001                                  | Black Shimmer                                                                                                                                                                                                                                                                              | <input type="checkbox"/> 1005015 Sky Blue           |                                                                              |
| <input type="checkbox"/>                                          | Other Standard PDG Premium Colour                                                                                                                                                                                                                                                          | Specify Colour _____                                | <a href="https://pdgmobility.com/colours">pdgmobility.com/colours</a> \$ 125 |

| Seat Width & Depth                                               |                                | Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW |                                                             |
|------------------------------------------------------------------|--------------------------------|------------------------------------------------------------------------|-------------------------------------------------------------|
| <b>Standard Seat Widths</b>                                      |                                |                                                                        |                                                             |
| <input type="checkbox"/> 1003014                                 | 14" Seat Width                 | <input type="checkbox"/> 1003015 15" Seat Width                        | <input type="checkbox"/> 1003016 16" Seat Width             |
| <input type="checkbox"/> 1003017                                 | 17" Seat Width                 | <input type="checkbox"/> 1003018 18" Seat Width                        | <input type="checkbox"/> 1003019 19" Seat Width             |
| <input type="checkbox"/> 1003020                                 | 20" Seat Width                 | <input type="checkbox"/> 1003021 21" Seat Width                        | <input type="checkbox"/> 1003022 22" Seat Width W111 \$ 250 |
| <b>Custom Seat Widths</b> Require 450 lbs Weight Capacity Frame. |                                |                                                                        |                                                             |
| <input type="checkbox"/> 1003023                                 | 23" Seat Width                 | <input type="checkbox"/> 1003024 24" Seat Width                        | <input type="checkbox"/> 1003025 25" Seat Width W111 \$ 500 |
| <input type="checkbox"/> 1003026                                 | 26" Seat Width                 | <input type="checkbox"/> 1003027 27" Seat Width                        | <input type="checkbox"/> 1003028 28" Seat Width W111 \$ 500 |
| <input type="checkbox"/> 1003029                                 | 29" Seat Width                 | <input type="checkbox"/> 1003030 30" Seat Width                        | <input type="checkbox"/> 1003031 31" Seat Width W111 \$ 750 |
| <input type="checkbox"/> 1003032                                 | 32" Seat Width                 |                                                                        | W111 \$ 750                                                 |
| <b>Seat Depths</b>                                               |                                |                                                                        |                                                             |
| <input type="checkbox"/> 1004016                                 | 16" Seat Depth                 | <input type="checkbox"/> 1004017 17" Seat Depth                        | <input type="checkbox"/> 1004018 18" Seat Depth             |
| <input type="checkbox"/> 1004019                                 | 19" Seat Depth                 |                                                                        |                                                             |
| <input type="checkbox"/> 1004020                                 | 20" Seat Depth                 | <input type="checkbox"/> 1004021 21" Seath Depth                       | <input type="checkbox"/> 1004022 22" Seat Depth W111 \$ 250 |
| <input checked="" type="checkbox"/> 1008001                      | <b>Solid Aluminum Seat Pan</b> |                                                                        |                                                             |

| Arm rests                                     |                        |                                                         |                                                |
|-----------------------------------------------|------------------------|---------------------------------------------------------|------------------------------------------------|
| <b>Single Post Height Adjustable Armrests</b> |                        |                                                         |                                                |
| <input type="checkbox"/> 1006111              | 8"-11" Range           | <input type="checkbox"/> 1006114 10"-14" Range          | <input type="checkbox"/> 1006116 12"-16" Range |
| <b>Dual Post Height Adjustable Armrests</b>   |                        |                                                         |                                                |
| <input type="checkbox"/> 1006211              | 8"-11" Range           | <input type="checkbox"/> 1006214 10"-14" Range          | <input type="checkbox"/> 1006216 12"-16" Range |
| <b>Regular Armpads</b>                        |                        |                                                         |                                                |
| <input type="checkbox"/> 1007100              | Desk Length - 2" x 10" | <input type="checkbox"/> 1007101 Full Length - 2" x 14" |                                                |
| <input type="checkbox"/> 1007310              | <b>Clothing Guards</b> |                                                         | \$ 125                                         |

Seat Height

➤ Don't forget

Select desired Seat-to-Floor Height    STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

| 6" Casters | Low |                          |                          |                          |                          |                          | High |                          |                          |                          |                          |                          |                          |
|------------|-----|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
|            | RW  | 13"                      | 14"                      | 15"                      | 16"                      | 17"                      | RW   | 15"                      | 16"                      | 17"                      | 18"                      | 19"                      | 20"                      |
|            | 12" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |                          |                          | 12"  | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |                          |                          |
|            | 16" |                          |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 16"  |                          |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
|            | 20" | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |                          |                          | 20"  | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |                          |                          |
|            | 22" |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |                          | 22"  |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |                          |
|            | 24" |                          |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 24"  |                          |                          | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

| 8" Casters | Low |     |     |                          |                          |                          | High |     |     |                          |                          |                          |                          |
|------------|-----|-----|-----|--------------------------|--------------------------|--------------------------|------|-----|-----|--------------------------|--------------------------|--------------------------|--------------------------|
|            | RW  | 13" | 14" | 15"                      | 16"                      | 17"                      | RW   | 15" | 16" | 17"                      | 18"                      | 19"                      | 20"                      |
|            | 16" |     |     | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 16"  |     |     | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
|            | 24" |     |     | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | 24"  |     |     | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Casters & Rear Wheels

Urethane Casters

- ☐ 1011061 6" x 1.25" Caster size
- ☐ 1011081 8" x 1.25" Caster size

Composite Wheels with Urethane Tires

- ☐ 1012012 12" Wheel size
- ☐ 1012016 16" Wheel size
- ☐ 1012020 20" Wheel size
- ☐ 1012022 22" Wheel size
- ☐ 1012024 24" Wheel size

☒ 1013000 Quick Release Rear Axles

Handrims    Handrims are only available with 20", 22", and 24" Rear Wheels.

- ☐ 1015010 Standard Aluminum
- ☐ 1015000 No Handrims
- ☐ 1015020 Plastic Coated

\$ 175

Wheel Locks

☒ 1016010 Push to Lock Wheel Locks

☐ 1016100 Wheel Lock Extensions

Anti-Tippers

☒ 1017010 Rear Anti-Tippers    Required Safety Option.

Back Post

☒ 1018010 1019010 Angle/Height Adjustable Stroller Handle    Button on Outside.

Back Post Height    Measured from the Seat Pan to the Handle.

- ☐ 1018425 Short (22")
- ☐ 1018431 Tall (25")

☒ 1019090 90° Seat to Back Angle

☒ 1019220 Fold Down Back

☒ 1019230 Spring Back Suspension

Tilt/Recline Options

☒ 1023011 Attendant Tilt only    Trigger on Right Back Post.

Front Rigging

Plug In Front Rigging

- ☐ 1025410 Easy-Out Fixed 70°    Effortless insertion and removal. No Hanger Bracket remains on frame.
- ☐ 1025110 Standard Plug In 70°    Legacy Design. Not compatible with Easy-Out.

\$ 95

Swing Away Front Rigging

- ☐ 1025140 Swing Away 70°
- ☐ 1025200 Swing Away Elevating Legrests    Includes Calf Pads.
- ☐ 1025300 1" Offset Hangers for Swing Away

W106 \$ 375

\$ 95

| Footplates & Options                 |         |                                                                                                          |             |        |
|--------------------------------------|---------|----------------------------------------------------------------------------------------------------------|-------------|--------|
| <b>Footplates</b>                    |         |                                                                                                          |             |        |
| <input type="checkbox"/>             | 1026010 | Composite                                                                                                |             |        |
| <input type="checkbox"/>             | 1026020 | Aluminum Angle Adjustable                                                                                | <i>W107</i> | \$ 165 |
| <input checked="" type="checkbox"/>  | 1027040 | <b>Length Adjustable Heel Loops</b>                                                                      |             |        |
| <b>Front Rigging Extension Tubes</b> |         |                                                                                                          |             |        |
| <input type="checkbox"/>             | 1028010 | Short    14"-18" for Easy-Out Front Rigging    13"-17" for Swing Away and Standard Plug In Front Rigging |             |        |
| <input type="checkbox"/>             | 1028020 | Long    18"-24" for Easy-Out Front Rigging    15"-19" for Swing Away and Standard Plug In Front Rigging  |             | \$ 75  |

Please send your completed order form by email to [sales.pdg@permobil.com](mailto:sales.pdg@permobil.com).