



Retail order form - January 2025

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.



The **conversion kit** converts the K110 mini stander to the K130 small. It includes the deck and complete body support of the K130. Request from customer service.

Conversion kit K116



*** Required**

*** Mobile Stander** frame includes body support with padding, seat pad, straps and fleece.

Specify color:
pink, blue or lime

CHOOSE

User dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Height	25–32	30–40	37–52	48–65
 Key user dimension: height Select the appropriate Mobile Stander by the user's overall height. Choose the model that allows for growth. Important: User's weight must not exceed the maximum working load.				
Item dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Overall length and width	25 x 25	25 x 25½	32 x 29	39 x 33½
Overall height	22½–28	23½–30	31–39	37½–49
Width without large wheels	19½	19½	24½	28¾
Wheel diameter	20	20	27	35
Width of body support	7–10	7–10	9–12	10–14
Height of body support	16½–22	22–29	29–38	36–48
Item weight (lb)	28	30	40	65
Max. working load (lb)	50	50	120	175

		K110 mini		K130 small		K150 medium		K170 large	
	CHOOSE	K101	\$3,609 ☐	K103	\$4,027 ☐	K105	\$4,894 ☐	K107	\$5,963 ☐
		Pink ☐		Pink ☐		Pink ☐		Pink ☐	
		Blue ☐		Blue ☐		Blue ☐		Blue ☐	
		Lime ☐		Lime ☐		Lime ☐		Lime ☐	
Wheels (pair) (include brakes)		K113	\$2,136 ☐	K113	\$2,136 ☐	K115	\$2,399 ☐	K117	\$2,894 ☐
Mounting bracket		K133	\$113 ☐	K133	\$113 ☐	K173	\$127 ☐	K173	\$127 ☐
Arm prompts (pair) (require mounting bracket)		K512	\$949 ☐	K512	\$949 ☐	K512	\$949 ☐	K512	\$949 ☐
Hand loops (require mounting bracket)	CHOOSE	Single	K578 \$178 ☐	K578	\$178 ☐	K578	\$178 ☐	K578	\$178 ☐
		Pair	K579 \$372 ☐	K579	\$372 ☐	K579	\$372 ☐	K579	\$372 ☐
Communication tray (requires mounting bracket)		K518	\$727 ☐	K518	\$727 ☐	K518	\$727 ☐	K518	\$727 ☐
Additional seat pad	CHOOSE	Small	N/A	N/A		K135	\$256 ☐	K135	\$256 ☐
		Large	N/A	N/A		K155	\$282 ☐	K155	\$282 ☐
Sandals with wedges (pair)	CHOOSE	Small	K114 \$198 ☐	K114	\$198 ☐	N/A		N/A	
		Medium	N/A	N/A		K154	\$210 ☐	N/A	
		Large	N/A	N/A		N/A		K174	\$224 ☐
Total		K110	\$ _____	K130	\$ _____	K150	\$ _____	K170	\$ _____