



Note: cover and foam insert required for use.

ROHO® Hybrid Elite Dual Compartment Cushion

This innovative cushion combines the skin and soft tissue protection and Shape Fitting Technology® of a ROHO DRY FLOATATION® cushion with stability from a contoured foam base.

The shape-fitting air insert provides critical pressure distribution under the pelvis for those with an obliquity. And the two independent air compartments let you set inflation levels to your shape and ideal position.



Combination of air and foam for aggressive positioning



Firm base for ease of transfers



Cushioned channels hold your legs in place



Two independent air compartments



Stability



Pressure injury prevention and treatment

Dimensions

Width: 13.75 in. to 23.75 in. / 35 cm to 60 cm

Height: 4.25 in. / 11 cm

Weight: Approx. 3 lbs. / 1.4 kg

Weight limit

500 lbs. / 227 kg

Limited warranty

Cushion: 24 months

Cover: 6 months



Our evidence shows the Hybrid Elite Dual Compartment is ideal for people who:

- Currently have, have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Require accommodation of pelvic asymmetry up to 2 in. (5 cm)
- Require a more stable environment for transfers
- Require pelvic and lower extremity positioning
- Require increased sitting stability