

permobil

ZRA

**The ultimate
ride quality.**

Experience the cutting edge technology behind our 1.25 inch Mono-Tube Titanium for unparalleled ride quality that will last a lifetime. Specifically selected for its durability and performance characteristics, our titanium puts you miles ahead of the rest.





TiFit. Personal Ergonomics.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. A TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do. TiFit...not standard.

Custom Configuration.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimise your riding experience. This includes choices of wheels, handrims, tyres, wheel locks, colour, tattoos, and more. TiLite offers more options for the individual.

Technical Specifications

Frame Features	Titanium Mono-Tube Rigid TiFit Limited Adjustability
Total Chair Width*	23"
Seat Width	10" – 20"
Seat Depth	10" – 20"
Front Seat Height	16" – 21"
Rear Seat Height	12.5" – 20"
Front Frame Angle	70°, 75°, 80°, 85°, or 90°
Seat Back Type	Folding Seat Back Adj. Height 8.5" – 21" Adj. Angle 80° – 101°

Camber	0°, 2°, 4°, 6° or 8°
Centre of Gravity Adjustment	6"
Front Seat Width (Seat Taper)	0", 1.5", 2.5"
Transit Options	Transport Bracket Transit Tie-Down
Transport Weight*	4.8 kg
Weight Limit	120.2 kg Standard 158.7 kg Heavy Duty
Warranty	Five Years on Frame One Year on TiLite Components

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.27 kg.


Colour Options

High Gloss Paint



TiLite Red Burgundy Neon Pink Super Wet Black

Matte Paint



Deep Blue Ice Titanium Grey Black Lava

Candy Paint




Ruby Red Azure Blue Electric Plum

Natural



Satin Titanium

Metallic Paint



Salsa Red Metallic Tangerine Metallic Acid Green Metallic Ocean Blue Metallic

Anodised Colours



Red Pepper Burnt Orange Sapphire Ultraviolet Fuchsia Black

Pearlescent Paint



Midnight Blue Pearl White Pearl