

Elevate with ActiveReach™



Benefits of seat elevation

- Facilitate vertical reach to promote independence with daily activities
- Protect upper extremities to reduce risk of repetitive strain injury related to overhead reaching
- Help improve safety and independence with transfers
- Allow eye contact with others reducing risk of cervical strain
- Provide psycho-social benefits of being on peer level and speaking eye-to-eye
- Promote enhanced safety with increased visibility in the community





Benefits of ActiveReach™

Your world isn't one dimensional. Your seat elevator shouldn't be either. ActiveReach™ technology is the ultimate advantage in accessing your environment. By combining Corpus seating functions with up to 14" of seat elevation and 45° of anterior tilt, Permobil's ActiveReach™ technology enables you to reach for more.

- Facilitate forward weight shift for transfers
- Assist with managing tone/abnormal reflexes
- Increase functional reach to promote greater independence with daily activities
- Promote improved postural control and ergonomics during forward reach activities
- Reduce risk of injury associated with lateral and overhead reaching
- Help aid digestion and speech production

Elevate & ActiveReach™ Packages

M3 Corpus

5°, 10°, 20°
+12" of elevation

M5 Corpus

5°, 10°, 20°
+14" of elevation

F3 Corpus

5°, 10°, 20°, 30°
+12" of elevation

F5 Corpus

5°, 10°, 20°, 45°
+14" of elevation

Please refer to your owner's manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions.