

Customer and Order Information			
Purchase Order *		Tag	v.0910
Dealer/Provider Name *			
Purchaser Name and Contact Information *			
Shipping Address *			
User Information (if available)		Weight _____ lbs	Height _____ ft/in
			*required
Order Notes/Instructions			

Base Chair & Frame Options			
☒	2601000	Fuze T50n 50° Manual Tilt Wheelchair	W654 \$ 2,685
☒	2601008	Narrow Frame Upgrade Reduces overall chair width to seat width plus only 4".	\$ 200
Frame Type			
☐	2602010	Fixed Width Frame 300 lbs Weight Capacity. Select a Fixed Seat Width below.	
☐	2602011	Fixed Width Frame 400 lbs Weight Capacity. Not compatible with Transit Tie Down System. Select a Fixed Seat Width below.	Prior Approval \$ 300
☐	2602025	Adjustable Width Frame 250 lbs Weight Capacity. Select and Adjustable Seat Width below.	\$ 200
☐	2602300	Transit Tie Down System (TTDS)	
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.			
PDG Premium Colours			
☐	2605001	Black Shimmer	
☐	2605015	Sky Blue	
☐		Other Standard PDG Premium Colour Specify Colour _____	pdgmobility.com/colours \$ 125

Seat Width & Depth				Overall chair width is +4.0"	
Fixed Seat Widths Order with Fixed Frame.					
☐	2603016	16" Seat Width	☐	2603017	17" Seat Width
☐	2603019	19" Seat Width	☐	2603018	18" Seat Width
☐	2603020	20" Seat Width	☐	2603021	21" Seat Width
☐	2603023	23" Seat Width	☐	2603022	22" Seat Width
					W111 \$ 250
					W111 \$ 500
Adjustable Seat Widths Order with Adjustable Frame.					
☐	2603216	16"-20" Range - Set at 16"	☐	2603217	16"-20" Range - Set at 17"
☐	2603219	16"-20" Range - Set at 19"	☐	2603218	16"-20" Range - Set at 18"
☐	2603220	16"-20" Range - Set at 20"			W111 \$ 250
Seat Depths					
☐	2604016	16" Seat Depth	☐	2604017	17" Seat Depth
☐	2604019	19" Seat Depth	☐	2604018	18" Seat Depth
☐	2604020	20" Seat Depth	☐	2604021	21" Seath Depth
☐			☐	2604022	22" Seat Depth
					W111 \$ 250
☒	2608001	Solid Aluminum Seat Pan			

Armrests	
Single Post Height Adjustable Armrests	
☐	2606114 10"-14" Range
☐	2606116 12"-16" Range
Regular Armpads	
☐	2607100 Desk Length - 2" x 10"
☐	2607101 Full Length - 2" x 14"

Seat Height



Select desired Seat-to-Floor Height STF Height is adjustable within the range of the chosen Axle Plate (Low/Standard/High).

Chairs with Swing-Away Front Rigging are set up with 3° of pre-tilt to avoid interference. This increases the actual STF Height at the front of the seat by up to 1" compared to the selected value. Easy-Out and Standard Plug In are not affected.

12" Rear Wheels					
Caster size	16"	17"	18"	19"	20"
5"	☐	☐			
6"	☐	☐	☐	☐	
8"	☐	☐	☐	☐	☐

Casters & Rear Wheels

Urethane Casters

☐ 2611051 5" x 1.25" Caster size ☐ 2611061 6" x 1.25" Caster size ☐ 2611081 8" x 1.25" Caster size

☒ 2614000 2614000 **12" Composite Wheels with Urethane Tires**

☒ 2613000 **Quick Release Rear Axles**

Wheel Locks

☒ 2616010 **Attendant Wheel Locks (Push to Lock)**

☐ 2616100 **Wheel Lock Extensions**

Anti-Tippers

☒ 2617010 **Rear Anti-Tippers** Required Safety Option.

Back Options

☒ 2618010 2619010 **Angle/Height Adjustable Stroller Handle** Button on Outside.

Back Post Height

 Measured from the Seat Pan to the Handle.

☐ 2618420 Short (17") ☐ 2618425 Tall (23")

☒ 2619090 **90° Seat to Back Angle**

☒ 2619220 **Fold Down Back**

Tilt/Recline Options

☒ 2623011 **Attendant Tilt only** Trigger on Right Back Post.

Front Rigging

Plug In Front Rigging

☐ 2625410 Easy-Out Fixed 70° Effortless insertion and removal. No Hanger Bracket remains on frame. \$ 95

☐ 2625110 Standard Plug In 70° Legacy Design. Not compatible with Easy-Out.

Swing Away Front Rigging

☐ 2625140 Swing Away 70°

☐ 2625200 Swing Away Elevating Legrests Includes Calf Pads. *W106* \$ 375

☐ 2625300 1" Offset Hangers for Swing Away \$ 95

Footplates & Options

Footplates

☐ 2626010 Composite

☐ 2626020 Aluminum Angle Adjustable *W107* \$ 165

☒ 2627040 **Length Adjustable Heel Loops**

Front Rigging Extension Tubes

☐ 2628010 Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging

☐ 2628020 Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging \$ 75

Please send your completed order form by email to sales.pdg@permobil.com.