

Permobil

M4 Nexa[®] LH



permobil

Closer to everyday life

With its compact base and a seat height of just 39 cm, the M4 Nexa LH (low height) helps users navigate everyday environments more easily, without compromising the power seating, performance, or durability expected from Permobil.

Made for everyday access

M4 Nexa LH is designed for users who need a lower seating position to get closer to daily life. The low seat-to-floor height of just 39 cm supports access to tables, vehicles, social environments, and transfers, while the compact base improves maneuverability in tight indoor spaces such as doorways, corridors, and bathrooms. Quieter driving and braking help create a smoother, more discreet mobility experience throughout the day.

Sit better, feel better, live better

M4 Nexa comes with a new biomechanical seating system designed to maintain a supported position with minimal displacement as power seat functions are used. The seating system allows for easy adjustment of seat width and depth, accommodating changes in body size, posture, function, or lifestyle.

Designed to reduce downtime

Service and reconfiguration are designed to be simple, with improved battery access, protected cabling, fewer components, and efficient seat adjustment using only two Allen keys.

Connectivity for more control

With Permobil +Connect, M4 Nexa supports remote insights, proactive service, and wireless programming. The MyPermobil app gives users practical information such as battery range, distance, and seating behaviour.





Specifications

Max. user weight	136 kg
Max. speed	10 km/h
Battery	52 Ah
Range (52 Ah battery)	33 km*
Weight incl. batteries 52 Ah/68 Ah	144 kg (52 Ah)
Obstacle climbing	70 mm
Ground clearance	74 mm
Base width	540 mm
Base length	850 mm
Seat depth	250–560 mm
Seat width	300–560 mm
Seat-to-floor height	390 mm
Back frame height	Low 300–410 mm** Standard 460–560 mm
R-net electronics Power Platform	120 A
Crash test ISO 7176-19:2022 tie-down	Yes—max 136 kg user weight
Tested according to EN12184:2022 CE	Yes

*Driving range (Actual driving range will vary based on driving, battery conditions and terrain)

**To fit back support sizes

Measurements are approximates and subject to change without notice.

Accent colours

BLACK	DARK BLUE	RED	GRAY	LIGHT BLUE	ORANGE
WHITE	DARK PURPLE	LIGHT PURPLE	YELLOW	PINK	GREEN

Actual product colors may vary.

Power functions



Anterior tilt
Active Reach 10°



Seat elevator
Active Height
300 mm



Power leg rest
85–175°



Posterior tilt
50°



Recline
150°

