



TiLite 2GX - Order Form



Dealer Name*:		Account #*:
City, State, Zip:		
Phone:	Email:	ATP:
Rep Name:	P.O. Number:	Quote Number:
Client Reference*:	Client Weight:	
Special considerations that need to be addressed here (e.g., diagnosis):		

*Required field

Frame Style

G2FS1	2GX Titanium Folding Manual Wheelchair Full adjustability, 265lb weight limit. Light weight and high strength. Not available with frame transit tie-down option.	K0005	\$ 4,692.00
-------	---	-------	-------------

Frame Options

XFS1	Swing Away Front Frame	Standard
XFS2	Fixed Front Frame	No charge
XFS3	Hemi Height Frame Available with a front seat height of 13.5"-20". Not with depth adjustable back/frame, 12"-13" seat depth or heavy duty frame.	No charge
XHDR1	Heavy Duty Frame 350lb weight limit. Includes tension adjustable by straps back upholstery. Available with 18"-22" seat widths and 16"-20" seat depths. Requires LiteSpeed billet aluminum front wheels and 0° camber. Not available with 90° front end or single sided forks.	No charge
XPWR1	Power Adaptable Frame Requires LiteSpeed billet aluminum front wheels and 0° camber. Not available with single-sided forks.	No charge
XRFR1	Reinforced Frame	No charge

Transit Options

XWTB	Transport Bracket Option Consists of four tie-down brackets. Unoccupied use only. Not available when swing away frame and heavy duty frame are both selected. Not with swing away tubular armrests or swing away flip back tubular armrests. Requires a minimum seat depth of 14", or 15" with hemi height frame.	E1022	\$ 379.00
------	--	-------	-----------

SmartDrive Power Assist System

XMAXSD1	Permobil SmartDrive MX2+ <i>Not available with rear seat heights less than 13.5". SmartDrive apps available on Apple App Store and Google Play Store for Apple Watch Series 5 or newer, and Samsung Galaxy Watch 4 or newer. Some exlusions may apply.</i>	E0986	\$ 7,529.00
XMAXSD18	Switch Control Buttons		No charge
XMAXSD8	SwitchControl Button with Mono Jack Option		No charge
XMAXSD19	Buddy Button <i>For use with SwitchControl button with mono jack option.</i>		No charge
XMAXSD4	Carrying Bag		\$ 84.00
XMAXSD10	Universal Padded Handle		\$ 22.00

Frame Finishes

XFF1	Satin Finish				No charge
	XSCBT1 Buffing Pads for Satin Frame				\$ 22.00
XFF	Powder Coated Frame <i>Select below.</i>				No charge
Select color:	XFF8 High Gloss Black	XFF10 Red Pepper	XFF29 Neon Pink	XFF14 Ocean Blue Metallic	
	XFF15 Tangerine Metallic	XFF30 Acid Green Metallic	XFF17 Midnight Blue Pearl	XFF18 White Pearl	
	XFF19 Black Lava	XFF27 Titanium Grey	XFF59 Pomegranate Red	XFF62 Blue Mist	
	XFF57 Indigo Purple	XFF40 Hyper Blue	XFF41 Flame Orange	XFF42 OD Green	
	XFF60 Silver Streak	XFF63 Lavendar Dusk	XFF58 Black Galaxy	XFF61 Seafoam Green	
	XFF64 Daffodil Yellow				
XFF3	Tattoos <i>Tattoos on front frame tubes. Finished in satin. Tattoos will be trimmed and adjusted to fit your frame.</i>				\$ 363.00
	XFF45 Need for Speed	XFF44 Flower Garden	XFF43 Mosaic		



Need for Speed



Flower Garden



Mosaic

Frame Measurements

Measurements shown below are our most popular selections. Not seeing what you want? Contact customer service or view our Creative Mobility Division (CMD) catalog [here](#) to view additional options. Please refer to Reference Diagrams at end of form for dimensional information.

Seat Width 'A' | Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

XSW1:	12"	13"	14"	15"	16"	17"	18"	19"	Standard
XSW2:	20"	21"	22"						E2201 \$ 536.00
21"-22" seat widths only available with heavy duty frame.									

Seat Depth 'B' | Measure from 0.5" in front of back posts to the front edge of seat sling. 12"-14" seat depths not with 25"-26" rear wheels.

XSD1:	12"	13"	14"	15"	16"	17"	18"	19"		Standard
XSD2:	20"								E2203	\$ 568.00
XSD3:	Custom Frame Depth. Specify:									No charge
<i>Frame depth longer than seat depth. Use 1" increments with a maximum of 2". A custom frame depth with a maximum of 20" is \$546.</i>										

Front Seat Height 'C' | Measure from floor to top of seat tube at front edge of seat sling. 13.5"-14" front seat height requires hemi height frame. 20.5"-21.5" front seat height not with hemi height frame.

XFSH1:	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"
	19.5"	20"	20.5"	21"								

Front Wheel	8" Front Wheel	19.5"-21"	18.5"-19.5"
	6" Front Wheel	17.5"-20"	16.5"-19"
	6" Front Wheel w/Glide	19"	17.5"-18"
	5" Front Wheel	16.5"-19.5"	16"-18.5"
	5" Front Wheel w/Glide	17.5"-18.5"	16.5"-18"
	4" Front Wheel	15.5"-19"	14.5"-18"
	4" Front Wheel w/Glide	17"-18"	15.5"-17"
	3" Front Wheel	14.5"-18.5"	13.5"-17.5"
	3" Front Wheel w/Glide	16.5"-17.5"	15"-17"

Std Front Seat Height

Hemi Front Seat Height

Front Seat Height Range

Rear Seat Height 'D' | Measure from floor to top of seat tube at back post.

XRSH1:	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"
	19"	19.5"	20"									

Rear Wheel	26"	15"-20"
	25"	14"-19"
	24"	14"-19"
	22"	13"-18"

Rear Seat Height Range

Back Supports

Comfort Back Supports | Backs over 20" wide or tall subject to increased pricing.

SBRS19	Comfort Unite 1 Back <i>Replaces Comfort Acta-Back. Inc 2-point quick release Compass hardware. Select shell height and cover.</i>				E2613	\$ 712
	12"	14"	16"	18"		No charge
	20"					\$ 100
	Comfort-Tek Cover		Diamond Spacer Cover			
	Privacy Shield					\$ 52
SBRS21	Comfort Unite 1 5" Wing Back <i>Replaces Comfort Acta-Back Deep w/5" Wing. Inc 2-point quick release Compass hardware. Select shell height and cover.</i>				E2620	\$ 845
	12"	14"	16"	18"		No charge
	20"					\$ 100
	Comfort-Tek Cover		Diamond Spacer Cover			
	Privacy Shield					\$ 52
SBRS20	Comfort Unite 2 Back <i>Replaces Comfort Acta-Back Deep w/3" Wing. Inc 2-point quick release Compass hardware. Select shell height and cover.</i>				E2615	\$ 820
	12"	14"	16"	18"		No charge
	20"					\$ 100
	Comfort-Tek Cover		Diamond Spacer Cover			
	Privacy Shield					\$ 52
SBRS22	Comfort Embrace Back <i>Replaces Comfort Acta-Embrace Back. Inc 2-point quick release Compass hardware. Select shell height.</i>				E2613	\$ 700
	10"	12"	14"	16"		
	Comfort-Tek Cover		Diamond Spacer Cover			
SBRS23	Comfort Relief Back <i>Replaces Comfort Acta-Relief Back. Inc 2-point quick release Compass hardware. Select shell height.</i>				E2615	\$ 886
	10"	12"	14"	16"		
	Comfort-Tek Cover		Diamond Spacer Cover			
	Privacy Shield					\$ 52

Backrest

Back Type

XSB13	Fixed Angle, Height Adjustable Back	Standard
XLB1	Folding Stabilizer Bar <i>Required with transit tie-down option unless depth adjustable backrest is selected. Includes tension adjustable by strap back upholstery.</i>	\$ 423.00
XSB14	Fixed Angle, Height Adjustable Back w/8° Bend <i>Not available with heavy duty frame.</i>	No charge
XSB18	Folding, Angle and Height Adjustable Back <i>Not available with heavy duty frame.</i>	\$ 250.00
XSB16	Depth Adjustable Back Frame <i>Not available with hemi height frame, heavy duty frame, 19"-20" seat depths, swing away tubular or swing away flip back armrests. Seat depths 12"-17" allow 3" of rearward adjustment. Seat depth of 18" allows 2" of rearward adjustment. Requires tension adjustable seat upholstery.</i>	\$ 244.00

Push Handles

XSB6	Integrated Push Handles	No charge
XSB7	Integrated Fold-Down Push Handles <i>Requires tension adjustable by straps back upholstery. Not available with heavy duty frame.</i>	\$ 125.00
XSB9	Bolt-On Height Adjustable Push Handles <i>Requires tension adjustable by straps back upholstery. Not available with heavy duty frame, folding angle and height adjustable back, swing away tubular or swing away flip back tubular armrests.</i>	\$ 363.00

Fixed Angle Adjustable Back Height 'E' | 11"-15" back height not with heavy duty frame, swing away tubular armrests 12"-15", bolt-on height adjustable push handles or folding stabilizer bar. Set at 1" increments.

Push Handles	None	11"-15"	16"-20"	Set at:
	Integrated	11"-15"	16"-20"	Set at:
	Fold-Down	11"-15"	16"-20"	Set at:
	Bolt-On Height Adj	N/A	16"-20"	Set at:
		Short	Tall	Specify
Fixed Angle Adj Back Height				

Folding Angle Adjustable Back Height 'E' | Not with folding stabilizer bar or heavy duty frame. Set at 1" increments.

Push Handles	None	11"-13"	13.5"-16.5"	16.5"-19.5"	Set at:
	Integrated	12.5"-14.5"	13.5"-16.5"	17"-19.5"	Set at:
	Fold-Down	15"-17"	16"-19"	18.5"-20.5"	Set at:
		Short	Medium	Tall	Specify
Folding Angle Adj Back Height					

Front Rigging

Front Angle 'H' | Front angle is measured from inside of front frame tube to top of bottom frame tube.

XFA3	Swing Away Front End	85° not with 6" front wheels. 90° requires aluminum footplate and not with heavy duty frame or 5"-6" front wheels.	No charge
	70° 80° 85° 90°	XFA0 Omit Hangers	
XFA6	Fixed Front End	85° not with 6" front wheels.	Standard
	70° 80° 85°		

Select one seat to footrest length and one footrest type.

Standard Footrests | For front angles 70°-85°. Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than front seat height. 70° fixed front has a minimum seat to footrest of 14.5".

Seat to Footrest 'E':	14" 20"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	
XFTR5	Flip-Up Footrest												No charge
XFTR9	Depth Adjustable Flip-Up Footrest												K0045 \$ 120.00
XFTR10	Depth and Angle Adjustable Flip-Up Footrest Available for 15"-20" seat to footrest only.												K0040, K0045 \$ 293.00
XFTR11	Angle Adjustable Flip-Up Footrest Available for 15"-20" seat to footrest only.												K0040, K0045 \$ 185.00
XHLP1	Heel Loops Not with seat width ≤ 13" or one piece footrests.												E0951 \$ 55.00

High Mount Footrests | Not with luggage carriers or impact guards.

Seat to Footrest 'E':	7" 13"	7.5" 13.5"	8" 14"	8.5" 14.5"	9"	9.5"	10"	10.5"	11"	11.5"	12"	12.5"	
XFTR20	High Mount Flip-Up												K0037 \$ 202.00
XFTR15	Angle Adjustable High Mount Flip-Up Footrest												K0040, K0037 \$ 352.00
XFTR16	Depth and Angle Adjustable High Mount Flip-Up Footrest												K0040, K0037 \$ 406.00
XHLP1	Heel Loops Not with seat width ≤ 13" or one piece footrests.												E0951 \$ 55.00

90° Footrests

Seat to Footrest 'E':	7"	8"	9"	10"	11"	12"	13"	
XFTR5	Flip-Up Footrest							No charge
XFTR9	Depth Adjustable Flip-Up Footrest							K0045 \$ 120.00
XFTR10	Depth and Angle Adjustable Flip-Up Footrest Available for 15"-20" seat to footrest only.							K0040, K0045 \$ 293.00
XHLP1	Heel Loops Not with seat width ≤ 13" or one piece footrests.							E0951 \$ 55.00

Elevating Legrests

Seat to Footrest 'E':	14"	15"	16"	17"	18"	19"	20"	21"		
XFA5	Depth Adjustable Flip-Up Footrest, Single								E0990	\$ 185.00
	Left	Right								
XFTR17	Depth Adjustable Flip-Up Footrest, Pair									\$ 363.00
XHLP1	Heel Loops <i>Not with seat width ≤ 13" or one piece footrests.</i>								E0951	\$ 55.00

Camber Options

Camber 'K'

XCBR8:	0°	2°	4°	6°	E2601	Standard
--------	----	----	----	----	-------	----------

Front Wheels

See front seat height ranges chart for front seat height 'C' limitations.

XFW24	3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	\$ 206.00
	Silver Hub Black Hub	
XFW14	4"x0.75" Plastic Wheel w/Poly Tire <i>Not available with power adaptable frame or heavy duty frame.</i>	No charge
XFW3	4"x0.75" Light-Up Micro Caster <i>Not available with power adaptable frame or heavy duty frame.</i>	\$ 88.00
XFW26	4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with power assist, power add-on system or heavy duty frame.</i>	\$ 112.00
XFW27	4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	\$ 217.00
	Silver Hub Black Hub	
XFW6	5"x1" Plastic Wheel w/Poly Tire <i>Not available with power adaptable frame, heavy duty frame or with 90° front angle.</i>	No charge
XFW28	5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with power adaptable frame, heavy duty frame or with 90° front angle.</i>	\$ 103.00
XFW30	5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not available with 90° front angle.</i>	\$ 233.00
	Silver Hub Black Hub	
XFW9	6"x1" Plastic Wheel w/Poly Tire <i>Not available with power adaptable frame, heavy duty frame or 85°-90° front angle.</i>	No charge
XFW32	6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with power adaptable frame, heavy duty frame or 85°-90° front angle.</i>	\$ 119.00
XFW33	6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not available with 85°-90° front angle.</i>	\$ 244.00
	Silver Hub Black Hub	
XFW16	8"x1" Plastic Wheel w/Poly Tire <i>Not available with power assist frame, power add-on frame, heavy duty frame or 12"-14" seat depths. Only with 70° front angle.</i>	No charge

Front Forks

XFK1	TiLite Standard Fork	No charge
XGL1	Glide Suspension Fork <i>Black anodize only. Not available with 8" front wheels or heavy duty frame. Select rider weight.</i>	E1015 \$ 526.00
	125lbs 125-180lbs 180-265lbs	

Rear Wheels

Rear Wheels

See rear seat height ranges chart for rear seat height 'D' limitations.

XOAD1	One Arm Drive System <i>Available with 24" mag wheels and standard aluminum handrim only. Rear seat height, camber or axle restrictions may apply.</i>	E0958 \$ 1,298.00
	Install on Left Install on Right	
XRW14	TiLite SHADOW Rear Wheel <i>Select size. 25" and 26" wheels require 15" minimum seat depth.</i>	Standard
	22" (501) 24" (540) 25" (559) 26" (590)	
XRW7	Spinergy SPOX Rear Wheel <i>Select size and spoke color. 25" and 26" wheels require 15" minimum seat depth.</i>	\$ 1,195
	22" (501) 24" (540) 25" (559) 26" (590)	
	Black Red Blue White Green	
XRW8	Spinergy LX Rear Wheel <i>Select size and spoke color. 25" and 26" wheels require 15" minimum seat depth. Only available with heavy duty frame when 0° or 2° camber is selected.</i>	\$ 1,295
	22" (501) 24" (540) 25" (559) 26" (590)	
	Black Red Blue White Orange Pink Purple Green	
ZRW9	Spinergy Design Lab LX Rear Wheel <i>Black spokes only. Select size and anodizing color. 25" and 26" wheels require 15" minimum seat depth.</i>	\$ 1,995
	22" (501) 24" (540) 25" (559) 26" (590)	
	Gold Green Blue Red Purple Orange	
XRW16	Spinergy Carbon Blade Rear Wheel <i>Select size. Not available with heavy duty frame. 25" wheels require 15" minimum seat depth.</i>	\$ 1,795
	24" (540) 25" (559)	
XRW17	X-Core Mag Rear Wheel <i>Select size. Not available with heavy duty frame, short tabs, quad release axles, spoke guards or omit handrims. Requires aluminum silver anodized handrims.</i>	\$ 225
	22" (501) 24" (540)	

Axles

XAXL1	Stainless Quick Release Axle		Standard
XAXL2	Titanium Quick Release Axle		\$ 277.00
XAXL3	Quad Release Axle <i>Not available with Mag or 26" (590) Spinergy Spox or LX rear wheels.</i>		\$ 169.00

Pneumatic Rear Tires

XRTR1	Treaded Tire 22"(501), 24"(540), 26"(590) <i>Not available with 25" rear wheels.</i>		Standard
XRTR2	Airless Inserts <i>Not available with Spinergy Carbon Blade rear wheels.</i>	E2213	\$ 98.00
XRTR31	MBL SpeedLite Tire 22"(501), 24"(540), 25"(559), 26"(590)		No charge
XRTR30	Primo Sentinel Tire 24"(540), 25"(559), 26"(590)		\$ 136.00
XRTR8	Schwalbe Marathon Plus Evolution Tire 24"(540), 25"(559), 26"(590)		\$ 169
XRTR32	Primo C-764 Knobby Tire 24"(540)		\$ 104.00

Solid Rear Tires

XRTR22	Primo Orion Tire 22"(501), 24"(540), 25"(559) <i>Not available with Spinergy Carbon Blade rear wheels.</i>		\$ 120.00
XRTR23	Primo Active Tire 24"(540), 25"(559) <i>Not available with Spinergy Carbon Blade rear wheels.</i>	E2213	\$ 147.00

Handrims

XRIM1	Aluminum Handrim, Silver Anodized		Standard
XRIM21	Aluminum Handrim, Black Anodized		No charge
XRIM12	Titanium Handrim, Satin Finish <i>Not with 22"(501) rear wheels or 24"(540) knobby tires.</i>		\$ 509.00
XRIM20	Q-Grip Handrim <i>Exclusive high-friction, anti-slip coating.</i>		\$ 358.00
XRIM5	Plastic Coated Handrim, Black		\$ 174.00
XRIM10	The Surge Handrim <i>Large profile. All-In-One oval with Gription strip.</i>		\$ 433.00
XRIM11	The Surge LT Handrim <i>Small profile. All-In-One oval with Gription strip.</i>		\$ 466.00
XRIM14	The Natural-Fit Handrim, Standard Grip <i>Original, large profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>		\$ 428.00
XRIM15	The Natural-Fit Handrim, Super Grip <i>Original, large profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>		\$ 466.00
XRIM16	The Natural-Fit Handrim, Without Thumb Grip <i>Original, large profile.</i>		\$ 390.00
XRIM17	The Natural-Fit LT Handrim, Standard Grip <i>Reduced weight, small profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>		\$ 440.00
XRIM18	The Natural-Fit LT Handrim, Super Grip <i>Reduced weight, small profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>		\$ 514.00
XRIM19	The Natural-Fit LT Handrim, Without Thumb Grip <i>Reduced weight, small profile.</i>		\$ 423.00
XRIM22	Omit Handrims		No charge

Handrim Tab Length

XRIM7	Long Tabs		Standard
XRIM8	Short Tabs <i>Not available with Mag rear wheels, The Natural-Fit handrims with thumb grip or The Natural-Fit handrims on Shadow wheels.</i>		\$ 28.00

Wheel Locks

XWLK7	Compact Composite Wheel Lock <i>Depending on configuration there may be compatibility issues with 12"-13" seat depths. Uni-Lock will be substituted at no charge.</i>		No charge
XWLK8	Compact Composite Pull to Lock <i>Depending on configuration there may be compatibility issues with 12"-13" seat depths. Uni-Lock will be substituted at no charge.</i>		No charge
XWLK10	Extension Handles for Compact Composite Push/Pull to Lock	E0961	\$ 93.00
XWLK11	MBL Composite Scissor Lock <i>Not available with swing away frame.</i>		\$ 158.00
XWLK12	MBL Performance Aluminum Scissor Lock <i>Not available with swing away frame.</i>		\$ 175
XWLK3	Uni-Lock		\$ 179.00
	Push to Lock Pull to Lock		
XWLK4	Extension Handles for Uni-Lock	E0961	\$ 93.00
XWLK15	Grade Aids, Push to Lock <i>Not available with solid rear tires.</i>	E0974	\$ 190.00

Back Upholstery

Back Upholstery

XBUP1	Padded Nylon Back Upholstery <i>Black only. Not available with heavy duty frame, folding back, folding stabilizer bar, depth adjustable frame, fold-down push handles or ultrasuede backrest stripe.</i>		Standard
XBUP3	Tension Adjustable Back Upholstery By Straps, Padded Nylon <i>Included with heavy duty frame. Required for aluminum fendered side guards and fold-down push handles.</i>	E2611	\$ 369.00
	XUPC3 Ultrasuede Backrest Stripe, Black		\$ 115.00
XBUP5	Tension Adjustable Back Upholstery By Straps, AIR Breathable <i>Included with heavy duty frame. Required for aluminum fendered side guards and fold-down push handles. Not with ultrasuede backrest stripe.</i>	E2611	\$ 369.00
XBUP0	Omit Back Upholstery		No charge

Seat Upholstery

XSUP4	Tension Adjustable Seat Upholstery, Black		Standard
XSUP0	Omit Seat Upholstery		No charge

Seat Cushions

XCSH2	Comfort Curve Cushion <i>Stretch-Air cover.</i>	E2601	\$ 122.00
XCSH23	Comfort Acta-Embrace Comfort-Tek Zero Elevation Cushion	E2607	\$ 428.00
XCSH15	ROHO Hybrid Elite SR Single Compartment Cushion <i>Fits seat widths 14"-19"</i>	E2622	\$ 552.00
XCSH14	ROHO Quadtro Select Cushion <i>Fits seat widths 14"-19".</i>	E2624	\$ 534.00
	Low Profile Mid Profile High Profile		
XCSH16	ROHO Hybrid Select Cushion <i>Fits seat widths 14"-18", and 20"</i>	E2624	\$ 571.00
XCSH17	Vicar Active O2 Low <i>Fits even seat widths 14"-20". Non-stock upcharges may apply.</i>	E2624	\$ 542.00

Side Guards

XSDG2	Aluminum Side Guards, Rigid Removable	E0956, E1034	\$ 217.00
	Small Standard		
XSDG3	Aluminum Side Guards, Quad Removable w/Hole	E0956, E1034	\$ 233.00
	Small Standard		
XSDG7	Carbon Side Guards, Rigid Removable <i>2x2 Twill Weave</i>	E0956, E1034	\$ 336.00
	Small Standard		
XSDG8	Carbon Side Guards, Quad Removable w/Hole <i>2x2 Twill Weave</i>	E0956, E1034	\$ 347.00
	Small Standard		
XSDG6	Aluminum Side Guards, Fendered Adjustable <i>Requires tension adjustable by straps back upholstery. Not available with angle adjustable fold-down back.</i>	E0956, E1034	\$ 363.00
CMD007	Custom Offset Profiled Side Guards (Removable or Fold-Down) <i>Please specify desired offset to nearest 0.5" and side guards type in the Notes section of this order form. Side guard material must be selected. Some limitations may apply.</i>		CMD

Armrests

XARM1	Swing Away Tubular Armrests <i>Select range. Not available with heavy duty frame, folding stabilizer bar with short or medium back height ranges, or depth adjustable back.</i>	E0973	\$ 244.00
	9"-12" 12"-15"		
XARM4	Swing Away Flip-Back Tubular Armrests <i>Not available with heavy duty frame, folding stabilizer bar or depth adjustable back.</i>	E0973	\$ 342.00
XARM2	Removable Height Adjustable Desk Arm with Rigid Side Guard <i>Select range and pad size. Requires 14" minimum seat depth.</i>	E0973	\$ 363.00
	9"-12" 12"-15"		
	Desk Length (10") Full Length (14")		
XARM3	Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop <i>Select range and pad size. Requires 16" minimum seat depth.</i>	E0973	\$ 406.00
	9"-12" 12"-15"		
	Desk Length (10") Full Length (14")		
XARM5	Armrests Covered in Black Ultrasuede <i>Not available with desk arms.</i>		\$ 120.00

Positioning Belts

XBLT1	Velcro Adjustable Belt, Black Nylon Webbing	E0978	\$ 61.00
XBLT2	Bodypoint 2" Belt, Push-Button Buckle (Auto Style)	E0978	\$ 103.00
XBLT4	Bodypoint 1" Padded Belt, Push-Button Buckle (Auto Style) <i>Not available with 19"-22" rear seat widths.</i>	E0978	\$ 196.00
XBLT5	Bodypoint 1.5" Padded Belt, Push-Button Buckle (Auto Style) <i>Not available with 19"-22" rear seat widths.</i>	E0978	\$ 185.00
XBLT10	Bodypoint Optional Band Clamp Mounting Hardware	K0108	\$ 22.00
XBLT16	Bodypoint EvoFlex Pelvic Stabilizer, Push Button	E0978	\$ 201.00

Calf Straps

XCLF2	Velcro Adjustable Small (9"-12" Footrest Widths) Medium (11"-15" Footrest Widths) Large (14"-17" Footrest Widths) Xtra Large (16"-19" Footrest Widths)	K0038	\$ 44.00
XCLF3	Bodypoint® Padded Velcro Adjustable Medium (12"-15" Footrest Widths) Large (16"-20" Footrest Widths)	K0038	\$ 109.00
XCLF1	Buckle Adjustable, Black Nylon Webbing	K0038	\$ 34.00

Residual Limb Supports

XRSL1	Residual Limb Support, 8"W x 10"L, Single Left Right	E1020	\$ 428.00
XRSL2	Residual Limb Support, 8"W x 10"L, Pair	E1020	\$ 855.00

Anti-Tips

XTIP1	Rear Anti-Tips, Aluminum Black	E0971	\$ 174.00
XTIP2	Rear Anti-Tips, User Friendly Flip-Up	E0971	\$ 363.00
XTIP9	Tip Assist	E0971	\$ 115.00

Accessories

XGRD1	Clear Spoke Guards <i>Not available with mag rear wheels.</i>	K0065	\$ 147.00
XBCK3	Permobil Medical Necessities Backpack		\$ 103.00
XPCH3	Permobil Medical Necessities Seat Pouch		\$ 66.00
XPCH4	Bodypoint Mobility Bag		\$ 109.00
XUEST1	Upper Extremity Support Tray <i>Requires armrests.</i> 12"-14" Seat Widths 13"-15" Seat Widths	E0950	\$ 315.00
XLUG1	Luggage Carriers <i>Not available with high mount swing away footrests or 1-piece flip-up footrest.</i>		\$ 158.00
XMPCT2	Black Neoprene Impact Guards with Ultrasuede Stripe <i>Not available with high mount footrests, 90° hangers or elevating legrests.</i>		\$ 115.00
XCH1	Crutch Holder	E2207	\$ 244.00
XAFP1	4" Military Service Patch Army Navy Air Force Marines Coast Guard		\$ 115.00

General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes

TiFit Reference Diagrams

