



Retail order form - January 2023

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.



User dimensions (inches)		E930 small	E940 med.	E950 large
Height		25 – 48	44 – 59	57 – 72
	Key user dimension: height			
	Select the appropriate Prone Stander by the user's overall height. Choose the model that allows for growth.			
		Important: User's weight must not exceed the maximum working load.		
Item dimensions (inches)		E930 small	E940 med.	E950 large
Board length	23½ (16 without kneeboard)		30	39
Board width (top/bottom)	12½ / 12½		11 / 16	14 / 19
Top of board to footboard	20 – 35		33 – 43½	44 – 57
Height when horizontal	29½		29½	31½
Height when vertical	38½		45½	61
Distance between laterals	8 – 16		8 – 16	12 – 22
Base width	23		26	30
Base length	35		39½	44
Tray inside width (widest point)	22½		22½	29½
Tray inside length	17½		17½	17½
Max. working load (lb)		100	150	200

*** Required**

* **Prone Stander** includes a pair of trunk laterals, hip strap, padding and footboard.

Specify color: red or blue

	E930 small		E940 medium		E950 large	
CHOOSE	E931	\$3,535 <input type="checkbox"/>	E941	\$4,084 <input type="checkbox"/>	E951	\$4,988 <input type="checkbox"/>
	Red <input type="checkbox"/>		Red <input type="checkbox"/>		Red <input type="checkbox"/>	
	Blue <input type="checkbox"/>		Blue <input type="checkbox"/>		Blue <input type="checkbox"/>	
	N/A		E614	\$352 <input type="checkbox"/>	E614	\$352 <input type="checkbox"/>
	E942	\$331 <input type="checkbox"/>	E942	\$331 <input type="checkbox"/>	E942	\$331 <input type="checkbox"/>
	N/A		E946	\$110 <input type="checkbox"/>	E956	\$121 <input type="checkbox"/>
	E947	\$1,007 <input type="checkbox"/>	E947	\$1,007 <input type="checkbox"/>	E957	\$1,140 <input type="checkbox"/>
	E933	\$142 <input type="checkbox"/>	E933	\$142 <input type="checkbox"/>	E933	\$142 <input type="checkbox"/>
	E934	\$121 <input type="checkbox"/>	E934	\$121 <input type="checkbox"/>	E934	\$121 <input type="checkbox"/>
	E655	\$472 <input type="checkbox"/>	E655	\$472 <input type="checkbox"/>	E655	\$472 <input type="checkbox"/>
	N/A		E945	\$404 <input type="checkbox"/>	E945	\$404 <input type="checkbox"/>
	N/A		E605	\$507 <input type="checkbox"/>	E605	\$507 <input type="checkbox"/>
	E422	\$151 <input type="checkbox"/>	E422	\$151 <input type="checkbox"/>	N/A	
	E433	\$164 <input type="checkbox"/>	E433	\$164 <input type="checkbox"/>	E433	\$164 <input type="checkbox"/>
	N/A		N/A		E432	\$174 <input type="checkbox"/>
	E424	\$110 <input type="checkbox"/>	E424	\$110 <input type="checkbox"/>	N/A	
	E435	\$117 <input type="checkbox"/>	E435	\$117 <input type="checkbox"/>	E435	\$117 <input type="checkbox"/>
	N/A		N/A		E434	\$125 <input type="checkbox"/>
	K821	\$258 <input type="checkbox"/>	K821	\$258 <input type="checkbox"/>	K821	\$258 <input type="checkbox"/>
	E930	\$ _____	E940	\$ _____	E950	\$ _____