



Frame colors



Component colors (TiLite Z only)



Specifications

	TiLite Z	TiLite X
Transport weight*	11.3 lb.	12.1 lb.
Frame material	Aluminum	Aluminum
Seat width	12" - 20"	14" - 22"
Seat depth	12" - 20"	14" - 20"
Front seat height	16" - 21"	13.5" - 20.5"
Rear seat height	13" - 20"	13" - 20"
Front frame angle	75°, 80°, 85°, or 90°	Swing Away 70°, 80°, 85°
Seat back height	9" - 20.5"	12" - 20.5"
Camber	0°, 2°, 4°, or 6°	0°, 2°, or 4°
Center of gravity	5"	4"
Backrest angle	80° - 101°	80° - 101°
Weight limit	275 lb.	Standard - 300 lb. Heavy Duty - 350 lb.
Warranty	Limited	Limited



Configure now at the Permobility Store

Transport weight is calculated based on a 16" x 16" seat dimension with 0 degrees of camber, 4" Poly casters, Comfort Company Compass 4 hardware, and standard seat upholstery. Wheel locks, and all items that can be removed without tools, such as anti-tippers, swing-away footrests, armrests, rear wheels, side guards, and seating products are excluded.

TiLite® X

Aluminum folding manual wheelchair

TiLite® Z

Aluminum rigid manual wheelchair



It's your TiLite

Introducing TiLite X and TiLite Z: truly personalized mobility with over 1 billion possible configurations and streamlined assembly, without compromising our signature quality and best-in-class performance.

With reduced transport weights and ultimate adjustability, these chairs are sure to adapt to your changing seating and mobility needs.



Performance comes standard



- New flip-back or removable T-style armrest
- Reinforced Aluminum Push/Pull wheel locks standard
- Re-designed frame and cross brace for best-in-class transport weight and performance
- GearLock caster adjustment
- Composite SoftRoll casters standard

Performance comes standard



- Made just for you
- Tension adjustable by strap seat sling
- Easier than ever to adjust
- Precision Lock caster adjustment
- Composite SoftRoll casters standard

Tailored, unmatched, dependable

- Result of 25 years of experience and 3 years of development through collaboration with industry experts.
- Individually built, just for you.
- Stronger, lighter, and faster than ever.

