



PDG Stellar GL - AADL W321



Dealer Name*: _____ Account #*: _____

City, Province, Post Code: _____

Phone: _____ Email: _____ ATP: _____

Rep Name: _____ P.O. Number: _____ Quote Number: _____

Client Reference*: _____ Client Weight: _____

Special considerations that need to be addressed here (e.g., diagnosis): _____ *Required field

October 2024, Revised January 2026. Prices are in Canadian Dollars (CAD)

Wheelchair & Frame

Stellar GL 45° Manual Tilt Wheelchair Also available: the lightweight Stellar GLT with aluminium frame. See separate order form. W321 \$ 2,685.00

Frame Type

- 250lbs Weight Capacity** Available for 14" - 22" seat widths. No charge
- 300lbs Weight Capacity** Available for 14" - 22" seat widths. Prior Approval \$ 200.00
- 450lbs Weight Capacity** Available for 20" - 32" seat widths. Prior Approval \$ 500.00

User Weight | Required for setup of tilt uplift force.

User Weight Range:	Up to 75lbs	75lbs - 100lbs	100lbs - 115lbs	115lbs - 150lbs	150lbs - 200lbs	200lbs - 250lbs
	250lbs - 300lbs	300lbs - 375lbs	375lbs - 450lbs			

Frame Options

Transit Tie Down System (TTDS) Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details. Not compatible with seat widths > 22". Tested to maximum weight capacity of 300lbs. No charge

Seat Measurements

Seat Width | Overall chair width is +9.5" with 12"/16" rear wheels, +10.5" with 20"/22"/24" rear wheels.

Select Width:	14"	15"	16"	17"	18"	19"		No charge
	20"	21"	22"				W111	\$ 250.00
	23"	24"	25"	26"	27"	28"	W111	\$ 500.00
	29"	30"	31"	32"			W111	\$ 750.00

Seat widths 23"-32" require 450lb Frame. Extended delivery lead times may apply.

Seat Depth

Select Depth:	16"	17"	18"	19"		No charge
	20"	21"	22"		W111	\$ 250.00

16" - 20" seat depths are adjustable. 21" and 22" seat depths require solid seat pan.

Frame Finishes

Select Color: High Gloss Black Hyper Blue

Armrests

Armrest Type

No Armrests	No charge
Single Post, Height Adjustable Armrests	No charge
8"-11" Range 10-14" Range 12"-16" Range	
Dual Post, Height Adjustable Armrests	No charge
8"-11" Range 10-14" Range 12"-16" Range	
Armrest Position: Rear Front	

If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.

Armpads

Regular Armpads	No charge
Desk Length 2" x 10" Full Length 2" x 14"	

Clothing Guards | *Not compatible with flip-up armrests.*

Clothing Guards, Single Post Armrests <i>Not compatible with 8"-11" height range.</i>	\$ 125.00
Clothing Guards, Dual Post Armrests <i>Not compatible with 8"-11" height range.</i>	\$ 125.00

Armrest Extensions | *Not compatible with flip-up armrests or desk length armpads.*

Armrest Extensions - 3" Forward	\$ 75.00
Armrest Extensions - 3" Rearward	\$ 75.00

Seat Height

Seat Height - 6" Casters | *Seat height is adjustable within the range of the chosen lower frame type (Low/High)*

RW		Low						High					
	24"												
	22"												
	20"												
	16"												
	12"												
		13"	14"	15"	16"	17"	18"	15"	16"	17"	18"	19"	20"
		Seat to Floor Height											

Seat Height - 8" Casters | *Seat height is adjustable within the range of the chosen lower frame type (Low/High)*

RW		Low				High			
	24"								
	16"								
		15"	16"	17"	18"	17"	18"	19"	20"
		Seat to Floor Height							

Seat Pan

No Seat Pan	No charge
Solid Aluminium Seat Pan <i>Required for 21" and 22" seat depths.</i>	No charge

Front Casters

Urethane Casters

6" x 1.25" Caster	No charge
8" x 1.25" Caster <i>Only compatible with 16" and 24" rear wheels.</i>	No charge

Rear Wheels

Composite Rear Wheel	No charge
12" 16" 20" 22" 24"	
Heavy Duty Spoke Wheel	\$ 350.00
24"	

Rear Wheel Position: Rear Front (Default)

If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

Axles

Quick Release Rear Axles	No charge
Fixed Rear Axles	No charge

Rear Tires

Urethane Tires	No charge
-----------------------	-----------

Handrims | Handrims are only available with 20", 22" and 24" rear wheels.

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims <i>Not compatible with drum brakes.</i>	\$ 175.00

Wheel Locks

User Wheel Locks	No charge
Push to Lock Pull to Lock	

Wheel Lock Options

Wheel Lock Extensions <i>Not compatible with attendant wheel lock.</i>	No charge
---	-----------

Anti-Tippers

Rear Anti-Tippers <i>Required safety option.</i>	No charge
---	-----------

Backrest

Back Post Style:	Stroller Handle Back
	<i>Seat widths greater than 22" require stroller handle back. Cane handle back may reduce tilt range. Cane handle back not compatible with power tilt with manual override, dynamic rocker back or endomorph back.</i>
Back Post Height:	Short (22") Tall (25")
	<i>Measured from seat pan to handle.</i>
Folding Style:	Fold-Down
	<i>Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.</i>
Seat to Back Angle:	Set at 90° Set at 98° Set at 106° Set at 114° Set at 122°

Backrest Options

Angle/Height Adjustable Stroller Handle: <i>Not compatible with cane handle back.</i>	No charge
Button on Outside Button on Inside	
Spring Back Suspension <i>Not compatible with reclining back pr dynamic rocker back.</i>	No charge

Tilt Activation

Attendant Tilt Only

Attendant Tilt Only, Right Back Post	No charge
Attendant Tilt Only, Left Back Post	No charge

Tilt Limit | *Optional. If no tilt limit is selected the full tilt range will be available.*

Tilt Limit:	Set at 20° Set at 30° Set at 40°
--------------------	----------------------------------

Anterior Pre-Tilt | *Add up to 5° of Anterior Tilt to the tilt range.*

Anterior Pre-Tilt	Set at 0° (Default) Set at 5°
--------------------------	-------------------------------

Front Rigging

No Front Rigging No charge

Easy-Out Front Rigging | *Effortless insertion and removal. No hanger brackets remain on the frame.*

Fixed 70° \$ 95.00

Swing Away Front Rigging | *Legacy design. Not compatible with Easy-Out.*

Swing Away 70° No charge

Swing Away Elevating Legrests (ELR) *Includes calf pads.* W106 \$ 375.00

1" Offset Hangers for Swing Away Front Rigging \$ 95.00

Plug In Front Rigging | *Legacy design. Not compatible with Easy-Out.*

Plug In 70° No charge

Footplates

Composite Footplates \$ 0.00

Aluminium Angle Adjustable Footplates W107 \$ 165.00

Extension Tubes | *Consider the height of your seat cushion.*

Standard Extension Tubes, Short No charge

Standard Extension Tubes, Long \$ 75.00

Footplate Height:

Specify:

Ext	Std Short	11" - 16"	12" - 16"	13" - 16"
	Std Long	14" - 22"	14" - 19"	14" - 17"
		Standard Easy-Out	60°, 70°, ELR Swing-Away	Standard Plug-In
		Legrest		

Front Rigging Accessories

Length Adjustable Heel Loops No charge

Standard Calf Strap Basic Seating \$ 50.00

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes