

Front wheel



Mid wheel



Rear wheel



| Paediatric Series    | Koala             | K300 PS Jr.                                      | M300 PS Jr.                              | K450 MX                                     |
|----------------------|-------------------|--------------------------------------------------|------------------------------------------|---------------------------------------------|
| Max user weight      | 35 kg             | 75 kg                                            | 75 kg                                    | 60 kg                                       |
| Max speed            | 6.5 km/h          | 8 km/h   10 km/h                                 | 10 km/h                                  | 8 km/h                                      |
| Range                | 24 km*            | 25 km*                                           | 20 km*                                   | 27 km*                                      |
| Base width           | 23"               | 24"                                              | 24"                                      | 25"                                         |
| Base length          | 32"               | 40"                                              | 36.5"                                    | 35"                                         |
| Min turning radius   | 21"               | 25"                                              | 20"                                      | 28"                                         |
| Ground clearance     | 2.5"              | 3"                                               | 3"                                       | 2.5"                                        |
| Weight w/ batteries  | 102 kg (M40)      | 123 kg (Gr 34)                                   | 142 kg (Gr 34)<br>151 kg (Gr 24)         | 153 kg (Gr 34)                              |
| Seat-to-floor height | 15" Elevator      | 16.5" Elevator<br>16.5", 17.25", 18.5" Fixed     | 17" Elevator<br>17", 17.75", 18.5" Fixed | 3" – 26" Elevator                           |
| Seat elevator        | 8"                | Optional 8"                                      | Optional 8"                              | Optional 8"                                 |
| Power tilt options   | 0° – 25° Power    | 0° – 45° Manual or Power                         | 0° – 45° Manual or Power                 | 0° – 45° Power                              |
| Recline options      | 85° – 120° Manual | 85° – 120° Manual                                | 85° – 120° Manual                        | 85° – 105° Manual                           |
| Legrest options      | 90° – 150° Manual | 70° or 90° Swing-Away<br>90° – 150° Centre Mount | 70° Swing-Away<br>90 – 150° Centre Mount | 90° – 125° Manual                           |
| Seat depth           | 10", 12", 14"     | 10" – 18"                                        | 10" – 18"                                | 10" – 18" (by 1" incr.)                     |
| Seat width           | 10", 12", 14"     | 11" – 16" (by 1" incr.)                          | 11" – 16" (by 1" incr.)                  | 10", 12", 14", 16"                          |
| Backrest height      | 14", 16", 20"     | 16.5" – 22.5" (by 1" incr.)                      | 16.5" – 22.5" (by 1" incr.)              | 14", 16", 20" (MX)<br>16.5" – 22.5" (PS Jr) |

All specifications may vary based upon configuration. | Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories. | \*Actual driving range will vary based on driving and battery conditions.

# Power Wheelchairs

Independence in motion.



Front wheel

Mid wheel

Rear wheel



| Rehab Series         | F3 Corpus                                         | F5 Corpus                                         | F5 Corpus VS            | M1                                                | M3 Corpus                                    | M5 Corpus                                    | M300 HD                               | X850                                  |
|----------------------|---------------------------------------------------|---------------------------------------------------|-------------------------|---------------------------------------------------|----------------------------------------------|----------------------------------------------|---------------------------------------|---------------------------------------|
| Max user weight      | 136 kg                                            | 136 kg                                            | 136 kg                  | 136 kg                                            | 136 kg                                       | 136 kg                                       | 204 kg                                | 136 kg                                |
| Max speed            | 8 km/h   10 km/h                                  | 10 km/h   12 km/h                                 | 10 km/h   12 km/h       | 8 km/h                                            | 10 km/h                                      | 10 km/h   12 km/h                            | 8 km/h                                | 10 km/h   15 km/h                     |
| Range                | 32.2 km*                                          | 32.2 km*                                          | 32.2 km*                | 25.8 km*                                          | 32.2 km*                                     | 32.2 km*                                     | 32.2 km*                              | 45 km*                                |
| Base width           | 24"                                               | 25.5"                                             | 25.5"                   | 24"                                               | 24"                                          | 25"                                          | 26"                                   | 28"                                   |
| Base length          | 40" with anti-tips<br>36" without anti-tips       | 42" with anti-tips<br>38" without anti-tips       | 43"                     | 35.5"                                             | 35.5"                                        | 37.5"                                        | 38"                                   | 49.5" – 57.5"                         |
| Min turning radius   | 26.5"                                             | 30"                                               | 30"                     | 20.5"                                             | 20.5"                                        | 21.5"                                        | 22"                                   | 66"                                   |
| Ground clearance     | 3"                                                | 3"                                                | 3"                      | 3"                                                | 3"                                           | 3"                                           | 3"                                    | 3.2"                                  |
| Weight w/ batteries  | 175 kg (Gr 34)   184 kg (Gr 24)                   | 186 kg (Gr 24)                                    | 190 kg (Gr 24)          | 145 kg (Gr 34)                                    | 175 kg (Gr 34)   184 kg (Gr 24)              | 192 kg (Gr 24)                               | 193 kg (Gr 34)   202 kg (Gr 24)       | 203 kg (M31)                          |
| Seat-to-floor height | 17.5", 18.5", 19.5"                               | 17.5", 18.5", 19.5"                               | 17.5", 18.5", 19.5"     | 16", 17", 18"                                     | 17.5", 18.5", 19.5"                          | 17.5", 18.5", 19.5"                          | 18", 19", 20"                         | 21"                                   |
| Seat elevator        | Optional 12"                                      | Optional 14"                                      | 14"                     | —                                                 | Optional 12"                                 | Optional 14"                                 | Optional 8"                           | Optional 8"                           |
| Power tilt options   | 0° – 50° Posterior<br>5°, 10°, 20° & 30° Anterior | 0° – 50° Posterior<br>5°, 10°, 20° & 45° Anterior | 0° – 50° Posterior      | 0° – 50° Posterior                                | 0° – 50° Posterior<br>5°, 10° & 20° Anterior | 0° – 50° Posterior<br>5°, 10° & 20° Anterior | 0° – 45° Posterior                    | 0° – 15° Posterior                    |
| Recline options      | 85° – 120° Manual<br>85° – 180° Power             | 85° – 120° Manual<br>85° – 180° Power             | 85° – 180° Power        | —                                                 | 85° – 120° Manual<br>85° – 180° Power        | 85° – 120° Manual<br>85° – 180° Power        | 85° – 120° Manual<br>85° – 150° Power | 90° – 130° Manual<br>90° – 120° Power |
| Legrest options      | 85° – 180° Manual or Power                        | 85° – 180° Manual or Power                        | 90° – 180° Power        | 80° – 135° Manual<br>65° – 80° Adjust. Swing-Away | 85° – 180° Manual or Power                   | 85° – 180° Manual or Power                   | 85° – 180° Manual or Power            | 110° – 160° Manual                    |
| Seat depth           | 14" – 22" (by 1" incr.)                           | 14" – 22" (by 1" incr.)                           | 14" – 22" (by 1" incr.) | 14" – 22" (by 1" incr.)                           | 14" – 22" (by 1" incr.)                      | 14" – 22" (by 1" incr.)                      | 19" – 25" (by 1" incr.)               | 14" – 22" (by 1" incr.)               |
| Seat width           | 17", 19", 21", 23"                                | 17", 19", 21", 23"                                | 17", 19", 21", 23"      | 14" – 22" (by 2" incr.)                           | 17", 19", 21", 23"                           | 17", 19", 21", 23"                           | 19", 21", 23", 25"                    | 17", 19", 21", 23"                    |
| Backrest height      | 20", 23" – 28" (by 1" incr.)                      | 20", 23" – 28" (by 1" incr.)                      | 23" – 28" (by 1" incr.) | 18" – 29"                                         | 20", 23" – 28" (by 1" incr.)                 | 20", 23" – 28" (by 1" incr.)                 | 21", 24" – 29" (by 1" incr.)          | 20", 23" – 28" (by 1" incr.)          |

All specifications may vary based upon configuration. | Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories. | \*Actual driving range will vary based on driving and battery conditions.