

## Fitting / Order form

Date: 2024-04-08

Dealer information	
Customer (name / code)	/
Contact	
Address	
Post code / City	/
Phone / Fax	/
E-mail	

## Promoted adaptations

Order information	
User reference	
Your reference	
Your order number	

Wheelchair serial number  
(for existing chair)

Fitting date	User weight
Fitting date (dd-mm-yyyy)	kg

### Chassis

#### Transportation tray

F3, F5, F5 VS, M3, M5, M VS


*This adaptation makes it possible to transport various objects on your wheelchair.*

*The tray is easy to detach.*

*If the wheelchair is equipped with a luggage box you are able to switch between the two, since the chair has the same attachment.*

*A strap/belt is included.*

*Max payload 15 kg.*



- A200001  Transportation tray W:320 mm D:250 mm H:50 mm
- A200002  Transportation tray W:400 mm D:300 mm H:50 mm
- A200003  Transportation tray W:460 mm D:400 mm H:50 mm

#### Hydrographic colour kit

F3, F5, F5 VS, M3, M5, M VS

*Individualise your colours to present your personality*



- | M3, M5  | MVS     | F3, F5  |                                       |
|---------|---------|---------|---------------------------------------|
| A200004 | A200007 | A200012 | <input type="checkbox"/> Carbon fibre |
| A200005 | A200008 | A200013 | <input type="checkbox"/> Fire ball    |
| A200006 | A200009 | A200014 | <input type="checkbox"/> Water lilly  |

### Seating

#### Extended seat depth 620mm

F3, F5, F5 VS, M5

*This extended seat depth of 620 mm is intended for tall users and provides improved biomechanics and optimises support for the legs.*


- A200062  Extended seat depth 620mm



**Seat width 370 mm**

F3, F5, F5 VS, M3, M5, M VS

*Reduction of seat width to accommodate smaller users*




A200068  **Seat width 370 mm**

**Ergo Air Seat cushions (special size)**

F3, F5, F5 VS, M3, M5, M VS

*Additional sizes of the Ergo Air seat cushion*



A200066  **Ergo air seat cushion W420 x D570**

A200067  **Ergo air seat cushion W520 x D570**

**Positioning**

**Chest support roll VS**

F5 Corpus VS, M VS

*This chest support with rolling padding gives you a softer feeling compared to the regular chest support.*

*The chest support padding will be mounted on the ordered chest support bracket.*

A200013  **Chest support roll (VS) Vinyl 340mm**

A200014  **Chest support roll (VS) Vinyl 390mm**

A200015  **Chest support roll (VS) Vinyl 440mm**

A200016  **Chest support roll (VS) Textile 340mm**

A200017  **Chest support roll (VS) Textile 390mm**

A200018  **Chest support roll (VS) Textile 440mm**



**Chest support roll Corpus**

Corpus

*This chest support with rolling padding gives you a softer feeling compared to the regular chest support.*

A200019  **Chest support roll Vinyl holder fixed 340mm**

A200020  **Chest support roll Vinyl holder fixed 390mm**

A200021  **Chest support roll Vinyl holder fixed 440mm**

A200022  **Chest support roll Textile holder fixed 340mm**

A200023  **Chest support roll Textile holder fixed 390mm**

A200024  **Chest support roll Textile holder fixed 440mm**

A200025  **Chest support roll Vinyl holder adjustable 340mm**

A200026  **Chest support roll Vinyl holder adjustable 390mm**

A200027  **Chest support roll Vinyl holder adjustable 440mm**

A200028  **Chest support roll Textile holder adjustable 340mm**

A200029  **Chest support roll Textile holder adjustable 390mm**

A200030  **Chest support roll Textile holder adjustable 440mm**




**Fixed armrest holders (Orthesis)**

F3, F5, M3, M5, X850

The fixed armrest holders are mounted on the seat frame and are adjustable in depth, height, and angle. These armrest DO NOT follow the backrest movement.

Note: The armrest needs to be ordered separately



A200069

Fixed armrest holders

A200070

Fixed armrest holders for seat width 57cm

Armrest height "H":  mm

**Power single flip-up armrest**

F3, F5, F5 VS, M3, M5, M VS

Provide you with the possibility to flip up armrest when entering the chair or when you need to have free space to the steering wheel in the car.

Specify operating buttons/switches and their position.



A200034

Pwr single flip-up armrest (Joystick)

A200035

Pwr single flip-up armrest (1 Button)

A200036

Pwr single flip-up armrest (2 Buttons)

left

right

Specify operating buttons/switches and their position in "Control buttons".

**Power independent dual flip-up armrests**

F3, F5, F5 VS, M3, M5, M VS

Provide you with the possibility to flip up armrests when entering the chair or when you need to have free space to the steering wheel in the car.

Specify operating buttons/switches and their position.



A200037

Pwr independent dual flip-up armrests (Joystick)

A200038

Pwr independent dual flip-up armrests (1 Button)

A200039

Pwr independent dual flip-up armrests (2 Buttons)

Specify operating buttons/switches and their position in "Control buttons".

**Double articulated armrest**

F3, F5, F5 VS, M3, M5, M VS

Sideways flip-up armrest reduces the distance between the armrest with 5 cm.

The armrest can still be flipped up due to the semi-sideways flip-up construction.



A200033

Double articulated armrest (left/right)

**Locking pin**

F3, F5, F5 VS, M3, M5, M VS

Locking pin is used to lock the swing away table.



A200040

Locking pin

Mounted

Distance: back - stomach cutout  cm

Not mounted

**Armrest left/right fixed (follows backrest)**

F3, F5, F5 VS, M3, M5, M VS

*Armrests stay perpendicular to the backrest during backrest recline and standing.  
You are still able to flip up the armrests.*



A200031  **Armrest left/right fixed (follows backrest)**

**Locked armrest**

F3, F5, F5 VS, M3, M5, M VS

*The armrest flip up function is locked. Pull the locking knob to flip up the armrests.  
The armrest keep the parallel motion when reclining backrest the backrest or going into stand position.*



A200032  **Locked armrest**


**Leg positioning**

**Swing-away legrests**

F3, F5, M3, M5

*This customization features detachable (tool less) swing away legrests for easy transfer in and out of the wheelchair.*

*Legrests angle (90°-115°) and footplate angles (Dorsi/Plantar Flexion up to 360° / Inversion/eversion 30°) are adjustable (tool required).*




A200041  **Swing-away legrests**

Lowerleg length (seatplate - footplate):  mm

**Extended knee support**

F3, F5, F5 VS, M3, M5, M VS

*This extends the range of the knee support in height and/or in depth.*



A200042  **Extended knee supports +50mm deep**

A200043  **Extended knee supports +100mm depth**

A200044  **Extended knee supports +50mm height**

A200045  **Extended knee supports +100mm height**

**Knee support with gel insert**

F3, F5, F5 VS, M3, M5, M VS

*This replaces the standard foam in the knee support with a Gel insert for a softer feeling.*



A200046  **Knee support with gel insert**

**Power flip-up single footplate**

F3, F5, F5 VS, M3, M5, M VS

*This adaptation gives you the option to control/flip up both footplates simultaneously through the joystick or buttons.*

*Minimum lower leg length 285 mm seatplate to footplate*

*Specify operating buttons/switches and their position.*

- A200047  **Pwr flip-up single footplate (Joystick)**
- A200048  **Pwr flip-up single footplate (1 Button)**
- A200049  **Pwr flip-up single footplate (2 Buttons)**



**Specify operating buttons/switches and their position in "Control buttons".**

Lowerleg length (seatplate - footplate):  mm

**Power independent flip-up footplates**

F3, F5, F5 VS, M3, M5, M VS

*This adaptation gives you the option to control/flip up your left and right footplate independently through the joystick or buttons.*

*Minimum lower leg length 285 mm seatplate to footplate*

*Specify operating buttons/switches and their position.*

- A200050  **Pwr independent flip-up footplates (Joystick)**
- A200051  **Pwr independent flip-up footplates (1 Button)**
- A200052  **Pwr independent flip-up footplates (2 Buttons)**



**Specify operating buttons/switches and their position in "Control buttons".**

Lowerleg length (seatplate - footplate):  mm

**Separate powered legrests**

F3, F5, M3, M5

*Left and right legrest are individually power adjustable, controlled via joystick or switch box.*

*They can be controlled independently or simultaneously*

To be operated with:

- A200053  **Separate powered legrests**
  - Joystick seating mode
  - Seat function switch box



**Short articulating legrest**

F3, F5, F5 VS, M3, M5, M VS

*For users that have short lowerlegs and who require more ground clearance and the possibility to lower the footplates to the ground to enable easy transfer.*

- A200064  **Short articulating legrest (-36mm)**





Alternative controls

**Rear view camera**

F3, F5, F5 VS, M3, M5, M VS

*Gives additional safety whilst maneuvering through tight spaces. Display is preferable mounted on one of the armrests with Permofix.*

*This Garmin system includes satellite navigation.*



A200054

**Rear view camera**

Display mounting position:  left  
 right

**Rain cover for joysticks**

All models

*Protects your joystick from the elements.*

A200055

**Rain cover CJSM2**

A200056

**Rain cover CJSM2+SB**

A200057

**Rain cover SB**

A200058

**Rain cover JS X850**

A200059

**Rain cover SB X850**

A200060

**Rain cover JSM Led**

A200061

**Rain cover JSM led + SB**



**panel holder power parallel swing away**

F3, F5, F5 VS, M3, M5, M VS

*Provides an easy lateral adjustment/movement of the parallel motion panel holder. The feature enables you to reposition the joystick electrically to your preference. The powered adjustable panelholder can be controlled via button/s or Joystick*

*ex:*

*1 button, with toggle function*

*2 buttons, one button for each direction*

*Via joystick in mode function.*



A200063

**panel holder power parallel swing away Left/right**

A200071

**panel holder power parallel swing away Left/right (2 buttons)**

**Specify operating buttons/switches and their position in "Control buttons".**

**Foot control**

F3, F5, F5 VS, M3, M5, M VS

*This foot control is a gas throttle with a left and right movement for steering. Switching between forwards and backwards driving is set in the OMNI parameters and is either a switch, auto toggle, or short command with the foot control.*



A200065

**Foot control (3 axis prop)**

**Control buttons**

All models

*A robust and heavy duty button. Available in three colours.*



	<i>Button 1</i>	<i>Button 2</i>	
599242-99-0	<input type="checkbox"/>	<input type="checkbox"/>	<b>Pico Button red</b>
599243-99-0	<input type="checkbox"/>	<input type="checkbox"/>	<b>Pico Button yellow</b>
599244-99-0	<input type="checkbox"/>	<input type="checkbox"/>	<b>Pico Button blue</b>
101766-99-0	<input type="checkbox"/>	<input type="checkbox"/>	<b>Buddy button black 63 mm</b>
101767-99-0	<input type="checkbox"/>	<input type="checkbox"/>	<b>Buddy button red 63 mm</b>
101768-99-0	<input type="checkbox"/>	<input type="checkbox"/>	<b>Buddy button yellow 63 mm</b>

Function Button 1

Position Button 1

Function Button 2

Position Button 2

