

Kombinationsavtal med drivaggregat



SmartDrive

Mototronik

	SmartDrive	Mototronik
Panthera X	●	●
Panthera X3	●	●
Panthera U3 Light	●	●
Panthera U3	●	●
Panthera S3*	●	●
Panthera S3 Swing	●	●
ETAC Cross	●	
Invacare Kuschall – K-series	●	
Sunrise Medical HELIUM	●	●
Sunrise Medical XENON	●	●
Sunrise Medical XENON HYBRID		●
Sunrise Medical NITRIUM		●
Sunrise Medical NITRIUM HYBRID		●
Sunrise Medical ARGON		●
Sunrise Medical EASY 200 / 300		●
Sunrise Medical LIFE		●
Sunrise Medical LIFE R/RT/T		●
Sunrise Medical NEON		●

*Mototronik - för fällbara fotplattor funkar enbart bredd 39, 42 och 45, då med assymetriska fotplattor