

New product

# Dreamline G7 Cushion

Learn more about the  
Dreamline G7 cushion



The Dreamline G7 cushion has a soft and comfortable feel, yet unyielding stability with a mild contoured base for easy transfers. Pressure redistribution is effective through its patented Cluster Gel Fluid packs, which are neatly hidden in a unique secure zip-up compartment.

## Key Features

- **Immediate Comfort.** The soft Comfort Wrap with the new FLO2™ fabric gives superior immersion and long-term sitting comfort.
- **Stability.** A firm, mildly contoured positioning base with a deep well, providing active users with a stable base of support.
- **Pressure Protection.** Large dual Cluster Gel Fluid packs nestled in 3cm of Enduro foam provides effective protection from skin breakdown, all secured within the FLO2™ cover for extra protection during transfers.



## FLO2™ 4-Way Stretch Fabric

FLO2™ is an advanced high-performance 4-way stretch fabric that utilises an engineered warp knit design to naturally facilitate moisture wicking and air exchange, helping to keep users cool and comfortable. The properties of the material enable the fabric to conform to the contours of asymmetrical configurations whilst providing skin protection through maximum immersion.



## Technical Specifications

<b>Weight Limit</b>	136kg (300lbs) - Up to 20"	
<b>Weight of Product</b>	3.46kg* (7.62lbs)	
<b>Width</b>	12" - 20" (30.5cm - 50.8cm)	
<b>Depth</b>	12" - 22" (30.5cm - 55.9cm)	
<b>Pelvic Well</b>	1.8" (4.5cm)	
<b>Positioning Base Height</b>	12"/14"	16" - 22"
	4.5cm	5.5cm
<b>Wrap Thickness</b>	2.5cm	3.0cm
<b>Gel Packs</b>	2 x Gel Fluid Packs	
<b>Base Foam</b>	EVA Closed Cell Foam	
<b>Comfort Wrap</b>	Enduro Foam with FLO2™ Cover	
<b>Covers</b>	FLO2™ Incontinence and Outer Cover	
	3D Spacer Side Panel	
	Neoprene Front Panel	
	Non-Slip Rubberised Urethane Base Panel	
<b>Adjustable Positioning</b>	Pelvic Obliquity / Pelvic rotation	
<b>Warranty</b>	18 Months	

\*Weight Based on 18" x 18" (45.8cm x 45.8cm)  
**Product Code D7C**

