

# M3 Corpus

# permobil



**\* Indicates a Required Field**

## Dealer Information

Order ID: \_\_\_\_\_

\*ATP: \_\_\_\_\_

Dealer Code: \_\_\_\_\_

\*Dealer Name: \_\_\_\_\_

Address: \_\_\_\_\_

\*City: \_\_\_\_\_

\* State/Zip: \_\_\_\_\_

\* Phone#: \_\_\_\_\_

PO#: \_\_\_\_\_

Email Address: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* First Name: \_\_\_\_\_

\* Last Name: \_\_\_\_\_

Diagnosis: \_\_\_\_\_

Funding Source: \_\_\_\_\_

Client Age: \_\_\_\_\_

## Client Measurements

\* Weight: \_\_\_\_\_

\*Height: \_\_\_\_\_

A) Top of Shoulders: \_\_\_\_\_

B) Chest Depth: \_\_\_\_\_

C) Chest Width: \_\_\_\_\_

D) Seat Depth: \_\_\_\_\_


E) Top of Head: \_\_\_\_\_

F) Elbow to Hand: \_\_\_\_\_

G) Seat Pan to Elbow: \_\_\_\_\_

H) Hip Width: \_\_\_\_\_

\*I) Knee to Foot: \_\_\_\_\_



Three column pricing on order/quote

Display HCPCS Codes on order/quote

\*\* Indicates HCPCS code is pending PDAC approval

Please send order/quote to fax#: **(800) 231-3256**

Email to: [sales@permobil.com](mailto:sales@permobil.com)

**Permobil Inc.**  
300 Duke Dr.  
Lebanon, TN 37090  
Tel: (800) 736-0925  
Fax: (800) 231-3256  
[www.permobil.com](http://www.permobil.com)

Prices effective **May 16, 2022**

# M3 Corpus

## Chair Model and Colors

<input type="checkbox"/> 11352	M3 Base Corpus - M3 SPO	8,457.00	K0856
	<i>Price Includes: MWD M3 Power Base Set Up For One Power Seat Function, Corpus Seat Frame, FlexLink Suspension, 7" Flat Free Casters, Tie-Down Hardware for Strap System, Charger, Flat Free or Pneumatic Tires and a Standard Positioning Belt. R-net 120amp Electronics. Max Speed = 6mph, Weight Capacity = 300lbs.</i>		
<input type="checkbox"/> 109052-99-0	M3 Base Corpus - M3 MPO	8,457.00	K0861
	<i>Price Includes: MWD M3 Power Base Set Up For Multiple Power Seat Functions, Corpus Seat Frame, FlexLink Suspension, 7" Flat Free Casters, Tie-Down Hardware for Strap System, Charger, Flat Free or Pneumatic Tires and a Standard Positioning Belt. R-net 120amp Electronics. Max Speed = 6mph, Weight Capacity = 300lbs.</i>		
<input type="checkbox"/> 11613	M3 Base Corpus - M3 SPO HD	9,957.00	K0858**
	<i>Price Includes: MWD M3 Power Base Set Up For One Power Seat Function, Corpus Seat Frame, FlexLink Suspension, 7" Flat Free Casters, Tie-Down Hardware for Strap System, Charger, Flat Free or Pneumatic Tires and a Standard Positioning Belt. R-net 120amp Electronics. Max Speed = 6mph, Weight Capacity = 330lbs.</i>		
<input type="checkbox"/> 11614	M3 Base Corpus - M3 MPO HD	9,957.00	K0862**
	<i>Price Includes: MWD M3 Power Base Set Up For Multiple Power Seat Functions, Corpus Seat Frame, FlexLink Suspension, 7" Flat Free Casters, Tie-Down Hardware for Strap System, Charger, Flat Free or Pneumatic Tires and a Standard Positioning Belt. R-net 120amp Electronics. Max Speed = 6mph, Weight Capacity = 330lbs.</i>		
<input type="checkbox"/> 108852-99-0	Permobil Connect	No Charge	
	<i>By selecting, user agrees to activate the Permobil Connect-ready wheelchair for connected services.</i>		

### CHOOSE DRIVE TIRES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 108414-99-0	Flat Free Tire 3"W x 14"T Black - Standard Tread	No Charge	
	<i>Designed for a smooth ride on indoor and outdoor surfaces. Includes flat free standard tread caster wheels.</i>		
<input type="checkbox"/> 108411-99-0	Pneumatic Tire 3"W x 14"T Black - Standard Tread	No Charge	
	<i>Includes flat free standard tread caster wheels.</i>		

### OPTIONAL DRIVE TIRES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 108412-99-0	Pneumatic/Snow Tire 3"W x 14"T Black - Aggr Tread	489.00	
	<i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Does not include caster wheels.</i>		

### CHOOSE SHROUD/RIM COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 109083-99-0	M3 Color Kit - Aqua	No Charge	
<input type="checkbox"/> 109084-99-0	M3 Color Kit - Sky	No Charge	
<input type="checkbox"/> 109085-99-0	M3 Color Kit - Violet	No Charge	
<input type="checkbox"/> 109086-99-0	M3 Color Kit - Flamingo	No Charge	
<input type="checkbox"/> 109087-99-0	M3 Color Kit - Chili	No Charge	
<input type="checkbox"/> 109088-99-0	M3 Color Kit - Canary	No Charge	
<input type="checkbox"/> 109089-99-0	M3 Color Kit - Spring	No Charge	
<input type="checkbox"/> 109090-99-0	M3 Color Kit - Arctic	No Charge	
<input type="checkbox"/> 109091-99-0	M3 Color Kit - Sunset	No Charge	
<input type="checkbox"/> 109176-99-0	M3 Color Kit - Midnight	No Charge	
<input type="checkbox"/> 109092-99-0	M3 Color Kit - Gold	No Charge	
<input type="checkbox"/> 109093-99-0	M3 Color Kit - Bronze	No Charge	
<input type="checkbox"/> 109094-99-0	M3 Color Kit - Platinum	No Charge	
<input type="checkbox"/> 109096-99-0	M3 Color Kit - Carbon Fiber	212.00	
<input type="checkbox"/> 109097-99-0	M3 Color Kit - Mossy Oak Break-Up Country	212.00	
<input type="checkbox"/> 109098-99-0	M3 Color Kit - CDG Poly	212.00	

# M3 Corpus

## Chair Model and Colors

### CHOOSE SHROUD/RIM COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 109100-99-0	M3 Color Kit - Fireball	212.00	
<input type="checkbox"/> 109101-99-0	M3 Color Kit - Water Lily	212.00	
<input type="checkbox"/> 109102-99-0	M3 Color Kit - Midnight Patriot	212.00	
<input type="checkbox"/> 109103-99-0	M3 Color Kit - American Liberty	212.00	

### CHOOSE HUB CAP KIT:

Part Number	Description	Price	HCPCS Code
	Hub Cap Color: <i>Please write in Hub Cap Color below. Hub Cap Color will be matched to the Shroud Color if no selection is made.</i> _____		
<input type="checkbox"/> 11343	Hydrographic Hub Cap Color: <i>Please write in Hub Cap Color below. Hub Cap Color will be matched to the Shroud Color if no selection is made.</i> _____	52.00	

## Base Options

### MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 11110	Active Height Power Adj. Seat Height - 12" Travel <i>Must also select Power Tilt when selecting this option.</i>	3,605.00	E2300
<input type="checkbox"/> 11117	Seat Tube Fixed	No Charge	

### SET SEAT HEIGHT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 107435-99-0	Set Seat Height to 17.5"	No Charge	
<input type="checkbox"/> 107433-99-0	Set Seat Height at 18.5" <i>If ordering Power Adjustable Seat Height, 1" of elevator travel is used to set seat at 18.5".</i>	No Charge	
<input type="checkbox"/> 107434-99-0	Set Seat Height at 19.5" <i>If ordering Power Adjustable Seat Height, 2" of elevator travel is used to set seat at 19.5".</i>	No Charge	

### MUST CHOOSE BATTERIES:

Part Number	Description	Price	HCPCS Code
	BATTERY OPTIONS:		
<input type="checkbox"/> 108999-99-0	Permobil Batteries, Grp 24 AGM GEL, Pair <i>Includes Permobil VoltPro charger optimized for Permobil batteries.</i>	1,042.00	E2363
<input type="checkbox"/> 109183-99-0	Permobil Batteries, Grp 34 AGM GEL, Pair <i>Includes Permobil VoltPro charger optimized for Permobil batteries.</i>	998.00	E2359

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 107998-99-0	Enhanced Steering Perform Unit (ESP), R-net <i>Provided at no charge when ordered with an Omni2 Control Module.</i>	1,580.00	
<input type="checkbox"/> 109170-99-0	Essentials Box with Reflectors <i>Please Note: The total backrest angle will be limited to 170°, or further depending on other backrest-mounted items, to prevent contact with the Essentials Box.</i>	509.00	
<input type="checkbox"/> 109169-99-0	Essentials Box with Lights <i>Please Note: The total backrest angle will be limited to 170°, or further depending on other backrest-mounted items, to prevent contact with the Essentials Box. Must also select Light Kit.</i>	830.00	
<input type="checkbox"/> 109375-99-0	LED Light Kit for M3/M5 <i>No charge with 4 power seat function configurations.</i>	1,035.00	
<input type="checkbox"/> 10715	Compact Laptop Ventilator Tray - Corpus Seat <i>Mounts to rear of Corpus backrest hinge. Accommodates Pulmonetic Systems / CareFusion LTV Series and similar sized Laptop Ventilators. The tray is 13"T, adjustable in width from 8" - 13", and adjustable in depth from 2" - 6.5". Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option. Additional drive restrictions may apply.</i>	1,671.00	E1029

# M3 Corpus

## Base Options

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10787	Battery Tray for Laptop Ventilator Tray <i>Measures 11 3/8" W x 6 1/2" H x 1 5/8" D. Accomodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	280.00	
<input type="checkbox"/> 10788	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 7" Tall. Width infinitely adjustable 3.5" - 7". Depth infinitely adjustable 3.75" - 5". Custom Option.</i>	355.00	E2208
<input type="checkbox"/> 109620-99-0	Programming Key for R-net - Permobil Version	547.00	

## Joystick Options

### CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
PLEASE NOTE: <i>For additional input devices, please refer to the Alternative Drive Systems Section.</i>			
<input type="checkbox"/> 107779-99-0	Color Joystick Module (CJSM) <i>This R-net color joystick features: A Large High-Resolution 3.5" VGA color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	1,232.00	
<input type="checkbox"/> 11186	Permobil Joystick Module, NC <i>Designed exclusively for Permobil, this R-net color joystick features: a 2.8" color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	No Charge	
<input type="checkbox"/> 107829-99-0	3.5mm Y-Splitter Cable for PJSM/CJSM <i>For use with the PJSM or CJSM Mode/Profile jack. Allows 2 mechanical switches to control 4 separate functions. Requires PC Programming to assign the short hold and long hold commands to each switch. Switches sold separately.</i>	37.00	
<input type="checkbox"/> 105246-99-0	R-net LED Remote Joystick - No Display <i>This Joystick does NOT have a Display Screen, does NOT contain mono jack ports, and does NOT include Bluetooth or Assignable Buttons. PC Programming required and profiles limited to five.</i>	No Charge	
<input type="checkbox"/> 105700-99-0	R-net Remote Color JS w/ Lights, Heavy Duty Gimbal <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Color LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	1,281.00	

### CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10045	Cone Shaped Joystick Knob - Standard	No Charge	
<input type="checkbox"/> 100960-99-0	Chin Cup For Joystick Knob	132.00	E2324
<input type="checkbox"/> 100961-99-0	"T" Handle For Joystick Knob	186.00	E2323
<input type="checkbox"/> 100962-99-0	"Large Ball" For Joystick Knob	132.00	E2323
<input type="checkbox"/> 100963-99-0	"Softball" For Joystick Knob	174.00	E2323
<input type="checkbox"/> 100964-99-0	"Mushroom" For Joystick Knob	174.00	E2323
<input type="checkbox"/> 100965-99-0	"Stick" For Joystick Knob	174.00	E2323
<input type="checkbox"/> 105442-99-0	1.75" Black Foam Ball For Joystick Knob	132.00	E2323
<input type="checkbox"/> 105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	174.00	E2323
<input type="checkbox"/> 105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	174.00	E2323
<input type="checkbox"/> 105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	174.00	E2323

### CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10664	Joystick Mount Right - Fixed	No Charge	
<input type="checkbox"/> 10665	Joystick Mount Left - Fixed	No Charge	
<input type="checkbox"/> 102537-99-0	Joystick Mount Right - Swing Away	389.00	E1028
<input type="checkbox"/> 102538-99-0	Joystick Mount Left - Swing Away	389.00	E1028

# M3 Corpus

## Joystick Options

### CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 108037-99-0	Slimline Retractable Joystick Mount - Right	523.00	E1028
<input type="checkbox"/> 108038-99-0	Slimline Retractable Joystick Mount - Left	523.00	E1028
<input type="checkbox"/> 11232	Adjustable Height Panel Bracket <i>Provides 1.77" of vertical adjustment below the armrest and 25° of angle adjustment. Compatible with all Permobil joystick mounts except for the M1 Fixed Joystick Mount.</i>	284.00	E1028
<input type="checkbox"/> 105210-99-0	Mount Standard Joystick on Corpus Back Hinge <i>Standard R-net Joystick is mounted to backrest hinge. Can be chosen instead of PN 1105408-99-0. Mount on Left: _____ Mount on Right: _____</i>	358.00	
<input type="checkbox"/> 108268-99-0	Bodypoint Tri-Lock Midline Mount Assembly <i>Includes L-shaped joystick arm and height adjustable mounting clamp. Joystick mount sold separately (select from the mounts listed below). Please indicate if you want midline mounted on left or right armrest. Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	668.00	
<input type="checkbox"/> 108280-99-0	Bodypoint Tri-Lock Mount for R-net JS	278.00	
<input type="checkbox"/> 108279-99-0	Bodypoint Tri-Lock Mount for Compact Joystick	278.00	
<input type="checkbox"/> 10842	PERMOFix Mount for Mini/Compact Joystick <i>Includes a shim washer that allows the PERMOFix joint to be mounted on the joystick arm. Provides an adjustable mount for a Mini Joystick or Compact Joystick.</i>	278.00	

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10411	Expandable Controller - R-net <i>This item is required when ordering power tilt, power recline, and power elevating legrests together OR when ordering an Alternative Drive Control.</i>	727.00	E2377
<input type="checkbox"/> 10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller.</i>	645.00	E2313
<input type="checkbox"/> 105225-99-0	Co-Pilot Attendant Control R-net - Corpus	5,320.00	E2331
<input type="checkbox"/> 105408-99-0	Attendant Control R-net	1,427.00	E2331
<input type="checkbox"/> 10582	Output Module (IOM) for R-net Control System <i>Enables a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	1,344.00	

## Seating System

### CHOOSE POSTERIOR TILT/RECLINE PACKAGE:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>Certain seat positions trigger speed limitations and/or drive inhibits based on product configuration and/or user weight.</i>		
<input type="checkbox"/> 11193	Corpus Seat w/ No Tilt & 180° Power Recline	9,100.00	E1004
<input type="checkbox"/> 11194	Corpus Seat w/ 50° Tilt <i>Backrest angle can be manually adjusted up to 120°.</i>	8,119.00	E1002
<input type="checkbox"/> 11195	Corpus Seat w/ 50° Tilt & 180° Power Recline	12,849.00	E1007

### ACTIVE REACH PACKAGES:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>Availability of ActiveReach and drive capabilities may be limited in certain configurations.</i>		
<input type="checkbox"/> 11141	5° Active Reach	No Charge	
	<i>When ordered with a Power Adjustable Seat Height, tilt range will be -5° to +50°. When ordered without a Power Adjustable Seat Height, tilt range will be -5° to +45°. If this option is selected with a fixed seat, the minimum seat to floor height at 0° of tilt will be approximately 19.5". Power Legrest recommended for full functionality. Minimum legrest angle will be ~ 97° when manual legrest is selected.</i>		
<input type="checkbox"/> 107090-99-0	10° Active Reach	No Charge	

# M3 Corpus

## Seating System

### ACTIVE REACH PACKAGES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 107090-99-0	10° Active Reach <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height and Power Tilt. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Power Legrest recommended for full functionality. Minimum legrest angle will be ~ 103° when manual legrest is selected.</i>	No Charge	
<input type="checkbox"/> 107091-99-0	20° Active Reach <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height, Power Tilt, Power Recline and Power Legrest Elevation. Must select Corpus VS Power Elevating Legrest, Corpus VS Footplates and Corpus VS Knee Support from VS Legrest Options section. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Not compatible with retractable positioning belts.</i>	2,825.00	

## Seating System Options

### CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10418	Seat Function Control Switch for R-net <i>Does not include power seat function control through the drive control.</i>	No Charge	
<input type="checkbox"/> 10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control. Order with R-net electronics or remote joystick options.</i>	2,009.00	E2310
<input type="checkbox"/> 10585	Multiple Seat Function Control Kit for R-net <i>For operation of two or more seat functions through the wheelchair electronics via the drive control. Order with R-net electronics or remote joystick options.</i>	3,803.00	E2311

### CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price	HCPCS Code
PLEASE NOTE: <i>ICS switchbox availability may be limited to certain configurations.</i>			
<input type="checkbox"/> 102798-99-0	Push Buttons on ICS Switchbox (1.4m)	No Charge	
<input type="checkbox"/> 102799-99-0	Toggle Switches on ICS Switchbox (1.4m)	340.00	
<input type="checkbox"/> 105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge	

### CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> 104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge	
<input type="checkbox"/> 106256-99-0	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	163.00	

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,329.00	
<input type="checkbox"/> 108416-99-0	Remote Stop System <i>Consists of a 315MHz Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair.</i>	1,315.00	
<input type="checkbox"/> 10750	Memory Seat Program, Corpus w/ R-net	137.00	
<input type="checkbox"/> 106209-99-0	Legrest Movement Linked to Recline Electronically <i>Power Legrest will raise when Power Recline moves backward and will lower when back moves forward. If Power Legrest is operated, Power Recline will NOT automatically reposition. Feature works through ICS Switchbox and Driver Control.</i>	52.00	

# M3 Corpus

## Seating System Options

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10898	<b>Enable Independent Repositioning Mode</b> <i>IRM allows for sequential movement of power seat functions into a posterior position (35° tilt, 115° recline, 135° legrest) and returns to upright position (5° tilt, 100° recline, 110° legrest). This can be accessed using a standard drive control or a single switch. On the Permobil Joystick Module, IRM is accessible through the mode/profile jack (long switch hold). Please see the Alt Drive Accessories section for switch options.</i>	No Charge	

## Backrest Options

### CHOOSE FIXED OR SLIDING BACKREST KIT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105067-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No Charge	
<input type="checkbox"/> 105066-99-0	Fixed Backrest Kit - No Mechanical Shear	41.00	

### CORPUS BACKREST OPTIONS:

Part Number	Description	Price	HCPCS Code
<b>BACKREST HEIGHT NOTE:</b>			
<i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Corpus Ergo Backs include lateral wedges and foam lumbar pads. Solid Backs do not include foam lumbar pads.</i>			
<input type="checkbox"/> 106545-99-0	Corpus Ergo Back - 14"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,135.00	E2620
<input type="checkbox"/> 106546-99-0	Corpus Ergo Back - 16"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,135.00	E2620
<input type="checkbox"/> 106547-99-0	Corpus Ergo Back - 18"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,135.00	E2620
<input type="checkbox"/> 106548-99-0	Corpus Ergo Back - 20"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	1,135.00	E2620
<input type="checkbox"/> 11104	Solid Corpus Ergo Back - 16"W x 20"T <i>Power Recline will be limited to 150° when ordering this backrest.</i>	No Charge	
<input type="checkbox"/> 11105	Solid Corpus Ergo Back - 18"W x 20"T <i>Power Recline will be limited to 150° when ordering this backrest.</i>	No Charge	
<input type="checkbox"/> 11106	Solid Corpus Ergo Back - 20"W x 20"T <i>Power Recline will be limited to 150° when ordering this backrest.</i>	No Charge	

### PLEASE NOTE:

*Taller backrests listed below include: lateral wedges, foam lumbar pads\*, UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided. \* Foam lumbar pads not included with Solid Backs.*

<input type="checkbox"/> 106532-99-0	Corpus Ergo Back - 14"W x 23" - 28"T (C3GERGBCA)	1,139.00	E2620
<input type="checkbox"/> 106533-99-0	Corpus Ergo Back - 16"W x 23" - 28"T (C3GERGBCA)	1,139.00	E2620
<input type="checkbox"/> 106534-99-0	Corpus Ergo Back - 18"W x 23" - 28"T (C3GERGBCA)	1,139.00	E2620
<input type="checkbox"/> 106535-99-0	Corpus Ergo Back - 20"W x 23" - 28"T (C3GERGBCA)	1,139.00	E2620
<input type="checkbox"/> 11107	Solid Corpus Ergo Back - 16"W x 23" - 28"T	No Charge	
<input type="checkbox"/> 11108	Solid Corpus Ergo Back - 18"W x 23" - 28"T	No Charge	
<input type="checkbox"/> 11109	Solid Corpus Ergo Back - 20"W x 23" - 28"T	No Charge	

### MUST ALSO SELECT BACK HEIGHT:

- Set Corpus Backrest at 23" \_\_\_\_\_
- Set Corpus Backrest at 24" \_\_\_\_\_
- Set Corpus Backrest at 25" \_\_\_\_\_
- Set Corpus Backrest at 26" \_\_\_\_\_
- Set Corpus Backrest at 27" \_\_\_\_\_

# M3 Corpus

## Backrest Options

### CORPUS BACKREST OPTIONS:

Part Number	Description	Price	HCPCS Code
	Set Corpus Backrest at 28" _____		

### CORPUS BACKREST CUSHION OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 11335	Leatherette Cushion - Corpus Ergo Back <i>Includes Dual Density Backrest Foam.</i>	No Charge	
<input type="checkbox"/> 11336	Stretch-Air Cushion - Corpus Ergo Back <i>Includes Dual Density Backrest Foam with Stretch-Air breathable fabric cover.</i>	No Charge	
<input type="checkbox"/> 11339	Foam Lumbar Pads for Solid Back - Corpus	42.00	

### AGILITY BACKREST OPTIONS:

Part Number	Description	Price	HCPCS Code
	<b>ROHO AGILITY NOTE:</b> <i>The ROHO AGILITY Back Systems are available in widths of 14", 16", 18" and 20" and shell heights of 16", 18" and 20". Includes backrest shell, cushion, air inserts, UT rail for headrest and Direct Mount adapter. If ordered with a headrest, extended headrest receiver will be provided for backrests heights less than 25".</i>		
<input type="checkbox"/> 11258	ROHO AGILITY Max Contour Direct Mount - 16" Tall <i>Backrest can be installed at a height (seat pan to top of AGILITY) between 22" and 25". Specify Shell Width (14", 16", 18", 20"W): _____</i>	832.00	E2620
<input type="checkbox"/> 11259	ROHO AGILITY Max Contour Direct Mount - 18" Tall <i>Backrest can be installed at a height (seat pan to top of AGILITY) between 22" and 26". Specify Shell Width (14", 16", 18", 20"W): _____</i>	832.00	E2620
<input type="checkbox"/> 11260	ROHO AGILITY Max Contour Direct Mount - 20" Tall <i>Backrest can be installed at a height (seat pan to top of AGILITY) between 24" and 27". Specify Shell Width (14", 16", 18", 20"W): _____</i>	832.00	E2620

### MUST ALSO SELECT BACK HEIGHT:

- Set AGILITY Backrest at 22" \_\_\_\_\_
- Set AGILITY Backrest at 23" \_\_\_\_\_
- Set AGILITY Backrest at 24" \_\_\_\_\_
- Set AGILITY Backrest at 25" \_\_\_\_\_
- Set AGILITY Backrest at 26" \_\_\_\_\_
- Set AGILITY Backrest at 27" \_\_\_\_\_

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 109029-99-0	BodiLink Premium 3" x 5" Lateral Supports - Pair <i>Must also select I109031-99-0. Includes foam wedge for additional support. Only compatible with Corpus Ergo backrest shells.</i>	446.00	E0956
<input type="checkbox"/> 109030-99-0	BodiLink Premium 4" x 6" Lateral Supports - Pair <i>Must also select I109031-99-0. Includes foam wedge for additional support. Only compatible with Corpus Ergo backrest shells.</i>	446.00	E0956
<input type="checkbox"/> 109031-99-0	Swing-Away Adj BodiLink Lateral Hdwr Corpus Pair	625.00	E1028

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105782-99-0	Tubular Backrest Frame - Corpus (1" Canes) <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Canes are attached with 2" long UniTrack linkage disks that can rotate to provide space for backrest hardware. The tops of the canes can be adjusted up to a height of 27" above the seat pan. Depending on the aftermarket backrest, several inches of seat depth may be lost when mounting this system.</i>	530.00	
<input type="checkbox"/> 106072-99-0	Direct Backrest Frame - Corpus <i>Used to mount ROHO Agility, Jay 3, Matrx Elite, Varilite Icon and Ride Java backrests ranging in width from 14" to 20". Slotted brackets and 1" UniTrack linkage disks provide a close, adjustable mount. Backrest pitch can be tuned by rotating the linkage disks. Some seat depth may be lost due to the position and orientation of the backrest.</i>	530.00	



# M3 Corpus

## Backrest Options

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10725	Split Backrest Cushion - Corpus <i>Custom Option.</i>	359.00	
<input type="checkbox"/> 10786	Tailored Backrest Cushion - Corpus <i>Sides of cushion are trimmed to provide additional clearance for lateral trunk supports. Must also select a backrest cushion (110659 or 110660) when ordering this custom option.</i>	134.00	

## Seat Sizing

### CHOOSE SEAT WIDTH:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 107023-99-0	Set Seat Width at 17": _____	No Charge	
<input type="checkbox"/> 107024-99-0	Set Seat Width at 19": _____	No Charge	
<input type="checkbox"/> 107025-99-0	Set Seat Width at 21": _____	No Charge	
<input type="checkbox"/> 107026-99-0	Set Seat Width at 23": _____	No Charge	

### CHOOSE SEAT DEPTH:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>There are three different seat depth kits, each providing three seat depth settings: (1) 14" - 16", (2) 17" - 19" and (3) 20" - 22". Please select the required seat depth and the appropriate kit will be provided.</i>		
<input type="checkbox"/> 107013-99-0	Seat Depth 14"	No Charge	
<input type="checkbox"/> 107014-99-0	Seat Depth 15"	No Charge	
<input type="checkbox"/> 107015-99-0	Seat Depth 16"	No Charge	
<input type="checkbox"/> 107016-99-0	Seat Depth 17"	No Charge	
<input type="checkbox"/> 107017-99-0	Seat Depth 18"	No Charge	
<input type="checkbox"/> 107018-99-0	Seat Depth 19"	No Charge	
<input type="checkbox"/> 107019-99-0	Seat Depth 20"	No Charge	
<input type="checkbox"/> 107020-99-0	Seat Depth 21"	No Charge	
<input type="checkbox"/> 107021-99-0	Seat Depth 22"	No Charge	

### CHOOSE SEAT CUSHION:

Part Number	Description	Price	HCPCS Code
	NOTE: <i>Corpus Ergo seat cushions are not available in all seat sizes.</i>		
<input type="checkbox"/> 11337	Corpus Ergo Cushion - Leatherette Cover	601.00	E2605
<input type="checkbox"/> 11338	Corpus Ergo Cushion - Stretch-Air Cover	601.00	E2605
<input type="checkbox"/> 10965	Corpus Ergo Air Cushion (C3GRHOSCA) <i>Features ROHO Dry Floatation technology. The removable ROHO cell pad provides adjustable skin protection while the contoured foam and thigh support wedges help maintain proper positioning and stability. Includes removable stretch fabric cover. Available in widths of 17", 19" and 21" and depths of 17", 19" and 21", as well as 17" x 23". Not Available in 21" x 17"</i>	635.00	E2622
<input type="checkbox"/> 11113	Incontinence Cover for Ergo Cushion <i>Made from a stretch Dartex fabric and suitable for all Ergo cushions including the Corpus Ergo Air. The cover is cut to the contours of the Ergo cushion and is designed to be used as an outer layer of protection. Cover size will be matched to cushion size. Sent as a part.</i>	193.00	
<input type="checkbox"/> 11166	ROHO QUADTRO SELECT HIGH PROFILE < 22"W <i>Registered in PDAC under ROHO model numbers QS98C - QS1114C.</i>	506.00	E2624
<input type="checkbox"/> 11167	ROHO QUADTRO SELECT HIGH PROFILE > or = 22"W <i>Registered in PDAC under ROHO model numbers QS128C - QS1413C.</i>	563.00	E2625
<input type="checkbox"/> 11168	ROHO HIGH PROFILE Single Comp < 22"W <i>Registered in PDAC under ROHO model numbers 1R98C - 1R1114C.</i>	506.00	E2622
<input type="checkbox"/> 11169	ROHO HIGH PROFILE Single Comp > or = 22"W <i>Registered in PDAC under ROHO model numbers 1R138C - 1R1414C.</i>	563.00	E2623

## M3 Corpus

### Seat Sizing

#### CHOOSE SEAT CUSHION:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 11173	Incontinence Cover for ROHO Cushion <i>Made from a stretch Dartex fabric and suitable for all ROHO seat cushions. Cover size will be matched to the cushion size. Sent as a part.</i>	78.00	
<input type="checkbox"/> 105118-99-0	Omit Seat Cushion	No Charge	

### Armrest Options

#### CHOOSE ARMBAR KIT AND MOUNTING POSITION:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>Narrow Armbar Kit can be adjusted from 15"W to 19"W between the armrest pads without any additional parts. Wide Armbar Kit can be adjusted from 19"W to 23"W between the armrest pads without any additional parts. Please select both Armbar Kit and Width Setting below.</i>		
<input type="checkbox"/> 106543-99-0	Narrow Armbar Kit, Inf Adjustable 15"W - 19"W	No Charge	
<input type="checkbox"/> 106544-99-0	Wide Armrest Kit, Inf Adjustable 19"W - 23"W <i>Required when selecting a 20" wide backrest.</i> Set the Width Between the Armrests to: _____ <i>Must be wider than backrest width to ensure proper fit. Please round to nearest 1".</i>	No Charge	

#### CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105069-99-0	Height Adjustable Armrests - NC <i>No charge when ordering power tilt and/or recline.</i>	No Charge	
<input type="checkbox"/> 105081-99-0	Height Adjustable Armrests <i>Must be selected when the NPO base is selected.</i> SET ARMREST HEIGHT AT: _____ <i>Corpus Armrests are infinitely adjustable from 8" - 15", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.</i>	372.00	

#### CHOOSE ARMREST SWIVEL FIXED OR ENABLED:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105071-99-0	Left Armrest Bracket - Fixed <i>Medial swivel cannot be enabled when selecting this option. Must be selected when ordering arm trough adapters.</i>	No Charge	
<input type="checkbox"/> 105070-99-0	Right Armrest Bracket - Fixed <i>Medial swivel cannot be enabled when selecting this option. Must be selected when ordering arm trough adapters.</i>	No Charge	
<input type="checkbox"/> 105073-99-0	Left Armrest Bracket - Medial Swivel Enabled	No Charge	
<input type="checkbox"/> 105072-99-0	Right Armrest Bracket - Medial Swivel Enabled	No Charge	

#### CHOOSE ARMREST PAD RIGHT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> 105105-99-0	4"x16" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> 105999-99-0	4"x18" Corpus HD Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> 106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> 107214-99-0	4"x16" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> 105995-99-0	4"x18" Corpus HD Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> 107593-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Right	187.00	
<input type="checkbox"/> 107595-99-0	4"x16" Corpus Arm w/ Permobil Gel Pad - Right	192.00	
<input type="checkbox"/> 107597-99-0	4"x18" Corpus HD Arm w/ Permobil Gel Pad - Right	198.00	
<input type="checkbox"/> 108861-99-0	Moldable Forearm Support - Right <i>11.3"L x 5"W. Not compatible with armrest pouch, trays (except Small Tray) or the VS Chest Support.</i>	268.00	

# M3 Corpus

## Armrest Options

### CHOOSE ARMREST PAD RIGHT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 108863-99-0	Contoured Hand Pad Medium - Right <i>6.4"L x 5.5" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	133.00	
<input type="checkbox"/> 108865-99-0	Contoured Hand Pad Large - Right <i>7.7"L x 6.1" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	140.00	

### CHOOSE ARMREST PAD LEFT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> 105107-99-0	4"x16" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> 105998-99-0	4"x18" Corpus HD Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> 106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> 107215-99-0	4"x16" Corpus Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> 105994-99-0	4"x18" Corpus HD Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> 107592-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Left	187.00	
<input type="checkbox"/> 107594-99-0	4"x16" Corpus Arm w/ Permobil Gel Pad - Left	192.00	
<input type="checkbox"/> 107596-99-0	4"x18" Corpus HD Arm w/ Permobil Gel Pad - Left	198.00	
<input type="checkbox"/> 108860-99-0	Moldable Forearm Support - Left <i>11.3"L x 5"W. Not compatible with armrest pouch, trays (except Small Tray) or the VS Chest Support.</i>	268.00	
<input type="checkbox"/> 108862-99-0	Contoured Hand Pad Medium - Left <i>6.4"L x 5.5" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	133.00	
<input type="checkbox"/> 108864-99-0	Contoured Hand Pad Large - Left <i>7.7"L x 6.1" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	140.00	

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105463-99-0	Corpus Arm Trough Adapter 1" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked (I10645 &amp; I10646) when selecting this option.</i>	359.00	
<input type="checkbox"/> 105462-99-0	Corpus Arm Trough Adapter 7/8" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked (I10645 &amp; I10646) when selecting this option.</i> Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____	359.00	
<input type="checkbox"/> 107339-99-0	Armrest Pouch - Right	254.00	
<input type="checkbox"/> 107338-99-0	Armrest Pouch - Left	254.00	
<input type="checkbox"/> 107840-99-0	UT Armrest Side Support - Right Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	366.00	
<input type="checkbox"/> 107863-99-0	UT Armrest Side Support - Left Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	366.00	
<input type="checkbox"/> 107864-99-0	UT Armrest Side Support - Right Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i>	366.00	
<input type="checkbox"/> 107865-99-0	UT Armrest Side Support - Left Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i>	366.00	
<input type="checkbox"/> 107839-99-0	UT Elbow Support - Right <i>Must also select PN I10881.</i>	181.00	

# M3 Corpus

## Armrest Options

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 107838-99-0	UT Elbow Support - Left <i>Must also select PN I10882.</i>	181.00	
<input type="checkbox"/> 10881	Adjustable Elbow Support Hardware - Right	266.00	
<input type="checkbox"/> 10882	Adjustable Elbow Support Hardware - Left	266.00	

## Legrest Options

### CORPUS LEGREST OPTIONS:

Part Number	Description	Price	HCPCS Code
CHOOSE LEGREST TYPE:			
<input type="checkbox"/> 107763-99-0	Corpus Manual Elevating Legrest	No Charge	
<input type="checkbox"/> 107757-99-0	Corpus Power Elevating Legrest - 170° <i>Includes mechanical articulation.</i>	3,906.00	E1012
CHOOSE FOOTPLATE STYLE:			
<input type="checkbox"/> 11197	2pc Narrow Corpus Footplates - 9"D x 6"W	No Charge	
<input type="checkbox"/> 11196	2pc Wide Corpus Footplates - 9"D x 7.5"W	No Charge	
<input type="checkbox"/> 106760-99-0	1pc Narrow Corpus Footpl - 9.75"D x 12.5"W <i>Footplate is 12.5"W in front, 8"W in back.</i>	280.00	
CALF SUPPORT OPTIONS:			
<input type="checkbox"/> 105252-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad.</i>	No Charge	
<input type="checkbox"/> 105258-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	No Charge	
<input type="checkbox"/> 105250-99-0	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad.</i>	No Charge	
<input type="checkbox"/> 105256-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i>	399.00	
<input type="checkbox"/> 105254-99-0	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8"W x 8"T Pad.</i>	399.00	

### CORPUS VS LEGREST OPTIONS (FOR SELECT ACTIVE REACH PACKAGES):

Part Number	Description	Price	HCPCS Code
CORPUS VS LEGREST AND KNEE SUPPORT:			
<i>NOTE: The Corpus VS Legrest is only available with 20°, 30° and 45° of Active Reach.</i>			
<input type="checkbox"/> 107758-99-0	Corpus VS Power Elevating Legrest - 180° <i>Includes power articulation.</i>	3,906.00	E1012
<input type="checkbox"/> 107038-99-0	Knee Support/Abductor - Corpus VS <i>Select when ordering Corpus VS Legrest. Must also select I11129.</i>	278.00	E0957
<input type="checkbox"/> 11129	Adj Removable Knee Support Hardware (1 Each)	291.00	E1028
CHOOSE FOOTPLATE STYLE:			
<input type="checkbox"/> 11198	2pc Narrow Corpus VS Footplates - 9"D x 6"W	No Charge	
CALF SUPPORT OPTIONS:			
<input type="checkbox"/> 107211-99-0	VS Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad. Compatible with footplate heights of 10" - 18".</i>	No Charge	
<input type="checkbox"/> 106190-99-0	VS Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad. Compatible with footplate heights of 10" - 18".</i>	No Charge	
<input type="checkbox"/> 106188-99-0	VS Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad. Compatible with footplate heights of 11" - 18".</i>	393.00	

# M3 Corpus

## Legrest Options

### CORPUS VS LEGREST OPTIONS (FOR SELECT ACTIVE REACH PACKAGES):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 106186-99-0	VS Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8"W x 8"T Pad. Compatible with footplate heights of 12" - 18".</i>	393.00	

### FOOTPLATE HEIGHT POSITION:

Part Number	Description	Price	HCPCS Code
	SET FOOTPLATE HEIGHT AT: _____ <i>This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90°. DO NOT INCLUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments, please indicate desired footplate height. NOTE: Corpus footplate height range is 6"-17" and Corpus VS height range is 9"-17". Height may need to be reduced for shipping. Calf Pads are not available at heights less than 9".</i>		

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 106530-99-0	Stealth Gel Pads - Narrow Footplates (Pair)	447.00	
<input type="checkbox"/> 106225-99-0	Stealth Gel Pads - Wide Footplates (Pair)	458.00	

## Accessories

<input type="checkbox"/> 104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge	
<input type="checkbox"/> 104520-99-0	Bodypoint Padded Hip Belt - Medium	159.00	
<input type="checkbox"/> 104514-99-0	Bodypoint Non-Padded Hip Belt - Large	135.00	
<input type="checkbox"/> 104516-99-0	Bodypoint Padded Hip Belt - Large	212.00	
<input type="checkbox"/> 100519-99-0	Positioning Belt Retractable	316.00	
<input type="checkbox"/> 100600-99-0	Positioning Belt Retractable Long	328.00	
<input type="checkbox"/> 108267-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Med	213.00	
<input type="checkbox"/> 108266-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Large	224.00	

### HEADREST OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 107470-99-0	Permobil 7"W Headrest <i>Must also select headrest mounting hardware.</i>	413.00	E0955
<input type="checkbox"/> 107469-99-0	Permobil 10"W Headrest <i>Must also select headrest mounting hardware.</i>	413.00	E0955
<input type="checkbox"/> 107468-99-0	Permobil 14"W Headrest <i>Must also select headrest mounting hardware.</i>	413.00	E0955
<input type="checkbox"/> 107480-99-0	Adj/Removable Headrest Link Hardware - Standard	284.00	E1028
<input type="checkbox"/> 109178-99-0	Extra Link for Link Headrest Hardware <i>Used to extend range of standard link hardware. Sent as part.</i>	52.00	
<input type="checkbox"/> 108721-99-0	Adj/Removable BodiLink Head Support Hdwr - MAX 8 <i>Must be selected when ordering more than 2 swing-away headrest accessories or chin boom assemblies.</i>	347.00	E1028
<input type="checkbox"/> 109624-99-0	Extra 4" Link for BodiLink MAX 8 Hdwr <i>Used to extend range of MAX 8 hardware. Sent as part.</i>	76.00	
<input type="checkbox"/> 11305	Swing Away Egg Switch Assembly <i>Mount on Upper Left:_____ Mount on Upper Right:_____</i>	609.00	E1028
<input type="checkbox"/> 11306	Swing Away Microlight Assembly <i>Mount on Upper Left:_____ Mount on Upper Right:_____</i>	695.00	E1028
<input type="checkbox"/> 11307	Swing Away Facial Pad Assembly 4"L x 2"T <i>Mount on Upper Left:_____ Mount on Upper Right:_____</i>	494.00	E1028
<input type="checkbox"/> 11308	Swing Away Facial Pad Assembly 6"L x 2"T <i>Mount on Upper Left:_____ Mount on Upper Right:_____</i>	494.00	E1028
<input type="checkbox"/> 11276	Provide Headrest Receiver for AGILITY Backrest <i>Select for AGILITY backs and third party backs ordered off chair.</i>	No Charge	

# M3 Corpus

## Accessories

### HEADREST OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 106331-99-0	Permobil 10"W Headrest, Leatherette <i>Not available with ROHO AGILITY Back Systems. Must also select PN I10552.</i>	378.00	E0955
<input type="checkbox"/> 10552	Adj/Removable Headrest Hardware	291.00	E1028
<input type="checkbox"/> 107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	287.00	

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	172.00	E1020
<input type="checkbox"/> 103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	172.00	E1020
<input type="checkbox"/> 10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left: _____ Mount Amp Suppt on Right: _____</i>	291.00	
<input type="checkbox"/> 109040-99-0	BodiLink Premium LPTS Pad (B), Stretch-Air - Pair <i>4.5"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	341.00	E0953
<input type="checkbox"/> 109043-99-0	BodiLink Premium LPTS Pad (A), Stretch-Air - Pair <i>6"W x 5.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	341.00	E0953
<input type="checkbox"/> 109041-99-0	BodiLink Premium LPTS Pad (C), Stretch-Air - Pair <i>8"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	341.00	E0953
<input type="checkbox"/> 109042-99-0	BodiLink Premium LPTS Pad (G), Stretch-Air - Pair <i>12"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	341.00	E0953
<input type="checkbox"/> 109044-99-0	Adj/Removable BodiLink LPTS Hdwr Med Corpus Pair <i>Medium TT Power Mount hardware can be extended upward to clear cushions of 3-5.5" thick.</i>	656.00	E1028
<input type="checkbox"/> 109194-99-0	Adj/Removable Lateral Thigh/Hip Hardware - Pair	573.00	E1028
<input type="checkbox"/> 106762-99-0	Abductor Kit Corpus - Leatherette <i>Not compatible with Corpus VS Legrest. Must also select PN I10555.</i>	402.00	E0957
<input type="checkbox"/> 10555	Adj Removable - Abductor Hardware (1 Each)	291.00	E1028
<input type="checkbox"/> 105260-99-0	Transfer Handles, Corpus Seat	523.00	
<input type="checkbox"/> 105292-99-0	Push Handles, Corpus Seat	465.00	
<input type="checkbox"/> 106210-99-0	RAM Self-Leveling Cup Holder	88.00	
<input type="checkbox"/> 106806-99-0	RAM X-Grip Phone Holder	88.00	
<input type="checkbox"/> 107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount (I103476-99-0). Customization may be required if ordered with an Omni.</i>	444.00	
<input type="checkbox"/> 107189-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. No charge with 4 power seat function configurations. Mount on left armrest: _____ Mount on right armrest: _____</i>	220.00	
<input type="checkbox"/> 105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordered with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest. May ship as a part at shorter backrest heights.</i>	105.00	

# M3 Corpus

## Accessories

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 102059-99-0	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest.</i>	305.00	
<input type="checkbox"/> 105840-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as a part.</i>	194.00	E0960
<input type="checkbox"/> 105839-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as a part.</i>	194.00	E0960
<input type="checkbox"/> 105850-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as a part.</i>	194.00	E0960
<input type="checkbox"/> 105262-99-0	Joystick Cutout Tray - 16" W <i>Must also select PN I10556.</i>	417.00	E0950
<input type="checkbox"/> 105263-99-0	Joystick Cutout Tray - 18" W <i>Must also select PN I10556.</i>	417.00	E0950
<input type="checkbox"/> 10843	Small Angle Adjustable Tray <i>Must also select PN I10556.</i> <i>Mount on Left: _____ Mount on Right: _____</i>	321.00	E0950
<input type="checkbox"/> 10556	Adj Removable - Tray Hardware (1 Each)	291.00	
<input type="checkbox"/> 108230-99-0	Permobil Adjustable Tray - 13"W	270.00	E0950
<input type="checkbox"/> 108228-99-0	Permobil Adjustable Tray - 15"W	270.00	E0950
<input type="checkbox"/> 108227-99-0	Permobil Adjustable Tray - 17"W	270.00	E0950
<input type="checkbox"/> 108229-99-0	Permobil Adjustable Tray - 19"W	270.00	E0950
<input type="checkbox"/> 108232-99-0	Permobil Tray Hardware Right Side Adjustable	198.00	
<input type="checkbox"/> 108233-99-0	Permobil Tray Hardware Left Side Adjustable	198.00	
<input type="checkbox"/> 105224-99-0	Oxygen Holder, Corpus	286.00	E2208
<input type="checkbox"/> 105223-99-0	Crutch Holder, Corpus	286.00	E2207

## Alternative Drive Systems

### ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10751	Mini-Joystick for Omni - (9-pin connector) <i>Operating Force = 0.5 oz (13 g) Neutral to Full Deflection = 0.2" (5mm)</i>	3,447.00	E2312
<input type="checkbox"/> 103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 7.8 oz (220 g) Neutral to Full Deflection = 1.1" (28mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,686.00	E2373
<input type="checkbox"/> 103894-99-0	Compact Joystick Lite for Omni - (9-pin connector) <i>Operating Force = 1.8 oz (50 g) Neutral to Full Deflection = 0.7" (17mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,907.00	E2373

### ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
NOTE: <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options. The Swing Away Chin Boom Assemblies are Compatible with Permobil Headrests I107470-99-0, I107469-99-0 or I107468-99-0 and must be selected with I108721-99-0.</i>			
<input type="checkbox"/> 107946-99-0	Swing Away Chin Boom Assm for Compact Lite - Left	442.00	E1028
<input type="checkbox"/> 107948-99-0	Swing Away Chin Boom Assm for Compact Lite - Right	442.00	E1028
<input type="checkbox"/> 107857-99-0	Swing Away Chin Boom Assm for Mini Joystick - Left	442.00	E1028
<input type="checkbox"/> 107947-99-0	Swing Away Chin Boom Assm for Mini Joystick - Right	442.00	E1028

## M3 Corpus

### Alternative Drive Systems

#### ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 104690-99-0	Alt JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	295.00	E1028
<input type="checkbox"/> 107135-99-0	Alt JS / Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod.</i> <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	558.00	E1028
<input type="checkbox"/> 104143-99-0	Alt JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	920.00	E1028
<input type="checkbox"/> 10752	Chin Control Harness, Large <i>For neck width up to 7½".</i>	654.00	E1028
<input type="checkbox"/> 10753	Chin Control Harness, Small <i>For neck width up to 5½".</i>	654.00	E1028

#### ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches set-up for Mode &amp; On/Off function. Switch operation force = 4.5oz (130 g)</i> <i>Mount at front of Joystick: _____ Mount at rear of Joystick: _____</i> <i>Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	274.00	
<input type="checkbox"/> 104641-99-0	Dual Switches for Compact Joystick -Gooseneck Mount <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches set-up for Mode &amp; On/Off function. Switch operation force = 4.5oz (130 g)</i>	493.00	

#### SIP & PUFF SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Includes breath tube kit.</i>	1,844.00	E2325,E2326
<input type="checkbox"/> 10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Includes breath tube kit.</i>	2,454.00	E2325,E2326
<input type="checkbox"/> 10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573 or I10574.</i>	295.00	

#### PROPORTIONAL FOOT CONTROL:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10767	Proportional Foot Control <i>Not available with Active Reach Packages.</i> <i>Mount on Left: _____ Mount on Right: _____</i>	3,514.00	E2328
<input type="checkbox"/> 10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	295.00	

### Head Array Options

#### TOTAL CONTROL HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 109558-99-0	Permobil Total Control Head Array <i>Features an occipital pad sensor and two lateral arms configured with proximity sensors. Must select swing away hardware and mounting hardware below to complete assembly. Choose additional sensors/switches, including mode switch, from options below. ESP recommended. Choose Omni Control Module separately.</i>	4,589.00	E2330
<input type="checkbox"/> 11575	Swing Away Hardware for Side Sensor - Left	250.00	E1028
<input type="checkbox"/> 11576	Swing Away Hardware for Side Sensor - Right	250.00	E1028
<input type="checkbox"/> 11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes &amp; Comments section. Choose Omni Control Module separately.</i>	No Charge	



# M3 Corpus

## Head Array Options

### TOTAL CONTROL HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):			
<input type="checkbox"/> 109563-99-0	Adj/Removable BodiLink HA Hdwr - MAX 8	347.00	E1028
<input type="checkbox"/> 109564-99-0	Adj/Removable BodiLink HA Hdwr - MAX 8 Flip Down <i>Features flip-down joint that allows head array to swing back 120° away from user's head.</i>	477.00	E1028
SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> 11601	Add Left Lateral Arm w Proximity Sensor	655.00	E1028
<input type="checkbox"/> 11602	Add Left Lateral Arm w Egg Switch	536.00	E1028
<input type="checkbox"/> 11603	Add Left Lateral Arm w Pad (No Sensor or Switch)	437.00	E1028
SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> 11604	Add Right Lateral Arm w Proximity Sensor	655.00	E1028
<input type="checkbox"/> 11605	Add Right Lateral Arm w Egg Switch	536.00	E1028
<input type="checkbox"/> 11606	Add Right Lateral Arm w Pad (No Sensor or Switch)	437.00	E1028

### ASL FUSION HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 109582-99-0	ASL FUSION Proportional and Digital Head Array <i>Includes occipital (ASL 935) and lateral (ASL 933) pads, and pulse mechanical switch (ASL 320BK) for reset/mode changes. Must select hardware listed below to complete assembly.</i>	6,395.00	E2328
<input type="checkbox"/> 11580	Swing Away Hardware for Lateral Sensor - Left	325.00	E1028
<input type="checkbox"/> 11581	Swing Away Hardware for Lateral Sensor - Right	325.00	E1028
<input type="checkbox"/> 11582	Swing Away Hardware for Reset/Mode Switch - Right	325.00	E1028
<input type="checkbox"/> 11583	Adj/Removable Headrest Mounting Hardware <i>Includes ASL 606 and BodiLink receiver.</i>	325.00	E1028
<input type="checkbox"/> 109583-99-0	ASL FUSION Display, Prgmr and Att Control w/Mount <i>Includes ASL 165 and quick release mount.</i>	695.00	K0108

## Alt Drive Systems Setup

### CHOOSE OMNI (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 108500-99-0	Omni2 Control Module for R-net - Corpus <i>Two-section design features 3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology for control of up to 4 Bluetooth devices. Must also select mounting hardware below. Please specify which type of alternative drive input device you will be using. Input Device: _____</i>	1,962.00	
<input type="checkbox"/> 103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	295.00	

### CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10080	Mount Seat Control Switchbox on Backrest	No Charge	

## Alt Drive Accessories

### PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	297.00	
<input type="checkbox"/> 106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	522.00	
<input type="checkbox"/> 106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	748.00	

# M3 Corpus

## Alt Drive Accessories

### PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105907-99-0	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	378.00	E1028
<input type="checkbox"/> 105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	255.00	
<input type="checkbox"/> 105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	154.00	
<input type="checkbox"/> 106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	135.00	

### OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 101769-99-0	Microlight Switch Qty: _____	206.00	
<input type="checkbox"/> 101767-99-0	Buddy Button Qty: _____	174.00	
<input type="checkbox"/> 106293-99-0	Egg Switch Black Qty: _____	174.00	
<input type="checkbox"/> 101771-99-0	Mini Cup Switch Qty: _____	174.00	
<input type="checkbox"/> 107500-99-0	Proximity Switch Qty: _____	688.00	
<input type="checkbox"/> 107495-99-0	Mini Joystick 5 Switch OC for Omni - (9-pin)	747.00	
<input type="checkbox"/> 10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 6" rod, one 8" rod and one 9" gooseneck switch with 50" cable and mounting clamp.</i>	545.00	
<input type="checkbox"/> 10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.</i>	806.00	
<input type="checkbox"/> 105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	34.00	
<input type="checkbox"/> 105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	37.00	
<input type="checkbox"/> 105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	44.00	
<input type="checkbox"/> 105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	48.00	
<input type="checkbox"/> 105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	37.00	

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgement:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer. \*\* Indicates HCPCS code is pending PDAC approval.*