

permobil

TiLite Pilot

Ready for take off.





Optimal rear wheel access

Built-in forward growth

Ready to grow.

The TiLite Pilot has built-in forward growth adjustment to grow with your child and provide optimal rear wheel access from the very beginning. This gives kids the most efficient propulsion possible without needing to adjust or relearn their push mechanics. The TiLite Pilot lets kids cruise their world like never before.

Designed to move.

The clean frame design features a seat pan and side panels that provides easy access for adjustments, making changes simple and intuitive. Unlike on other chairs, all of these adjustments can be done with little or no compromise to rear wheel access.



Easy to use wheel locks are within reach

Frame designed for maximum stability



Lightweight for easier transport

Cleared for landing.

The ultra-lightweight design and carbon fibre frame panels make it 50% lighter and a breeze for kids to navigate all the adventures of their great big world. Quick release wheels and transit tie-downs mean that bringing the pilot everywhere is also simple for mum and dad.

Specifications


Frame features	Carbon frame with 1" lower aluminium tubing
Seat width	8" – 14" (20 - 35 cm)
Seat depth	9" – 15" (23 - 38 cm)
Front seat height	12" – 16.5" (30 - 42 cm)
Rear seat height	12" – 16" (30 - 41 cm)
Front frame angle	Unlimited
Seat back height	11" to 13" (8 - 33 cm) with standard back

Centre of gravity adjustment	0.5" – 3.5" (1 - 9 cm)
Transit options	Transit Tie-Down, Transport
Transport weight	4.9 kg [*]
Weight limit	43 kg transit, 54 kg standard
Seat to footrest	3" – 14.5" (8 - 37 cm) without high mount parts

[†]Specifications may vary based upon configuration.

^{*}Weight varies by configuration and does not include rear wheels.

Frame Colours



Acid Green Metallic

Neon Pink

Electric Plum

Ocean Blue Metallic

Anodised Colours



Sapphire

Black

Red Pepper

Ultraviolet