



Order form

# K300 PS Jr



Dealer Name\*:

Account #\*:

City, State, Zip\*:

Phone:

Email:

ATP:

Rep Name:

P.O. Number:

Quote Number:

Client Reference\*:

Client Weight\*:

*\*Required field*

Special considerations that need to be addressed here (e.g., diagnosis):

## Chair Model and Wheels

Front-wheel drive C300 power base includes: PS Jr seat frame, adjustable suspension, 4-point tie-down hardware, support wheels, VoltPro charger, and R-net 120 amp electronics. Pre-selected items are required for this configuration.

### Choose Base

103573-99-0	<b>K300 Base PS - C300 PS0</b> <i>FWD C300 power base set up for no power seat functions, max speed: 5 mph, weight capacity: 165 lb.</i>	\$8,616.00	K0848
10431	<b>K300 Base PS - C300 PS1</b> <i>FWD C300 power base set up for one power seat function, max speed: 5 mph, weight capacity: 165 lb.</i>	\$8,616.00	K0856
10520	<b>K300 Base PS - C300 PS2</b> <i>FWD C300 power base set up for two power seat functions, max speed: 5 mph, weight capacity: 165 lb.</i>	\$8,616.00	K0861

### Choose Drive Tires

103499-99-0	<b>Tri-Spoke Split Rim w/ Flat Free Tire 3x8" Gray</b>	No charge
103498-99-0	<b>Tri-Spoke Split Rim w/ Pneumatic Tire 3x8" Gray</b>	No charge
10782	<b>Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8"</b> <i>Not recommended for year-long use and must be ordered along with standard drive wheels. Does not include caster wheels.</i>	\$514.00

## Base Options

Choose Base Color | *Choose both shroud and rim color*

<b>Standard Shroud Colors</b>	Cobalt Blue	Purple Rain	Onyx Black	Phantom Gray	No charge
	Popstar Pink	Volcano Red			
<b>Standard Rim Colors</b>	Cobalt Blue	Purple Rain	Onyx Black	Phantom Gray	No charge
	Popstar Pink	Volcano Red	Platinum Sliver	Sublime Green	
	Extreme Yellow				

## Batteries

107032-99-0	<b>Batteries, Grp 34 Gel - Pair</b>	\$1,048.00	E2359
-------------	-------------------------------------	------------	-------

## Optional Items

ESP UNIT	<b>Enhanced Steering Perform Unit (ESP), R-net</b> <i>Required when ordering any ventilator tray. Recommended for switched input devices.</i>	\$1,659.00	
105803-99-0	<b>LED Light Kit for C/K300</b>	\$1,344.00	
102768-99-0	<b>Fixed Vent Shelf for C300</b> <i>Must also select ESP and expandable controller.</i>	\$523.00	E1029
109620-99-0	<b>Programming Key for R-net - Permobil Version</b>	\$575.00	

## Joystick Options

Choose joystick type, handle and mount. **NOTE:** For additional input devices, please refer to the Alternative Drive Systems section.

### Joystick Type

105246-99-0	<b>R-net LED Remote Joystick - No Display</b> <i>No display, no monojacks, limited to 5 profiles.</i>	No charge	
108319-99-0	<b>Permobil Joystick Module (PJSM)</b> <i>2.8" color display.</i>	\$1,182.00	
105700-99-0	<b>R-net Remote Color Joystick w/ Lights, Heavy Duty Gimbal</b> <i>Color LCD joystick, includes 3" knob with 3/4" diameter. Not compatible with alternative joystick handles.</i>	\$1,346.00	

### Joystick Handle

10045	<b>Cone Shaped Handle - Standard</b>	No charge	
100960-99-0	<b>Chin Cup Handle</b>	\$139.00	E2324
100961-99-0	<b>T Handle</b>	\$196.00	E2323
100962-99-0	<b>Large Ball</b>	\$139.00	E2323
100963-99-0	<b>Softball</b>	\$183.00	E2323
100964-99-0	<b>Mushroom Handle</b>	\$183.00	E2323
100965-99-0	<b>Stick Handle</b>	\$183.00	E2323
105442-99-0	<b>1.75" Black Foam Ball</b>	\$139.00	E2323
105846-99-0	<b>Bodypoint 3" U-Shaped Handle w/ Flex-Shaft</b>	\$183.00	E2323
105847-99-0	<b>Bodypoint 4" U-Shaped Handle w/ Flex-Shaft</b>	\$183.00	E2323
105848-99-0	<b>Bodypoint Rubber Dome Handle</b>	\$183.00	E2323

## Joystick Options

### Joystick Mount

<b>Fixed</b>	Mount on Left	Mount on Right	No charge
<b>Swing-Away</b>	Mount on Left	Mount on Right	\$409.00 E1028
<b>Slimline Retractable</b>	Mount on Left	Mount on Right	\$550.00 E1028
11232	<b>Adjustable Height Panel Bracket</b> <i>Provides 1.77" of vertical adjustment and 25° of angle adjustment.</i>		\$299.00 E1028
103472-99-0	<b>Mount Standard Joystick On Backrest, PS, PSJR, KO</b> <i>Can be chosen instead of attendant control.</i>		\$340.00
	Mount on Left	Mount on Right	
108268-99-0	<b>Bodypoint Tri-Lock Midline Mount Assembly</b> <i>Must select specific joystick attachment separately.</i>		\$702.00
	Mount on Left	Mount on Right	
108280-99-0	<b>Bodypoint Tri-Lock Mount for R-net Joystick</b>		\$292.00
108279-99-0	<b>Bodypoint Tri-Lock Mount for Compact Joystick</b>		\$292.00
10842	<b>PERMOFix Mount for mo-vis/Compact Joystick</b> <i>PermoFix mount for attaching Compact and mo-vis joysticks to Tri-Lock midline mount.</i>		\$292.00

### Other Options

10411	<b>Expandable Controller - R-net</b> <i>This item is required when ordering an Alternative Drive Control.</i>	\$764.00 E2377
10513	<b>Harness for Expandable Controller</b> <i>Must be ordered with any expandable controller.</i>	\$678.00 E2313
109779-99-0	<b>R-net Attendant Control</b>	\$1,499.00 E2331
103811-99-0	<b>Stroller Handle for PS Jr Seating System</b> <i>Angle adjustable.</i>	\$264.00
10582	<b>Output Module (IOM) for R-net Control System</b> <i>Enables user to control an external device (e.g., communication device, door opener, etc.) through the drive input (e.g., JSM, head array, etc.). May require specialty cable or wiring to be ordered separately.</i>	\$1,412.00

## Power Seating

### Adjustable Seat Height

105116-99-0	<b>Seat Tube Fixed</b>	No charge
106748-99-0	<b>Power Adjustable Seat Height - 8" Travel</b>	\$3,429.00 E2300

### Power Seat Functions

103712-99-0	<b>PS Jr Seat Frame - No Tilt</b>	No charge
103716-99-0	<b>45° Power Tilt</b>	\$7,297.00 E1002

### Seat Function Control Options | ICS switchbox availability may be limited to certain configurations.

10418	<b>Seat Function Control Switch for R-net</b> <i>Does not include operation of power seat functions through drive control.</i>	No charge
10419	<b>Single Seat Function Control Kit for R-net</b> <i>Allows operation of a single power seat function through drive control.</i>	\$2,110.00 E2310
10585	<b>Multiple Seat Function Control Kit for R-net</b> <i>Allows operation of a two or more power seat functions through drive control.</i>	\$3,994.00 E2311
104135-99-0	<b>Push Buttons on ICS Switchbox</b>	No charge
104136-99-0	<b>Toggle Switches on ICS Switchbox</b>	\$357.00

## Power Seating

### Switchbox Mounting

10626	<b>Mount ICS Switchbox onto R-net Joystick Bracket</b>	No charge
104651-99-0	<b>Mount ICS Switchbox onto Armrest</b> Mount on Left      Mount on Right	No charge
106256-99-0	<b>Mount ICS Switchbox on Back</b> <i>No charge when ordering Omni2 Control Module.</i> Mount on Left      Mount on Right	\$172.00
105409-99-0	<b>Do Not Install ICS Switchbox, Send as Part</b>	No charge

### Power Seating Options

104196-99-0	<b>ICS Alternative Switchbox for R-net Electronics</b> <i>Accommodates up to 8 momentary switches.</i>	\$1,396.00
106743-99-0	<b>Remote Stop System - All Electronics</b> <i>Consists of a 315MHz radio frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair.</i>	\$1,381.00 E2321

## Backrest Options

### PS Jr Backrest Options

103717-99-0	<b>PS Jr Back Canes, Height &amp; Angle Adjustable</b> <i>Must select PS Jr style armrests with this option. Backrest width should match seat width.</i>	No charge
<b>PS Jr Solid Backrest</b>	11" W    12" W    13" W    14" W    15" W    16" W	No charge
10449	<b>Omit Backrest</b>	No charge
104616-99-0	<b>Backrest Adjustment Quick Lock, PS / PS Jr</b> <i>Provides manual back angle adjustment.</i>	\$217.00

## Seat Sizing

### Seat Pan | Choose width and depth

<b>Seat Width</b>	11" W    12" W    13" W    14" W    15" W    16" W	No charge
<b>Seat Depth</b>	5" D    6" D    7" D    8" D    9" D    10" D    11" D    12" D 13" D    14" D    15" D    16" D    17" D    18" D	No charge

## Seat Cushion Options

### Mini-Flex Seat Cushions

<b>Leatherette</b>	10" W x 10" D    12" W x 12" D    14" W x 14" D    16" W x 16" D	\$653.00
<b>Black Mesh</b>	10" W x 10" D    12" W x 12" D    14" W x 14" D    16" W x 16" D	\$657.00

### Cushion Options

105118-99-0	<b>Omit Seat Cushion</b>	No charge E2607
-------------	--------------------------	-----------------

## Armrest Options

### Armrest Set-Up

10825	<b>Flip Back, Height Adjustable Armrests - NC</b> <i>No charge with power tilt. Set Armrest Height to:</i>										No charge
		6"	7"	8"	9"	10"	11"	12"	13"	14"	
10824	<b>Flip Back, Height Adjustable Armrests</b> <i>Set Armrest Height to:</i>										\$426.00 E0973
		6"	7"	8"	9"	10"	11"	12"	13"	14"	

### Choose Armrest Pads | *Must select left and right. Priced per side.*

CL Pads	<b>Corpus Leatherette</b>										No charge
		10" Left	13" Left	10" Right	13" Right						
CF Pads	<b>Corpus Black Fabric</b>										No charge
		10" Left	13" Left	10" Right	13" Right						
UT Pads	<b>UniTrack Armrest Pad</b>										No charge
		8" Left	12.5" Left	8" Right	12.5" Right						
CG Pads	<b>Corpus Gel</b>										\$192.00
		10" Left	13" Left	10" Right	13" Right						

### Optional Items | *Only available with Corpus style armrest pads. Priced per side/support.*

<b>Armrest Pouch</b>	Left	Right									\$267.00
<b>Armrest Supports</b>	Left Outer	Left Inner	Right Outer	Right Inner							\$385.00

## Legrest Options

### Legrest Type

<b>Swing-Away</b>	70° w/ Angle Adjustable Footplate	90° w/ Angle Adjustable Footplate				No charge
<b>Center Mount (Aeron)</b>	2pc Footplate	1pc Footplate				\$403.00

### Optional Items

104684-99-0	<b>UT Calf Support Kit (B), Mesh Cover &amp; Adj Hardware</b> <i>3.5"W x 4.5"T Pad. Cannot be used with swing-away legrests.</i>	No charge
104686-99-0	<b>UT Calf Support Kit (J), Mesh Cover &amp; Adj Hardware</b> <i>5"W x 6.5"T Pad. Cannot be used with swing-away legrests.</i>	No charge
107503-99-0	<b>BP Ankle Huggers, Med Side Release, Adj - Pair</b>	\$315.00

## Accessories

### Belt Options

104518-99-0	<b>Bodypoint Non-Padded Hip Belt - Medium</b>	No charge
105849-99-0	<b>Bodypoint Padded Hip Belt - Extra Small</b>	\$167.00
104520-99-0	<b>Bodypoint Padded Hip Belt - Medium</b>	\$167.00
100519-99-0	<b>Positioning Belt Retractable</b>	\$332.00
108267-99-0	<b>Bodypoint Padded Hip Belt "Airplane" Latch - Medium</b>	\$224.00

## Accessories

### Head Support Options

100749-99-0	<b>MF Headrest for Koala/PS Jr - Leatherette</b> <i>Must also select hardware.</i>	\$374.00	E0955
10552	<b>Adj/Removable Headrest Hardware</b>	\$306.00	E1028
107499-99-0	<b>Universal Headrest Adapter</b>	\$302.00	

### Leg Support Options

103077-99-0	<b>UT Hip Suppt (B), Mesh Cover - Low - Pair</b> <i>4.5"W x 3.5"T pad. Must also select hardware.</i>	\$315.00	E0953
103079-99-0	<b>UT Hip Suppt (A), Mesh Cover - Low - Pair</b> <i>6"W x 5.5"T pad. Must also select hardware.</i>	\$315.00	E0953
103095-99-0	<b>UT Thigh/Hip Suppt (C), Mesh Cover - High - Pair</b> <i>8"W x 3.5"T pad. Must also select hardware.</i>	\$315.00	E0953
10550	<b>Adj Removable - Thigh Support Hardware - Pair</b>	\$602.00	E1028

### Other Options

106210-99-0	<b>RAM Self-Leveling Cup Holder</b>	\$93.00	
106806-99-0	<b>RAM X-Grip Phone Holder</b>	\$93.00	
107100-99-0	<b>RAM X-Grip 10" Tablet Holder w/ Swing-Away Mount</b>	\$467.00	
107189-99-0	<b>USB Charger</b> <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets.</i> Mount on Left Armrest      Mount on Right Armrest	\$231.00	
105805-99-0	<b>Medical Necessities Bag Hooks</b>	\$111.00	E0950
102059-99-0	<b>Medical Necessities Bag</b>	\$321.00	
108276-99-0	<b>Bodypoint PivotFit Shoulder Harness - Medium</b>	\$289.00	
100578-99-0	<b>Small Angle Adjustable Tray</b> <i>Must also select hardware.</i>	\$339.00	
10556	<b>Adj Removable - Tray Hardware - 1 Each</b> Mount on Left      Mount on Right	\$306.00	
108231-99-0	<b>Permobil Adjustable Tray - 11"W</b> <i>Must also select hardware.</i>	\$284.00	E0960
108230-99-0	<b>Permobil Adjustable Tray - 13"W</b> <i>Must also select hardware.</i>	\$284.00	E0950
	<b>Permobil Adjustable Tray Hardware</b> Left Side Adjustable      Right Side Adjustable	\$208.00	
103798-99-0	<b>After Market Seating System Adapter</b> <i>For UniTrack seat rail.</i>	\$109.00	

## Alternative Drive Systems

*Omni2 Control Module required for all alternative drive systems. Enhanced Steering Performance (ESP) recommended for switched input devices.*

### Alt Proportional Joysticks

109525-99-0	<b>mo-vis Micro Joystick for Omni</b> <i>Operating Force = 8 g (0.28 oz). Neutral to Full Deflection = 0.13" (3.3 mm).</i>	\$3,200.00	E2312
109527-99-0	<b>mo-vis Multi Joystick for Omni</b> <i>Operating Force = 50 g (1.8 oz). Neutral to Full Deflection = 0.31" (8 mm).</i>	\$3,200.00	E2312
103894-99-0	<b>Compact Joystick Lite for Omni - (9-pin connector)</b> <i>Operating Force = 50 g (1.8 oz). Neutral to Full Deflection = 0.7" (17mm).</i>	\$2,003.00	E2373
103893-99-0	<b>Compact Joystick for Omni - (9-pin connector)</b> <i>Operating Force = 220 g (7.8 oz). Neutral to Full Deflection = 1.1" (28mm).</i>	\$1,771.00	E2373

## Alternative Drive Systems

**Alt Proportional Joystick Mounts** | Please see *Alt Drive Accessories* section for additional mounting options.

104690-99-0	<b>Alt Joystick Swing-Away Mount Kit - Armrest</b> Mount on Left      Mount on Right	\$310.00	E1028
107135-99-0	<b>Alt Joystick / Omni Combo Swing-Away Mount Kit - Armrest</b> Mount on Left      Mount on Right	\$586.00	E1028
104143-99-0	<b>Alt Joystick Adjustable Mount Kit - Chin Boom</b> <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	\$966.00	E1028
10752	<b>Chin Control Harness - Large</b> <i>For neck width up to 7½".</i>	\$687.00	E1028
10753	<b>Chin Control Harness - Small</b> <i>For neck width up to 5½".</i>	\$687.00	E1028

### Alt Proportional Joystick Accessories

104904-99-0	<b>Dual Switches for Compact Joystick - Fixed Mount</b> Mount Left of JS      Mount Right of JS      Mount front of JS      Mount back of JS	\$288.00	
104641-99-0	<b>Dual Switches for Compact Joystick - Gooseneck Mount</b>	\$518.00	

### Total Control Head Array

109558-99-0	<b>Permobil Total Control Head Array</b> <i>Features an occipital pad sensor and two lateral arms configured with proximity sensors. Must select swing away hardware and mounting hardware below to complete assembly.</i>	\$4,819.00	E2330
11575	<b>Swing-Away Hardware for Side Sensor - Left</b>	\$263.00	E1028
11576	<b>Swing-Away Hardware for Side Sensor - Right</b>	\$263.00	E1028
11072	<b>Program Head Array as a 5-Switch System</b> <i>Select additional arms/sensors below and indicate the intended use of each sensor in the notes section.</i>	No charge	
109563-99-0	<b>Adj/Removable BodiLink HA Hardware - MAX 8</b>	\$365.00	E1028
109564-99-0	<b>Adj/Removable BodiLink HA Hardware - MAX 8 Flip Down</b> <i>Flip-down joint that allows head array to swing back 120°.</i>	\$501.00	E1028
11601	<b>(Optional) Add Left Lateral Arm w/ Proximity Sensor</b>	\$688.00	E1028
11602	<b>(Optional) Add Left Lateral Arm w/ Egg Switch</b>	\$563.00	E1028
11603	<b>(Optional) Add Left Lateral Arm w/ Pad (No Sensor or Switch)</b>	\$459.00	E1028
11604	<b>(Optional) Add Right Lateral Arm w/ Proximity Sensor</b>	\$688.00	E1028
11605	<b>(Optional) Add Right Lateral Arm w/ Egg Switch</b>	\$563.00	E1028
11606	<b>(Optional) Add Right Lateral Arm w/ Pad (No Sensor or Switch)</b>	\$459.00	E1028

### ASL FUSION Head Array

109582-99-0	<b>ASL FUSION Proportional and Digital Head Array</b> <i>Includes occipital (ASL 935) and lateral (ASL 933) pads, and pulse mechanical switch (ASL 320BK) for reset/mode changes. Must select hardware listed below to complete assembly.</i>	\$6,715.00	E2328
11580	<b>Swing-Away Hardware for Lateral Sensor - Left</b>	\$342.00	E1028
11581	<b>Swing-Away Hardware for Lateral Sensor - Right</b>	\$342.00	E1028
11582	<b>Swing-Away Hardware for Reset/Mode Switch - Right</b>	\$342.00	E1028
11583	<b>Adj/Removable Headrest Mounting Hardware</b> <i>Includes ASL 606 and BodiLink receiver.</i>	\$342.00	E1028
109583-99-0	<b>ASL FUSION Display, Programmer and Attendant Control w/ Mount</b> <i>Includes ASL 165 and quick release mount.</i>	\$730.00	K0108

## Alternative Drive Systems

### Sip & Puff Systems

10573	<b>Sip &amp; Puff System for Omni</b>	\$1,937.00	E2325, E2326
10574	<b>Sip &amp; Puff System with Fiber Optic Opti-Stop</b>	\$2,577.00	E2325, E2326
10762	<b>Adjustable Mounting Hardware - Sip &amp; Puff</b> <i>Required with any sip &amp; puff system.</i>	\$310.00	

## Alt Drive System Setup

### Omni2 Options | *Required for all alternative drive systems.*

108523-99-0	<b>Omni2 Control Module for R-net - PS Jr</b> <i>3.5" VGA high resolution display, USB charger port, built-in IR control and integrated Bluetooth technology. Must select mounting hardware below. Specify alternative input device being used:</i>	\$2,061.00	
103476-99-0	<b>Adjustable Swing-Away Hardware - Omni/Tablet</b> Mount on Left      Mount on Right	\$310.00	

## Alt Drive Accessories

### PERMOfix Mounting Solutions

106655-99-0	<b>PERMOfix UniTrack Fixed Mount Kit</b>	\$312.00	
106656-99-0	<b>PERMOfix UniTrack Removable Mount Kit</b>	\$549.00	
106657-99-0	<b>PERMOfix UniTrack Swing-Away Mount Kit</b>	\$786.00	
105907-99-0	<b>Proximity Sensor + Swivel Assembly</b> <i>Swivel assembly allows sensor to be mounted with PERMOfix. Since Proximity Sensors require power, CANNOT be used as On/Off switches.</i>	\$397.00	E1028
105919-99-0	<b>Egg Switch + Swivel Mount Assembly</b> <i>Swivel assembly allows switch to be mounted with PERMOfix.</i>	\$268.00	
105912-99-0	<b>Solid Oval Pad + Swivel Assembly</b> <i>Swivel assembly allows switch to be mounted with PERMOfix.</i>	\$162.00	
106675-99-0	<b>Oval Plate + Swivel Assembly - No Switch</b> <i>Blank plate provides mounting option for Microlight, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix.</i>	\$142.00	

### Switches and Cables

101769-99-0	<b>Microlight Switch   Qty:</b>	\$217.00	
101767-99-0	<b>Buddy Button   Qty:</b>	\$183.00	
106293-99-0	<b>Egg Switch Black   Qty:</b>	\$183.00	
101771-99-0	<b>Mini Cup Switch   Qty:</b>	\$183.00	
107500-99-0	<b>Proximity Switch   Qty:</b>	\$723.00	
105935-99-0	<b>3.5mm Ext. Cable for Sensor/Switch, 0.3m Long</b>	\$36.00	
105931-99-0	<b>3.5mm Ext. Cable for Sensor/Switch, 1m Long</b>	\$39.00	
105932-99-0	<b>3.5mm Ext. Cable for Sensor/Switch, 2m Long</b>	\$47.00	
105933-99-0	<b>3.5mm Ext. Cable for Sensor/Switch, 3m Long</b>	\$51.00	
105934-99-0	<b>3.5mm Y-splitter Cable for Sensor/Switch</b> <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	\$39.00	

## Additional Order Instructions (for Permobil):

## Notes & Comments (not for Permobil):

## Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- The client's name or code and the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

## Order Acknowledgement:

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

*NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*