



permobil

# TRA

**Roll better,  
live better.**

With the power and style of its titanium dual-tube frame, TiFit precision, and full adjustability, the TRA just rolls better. The result is better performance and comfort built around you.



## TiFit. Personal Ergonomics.

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. A TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do. TiFit...not standard.

## Custom Configuration.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimise your riding experience. This includes choices of wheels, handrims, tyres, wheel locks, colour, tattoos, and more. TiLite offers more options for the individual.

## Technical Specifications

Frame Features	Titanium Dual-Tube Rigid TiFit Limited Adjustability	Camber	0°, 2°, 4°, 6° or 8°
Total Chair Width*	23"	Centre of Gravity Adjustment	5"
Seat Width	12" – 20"	Front Seat Width (Seat Taper)	0", 1", 2", 3"
Seat Depth	12" – 20"	Transit Options	Transport Bracket
Front Seat Height	15" – 21"	Transport Weight*	5 kg.
Rear Seat Height	13" – 21"	Weight Limit	120.2 kg Standard 158.7 kg Heavy Duty
Front Frame Angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five Years on Frame One Year on TiLite Components
Seat Back Type	Folding Seat Back Adj. Height 8.5" – 21" Adj. Angle 80° – 101°	* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.27 kg.	


## Colour Options

### High Gloss Paint



TiLite Red    Burgundy    Neon Pink    Super Wet Black

### Matte Paint



Deep Blue Ice    Titanium Grey    Black Lava

### Candy Paint



Ruby Red    Azure Blue    Electric Plum

### Natural



Satin Titanium

### Metallic Paint



Salsa Red Metallic    Tangerine Metallic    Acid Green Metallic    Ocean Blue Metallic

### Anodised Colours



Red Pepper    Burnt Orange    Sapphire    Ultraviolet    Fuchsia    Black

### Pearlescent Paint



Midnight Blue Pearl    White Pearl