

M300 Corpus HD

permobil



*** Indicates a Required Field**

Dealer Information

Order ID: _____

***ATP:** _____

Dealer Code: _____

***Dealer Name:** _____

Address: _____

***City:** _____

*** State/Zip:** _____

*** Phone#:** _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

*** First Name:** _____

*** Last Name:** _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

*** Weight:** _____

***Height:** _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____


E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

***I) Knee to Foot:** _____



Three column pricing on order/quote

Display HCPCS Codes on order/quote

** Indicates HCPCS code is pending PDAC approval

Please send order/quote to fax#: **(800) 231-3256**

Email to: sales@permobil.com

Permobil Inc.

300 Duke Dr.

Lebanon, TN 37090

Tel: (800) 736-0925

Fax: (800) 231-3256

www.permobil.com

Prices effective **May 16, 2022**

M300 Corpus HD

Chair Model and Colors

CHOOSE BASE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10908	M300 Base Corpus HD - SPO <i>Price Includes: M300 Power HD Base Set Up For One Power Seat Function, Corpus HD Seat Frame, Shock Absorbing Suspension System, 8" Flat-Free Casters, Tie-Down Hardware for Strap Systems, Charger, Flat Free or Pneumatic Tires and a Standard Positioning Belt. R-net 120amp Electronics. Max Speed = 5.0mph, Weight Capacity = 450lbs.</i>	10,899.00	K0858
<input type="checkbox"/> 105799-99-0	M300 Base Corpus HD - MPO <i>Price Includes: M300 HD Power Base Set Up For Multiple Power Seat Functions, Corpus HD Seat Frame, Shock Absorbing Suspension System, 8" Flat-Free Casters, Tie-Down Hardware for Strap Systems, Charger, Flat Free or Pneumatic Tires and a Standard Positioning Belt. R-net 120amp Electronics. Max Speed = 5.0mph, Weight Capacity = 450lbs.</i>	10,899.00	K0862

CHOOSE WHEEL PACKAGE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 11081	Gray Wheel Package - MWD <i>Includes grey tires on 8" silver casters and drive wheel rims. Hub cap color is selected below. Please select drive wheel type: Flat free drive tires: _____ Pneumatic drive tires: _____</i>	No Charge	
<input type="checkbox"/> 11082	Black Wheel Package - MWD <i>Includes black tires on 8" silver casters and drive wheel rims. Hub cap color is selected below. Please select drive wheel type: Flat free drive tires: _____ Pneumatic drive tires: _____</i>	No Charge	
<input type="checkbox"/> 11083	Blackout Wheel Package - MWD <i>Custom Option. Includes black tires and custom matte black 8" casters, drive wheels and hub caps. Please select drive wheel type: Flat free drive tires: _____ Pneumatic drive tires: _____</i>	353.00	

OPTIONAL DRIVE TIRES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10782	Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8" <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Not available with black tires or black rims.</i>	489.00	

CHOOSE SHROUD COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105885-99-0	M300 HD Shroud - Onyx Black	No Charge	
<input type="checkbox"/> 105881-99-0	M300 HD Shroud - Platinum Silver	No Charge	
<input type="checkbox"/> 105887-99-0	M300 HD Shroud - Sublime Green	No Charge	
<input type="checkbox"/> 106169-99-0	M300 HD Shroud - Cobalt Blue	No Charge	
<input type="checkbox"/> 106171-99-0	M300 HD Shroud - Phantom Gray	No Charge	
<input type="checkbox"/> 106159-99-0	M300 HD Shroud - Pop Star Pink	No Charge	
<input type="checkbox"/> 106170-99-0	M300 HD Shroud - Purple Rain	No Charge	
<input type="checkbox"/> 106160-99-0	M300 HD Shroud - Volcano Red	No Charge	

CHOOSE HUBCAP COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10871	Match Hub Cap Color to Shroud: _____ <i>Note: Hub caps will be matched to shroud color if no selection is made.</i> Solid Color Hub Caps - Pair <i>Platinum Silver: _____ Onyx Black: _____ Cobalt Blue: _____ Phantom Grey: _____ Pop Star Pink: _____ Volcano Red: _____ Purple Rain: _____ Titan Blue: _____ Sublime Green: _____ Crystal Blue: _____ Vintage Wine: _____ Extreme Yellow: _____ Chili Pepper Red: _____</i>	No Charge	

Base Options

M300 Corpus HD

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105875-99-0	HD Power Adjustable Seat Height - 8" Travel <i>Please Note: Drive is restricted when seat is in an elevated position.</i>	3,605.00	E2300
<input type="checkbox"/> 105880-99-0	HD Seat Tube Fixed	No Charge	

MUST CHOOSE BATTERIES:

Part Number	Description	Price	HCPCS Code
BATTERY OPTIONS:			
<input type="checkbox"/> 108999-99-0	Permobil Batteries, Grp 24 AGM GEL, Pair <i>Includes Permobil VoltPro charger optimized for Permobil batteries.</i>	1,042.00	E2363
<input type="checkbox"/> 107032-99-0	Batteries, Grp 34 Gel, Pair <i>Includes Permobil VoltPro charger for non-Permobil gel batteries.</i>	998.00	E2359

CHOOSE SEAT PAN TO FLOOR HEIGHT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105878-99-0	Set Seat to Floor Height at 18"	No Charge	
<input type="checkbox"/> 105877-99-0	Set Seat to Floor Height at 19"	No Charge	
<input type="checkbox"/> 105876-99-0	Set Seat to Floor Height at 20"	No Charge	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105428-99-0	Enhanced Steering Perform Unit (ESP), R-net <i>Recommended item with switched drive inputs for improved tracking.</i>	1,580.00	K0108
<input type="checkbox"/> 105297-99-0	LED Light Kit for M300 / M400 w/ R-net Electronics	1,282.00	
<input type="checkbox"/> 10715	Compact Laptop Ventilator Tray - Corpus Seat <i>Mounts to rear of Corpus backrest hinge. Accommodates Pulmonetic Systems / CareFusion LTV Series and similar sized Laptop Ventilators. The tray is 13"T, adjustable in width from 8" - 13", and adjustable in depth from 2" - 6.5". Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option. Additional drive restrictions may apply.</i>	1,671.00	E1029
<input type="checkbox"/> 10787	Battery Tray for Laptop Ventilator Tray <i>Measures 11 3/8" W x 6 1/2" H x 1 5/8" D. Accommodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	280.00	
<input type="checkbox"/> 10788	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 7" Tall. Width infinitely adjustable 3.5" - 7". Depth infinitely adjustable 3.75" - 5". Custom Option.</i>	355.00	E2208
<input type="checkbox"/> 109620-99-0	Programming Key for R-net - Permobil Version	547.00	

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
PLEASE NOTE: <i>For additional input devices, please refer to the Alternative Drive Systems Section.</i>			
<input type="checkbox"/> 107779-99-0	Color Joystick Module (CJSM) <i>This R-net color joystick features: A Large High-Resolution 3.5" VGA color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	1,232.00	
<input type="checkbox"/> 11186	Permobil Joystick Module, NC <i>Designed exclusively for Permobil, this R-net color joystick features: a 2.8" color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	No Charge	
<input type="checkbox"/> 107829-99-0	3.5mm Y-Splitter Cable for PJSM/CJSM <i>For use with the PJSM or CJSM Mode/Profile jack. Allows 2 mechanical switches to control 4 separate functions. Requires PC Programming to assign the short hold and long hold commands to each switch. Switches sold separately.</i>	37.00	

M300 Corpus HD

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105246-99-0	R-net LED Remote Joystick - No Display <i>This Joystick does NOT have a Display Screen, does NOT contain mono jack ports, and does NOT include Bluetooth or Assignable Buttons. PC Programming required and profiles limited to five.</i>	No Charge	
<input type="checkbox"/> 105700-99-0	R-net Remote Color JS w/ Lights, Heavy Duty Gimbal <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Color LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	1,281.00	

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10045	Cone Shaped Joystick Knob - Standard	No Charge	
<input type="checkbox"/> 100960-99-0	Chin Cup For Joystick Knob	132.00	E2324
<input type="checkbox"/> 100961-99-0	"T" Handle For Joystick Knob	186.00	E2323
<input type="checkbox"/> 100962-99-0	"Large Ball" For Joystick Knob	132.00	E2323
<input type="checkbox"/> 100963-99-0	"Softball" For Joystick Knob	174.00	E2323
<input type="checkbox"/> 100964-99-0	"Mushroom" For Joystick Knob	174.00	E2323
<input type="checkbox"/> 100965-99-0	"Stick" For Joystick Knob	174.00	E2323
<input type="checkbox"/> 105442-99-0	1.75" Black Foam Ball For Joystick Knob	132.00	E2323
<input type="checkbox"/> 105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	174.00	E2323
<input type="checkbox"/> 105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	174.00	E2323
<input type="checkbox"/> 105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	174.00	E2323

CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10664	Joystick Mount Right - Fixed	No Charge	
<input type="checkbox"/> 10665	Joystick Mount Left - Fixed	No Charge	
<input type="checkbox"/> 102537-99-0	Joystick Mount Right - Swing Away	389.00	E1028
<input type="checkbox"/> 102538-99-0	Joystick Mount Left - Swing Away	389.00	E1028
<input type="checkbox"/> 108037-99-0	Slimline Retractable Joystick Mount - Right	523.00	E1028
<input type="checkbox"/> 108038-99-0	Slimline Retractable Joystick Mount - Left	523.00	E1028
<input type="checkbox"/> 11232	Adjustable Height Panel Bracket <i>Provides 1.77" of vertical adjustment below the armrest and 25° of angle adjustment. Compatible with all Permobil joystick mounts except for the M1 Fixed Joystick Mount.</i>	284.00	E1028
<input type="checkbox"/> 105210-99-0	Mount Standard Joystick on Corpus Back Hinge <i>Standard R-net Joystick is mounted to backrest hinge. Can be chosen instead of PN 1105408-99-0. Mount on Left: _____ Mount on Right: _____</i>	358.00	
<input type="checkbox"/> 108268-99-0	Bodypoint Tri-Lock Midline Mount Assembly <i>Includes L-shaped joystick arm and height adjustable mounting clamp. Joystick mount sold separately (select from the mounts listed below). Please indicate if you want midline mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	668.00	
<input type="checkbox"/> 108280-99-0	Bodypoint Tri-Lock Mount for R-net JS	278.00	
<input type="checkbox"/> 108279-99-0	Bodypoint Tri-Lock Mount for Compact Joystick	278.00	
<input type="checkbox"/> 10842	PERMOFix Mount for Mini/Compact Joystick <i>Includes a shim washer that allows the PERMOFix joint to be mounted on the joystick arm. Provides an adjustable mount for a Mini Joystick or Compact Joystick.</i>	278.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10411	Expandable Controller - R-net <i>This item is required when ordering power tilt, power recline, and power elevating legrests together OR when ordering an Alternative Drive Control.</i>	727.00	E2377

M300 Corpus HD

Joystick Options

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller.</i>	645.00	E2313
<input type="checkbox"/> 105225-99-0	Co-Pilot Attendant Control R-net - Corpus	5,320.00	E2331
<input type="checkbox"/> 105408-99-0	Attendant Control R-net	1,427.00	E2331
<input type="checkbox"/> 10582	Output Module (IOM) for R-net Control System <i>Enables a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	1,344.00	

Seating System

CHOOSE SEATING SYSTEM:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>Weight Limit 450lbs. For the Corpus HD Seat on the M300, there is a half speed drive limit when the Total Backrest Angle (TBA) > 135° and a drive inhibit at 150° of TBA. There is also a drive inhibit at 20° of tilt.</i>		
<input type="checkbox"/> 10911	HD Seat w/ No Tilt & 150° Power Recline	9,092.00	E1004
<input type="checkbox"/> 10912	Corpus HD Seat w/ 45° DuraCore Tilt <i>Backrest angle can be manually adjusted up to 120°</i>	9,262.00	E1002
<input type="checkbox"/> 10913	HD Seat w/ 45° DuraCore Tilt & 150° Power Recline	14,014.00	E1007

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10418	Seat Function Control Switch for R-net <i>Does not include power seat function control through the drive control.</i>	No Charge	
<input type="checkbox"/> 10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control. Order with R-net electronics or remote joystick options.</i>	2,009.00	E2310
<input type="checkbox"/> 10585	Multiple Seat Function Control Kit for R-net <i>For operation of two or more seat functions through the wheelchair electronics via the drive control. Order with R-net electronics or remote joystick options.</i>	3,803.00	E2311

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>ICS Switchbox is only available with chairs that include an R-Net LED Remote Joystick.</i>		
<input type="checkbox"/> 102798-99-0	Push Buttons on ICS Switchbox (1.4m)	No Charge	
<input type="checkbox"/> 102799-99-0	Toggle Switches on ICS Switchbox (1.4m)	340.00	
<input type="checkbox"/> 105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge	

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> 104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge	
<input type="checkbox"/> 106256-99-0	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	163.00	

M300 Corpus HD

Seating System Options

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,329.00	
<input type="checkbox"/> 106743-99-0	Remote Stop System - All Electronics <i>Consists of a 315MHz Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.</i>	1,315.00	E2321
<input type="checkbox"/> 10750	Memory Seat Program, Corpus w/ R-net	137.00	
<input type="checkbox"/> 106209-99-0	Legrest Movement Linked to Recline Electronically <i>Power Legrest will raise when Power Recline moves backward and will lower when back moves forward. If Power Legrest is operated, Power Recline will NOT automatically reposition. Feature works through ICS Switchbox and Driver Control.</i>	52.00	
<input type="checkbox"/> 10898	Enable Independent Repositioning Mode <i>IRM allows for sequential movement of power seat functions into a posterior position (35° tilt, 115° recline, 135° legrest) and returns to upright position (5° tilt, 100° recline, 110° legrest). This can be accessed using a standard drive control or a single switch. On the Permobil Joystick Module, IRM is accessible through the mode/profile jack (long switch hold). Please see the Alt Drive Accessories section for switch options.</i>	No Charge	

Backrest Options

CHOOSE FIXED OR SLIDING BACKREST KIT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105871-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No Charge	
<input type="checkbox"/> 105869-99-0	Fixed Backrest Kit - No Mechanical Shear	35.00	

CHOOSE BACKREST TYPE:

Part Number	Description	Price	HCPCS Code
BACKREST HEIGHT NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°.</i>			
<input type="checkbox"/> 106016-99-0	Corpus HD Ergo Backrest - 20"W x 21"T (C3GERGBCA)	1,114.00	E2620
<input type="checkbox"/> 106017-99-0	Corpus HD Ergo Backrest - 22"W x 21"T	1,089.00	E2621
<input type="checkbox"/> 106018-99-0	Corpus HD Ergo Backrest - 24"W x 21"T	1,089.00	E2621
<input type="checkbox"/> 10915	Solid Back 20"W X 21"T - Corpus HD Seat	No Charge	
<input type="checkbox"/> 10916	Solid Back 22"W x 21"T - Corpus HD Seat	No Charge	
<input type="checkbox"/> 10917	Solid Back 24"W x 21"T - Corpus HD Seat	No Charge	

PLEASE NOTE:

Taller backrests listed below include: foam lumbar pads, a UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided. * Foam lumbar pads not included with Solid Backs.*

<input type="checkbox"/> 106019-99-0	Corpus HD Ergo Back -20"W x 24 - 29"T (C3GERGBCA)	1,114.00	E2620
<input type="checkbox"/> 106020-99-0	Corpus HD Ergo Back -22"W x 24 - 29"T	1,089.00	E2621
<input type="checkbox"/> 106021-99-0	Corpus HD Ergo Back -24"W x 24 - 29"T	1,089.00	E2621
<input type="checkbox"/> 10919	Solid Back 20"W x 24 - 29"T - Corpus HD Seat	No Charge	
<input type="checkbox"/> 10920	Solid Back 22"W x 24 - 29"T - Corpus HD Seat	No Charge	
<input type="checkbox"/> 10921	Solid Back 24"W x 24 - 29"T - Corpus HD Seat	No Charge	
	Set Backrest at 24" _____		
	Set Backrest at 25" _____		
	Set Backrest at 26" _____		
	Set Backrest at 27" _____		
	Set Backrest at 28" _____		

M300 Corpus HD

Backrest Options

CHOOSE BACKREST TYPE:

Part Number	Description	Price	HCPCS Code
	Set Backrest at 29" _____		

CHOOSE BACKREST CUSHION:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10659	Leatherette Cushion - Corpus Ergo Back	No Charge	
<input type="checkbox"/> 10660	Black Mesh Cushion - Corpus Ergo Back	45.00	
<input type="checkbox"/> 11340	Foam Lumbar Pads for Solid Back - Corpus HD	42.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10725	Split Backrest Cushion - Corpus <i>Custom Option.</i>	359.00	

Seat Sizing

CHOOSE SEAT WIDTH SETTING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105859-99-0	Set Seat Width at 19"	No Charge	
<input type="checkbox"/> 105860-99-0	Set Seat Width at 21"	No Charge	
<input type="checkbox"/> 105861-99-0	Set Seat Width at 23"	No Charge	
<input type="checkbox"/> 105862-99-0	Set Seat Width at 25"	No Charge	

CHOOSE SEAT DEPTH:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>Odd numbered seat depths can be made 1" shorter without any additional parts. For example, a 21"D Seat can also be set at 20"D.</i>		
<input type="checkbox"/> 105863-99-0	Corpus HD Seat Depth 19"	No Charge	
<input type="checkbox"/> 105864-99-0	Corpus HD Seat Depth 20"	No Charge	
<input type="checkbox"/> 105865-99-0	Corpus HD Seat Depth 21"	No Charge	
<input type="checkbox"/> 105866-99-0	Corpus HD Seat Depth 22"	No Charge	
<input type="checkbox"/> 105868-99-0	Corpus HD Seat Depth 23"	No Charge	
<input type="checkbox"/> 105856-99-0	Corpus HD Seat Depth 24"	No Charge	
<input type="checkbox"/> 105857-99-0	Corpus HD Seat Depth 25"	No Charge	

CHOOSE SEAT CUSHION:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10651	Corpus Ergo Cushion - Leatherette (C3GERGSCA)	601.00	E2605
<input type="checkbox"/> 10652	Corpus Ergo Cushion - Black Fabric	601.00	E2605
<input type="checkbox"/> 11113	Incontinence Cover for Ergo Cushion <i>Made from a stretch Dartex fabric and suitable for all Ergo cushions including the Corpus Ergo Air. The cover is cut to the contours of the Ergo cushion and is designed to be used as an outer layer of protection. Cover size will be matched to cushion size. Sent as a part.</i>	193.00	
<input type="checkbox"/> 11166	ROHO QUADTRO SELECT HIGH PROFILE < 22"W <i>Registered in PDAC under ROHO model numbers QS98C - QS1114C.</i>	506.00	E2624
<input type="checkbox"/> 11167	ROHO QUADTRO SELECT HIGH PROFILE > or = 22"W <i>Registered in PDAC under ROHO model numbers QS128C - QS1413C.</i>	563.00	E2625
<input type="checkbox"/> 11168	ROHO HIGH PROFILE Single Comp < 22"W <i>Registered in PDAC under ROHO model numbers 1R98C - 1R1114C.</i>	506.00	E2622
<input type="checkbox"/> 11169	ROHO HIGH PROFILE Single Comp > or = 22"W <i>Registered in PDAC under ROHO model numbers 1R138C - 1R1414C.</i>	563.00	E2623
<input type="checkbox"/> 11173	Incontinence Cover for ROHO Cushion <i>Made from a stretch Dartex fabric and suitable for all ROHO seat cushions. Cover size will be matched to the cushion size. Sent as a part.</i>	78.00	
<input type="checkbox"/> 105118-99-0	Omit Seat Cushion	No Charge	

M300 Corpus HD

Armrest Options

CHOOSE ARMBAR MOUNTING POSITION:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 106022-99-0	Set Armbar at 20"W Between Armrests: _____	No Charge	
<input type="checkbox"/> 106023-99-0	Set Armbar at 21"W Between Armrests: _____	No Charge	
<input type="checkbox"/> 106024-99-0	Set Armbar at 22"W Between Armrests: _____	No Charge	
<input type="checkbox"/> 106025-99-0	Set Armbar at 23"W Between Armrests: _____	No Charge	
<input type="checkbox"/> 106026-99-0	Set Armbar at 24"W Between Armrests: _____	No Charge	
<input type="checkbox"/> 106027-99-0	Set Armbar at 25"W Between Armrests: _____	No Charge	

CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10923	Height Adjustable Armrests	371.00	
<input type="checkbox"/> 10924	Height Adjustable Armrests - NC	No Charge	

No charge when ordering power tilt and/or recline.

SET ARMREST HEIGHT AT: _____

Corpus HD Armrests are infinitely adjustable from 8" - 16", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.

CHOOSE ARMREST PAD RIGHT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> 105997-99-0	4"x16" Corpus HD Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> 105999-99-0	4"x18" Corpus HD Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> 106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> 105993-99-0	4"x16" Corpus HD Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> 105995-99-0	4"x18" Corpus HD Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> 107593-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Right	187.00	
<input type="checkbox"/> 108246-99-0	4"x16" Corpus HD Arm w/ Permobil Gel Pad - Right	192.00	
<input type="checkbox"/> 107597-99-0	4"x18" Corpus HD Arm w/ Permobil Gel Pad - Right	198.00	
<input type="checkbox"/> 108861-99-0	Moldable Forearm Support - Right	268.00	
	<i>11.3"L x 5"W. Not compatible with armrest pouch, trays (except Small Tray) or the VS Chest Support.</i>		
<input type="checkbox"/> 108863-99-0	Contoured Hand Pad Medium - Right	133.00	
	<i>6.4"L x 5.5" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>		
<input type="checkbox"/> 108865-99-0	Contoured Hand Pad Large - Right	140.00	
	<i>7.7"L x 6.1" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>		

CHOOSE ARMREST PAD LEFT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> 105996-99-0	4"x16" Corpus HD Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> 105998-99-0	4"x18" Corpus HD Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> 106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> 105992-99-0	4"x16" Corpus HD Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> 105994-99-0	4"x18" Corpus HD Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> 107592-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Left	187.00	
<input type="checkbox"/> 108247-99-0	4"x16" Corpus HD Arm w/ Permobil Gel Pad - Left	192.00	
<input type="checkbox"/> 107596-99-0	4"x18" Corpus HD Arm w/ Permobil Gel Pad - Left	198.00	
<input type="checkbox"/> 108860-99-0	Moldable Forearm Support - Left	268.00	
	<i>11.3"L x 5"W. Not compatible with armrest pouch, trays (except Small Tray) or the VS Chest Support.</i>		

M300 Corpus HD

Armrest Options

CHOOSE ARMREST PAD LEFT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 108862-99-0	Contoured Hand Pad Medium - Left <i>6.4"L x 5.5" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	133.00	
<input type="checkbox"/> 108864-99-0	Contoured Hand Pad Large - Left <i>7.7"L x 6.1" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	140.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105463-99-0	Corpus Arm Trough Adapter 1" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked (I10645 & I10646) when selecting this option.</i>	359.00	
<input type="checkbox"/> 105462-99-0	Corpus Arm Trough Adapter 7/8" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked (I10645 & I10646) when selecting this option.</i> Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____	359.00	
<input type="checkbox"/> 107339-99-0	Armrest Pouch - Right	254.00	
<input type="checkbox"/> 107338-99-0	Armrest Pouch - Left	254.00	
<input type="checkbox"/> 107840-99-0	UT Armrest Side Support - Right Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	366.00	
<input type="checkbox"/> 107863-99-0	UT Armrest Side Support - Left Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	366.00	
<input type="checkbox"/> 107864-99-0	UT Armrest Side Support - Right Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i>	366.00	
<input type="checkbox"/> 107865-99-0	UT Armrest Side Support - Left Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i>	366.00	
<input type="checkbox"/> 107839-99-0	UT Elbow Support - Right <i>Must also select PN I10881.</i>	181.00	
<input type="checkbox"/> 107838-99-0	UT Elbow Support - Left <i>Must also select PN I10882.</i>	181.00	
<input type="checkbox"/> 10881	Adjustable Elbow Support Hardware - Right	266.00	
<input type="checkbox"/> 10882	Adjustable Elbow Support Hardware - Left	266.00	

Legrest Options

CHOOSE LEGREST TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 106053-99-0	Corpus HD Manual Adj Legrest Elevation	No Charge	
<input type="checkbox"/> 106052-99-0	Power Elevating Center Mount Legrests - 175° <i>Includes mechanical articulation. HCPCS Determined by Pay Source.</i>	3,935.00	E1012

CHOOSE FOOTPLATE STYLE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 106008-99-0	2pc Wide Corpus Footplates - 9"D x 7.5"W	No Charge	
<input type="checkbox"/> 106009-99-0	2pc Extra Wide Corpus Footplates - 10"D x 9"W	279.00	

M300 Corpus HD

Legrest Options

FOOTPLATE HEIGHT POSITION:

Part Number	Description	Price	HCPCS Code
	SET FOOTPLATE HEIGHT AT: _____ <i>This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90° DO NOT INCLUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments, please indicate desired footplate height. NOTE: Standard footplate height range is 6" - 18"; however, height may need to be reduced for shipping. Calf Pads are not available at heights less than 10".</i>		

CHOOSE LATERAL FOOTPLATE ROTATION OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 106006-99-0	Lateral Footplate Rotation - Fixed	No Charge	
<input type="checkbox"/> 106007-99-0	Lateral Footplate Rotation - Enabled	34.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 106190-99-0	VS Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad. Compatible with footplate heights of 10" - 18".</i>	No Charge	
<input type="checkbox"/> 106188-99-0	VS Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad. Compatible with footplate heights of 11" - 18".</i>	393.00	
<input type="checkbox"/> 106186-99-0	VS Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8"W x 8"T Pad. Compatible with footplate heights of 12" - 18".</i>	393.00	
<input type="checkbox"/> 106225-99-0	Stealth Gel Pads - Wide Footplates (Pair)	458.00	
<input type="checkbox"/> 106224-99-0	Stealth Gel Pads - Extra Wide Footplates (Pair)	468.00	

Accessories

<input type="checkbox"/> 104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge
<input type="checkbox"/> 104520-99-0	Bodypoint Padded Hip Belt - Medium	159.00
<input type="checkbox"/> 10940	Bodypoint Non-Padded Hip Belt, UniTrack - Large	No Charge
<input type="checkbox"/> 104516-99-0	Bodypoint Padded Hip Belt - Large	212.00
<input type="checkbox"/> 106531-99-0	Bodypoint Padded Hip Belt - Extra Large	252.00
<input type="checkbox"/> 100519-99-0	Positioning Belt Retractable	316.00
<input type="checkbox"/> 100600-99-0	Positioning Belt Retractable Long	328.00
<input type="checkbox"/> 108267-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Med	213.00
<input type="checkbox"/> 108266-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Large	224.00

HEADREST OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 107470-99-0	Permobil 7"W Headrest <i>Must also select headrest mounting hardware.</i>	413.00	E0955
<input type="checkbox"/> 107469-99-0	Permobil 10"W Headrest <i>Must also select headrest mounting hardware.</i>	413.00	E0955
<input type="checkbox"/> 107468-99-0	Permobil 14"W Headrest <i>Must also select headrest mounting hardware.</i>	413.00	E0955
<input type="checkbox"/> 107480-99-0	Adj/Removable Headrest Link Hardware - Standard	284.00	E1028
<input type="checkbox"/> 109178-99-0	Extra Link for Link Headrest Hardware <i>Used to extend range of standard link hardware. Sent as part.</i>	52.00	
<input type="checkbox"/> 108721-99-0	Adj/Removable BodiLink Head Support Hdwr - MAX 8 <i>Must be selected when ordering more than 2 swing-away headrest accessories or chin boom assemblies.</i>	347.00	E1028
<input type="checkbox"/> 109624-99-0	Extra 4" Link for BodiLink MAX 8 Hdwr <i>Used to extend range of MAX 8 hardware. Sent as part.</i>	76.00	
<input type="checkbox"/> 11305	Swing Away Egg Switch Assembly <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	609.00	E1028

M300 Corpus HD

Accessories

HEADREST OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 11306	Swing Away Microlight Assembly <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	695.00	E1028
<input type="checkbox"/> 11307	Swing Away Facial Pad Assembly 4"L x 2"T <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	494.00	E1028
<input type="checkbox"/> 11308	Swing Away Facial Pad Assembly 6"L x 2"T <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	494.00	E1028
<input type="checkbox"/> 106331-99-0	Permobil 10"W Headrest, Leatherette <i>Not available with ROHO AGILITY Back Systems. Must also select PN I10552.</i>	378.00	E0955
<input type="checkbox"/> 10552	Adj/Removable Headrest Hardware	291.00	E1028
<input type="checkbox"/> 107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	287.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	172.00	E1020
<input type="checkbox"/> 103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	172.00	E1020
<input type="checkbox"/> 10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left: _____ Mount Amp Suppt on Right: _____</i>	291.00	
<input type="checkbox"/> 109040-99-0	BodiLink Premium LPTS Pad (B), Stretch-Air - Pair <i>4.5"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	341.00	E0953
<input type="checkbox"/> 109043-99-0	BodiLink Premium LPTS Pad (A), Stretch-Air - Pair <i>6"W x 5.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	341.00	E0953
<input type="checkbox"/> 109041-99-0	BodiLink Premium LPTS Pad (C), Stretch-Air - Pair <i>8"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	341.00	E0953
<input type="checkbox"/> 109042-99-0	BodiLink Premium LPTS Pad (G), Stretch-Air - Pair <i>12"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	341.00	E0953
<input type="checkbox"/> 109044-99-0	Adj/Removable BodiLink LPTS Hdwr Med Corpus Pair <i>Medium TT Power Mount hardware can be extended upward to clear cushions of 3-5.5" thick.</i>	656.00	E1028
<input type="checkbox"/> 109194-99-0	Adj/Removable Lateral Thigh/Hip Hardware - Pair	573.00	E1028
<input type="checkbox"/> 105260-99-0	Transfer Handles, Corpus Seat	523.00	
<input type="checkbox"/> 105292-99-0	Push Handles, Corpus Seat	465.00	
<input type="checkbox"/> 106210-99-0	RAM Self-Leveling Cup Holder	88.00	
<input type="checkbox"/> 106806-99-0	RAM X-Grip Phone Holder	88.00	
<input type="checkbox"/> 107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount (I103476-99-0). Customization may be required if ordered with an Omni.</i>	444.00	
<input type="checkbox"/> 107189-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. No charge with 4 power seat function configurations. Mount on left armrest: _____ Mount on right armrest: _____</i>	220.00	
<input type="checkbox"/> 105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordered with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest. May ship as a part at shorter backrest heights.</i>	105.00	

M300 Corpus HD

Accessories

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 102059-99-0	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest.</i>	305.00	
<input type="checkbox"/> 105840-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as a part.</i>	194.00	E0960
<input type="checkbox"/> 105839-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as a part.</i>	194.00	E0960
<input type="checkbox"/> 105850-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as a part.</i>	194.00	E0960
<input type="checkbox"/> 10556	Adj Removable - Tray Hardware (1 Each)	291.00	
<input type="checkbox"/> 108227-99-0	Permobil Adjustable Tray - 17"W	270.00	E0950
<input type="checkbox"/> 108229-99-0	Permobil Adjustable Tray - 19"W	270.00	E0950
<input type="checkbox"/> 108232-99-0	Permobil Tray Hardware Right Side Adjustable	198.00	
<input type="checkbox"/> 108233-99-0	Permobil Tray Hardware Left Side Adjustable	198.00	
<input type="checkbox"/> 105224-99-0	Oxygen Holder, Corpus	286.00	E2208
<input type="checkbox"/> 105223-99-0	Crutch Holder, Corpus	286.00	E2207

Stealth Accessories

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 106221-99-0	Stealth 4" x 6" Lateral Supports, HD <i>Must also select I10723. Kit contains thin, Stealth hardware with swing away brackets and two height adjustment rails. The rails are bolted to the backrest to provide mounts for the lateral hardware. PLEASE NOTE: The spacing between the armrests must be at least 1" wider than the backrest width when ordering these laterals.</i>	444.00	E0956
<input type="checkbox"/> 10723	Adjustable Stealth Lateral Hardware (Pair)	602.00	E1028

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10751	Mini-Joystick for Omni - (9-pin connector) <i>Operating Force = 0.5 oz (13 g) Neutral to Full Deflection = 0.2" (5mm)</i>	3,447.00	E2312
<input type="checkbox"/> 103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 7.8 oz (220 g) Neutral to Full Deflection = 1.1" (28mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,686.00	E2373
<input type="checkbox"/> 103894-99-0	Compact Joystick Lite for Omni - (9-pin connector) <i>Operating Force = 1.8 oz (50 g) Neutral to Full Deflection = 0.7" (17mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,907.00	E2373

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
NOTE: <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options. The Swing Away Chin Boom Assemblies are Compatible with Permobil Headrests I107470-99-0, I107469-99-0 or I107468-99-0 and must be selected with I108721-99-0.</i>			
<input type="checkbox"/> 107946-99-0	Swing Away Chin Boom Assm for Compact Lite - Left	442.00	E1028
<input type="checkbox"/> 107948-99-0	Swing Away Chin Boom Assm for Compact Lite - Right	442.00	E1028
<input type="checkbox"/> 107857-99-0	Swing Away Chin Boom Assm for Mini Joystick - Left	442.00	E1028
<input type="checkbox"/> 107947-99-0	Swing Away Chin Boom Assm for Mini Joystick - Right	442.00	E1028

M300 Corpus HD

Alternative Drive Systems

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 104690-99-0	Alt JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	295.00	E1028
<input type="checkbox"/> 107135-99-0	Alt JS / Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod.</i> <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	558.00	E1028
<input type="checkbox"/> 104143-99-0	Alt JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	920.00	E1028
<input type="checkbox"/> 10752	Chin Control Harness, Large <i>For neck width up to 7½".</i>	654.00	E1028
<input type="checkbox"/> 10753	Chin Control Harness, Small <i>For neck width up to 5½".</i>	654.00	E1028

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches set-up for Mode & On/Off function. Switch operation force = 4.5oz (130 g)</i> <i>Mount at front of Joystick: _____ Mount at rear of Joystick: _____</i> <i>Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	274.00	
<input type="checkbox"/> 104641-99-0	Dual Switches for Compact Joystick -Gooseneck Mount <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches set-up for Mode & On/Off function. Switch operation force = 4.5oz (130 g)</i>	493.00	

SIP & PUFF SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Includes breath tube kit.</i>	1,844.00	E2325,E2326
<input type="checkbox"/> 10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Includes breath tube kit.</i>	2,454.00	E2325,E2326
<input type="checkbox"/> 10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573 or I10574.</i>	295.00	

PROPORTIONAL FOOT CONTROL:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10767	Proportional Foot Control <i>Not available with Active Reach Packages.</i> <i>Mount on Left: _____ Mount on Right: _____</i>	3,514.00	E2328
<input type="checkbox"/> 10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	295.00	

Head Array Options

TOTAL CONTROL HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 109558-99-0	Permobil Total Control Head Array <i>Features an occipital pad sensor and two lateral arms configured with proximity sensors. Must select swing away hardware and mounting hardware below to complete assembly. Choose additional sensors/switches, including mode switch, from options below. ESP recommended. Choose Omni Control Module separately.</i>	4,589.00	E2330
<input type="checkbox"/> 11575	Swing Away Hardware for Side Sensor - Left	250.00	E1028
<input type="checkbox"/> 11576	Swing Away Hardware for Side Sensor - Right	250.00	E1028
<input type="checkbox"/> 11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module separately.</i>	No Charge	

M300 Corpus HD

Head Array Options

TOTAL CONTROL HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):			
<input type="checkbox"/> 109563-99-0	Adj/Removable BodiLink HA Hdwr - MAX 8	347.00	E1028
<input type="checkbox"/> 109564-99-0	Adj/Removable BodiLink HA Hdwr - MAX 8 Flip Down <i>Features flip-down joint that allows head array to swing back 120° away from user's head.</i>	477.00	E1028
SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> 11601	Add Left Lateral Arm w Proximity Sensor	655.00	E1028
<input type="checkbox"/> 11602	Add Left Lateral Arm w Egg Switch	536.00	E1028
<input type="checkbox"/> 11603	Add Left Lateral Arm w Pad (No Sensor or Switch)	437.00	E1028
SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> 11604	Add Right Lateral Arm w Proximity Sensor	655.00	E1028
<input type="checkbox"/> 11605	Add Right Lateral Arm w Egg Switch	536.00	E1028
<input type="checkbox"/> 11606	Add Right Lateral Arm w Pad (No Sensor or Switch)	437.00	E1028

ASL FUSION HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 109582-99-0	ASL FUSION Proportional and Digital Head Array <i>Includes occipital (ASL 935) and lateral (ASL 933) pads, and pulse mechanical switch (ASL 320BK) for reset/mode changes. Must select hardware listed below to complete assembly.</i>	6,395.00	E2328
<input type="checkbox"/> 11580	Swing Away Hardware for Lateral Sensor - Left	325.00	E1028
<input type="checkbox"/> 11581	Swing Away Hardware for Lateral Sensor - Right	325.00	E1028
<input type="checkbox"/> 11582	Swing Away Hardware for Reset/Mode Switch - Right	325.00	E1028
<input type="checkbox"/> 11583	Adj/Removable Headrest Mounting Hardware <i>Includes ASL 606 and BodiLink receiver.</i>	325.00	E1028
<input type="checkbox"/> 109583-99-0	ASL FUSION Display, Prgmr and Att Control w/Mount <i>Includes ASL 165 and quick release mount.</i>	695.00	K0108

Alt Drive Systems Setup

CHOOSE OMNI (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 108500-99-0	Omni2 Control Module for R-net - Corpus <i>Two-section design features 3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology for control of up to 4 Bluetooth devices. Must also select mounting hardware below. Please specify which type of alternative drive input device you will be using. Input Device: _____</i>	1,962.00	
<input type="checkbox"/> 103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	295.00	

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 10080	Mount Seat Control Switchbox on Backrest	No Charge	

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	297.00	
<input type="checkbox"/> 106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	522.00	
<input type="checkbox"/> 106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	748.00	

M300 Corpus HD

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 105907-99-0	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	378.00	E1028
<input type="checkbox"/> 105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	255.00	
<input type="checkbox"/> 105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	154.00	
<input type="checkbox"/> 106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	135.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> 101769-99-0	Microlight Switch Qty: _____	206.00	
<input type="checkbox"/> 101767-99-0	Buddy Button Qty: _____	174.00	
<input type="checkbox"/> 106293-99-0	Egg Switch Black Qty: _____	174.00	
<input type="checkbox"/> 101771-99-0	Mini Cup Switch Qty: _____	174.00	
<input type="checkbox"/> 107500-99-0	Proximity Switch Qty: _____	688.00	
<input type="checkbox"/> 107495-99-0	Mini Joystick 5 Switch OC for Omni - (9-pin)	747.00	
<input type="checkbox"/> 10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 6" rod, one 8" rod and one 9" gooseneck switch with 50" cable and mounting clamp.</i>	545.00	
<input type="checkbox"/> 10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.</i>	806.00	
<input type="checkbox"/> 105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	34.00	
<input type="checkbox"/> 105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	37.00	
<input type="checkbox"/> 105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	44.00	
<input type="checkbox"/> 105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	48.00	
<input type="checkbox"/> 105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	37.00	

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer. ** Indicates HCPCS code is pending PDAC approval.*