



TiLite CR1 - Order Form



Dealer Name*:

Account #*:

City, State, Zip:

Phone:

Email:

ATP:

Rep Name:

P.O. Number:

Quote Number:

Client Reference*:

Client Weight:

Special considerations that need to be addressed here (e.g., diagnosis):

*Required field

Frame Style

☐ CFS1 **CR1 Carbon Fiber Frame** Minimal adjustability. 231lb weight limit. \$ 8,320.00

SmartDrive Power Assist System

☐ CMAXSD1 **Permobil SmartDrive MX2+** Requires CCBR2 SmartDrive ready carbon fiber camber tube. E0986 \$ 7,529.00

☐ CMAXSD18 **Switch Control Buttons** No charge

☐ CMAXSD8 **SwitchControl Button with Mono Jack Option** No charge

☐ CMAXSD19 **Buddy Button** For use with SwitchControl button with mono jack option. No charge

☐ CMAXSD4 **Carrying Bag** \$ 84.00

☐ CMAXSD10 **Universal Padded Handle** \$ 22.00

Frame Finishes

Color Anodize Package

Anodized Package: ☐ Charcoal Grey ☐ Gold ☐ Green ☐ Purple ☐ Orange ☐ Blue

☐ Red

Includes rear strut band clamps, footrest clamps, caster centre caps and side guard caps. If no anodize colour is selected, default will be charcoal grey.

Trim Rings: ☐ Charcoal Grey ☐ Polished Aluminum

Frame and backrest.

TiFit: Tailored Individualized Frame

Please refer to Reference Diagrams at end of form for more information.

Front End

☒ CFE1 **2.1" Tapered Front End 'I'** Footrest width is 4.2" narrower than seat width.

☒ CFA1 **Front Frame Angle 85°**

Seat Width 'A' | Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

CSW1: ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18"

13" seat width requires CBUP3 padded tension adjustable by strap back upholstery.

Seat Depth 'B' | Measure from front of back post (at intersection with seat tube) to front edge of seat sling. Frame bend will start 1" in front of seat sling unless additional custom frame depth is specified or the combination of specifications affect the rolling dynamics or functionality of the chair. Please note that 1" is the equivalent to 0.5" in front of the seat sling on other TiLite frame styles.

CSD1: ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19"

Seat Heights 'C' & 'D' | 17.5-17.75" front seat to floor requires 3" casters and 24" rear wheels. Extra frame depth is measured from front of seat sling to beginning of frame bend, standard is 1".

Seat Depth	Max Extra Frame Depth	Extra Frame Depth CSD3 Set in 0.5" increments.	Front STF Range	Front STF 'C' Set in 0.25" increments.	Rear STF 'D' Amount lower than Front STF	Frame Type
13"	☐ 2"	☐ Set at	☐ 17.5"-17.75" ☐ 18"-20.75"	☐ Set at ☐ Set at	☐ 1.6" lower ☐ 1.6" lower ☐ 2.3" lower ☐ 2.9" lower	Frame 1 Frame 1 Frame 2 Frame 3
14"	☐ 2"	☐ Set at	☐ 17.5"-17.75" ☐ 18"-20.75"	☐ Set at ☐ Set at	☐ 1.7" lower ☐ 1.7" lower ☐ 2.4" lower ☐ 3.1" lower	Frame 1 Frame 1 Frame 2 Frame 3
15"	☐ 2"	☐ Set at	☐ 17.5"-17.75" ☐ 18"-20.75"	☐ Set at ☐ Set at	☐ 1.8" lower ☐ 1.8" lower ☐ 2.6" lower ☐ 3.4" lower	Frame 1 Frame 1 Frame 2 Frame 3
16"	☐ 2"	☐ Set at	☐ 17.5"-17.75" ☐ 18"-20.75"	☐ Set at ☐ Set at	☐ 1.9" lower ☐ 1.9" lower ☐ 2.8" lower ☐ 3.6" lower	Frame 1 Frame 1 Frame 2 Frame 3
17"	☐ 2"	☐ Set at	☐ 18"-20.75"	☐ Set at	☐ 2.1" lower ☐ 3.0" lower ☐ 3.8" lower	Frame 1 Frame 2 Frame 3
18"	☐ 1"	☐ Set at	☐ 18"-20.75"	☐ Set at	☐ 2.2" lower ☐ 3.1" lower ☐ 4.0" lower	Frame 1 Frame 2 Frame 3
19"	☐ N/A	☐ N/A	☐ 18"-20.75"	☐ Set at	☐ 2.3" lower ☐ 3.3" lower ☐ 4.3" lower	Frame 1 Frame 2 Frame 3

Footrest

Seat to Footrest 'E' | Measure from front edge of seat sling to top middle of footrest. Measurement must be 1"-3" shorter than front seat height.

Footrest Height: ☐ Set At

0.25" increments. CMD25 is available for individuals needing a shorter footrest height.

☐ CFTR13	Carbon Fiber Open Loop Footrest With center support tube.	Standard
☐ CFTR14	Carbon Fiber Footrest w/Flat Carbon Fiber Cover	\$ 200.00
☐ CFCT1	Grip Tape for Flat Carbon Fiber Cover	\$ 15.00

2" sided custom width/depth aluminum footrest cover CMD49 is available for individuals wanting to keep their feet centred on the footrest.

Footrest Trim Rings

☐ CFAR3	Stainless Steel	Standard
☐ CFAR4	Titanium	\$ 120.00

Back Supports

ROHO Back Supports

☐ SBRS14	Agility Active Aluminum, Min Contour with Foam Includes fixed hardware. Select shell height.	E2620	\$ 738.00
☐ SBRS15	Agility Active Aluminum, Mid Contour with Foam Includes fixed hardware. Select shell height.	E2620	\$ 738.00
	☐ 8" ☐ 10"		
☐ SBRS7	Agility Carbon Min Contour Back with Foam/Air Insert 2x2 twill weave carbon fiber. Includes fixed hardware. Select shell height.	E2620	\$ 738.00
☐ ACB1	Agility Carbon Min Upgrade Required with Agility Carbon Backs. Not available with bolt-on push handles.		\$ 464.00
☐ SBRS9	Agility Carbon Mid Contour Back with Foam/Air Insert 2x2 twill weave carbon fiber. Includes fixed hardware. Select shell height.	E2620	\$ 738.00
☐ ACB2	Agility Carbon Mid Upgrade Required with Agility Carbon Backs. Not available with bolt-on push handles.		\$ 464.00
	☐ 10" ☐ 13" ☐ 16"		

Back Upholstery

☐ CBUP3	Tension Adjustable By Straps, Padded Nylon Required for 13" seat width.	Standard
☐ CBUP0	Omit Back Upholstery	No charge

Backrest

Back Type

☐ CSB15	Carbon Folding Fixed Height Backrest 10"-17" height 'F', 0.25" increments.	Standard
	☐ Set At	
Back Angle 'G':	☐ 80° ☐ 81° ☐ 82° ☐ 83° ☐ 84° ☐ 85° ☐ 86° ☐ 87° ☐ 88° ☐ 89° ☐ 90° ☐ 91° ☐ 92° ☐ 93° ☐ 94° ☐ 95° ☐ 96° ☐ 97° ☐ 98° ☐ 99° ☐ 100° ☐ 101°	

Side Guards

☐ CSG1	Download Pivoting Fixed Lockout Carbon Fiber Side Guards Profiled 0.125" above rear wheels or as close as possible depending upon chair configuration.	Standard
☐ CSG2	Download Pivoting Fixed Lockout Fendered Carbon Fiber Side Guards Custom cut to profile of rear wheels.	\$ 460.00

Push Handles

☐ CSB10	Integrated Carbon Fiber Push Handles, Short	\$ 300.00
--------------------------------	--	-----------

Rear Wheel Positioning

Center of Gravity 'H' | Measure from front of back post at the seat tube to the center of the rear axle. Please note that due to its design and materials the CR1 has a lighter than normal front end compared to other TiLite models. This will affect the overall feel when selecting this measurement. There is 1" maximum adjustability in either direction from specified location.

CCOG1: ☐ 2.25" ☐ 2.5" ☐ 2.75" ☐ 3" ☐ 3.25" ☐ 3.5" ☐ 3.75" ☐ 4" ☐ 4.25" ☐ 4.5" ☐ 4.75" ☐ 5"

Camber 'K'

Camber Angle: ☐ 0° ☐ 2° ☐ 4°

☐ CCB1	Aero Full Carbon Fiber Camber Tube Includes integrated zero play carbon fiber camber plugs. Not available with 13" seat width or certain configuration combinations.	Standard
☐ CCB2	SmartDrive Ready (Round) Carbon Fiber Camber Tube Includes integrated zero play carbon fiber camber plugs.	No charge

Seat Upholstery

☒ CSUP2	Tension Adjustable By Straps, Black	Standard
☐ CMD110	Carbon Fiber Seat Pan Please specify size in Notes section at the end of this order form. Maximum size is 17"w x 19"d.	CMD

Seat Cushions

☐ CASH2	Comfort Curve Cushion Stretch-Air cover.	E2601	\$ 122.00
☐ CASH23	Comfort Acta-Embrace Comfort-Tek Zero Elevation Cushion	E2607	\$ 428.00
☐ CASH11	ROHO Single Compartment Fits seat widths 14"-19".	E2622	\$ 534.00
	☐ Low Profile ☐ Mid Profile ☐ High Profile		
☐ CASH14	ROHO Quadtro Select Cushion Fits seat widths 14"-19".	E2624	\$ 534.00
	☐ Low Profile ☐ Mid Profile ☐ High Profile		

Front Wheels

☐ CFW24	3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Required when front seat height is 17.5"-17.75"</i>	No charge
	☐ Silver Hub ☐ Black Hub	
☐ CFW25	4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	Standard
	☐ Silver Hub ☐ Black Hub	
☐ CFW26	4"x1.25" LiteSpeed-6 Carbon Fiber Wheel w/Soft Roll Tire	\$ 190.00

Front Forks

☐ CFK2	TiLite Slipstream CR1 Single-Sided Fork <i>Shorter length and width, stainless steel caster axle.</i>	Standard
☐ CTFS1	Titanium Fork Stem	Standard
☐ CTFA1	Titanium Caster Axle	\$ 130.00

Rear Wheels

☐ CRW14	Spinergy SPOX Rear Wheel, Black Spokes <i>Select size.</i>	Standard
	☐ 24" (540) ☐ 25" (559)	
☐ CRW16	Spinergy Carbon Blade Rear Wheel <i>2x2 twill weave carbon fiber. Select size. Not available with Natural-Fit or Natural-Fit LT handrims with thumb grip.</i>	\$ 1,795
	☐ 24" (540) ☐ 25" (559)	

Axles

☐ CAXL1	Stainless Quick Release Axle	Standard
☐ CAXL2	Titanium Quick Release Axle	\$ 277.00
☐ CAXL3	Quad Release Axle	\$ 169.00

Pneumatic Rear Tires

☐ CRTR1	MBL SpeedLite Tire 24"(540), 25"(559)	No charge
☐ CRTR10	MBL Gazelle Tire 24"(540), 25"(559) <i>Kevlar bead, lightest weight option.</i>	\$ 99.00
☐ CRTR8	Schwalbe Marathon Plus Evolution Tire 24"(540), 25"(559)	\$ 169
☐ CRIM21	Aluminum Handrim, Black Anodized	Standard

Handrims

☐ CRIM23	Spinergy Handrim, Aluminum Anodized <i>Anodized to match chair anodize package color.</i>	\$ 240.00
☐ CRIM12	Titanium Handrim, Satin Finish	\$ 509.00
☐ CRIM20	Q-Grip Handrim <i>Exclusive high-friction, anti-slip coating.</i>	\$ 358.00
☐ CRIM17	The Natural-Fit LT Handrim, Standard Grip <i>Reduced weight, small profile. Not with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$ 440.00
☐ CRIM18	The Natural-Fit LT Handrim, Super Grip <i>Reduced weight, small profile. Not with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$ 514.00
☐ CRIM19	The Natural-Fit LT Handrim, Without Thumb Grip <i>Reduced weight, small profile.</i>	\$ 423.00

Handrim Tab Length

☐ CRIM7	Long Tabs	Standard
☐ CRIM8	Short Tabs <i>Trimmed. Not available with The Natural-Fit handrims with thumb grip.</i>	\$ 30.00

Wheel Locks

☐ CWLK9	MBL Composite Scissor Lock	Standard
☐ CWLK12	MBL Performance Aluminum Scissor Lock	\$ 175
☐ CWLK10	Omit Wheel Locks <i>Selecting this option prevents adding wheel locks at a later date.</i>	No charge

Positioning Belts

☐ CBLT1	Velcro Adjustable Positioning Belt, Black Nylon Webbing	E0978	\$ 61.00
--------------------------------	--	-------	----------

Calf Straps

☐ CCLF3	Bodypoint Padded Velcro Adjustable Calf Strap	K0038	\$ 109.00
	☐ Medium (13"-16" footrest widths) ☐ Large (17"-18" footrest widths)		

Anti-Tips

☐ CTIP1	Rear Anti-Tips, Aluminum Black	E0971	\$ 174.00
--------------------------------	---------------------------------------	-------	-----------

Accessories

☐ CMPCT2	Black Neoprene Impact Guards with Ultrasuede Stripe	\$ 120.00
☐ CBACK3	Permobil Medical Necessities Backpack	\$ 103.00
☐ CPCH3	Permobil Medical Necessities Seat Pouch	\$ 66.00
☐ CPCH4	Bodypoint Mobility Bag	\$ 109.00
☐ CAFP1	4" Military Service Patch	\$ 115.00
☐ Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard		

General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes

TiFit Reference Diagrams

