

# TiLite Pilot Canadian Order Form

Effective May, 2022 (New Items in Red)

CARBON RIGID, GROWABLE FRAME

Page 1 of 5

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		<b>SHIP TO INFORMATION</b>	
Mark For:	User Weight:	Name:	Ship Via:
TiLite Sales Rep:		Address:	
<b>TiLite Information:</b>	Promo Code:	City:	Province: Postal Code:
Ship Date:	Order #:	Phone:	Fax:

**Please Email or Fax to the Permobil Manual Customer Service Team at [tilitesales@permobil.com](mailto:tilitesales@permobil.com) or 800.736.0925, option 2.**

## FRAME STYLE

- TPFS1 ☐ **TiLite Pilot INCLUDING Growth Kit** (Full Adjustability. Seat Width **A** grows 2" and Seat Depth **B** grows 3".) 120 Pound Weight Limit. Price includes one-time growth kit and Transit Tie-Down option. Additional growth kits sold separately.).....\$3833
- TPFS2 ☐ **TiLite Pilot OMIT Growth Kit**.....\$3544

## SMARTDRIVE POWER ASSIST SYSTEM



- TPMAXSD1 ☐ **Permobil Smart Drive MX2+** (N/A with rear seat heights less than 13.5", or 20" rear wheels).....\$7870
- TPMAXSD18 ☐ SwitchControl Buttons.....NCO
- TPMAXSD13 ☐ **NEW PushTracker E3** (Not required).....\$70
- TPMAXSD8 ☐ SwitchControl Button with Mono Jack Option.....\$204
- TPMAXSD19 ☐ **NEW Buddy Button**.....\$362
- TPMAXSD4 ☐ Carrying Bag.....\$108
- TPMAXSD10 ☐ Universal Padded Handle.....\$30

\*SmartDrive apps available on Watch OS (Apple Watch Series 5 or newer).

## TRANSIT TIE-DOWN OPTION (Refer to the Transit Tie-Down Option Brochure for important information regarding the Transit Tie-Down options.)

- TPWTB ☐ Transport Brackets (Consists of four Tie-Down Brackets.).....NCO
- TPTTDO ☐ Transit Tie-Down (WC-19 Tie-Down Brackets require the TiLite back and head support. 95 lb weight limit.).....\$410

## FRAME FINISHES

- ☐ TPFF29 Neon Pink ☐ TPFF14 Ocean Blue Metallic ☐ TPFF30 Acid Green Metallic ☐ TPFF25 Electric Plum.....STD

## COLOR ANODIZE PACKAGE

(Color Anodize Package includes: Caster Mount, Caster Fork, and Caster Caps. If no anodize color is selected, the default will be black.)


- ☐ ZANO1 Black.....STD
- ☐ ZAN12 Gold ☐ ZAN15 Green ☐ ZAN13 Purple ☐ ZAN14 Orange ☐ ZAN10 Blue ☐ ZAN11 Red.....NCO

## FRAME MEASUREMENTS

### FRONT SEAT TO FLOOR HEIGHT, REAR SEAT TO FLOOR HEIGHT, CENTER OF GRAVITY AND SEAT TO FOOTREST

Please see the table below for available configuration adjustment ranges. Center of gravity will be set at the factory. See documentation shipped with chair for detailed information on adjustments.

SEAT DEPTH	MEASUREMENT:	20" REAR WHEEL	22" REAR WHEEL	24" REAR WHEEL
		AVAILABLE RANGE (IN):	AVAILABLE RANGE (IN):	AVAILABLE RANGE (IN):
9"	FRONT SEAT HEIGHT	12.5-15	N/A	N/A
	REAR SEAT HEIGHT	12-13.5	N/A	N/A
	CENTER OF GRAVITY	0.5-1.5	N/A	N/A
10"	FRONT SEAT HEIGHT	12-15	N/A	N/A
	REAR SEAT HEIGHT	12-13.5	N/A	N/A
	CENTER OF GRAVITY	0.5-2.25	N/A	N/A
11"	FRONT SEAT HEIGHT	12-15	13.5-16.5	14-16.5
	REAR SEAT HEIGHT	12-13.5	13-15	13.5-16
	CENTER OF GRAVITY	0.5-3	0.5-2.25	0.5-2.25
12"	FRONT SEAT HEIGHT	12-15	13.5-16	14-16.5
	REAR SEAT HEIGHT	12-13.5	13-15	13.5-16
	CENTER OF GRAVITY	0.5-3.5	0.5-3.25	0.5-3.25
13"	FRONT SEAT HEIGHT	N/A	13-16	13.5-16.5
	REAR SEAT HEIGHT	N/A	13-15	13.5-16
	CENTER OF GRAVITY	N/A	0.5-3.5	0.5-3.5



Permobil.com  
2701 West Court Street / Pasco, Washington 99301  
Tel 800 736 0925 or Email [tilitesales@permobil.com](mailto:tilitesales@permobil.com)

Rev H JR 8/10/2018

# TiLite Pilot Canadian Order Form

Effective May 2022 (New Items in Red)

CARBON RIGID, GROWABLE FRAME

Page 2 of 5

14"	FRONT SEAT HEIGHT	N/A	13-16	13.5-16.5
	REAR SEAT HEIGHT	N/A	13-15	13.5-16
	CENTER OF GRAVITY	N/A	0.5-3.5	0.5-3.5
15"	FRONT SEAT HEIGHT	N/A	13-16	13.5-16.5
	REAR SEAT HEIGHT	N/A	13-15	13.5-15.5
	CENTER OF GRAVITY	N/A	1.5-3.5	1.5-3.5

## FRAME MEASUREMENTS

**SEAT WIDTH A** Measure inside width from side panel to side panel. Seat width grows 2" if selected width is between 8" and 12". The 13" width will grow a maximum of 1", and the 14" width will not grow any wider. Refer to chart for seat depth availability.

TPSW1 ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14".....STD

**SEAT DEPTH RANGE:** The seat depth can be grown within the selected range.

☐ 9"-12" ☐ 11"-15".....STD

**SEAT DEPTH B** Measure from front of seat pan to front of back support cushion.

TPSD1 ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15".....

**FRONT SEAT HEIGHT C** Measure from the floor to the top edge of the front of the seat pan. Adjustable

TPFSH1 ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5".....

**REAR SEAT HEIGHT D** Measure from the floor to the top edge of the of the seat pan at the intersection of support pad and the seat pan. Adjustable in 0.5" increments.

TPRSH1 ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5".....

## FOOTREST TYPE

TPFTR3 ☐ Angle Adjustable Footrest.....

TPFTR24 ☐ Flip-Back Depth and Angle Adjustable Footrest (Not available on 8", 9" or 14" seat widths).....

**SEAT TO FOOTREST E** Measure from the top of the footplate to the front of the seat pan. For an angle adjustable footrest, this measurement must be at least 2" less than front seat height. For a Flip-Back Depth and Angle Adjustable footrest, this measurement must be at least 3" less than front seat height.

☐ 3" ☐ 3.5" ☐ 4" ☐ 4.5" ☐ 5" ☐ 5.5" ☐ 6" ☐ 6.5" ☐ 7" ☐ 7.5" ☐ 8" ☐ 8.5" ☐ 9" ☐ 9.5" ☐ 10" ☐ 10.5" ☐ 11" ☐ 11.5" ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5".....STD

**COMFORT BACK SUPPORTS** (WC-20 certification for integrated hardware is pending testing. Cane mounted hardware is WC-20 certified.)

TPSB4 ☐ Acta-Back (Available on 12" to 14" seat widths)  
Back support length: ☐ 10" ☐ 12" ☐ 14".....\$723

TPSB5 ☐ Acta-Embrace (Available on 10" to 14" seat widths)  
Back support length: ☐ 10" ☐ 12" ☐ 14".....\$800

Comfort Back Options:

TPCHD1 Backrest Mounting Hardware Options: ☐ Integrated Compass Hardware ☐ Cane Mount Hardware.....NCO

TPFAB1 Cover Fabric: ☐ Comfort-Tek ☐ Stretch-Air.....NCO

TPLAT1 ☐ Bodilink 3"L x 4"D Lateral trunk support, each (Swing away, GlideWear cover, non-padded).....\$540  
Mounting Location: ☐ Left ☐ Right

**SEAT BACK TYPE** Measure from top of back support to seat pan. Back support shell height is 9.5" with scapular cutouts.

TPSB1 ☐ Standard Non-Folding Adjustable-Height Rigid Back with Pad (11" to 13" height range).....STD

TPSB2 ☐ Contoured Non-Folding Adjustable-Height Rigid Back with Pad (11" to 13" height range).....NCO

TPSB3 ☐ Aluminum Canes for Aftermarket back support, no Pad (Not available with height adjustable fenders).....NCO

**HEAD SUPPORT** (The Comfort head support is WC-20 certified on the Acta-Back and Acta-Embrace. The head support is WC-19 certified when used with the TiLite backrest.)

TPHDR1 Comfort BodiLink head support with GlideWear cover (Uses piano taper midline hardware. Available on Acta-Back, Acta-Embrace and TiLite rigid back supports).....\$604  
Width: ☐ 6" ☐ 8"

**PUSH HANDLES** (Height adjustable range is 33" - 40" as measured from the floor with 24" (540) wheels.)

TPSB6 ☐ Bolt-On Height-Adjustable Removable Single Push Handle.....\$242

TPSB7 ☐ Bolt-On Height-Adjustable Removable Dual Push Handles.....\$300

**CAMBER / CAMBER TUBES** (Aluminum Camber tube standard. Anodize color will match selection above.)

☐ 2° (Minimum 1.5" rear wheel spacing) ☐ 4° (Minimum 1.25" rear wheel spacing) ☐ 6° (Minimum 1" rear wheel spacing)  
☐ 8° ☐ 12".....STD

**REAR WHEEL SPACING** Measure from outside of side panel to inside of rear tire at top of wheel. (Not all rear wheel spacing is available with all accessories.)


TPRWS1 ☐ 0.75" ☐ 1" ☐ 1.25" ☐ 1.5" ☐ 1.75" ☐ 2".....STD

TPCBR12 ☐ Adjustable Rear Wheel Spacing Camber Plugs.....STD

TPCBR9 ☐ Additional Spacers to Increase Rear Wheel Spacing.....NCO

**FRONT WHEELS** (TiLite Standard Fork style will be provided.)

TPFW14 ☐ 4" X 0.75" Plastic Wheel w/Poly Tire.....STD



# TiLite Pilot Canadian Order Form

Effective May 2022 (New Items in Red)

CARBON RIGID, GROWABLE FRAME

Page 3 of 5

TPFW3	<input type="checkbox"/> 4" X 0.75" Light-Up Micro-Caster.....	\$84
TPFW26	<input type="checkbox"/> 4" X 1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire.....	\$100
TPFW27	<input type="checkbox"/> 4" X 1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire (Black Hub Only).....	\$200

## CUSTOM CONFIGURATION OPTIONS

### REAR WHEELS

TPRW14	TiLite SHADOW™ <input type="checkbox"/> 20"(451) Only available with seat depth range 9"-12". <input type="checkbox"/> 22"(501) Only available with seat depth range 11"-15". <input type="checkbox"/> 24"(540) Only available with seat depth range 11"-15".	STD
TPRW7	Spinergy SPOX <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White <input type="checkbox"/> Orange <input type="checkbox"/> Pink <input type="checkbox"/> Purple <input type="checkbox"/> Green (Orange, Pink, Purple And Green only available in 20" size.) <input type="checkbox"/> 20"(451) - Std. Hub <input type="checkbox"/> 22"(501) - Std. Hub Only available with seat depth range 11"-15". <input type="checkbox"/> 24"(540) - Std. Hub Only available with seat depth range 11"-15".	\$945
TPRW8	Spinergy LX <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White <input type="checkbox"/> Bright Orange <input type="checkbox"/> Flamingo Pink <input type="checkbox"/> Purple <input type="checkbox"/> Green.....	\$1108
TPRW9	New Spinergy Design Lab LX Wheels (If same color is selected, anodized hubs and rims match color anodized package on page 1. Available with black spokes only.) <input type="checkbox"/> Gold Anodize <input type="checkbox"/> Green Anodize <input type="checkbox"/> Blue Anodize <input type="checkbox"/> Red Anodize <input type="checkbox"/> Purple Anodize <input type="checkbox"/> Orange Anodize..... <input type="checkbox"/> 22"(501) - Std. Hub <input type="checkbox"/> 24"(540) - Std. Hub <input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" minimum Seat Depth B.) <input type="checkbox"/> 26"(590) - Std. Hub (Requires 15" minimum Seat Depth B.)	\$1853
TPRW16	New Spinergy® Carbon Blade with 2x2 Twill Weave Carbon Fiber (Not available with Heavy Duty Frame.) <input type="checkbox"/> 24"(540) <input type="checkbox"/> 25"(559) (Requires 15" minimum Seat Depth B.)	\$1617
TPRW17	X-Core Mag (Not available with Short Tabs, Quad Release Axles, or Omit Handrims Option.) <input type="checkbox"/> 20"(451) Only available with seat depth range 9"-12". <input type="checkbox"/> 22"(501) <input type="checkbox"/> 24"(540) Only available with seat depth range 11"-15".	\$111 \$116

### AXLES

TPAXL1	<input type="checkbox"/> Stainless Quick Release.....	STD
TPAXL3	<input type="checkbox"/> Quad Release (Not available with Mag Rear Wheels.).....	\$179

### PNEUMATIC REAR TIRES

TPRTR1	<input type="checkbox"/> Treaded Treaded 22"(501), 24"(540), and 26"(590).....	STD
TPRTR2	<input type="checkbox"/> Airless Inserts.....	\$121
ZRTR31	<input type="checkbox"/> MBL SpeedLite 20"(451), 22"(501), 24"(540).....	STD
TPRTR30	<input type="checkbox"/> Primo Sentinel (Available on all wheel sizes, 20" (451), 22" (501), and 24" (540)).....	\$168
TPRTR8	<input type="checkbox"/> Schwalbe Marathon (Only available on 24" (540)).....	\$163

### HANDRIMS

TPRIM1	<input type="checkbox"/> Aluminum – Silver Anodized.....	STD
TPRIM21	<input type="checkbox"/> Aluminum – Black Anodized.....	NCO

### Handrims by OUT-FRONT®

The Natural-Fit® LT (Reduced Weight, Smaller Profile) (Not available with 20"(451) Rear Wheels.)		
TPRIM17	<input type="checkbox"/> Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$525
TPRIM18	<input type="checkbox"/> Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$562
TPRIM19	<input type="checkbox"/> Without Thumb Grip.....	\$483
TPRIM11	<input type="checkbox"/> The Surge LT (All-In-One Oval with Gription™ Strip, Smaller Profile) (Not available with 20"(451) Rear Wheels.).....	\$494
TPRIM20	<input type="checkbox"/> Q-Grip (Exclusive High-Friction, Anti-Slip Coating) (Not available with 20"(451) Rear Wheels.).....	\$394

### Additional Handrim Options

TPRIM5	<input type="checkbox"/> Plastic Coated – Black (Not available with Mag Rear Wheels.).....	\$168
TPRIM22	<input type="checkbox"/> Omit Handrims.....	NCO

### Tab Length

TPRIM7	<input type="checkbox"/> Long Tabs.....	STD
TPRIM8	<input type="checkbox"/> Short Tabs (Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.).....	\$37

### WHEEL LOCKS


TPWLK30	<input type="checkbox"/> Compact Push to Lock.....	NCO
TPWLK31	<input type="checkbox"/> Compact Pull to Lock.....	NCO

Mounting Location:

☐ Standard Location F ☐ Forward Mount Location G (Mounted on forward front of side panel.)

### SOLID SEAT PAN

TPSP1	Solid Seat Pan (Required Option).....	\$205
-------	---------------------------------------	-------



# TiLite Pilot Canadian Order Form

Effective May 2022 (New Items in Red)

CARBON RIGID, GROWABLE FRAME

Page 4 of 5

## SEAT CUSHION

- TPCSH2 ☐ Comfort Curve Cushion (Stretch-Air cover.).....\$142  
 TPCSH23 ☐ Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion...\$497  
 SCSH11 ☐ ROHO Single Compartment Low Profile Cushion.....\$711

## CUSTOM CONFIGURATION OPTIONS

ROHO OPTIONS: (Made to order  
options have a 2-3 week lead time)

	SEAT DEPTH (B)						
	10"	11"	12"	13"	14"	15"	
SEAT WIDTH (A)	8"	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order
	9"	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order
	10"	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order
	11"	Made to Order	Stocked	Stocked	Stocked	Made to Order	Made to Order
	12"	Made to Order	Stocked	Stocked	Stocked	Made to Order	Made to Order
	13"	Made to Order	Made to Order	Made to Order	Stocked	Stocked	Stocked
	14"	Made to Order	Made to Order	Made to Order	Stocked	Stocked	Stocked

## SIDE GUARDS: (Side guards required)

- TPSDG30 ☐ Built in Side Guard with Height Adjustable Fender (Not available with 20" rear wheels, .).....\$609  
 TPSDG31 ☐ Built in Side Guard WITHOUT Fender.....STD

## ARMRESTS (Armrests attach to built in Side Guards.)

- TPARM2 ☐ Height Adjustable Desk Arm w/ 10" Pad (Works with Fendered Side Guard Attachment Option.).....\$452  
 TPARM31 ☐ Height Adjustable Armrest.....STD

## UPPER BODY POSITIONING

- TPUBP1 ☐ Bodypoint® PivotFit™ Dynamic shoulder harness with cam buckle mounts, rear-pull .....\$284  
☐ Extra Small ☐ Small  
 TPUBP2 ☐ Bodypoint® Monoflex™ chest support.....\$273  
☐ Extra Small ☐ Small

## PELVIC POSITIONING BELTS

- TPBLT1 ☐ Velcro Adjustable, black nylon webbing seat belt.....\$74  
 TPBLT8 ☐ Q'Straint® WC-19 Wheelchair-Anchored Pelvic Belt (Recommended for use with Transit Tie Down.).....\$210  
 TPBLT30 ☐ Bodypoint® 1" wide two-point padded hip belt with swivel buckle.....\$242  
 TPBLT16 ☐ Bodypoint® 1" wide Evoflex® belt with small push-button buckle.....\$252

## LOWER BODY POSITIONING

- TPCLF15 ☐ Bodypoint® Ankle Huggers®.....\$305  
☐ Extra Small ☐ Small

## ANTI-TIPS

- TPTIP1 ☐ Rear Aluminum - Black.....\$210  
 TPTIP7 ☐ Single Anti-tip.....\$168  
 Mounting Location: ☐ Right Side ☐ Left Side

## BACKPACK

- TPBCK3 ☐ Black Permobil Backpack.....\$100

## SPOKE GUARDS

- TPGRD30 ☐ Clear polycarbonate spoke guard with TiLite Pilot sticker pack (Not available with Mag Rear Wheels.).....\$142

## ACCESSORIES

- TPGS1 ☐ GRIP SolutionsLap Board (Lap board dimensions are 11" X 14").....\$58  
 TPCH1 ☐ Crutch Holder.....\$153

# TiLite Pilot Canadian Order Form

Effective May 2022 (New Items in Red)

CARBON RIGID, GROWABLE FRAME

Page 5 of 5

## GENERAL

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional.

**Get the latest order forms at permobilus.com.** Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. ROHO®, AGILITY®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q'Straint® is a registered trademark of Q'Straint Inc.

## NOTES

### QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

### ORDER THIS CHAIR