

## ROHO Back Systems

HCPCS	DESCRIPTION	ROHO® PRODUCTS
E2611	General use wheelchair back cushion, width less than 22 inches, any height, including any type mounting hardware.	ROHO® AGILITY® Minimum Contour Back System (foam only), ROHO Retroback®
E2617	Custom fabricated wheelchair back cushion, any size, including any type mounting hardware.	ROHO Custom Back, ROHO AGILITY Custom Back
E2613	Positioning wheelchair back cushion, posterior, width less than 22 inches, any height, including any type mounting hardware.	ROHO AGILITY Minimum Contour Back System (foam + lumbar)
E2615	Positioning wheelchair back cushion, posterior-lateral, width less than 22 inches, any height, including any type mounting hardware.	ROHO AGILITY Minimum Contour Back System, ROHO AGILITY Mid Contour Back System, ROHO AGILITY Max Contour Back System (foam with air insert + lumbar)
E2620	Positioning wheelchair back cushion, planar back with lateral supports, width less than 22 inches, any height, including any type mounting hardware.	ROHO AGILITY Minimum Contour Back System, ROHO AGILITY Mid Contour Back System, ROHO AGILITY Max Contour Back System (foam with air insert) ROHO AGILITY CARBON™ Minimum Contour Back System + Carbon Upgrade, ROHO AGILITY Mid Contour Back System + Carbon Upgrade



# AGILITY® Minimum Contour Back System Specifications

Depth, Width and Height	Depth: Foam 2.5" (6.5 cm) or Foam with Air 3" (7.5 cm) Width: 13" – 20" ( ) Height: 10", 13", 16" and 18"
Weight	Weight ranges do not include Back Cane Hardware Back Cane Hardware with Standard Bracket: 0.75 lbs. (0.3 kg) 14 in. x 10 in. (35.5 cm x 25.5 cm) 2½ lbs. (1.1 kg) 16 in. x 16 in. (40.5 cm x 40.5 cm) 4 lbs. (1.8 kg) 20 in. x 18 in. (51.0 cm x 45.5 cm) 5¼ lbs. (2.3 kg)
Weight Limit	Weight limit of 300 lbs. (136 kg) and must be properly sized to the individual
Limited Warranty	Back Shell and Hardware: 5 Years Air Insert, Cover and Foam Components: 18 months Custom AGILITY DRY FLOATATION Air Insert: 36 months
Available Options	F = Foam only FL = Foam Insert and Lumbar Support Pad A = Air Insert only AL = Air Insert and Lumbar Support Pad
Hardware Options	AGILITY Quick Release Hardware (1) or Fixed Hardware (2) (Fixed Hardware is NOT intended to be used on ROHO AGILITY Back Systems greater than 16 in. (40.5 cm) in height.)
Recommended for Individuals who	<ul style="list-style-type: none"> <li>• Require basic postural support</li> <li>• Have good trunk control</li> <li>• Need adjustability to enhance their range of motion and mobility</li> <li>• Have, or have potential for, shoulder or back pain</li> </ul>
U.S. Medicare Codes	E2611 (foam only, no lumbar support pad) E2613 (foam only, with 1" lumbar support pad) E2615 (air insert, with 1" lumbar support pad) E2620 (air insert, no lumbar support pad) E2617 Custom AGILITY (AG2CUSTC) – custom fabricated wheelchair back cushion, any size, including any type mounting hardware
Additional Product Features	<ul style="list-style-type: none"> <li>• Hardware designed for easy adjustment of angle, depth, and height – even with the client seated in the wheelchair to eliminate the need for transfers</li> <li>• Triangle bracket design allows 1 inch of rearward or forward adjustment and up to 4 inches in height adjustment</li> <li>• Height and width adjustment plate allows for client to adjust width without removing hardware from back shell</li> <li>• Single-point angle adjustment allows tilt adjustment up to 15° forward and 15° backward in 5° increments</li> </ul>



## Available Accessories:

Swing-Away Laterals and Fixed Laterals  
Headrest  
Headrest Adaptor Plate  
Privacy Shield

# AGILITY® Mid Contour Back System Specifications

Depth, Width and Height	Depth: Foam with Air 4" (10 cm) Width: 13" – 20" Height: 10", 13", 16", 18" and 20"
Weight	Weight ranges do not include Back Cane Hardware Back Cane Hardware with Standard Bracket: 0.75 lbs. (0.3 kg) 14 in. x 10 in. (35.5 cm x 25.5 cm) 2½ lbs. (1.1 kg) 16 in. x 16 in. (40.5 cm x 40.5 cm) 4 lbs. (1.9 kg) 20 in. x 20 in. (51.0 cm x 51.0 cm) 6 lbs. (2.8 kg)
Weight Limit	Weight limit of 300 lbs. (136 kg) and must be properly sized to the individual
Limited Warranty	Back Shell and Hardware: 5 Years Air Insert, Cover and Foam Components: 18 months Custom AGILITY DRY FLOATATION Air Insert: 36 months
Available Options	A = Air Insert only AL = Air Insert and Lumbar Support Pad
Hardware Options	Quick Release Hardware (1), Fixed Hardware (2) or Direct Mount Hardware (P) (Fixed Hardware is NOT intended to be used on ROHO AGILITY Back Systems greater than 16 in. (40.5 cm) in height.)
Recommended for Individuals who	<ul style="list-style-type: none"> <li>• Require postural support</li> <li>• Have minimal trunk support</li> <li>• Need adjustability to enhance their range of motion and mobility</li> <li>• Have, or have potential for, shoulder or back pain</li> </ul>
U.S. Medicare Codes	E2615 (air insert, with 1" lumbar support pad) E2620 (air insert, no lumbar support pad) E2617 Custom AGILITY (AG2CUSTC) – custom fabricated wheelchair back cushion, any size, including any type mounting hardware
Additional Product Features	<ul style="list-style-type: none"> <li>• Hardware designed for easy adjustment of angle, depth, and height – even with the client seated in the wheelchair to eliminate the need for transfers</li> <li>• Triangle bracket design allows 1 inch of rearward or forward adjustment and up to 4 inches in height adjustment</li> <li>• Height and width adjustment plate allows for client to adjust width without removing hardware from back shell</li> <li>• Single-point angle adjustment allows tilt adjustment up to 15° forward and 15° backward in 5° increments</li> </ul>



## Available Accessories:

Swing-Away Laterals and Fixed Laterals  
Headrest  
Headrest Adaptor Plate  
Privacy Shield

# AGILITY® Max Contour Back System Specifications

Depth, Width and Height	Depth: Foam with Air 7" (18 cm) Width: 13" – 20" Height: 16", 18" and 20"
Weight	Weight ranges do not include Back Cane Hardware Back Cane Hardware with Standard Bracket: 0.75 lbs. (0.3 kg) 14 in. x 16 in. (35.5 cm x 40.5 cm) 4¼ lbs. (2.0 kg) 20 in. x 20 in. (51.0 cm x 51.0 cm) 6 lbs. (2.8 kg)
Weight Limit	Weight limit of 300 lbs. (136 kg) and must be properly sized to the individual
Limited Warranty	Back Shell and Hardware: 5 Years Air Insert, Cover and Foam Components: 18 months Custom AGILITY DRY FLOATATION Air Insert: 36 months
Available Options	A = Air Insert only AL = Air Insert and Lumbar Support Pad
Hardware Options	Quick Release Hardware (1), Fixed Hardware (2) or Direct Mount Hardware (P) (Fixed Hardware is NOT available on ROHO AGILITY Back Systems greater than 16 in. (40.5 cm) in height.)
Recommended for Individuals who	<ul style="list-style-type: none"> <li>• Require significant postural support</li> <li>• Have minimal to maximum trunk support and control</li> <li>• Need adjustability to enhance their range of motion and mobility</li> <li>• Have, or have potential for, shoulder or back pain</li> </ul>
U.S. Medicare Codes	E2615 (air insert, with 1" lumbar support pad) E2620 (air insert, no lumbar support pad) E2617 Custom AGILITY (AG2CUSTC) – custom fabricated wheelchair back cushion, any size, including any type mounting hardware
Additional Product Features	<ul style="list-style-type: none"> <li>• Hardware designed for easy adjustment of angle, depth, and height – even with the client seated in the wheelchair to eliminate the need for transfers</li> <li>• Triangle bracket design allows 1 inch of rearward or forward adjustment and up to 4 inches in height adjustment</li> <li>• Height and width adjustment plate allows for client to adjust width without removing hardware from back shell</li> <li>• Single-point angle adjustment allows tilt adjustment up to 15° forward and 15° backward in 5° increments</li> </ul>



## Available Accessories:

Swing-Away Laterals and Fixed Laterals  
Headrest  
Headrest Adaptor Plate  
Privacy Shield