




## Retail order form - Sept. 2023

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.

### Range of adjustment



#### User dimensions (inches)

**S410 Size 1**

Height 71–107



#### Key user dimension: height

Select the appropriate stander by the user's overall height.  
Choose the model that allows for growth.

**Important:** User's weight must not exceed the maximum working load.

#### Item dimensions (inches)

**S410 Size 1**

Base length x width	32½ x 21½
Height of horizontal body support above floor	22–28
Sandal to knee height	6¾–11
Knee to hip height	6¾–10½
Sandal to top of headrest	31½–43
Width between trunk laterals	5½–10
Width between medial/lateral thigh supports	2½–4
Basic item weight PRONE (lb)	47
Basic item weight SUPINE (lb)	52
Max. working load (lb)	75

				S410 size 1		
* Frame				S411	\$2,995 <input type="checkbox"/>	
* Color				CHOOSE	Red <input type="checkbox"/>	
					Blue <input type="checkbox"/>	
* CHOOSE CONFIGURATION	Prone (with laterals, pelvic harness)			S412	\$1,795 <input type="checkbox"/>	
		Prone knee cuffs (pair)		S485	\$395 <input type="checkbox"/>	
	Supine (with laterals, knee strap)			S413	\$1,895 <input type="checkbox"/>	
		Head laterals (pair)		S497	\$395 <input type="checkbox"/>	
		Butterfly harness		S456	\$595 <input type="checkbox"/>	
		* Hip strap	CHOOSE	Standard	S472	\$595 <input type="checkbox"/>
				Rotation control	S476	\$595 <input type="checkbox"/>
	Multi-Position (with laterals, knee strap, pelvic harness)			S414	\$2,895 <input type="checkbox"/>	
		Head laterals		S497	\$395 <input type="checkbox"/>	
		Butterfly harness		S456	\$595 <input type="checkbox"/>	
		* Hip strap	CHOOSE	Standard	S472	\$595 <input type="checkbox"/>
				Rotation control	S476	\$595 <input type="checkbox"/>
		Prone knee cuffs (pair)		S485	\$395 <input type="checkbox"/>	
* Trunk strap			CHOOSE	Standard	S462	\$595 <input type="checkbox"/>
				Rotation control	S466	\$595 <input type="checkbox"/>
* Leg type			CHOOSE	Standard	S415	\$1,295 <input type="checkbox"/>
				Advanced	N/A	
Tray				S491	\$1,095 <input type="checkbox"/>	
Additional trunk laterals (pair)				S494	\$395 <input type="checkbox"/>	
Medial/lateral thigh supports (pair)				S482	\$395 <input type="checkbox"/>	
Additional medial/lateral thigh supports (pair)				S482	\$395 <input type="checkbox"/>	
Total				S410		

\* Required