



PDG Fuze T50n - AADL W654



Dealer Name*: _____ Account #*: _____

City, Province, Post Code: _____

Phone: _____ Email: _____ ATP: _____

Rep Name: _____ P.O. Number: _____ Quote Number: _____

Client Reference*: _____ Client Weight: _____

Special considerations that need to be addressed here (e.g., diagnosis): _____ *Required field

October 2024, Revised January 2026. Prices are in Canadian Dollars (CAD)

Wheelchair & Frame

Fuze T50 50° Manual Tilt Wheelchair Base Model	W654	\$ 2,685.00
Fuze T50n Narrow Frame Upgrade <i>Reduces overall width to seat width + 4".</i>		\$ 200.00

Frame Type

Fixed Width Frame, 300lbs Weight Capacity <i>Select a fixed seat width below.</i>		No charge
Fixed Width Frame, 400lbs Weight Capacity <i>Not compatible with transit tie down system. Select a fixed seat width below.</i>	Prior Approval	\$ 300.00
Adjustable Width Frame. 250lbs Weight Capacity <i>Select an adjustable seat width below.</i>		\$ 200.00

Frame Options

Transit Tie Down System (TTDS) <i>Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details. Tested to maximum weight capacity of 300lbs.</i>		No charge
--	--	-----------

Seat Measurements

Seat Width | Overall chair width is +4"

Fixed:	16"	17"	18"	19"		No charge
	20"	21"	22"		W111	\$ 250.00
	23"	24"			W111	\$ 500.00
Adjustable 16"-20":	16"	17"	18"	19"		No charge
	20"				W111	\$ 250.00

Seat Depth

Select Depth:	15"	16"	17"	18"	19"		No charge
	20"	21"	22"			W111	\$ 250.00

21"-22" seat depths are not adjustable and require straight stroller handle back. Not compatible with reclining or fold-down back.

Frame Finishes

Select Color:	High Gloss Black	Hyper Blue
---------------	------------------	------------

Armrests

Armrest Type

No Armrests No charge

Single Post, Height Adjustable Armrests No charge

10-14" Range 12"-16" Range

Armrest Position: Rear (Default) Front

If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.


Armpads

Regular Armpads No charge

Desk Length 2" x 10" Full Length 2" x 14"

Seat Height

Chairs with swing-away front rigging are set up with 3° of pre-tilt to avoid interference. This increases the actual seat to floor height at the front of the seat by up to 1" compared to the selected value. Easy-Out and standard plug-in front rigging are not affected.



Seat Pan

No Seat Pan No charge

Solid Aluminium Seat Pan No charge

Front Casters

Urethane Casters

5" x 1.25" Caster No charge

6" x 1.25" Caster No charge

8" x 1.25" Caster No charge

Rear Wheels

Composite Rear Wheel, 12" No charge

Rear Wheel Position: Rear (More Stable) Mid (Default) Front (Less Stable)

If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

Axles

Quick Release Rear Axles No charge

Fixed Rear Axles No charge

Rear Tires

Urethane Tires No charge

Handrims

No Handrims No charge

Wheel Locks

Attendant Wheel Locks No charge

Push to Lock Pull to Lock

Wheel Lock Options

Wheel Lock Extensions No charge

Anti-Tippers

Rear Anti-Tippers *Required safety option.*

No charge

Backrest

Back Post Style: Stroller Handle Back

Cane handle back not compatible with power tilt with manual override.

Back Post Type: Straight 8° Bend

Back Post Height: Short (17") Tall (23")

Measured from seat pan to handle.

Folding Style: Fold-Down

Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.

Seat to Back Angle: Set at 90° Set at 98° Set at 106° Set at 114° Set at 122°

Backrest Options

Angle/Height Adjustable Stroller Handle: *Not compatible with cane handle back.*

No charge

Button on Outside Button on Inside

Tilt Activation

Attendant Tilt Only

Attendant Tilt Only, Right Back Post

No charge

Attendant Tilt Only, Left Back Post

No charge

Tilt Limit | *Optional. If no tilt limit is selected the full tilt range will be available.*

Tilt Limit: Set at 20° Set at 30° Set at 40°

Front Rigging

No Front Rigging

No charge

Easy-Out Front Rigging | *Effortless insertion and removal. No hanger brackets remain on the frame.*

Fixed 70°

\$ 95.00

Swing Away Front Rigging | *Legacy design. Not compatible with Easy-Out.*

Swing Away 70°

No charge

Swing Away Elevating Legrests (ELR) *Includes calf pads.*

W106

\$ 375.00

1" Offset Hangers for Swing Away Front Rigging

\$ 95.00

Plug In Front Rigging | *Legacy design. Not compatible with Easy-Out.*

Plug In 70°

No charge

Footplates

Composite Footplates

\$ 0.00

Aluminium Angle Adjustable Footplates

W107

\$ 165.00

Extension Tubes | *Consider the height of your seat cushion.*

Standard Extension Tubes, Short

No charge

Standard Extension Tubes, Long

\$ 75.00

Footplate Height: Specify:

Ext	Short	11" - 16"	12" - 16"	4" - 9"	13" - 16"
	Long	14" - 22"	14" - 19"	12" - 16"	14" - 17"
		Standard Easy-Out	60°, 70°, ELR Swing-Away	90° Swing-Away	Standard Plug-In
Legrest					

Front Rigging Accessories

Length Adjustable Heel Loops

No charge

Standard Calf Strap

Basic
Seating

\$ 50.00

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes