

Manual wheelchairs

# TiLite CR1

permobil



### Key features

- Exceptionally lightweight frame
- $\pm 1$ " Centre of Gravity fine tuning adjustment
- Integrated, fold-down backrest for easy transport
- Individually made, just for you

### Custom profiled sideguard

- Supports backrest without need for backrest hinges
- Individually designed for tailored fit

[12.0096, 12.0116]




### Performance components

- Specially designed castor bearings for optimal manoeuvrability and weight savings
- PBO fibre and titanium reinforcements
- New carbon fibre castor wheels

### Dual-angle front end

- Designed to increase environmental access
- Optimised castor placement





## Component Colours

### Anodised Options



### COG Bracket & Trim Rings



## Specifications

<b>Transport weight</b>	3.9 kg
<b>Seat width</b>	13" – 18"
<b>Seat depth</b>	13" – 19"
<b>Front seat height</b>	17.5" – 21"
<b>Seat slope</b>	7°, 10°, 13°
<b>Front frame angle</b>	85°
<b>Seat back height</b>	10" – 17"

<b>Camber</b>	0°, 2°, 4°
<b>Centre of gravity adjustment</b>	2"
<b>Backrest angle</b>	80° - 101°
<b>Compatible with SmartDrive?</b>	Yes
<b>Weight limit</b>	105 kg
<b>Warranty</b>	Limited