

Document	Doc. Id	Version	Published Date	Page
Permobil Crash Test Certificare	PD0084	3	2020-07-03	1(3)
Issued by (Name)	Henrik Larsson		Date	2020-07-01
Reviewed by (Name)	Marita Brundin		Date	2020-07-03
Approved by (Name)	Johan Holmström		Date	2020-07-03

Crash test Certificate

To whom it may concern, Permobil AB hereby confirms that the wheelchairs listed below have been verified concerning crashworthiness in accordance with ISO 7176-19. This implies that it's safe to be transported in a vehicle when using the combination of listed wheelchair and fastening system.

Manufacturer: Permobil AB, Per Uddéns Väg 20, 861 23 Timrå, Sweden

The wheelchairs fulfill all applicable design and performance requirements according to ISO 7176-19 as described in the table below.

This certificate is valid assuming that no changes are made to the wheelchair after delivery from Permobil AB.

Wheelchair / Product	Fastening System	User weight Kg	ISO 7176-19 Wheelchairs for use as seats in motorvehicles
C350 Corpus 3G	Four point tie-down	136	Pass
F3 Corpus	Four point tie-down	150	Pass
F3 Corpus	Dahl docking station	136	Pass
F5 Corpus	Four point tie-down	150	Pass
F5 Corpus	Dahl docking station	136	Pass
F5 Corpus VS	Four point tie-down	136	Pass
F5 Corpus VS	Dahl docking station	136	Pass
K300 PS Jr	Four point tie-down	76	Pass
K450 MX	Four point tie-down	57	Pass
M1 Agility backrest	Four point tie-down	136	Pass
M1 Solid backrest	Four point tie-down	136	Pass
M300/M400 Corpus 3G HD	Four point tie-down	205	Pass
M300/M400 PS Junior	Four point tie-down	75	Pass
M3 Corpus	Four point tie-down	150	Pass
M3 Corpus	Dahl docking station	136	Pass
M5 Corpus	Four point tie-down	150	Pass
M5 Corpus	Dahl docking station	136	Pass
X850 Corpus 3G	Four point tie-down	136	Pass
F3 Corpus Agility Direct Mount	Dahl docking station	136	Pass
F3 Corpus Agility Direct Mount	Four point tie-down	136	Pass
F5 Corpus Agility Direct Mount	Dahl docking station	136	Pass

F5 Corpus Agility Direct Mount	Four point tie-down	136	Pass
F5 Corpus VS Agility Direct Mount	Dahl docking station	136	Pass
F5 Corpus VS Agility Direct Mount	Four point tie-down	136	Pass
M3 Corpus Agility Direct Mount	Dahl docking station	136	Pass
M3 Corpus Agility Direct Mount	Four point tie-down	136	Pass
M5 Corpus Agility Direct Mount	Dahl docking station	136	Pass
M5 Corpus Agility Direct Mount	Four point tie-down	136	Pass