

BODILINK® HEAD SUPPORT BASIC ORDER FORM



HCPCS CODING: E0955 - CUSHIONED HEADREST
HCPCS CODING: E1028 - GT REMOVABLE HARDWARE

Company Name/ACCT #: _____
P.O. Number: _____
Requested By: _____
Phone: _____ Fax: _____
Email: _____
Ship To: _____
Patient Reference: _____



GT Lateral Hardware
Shown on Acta-Back®

For best results, do not fill out in your browser. Interactive form should be completed using Adobe Reader after saving to your local drive. Then email or print and fax to Customer Support at orders.comfort@permobil.com or 1.406.522.8563.

DO NOT SEND PATIENT PROTECTED HEALTH INFORMATION. IT IS NOT NEEDED TO MAKE THE PRODUCT YOU ARE REQUESTING.

1 BODILINK® HEAD SUPPORT PAD

- You may make one selection per column.
- Choose pad size in 1.1.
- Choose a cover material in 1.2.

1.1 BASIC PAD SIZE

HCPCS CODING: E0955 MSRP \$136.00

Length (L) refers to the actual size dimension of the support from bottom to top edge.

	← W →	
L ↓	6W x 5L	BL-HSP3-6W5L
L ↓	8W x 5L	BL-HSP3-8W5L
L ↓	10W x 6L	BL-HSP3-10W6L
L ↓	12W x 6L	BL-HSP3-12W6L
L ↓	14W x 6L	BL-HSP3-14W6L
L ↓	18W x 6L	BL-HSP3-18W6L



Basic head support pads can be formed to any shape using the moldable steel insert inside the pad.

1.2 BASIC PAD COVER MATERIAL

COMFORT-TEK FOR FLUID PROTECTION & AN EASILY CLEANED SURFACE.

STRETCH-AIR FOR PATIENT COMFORT & HEAT DISSIPATION.

	Comfort-Tek (COMFORT-TEK)
	Stretch-Air (STRETCH-AIR)

2 BODILINK™ HEAD SUPPORT HARDWARE

- You may make one selection per column in either.
- Choose which mounting block option you require in 2.2.

2.1 GT HARDWARE, Removable

HCPCS E1028

RANGE OF MOTION

GT MODEL	TOTAL # OF LINKS	TOTAL LINK LENGTH	MAX. ANTERIOR REACH	MAX. POSTERIOR REACH	MAX. LATERAL REACH
LAT	5	10.5"	10"	2"	7.5"
MID	3	8.5"	8"	2"	N/A

STYLE	MSRP	PART NUMBER
LATERAL	\$240.00	(BL-HS-GT3RMLAT)
MID-LINE	\$201.00	(BL-HS-GT3RMMID)


GT HARDWARE, LATERAL



GT HARDWARE, MID-LINE



2.2 HEAD SUPPORT MOUNTING OPTION




MOUNTING BLOCK OR BRACKET WITH MOUNTING BLOCK	MSRP	PART NUMBER
POSTERIOR A MOUNTING BLOCK	\$0.00	(HS-MB-POSTA)
POSTERIOR B MOUNTING BLOCK	\$0.00	(HS-MB-POSTB)
ANTERIOR A MOUNTING BLOCK	\$0.00	(HS-MB-ANTA)
QUANTUM A MOUNTING BLOCK	\$90.00	(HS-MB-PW1A)
PERMOBIL A MOUNTING BLOCK	\$90.00	(HS-MB-PW2A)
PERMOBIL B MOUNTING BLOCK	\$90.00	(HS-MB-PW2B)
*ACTA-RELIEF™ 14W BRACKET WITH MOUNTING BLOCK	\$72.00	(HS-MB-AR14W)
*ACTA-RELIEF™ 16W BRACKET WITH MOUNTING BLOCK	\$72.00	(HS-MB-AR16W)
*ACTA-RELIEF™ 18W BRACKET WITH MOUNTING BLOCK	\$72.00	(HS-MB-AR18W)
*ACTA-RELIEF™ 20W BRACKET WITH MOUNTING BLOCK	\$72.00	(HS-MB-AR20W)

*For Acta Relief™, includes mounting bracket and POSTA block.
 If you already have an Acta-Relief™ head support bracket, order POSTA Block

3 HEAD SUPPORT ADDITIONAL MOUNTING HARDWARE

3.1 EXTRA MOUNTING HARDWARE



MOUNTING BLOCK OR BRACKET	MSRP	PART NUMBER
POSTERIOR A MOUNTING BLOCK	\$55.00	(BL-HS-MBPOSTA)
POSTERIOR B MOUNTING BLOCK	\$55.00	(BL-HS-MBPOSTB)
ANTERIOR A MOUNTING BLOCK	\$55.00	(BL-HS-MBANTA)
QUANTUM A MOUNTING BLOCK	\$95.00	(BL-HS-MBPW1A)
PERMOBIL A MOUNTING BLOCK	\$95.00	(BL-HS-MBPW2A)
PERMOBIL B MOUNTING BLOCK	\$95.00	(BL-HS-MBPW2B)
*ACTA-RELIEF™ 14W BRACKET	\$160.00	(BL-HS-AR14W)
*ACTA-RELIEF™ 16W BRACKET	\$160.00	(BL-HS-AR16W)
*ACTA-RELIEF™ 18W BRACKET	\$160.00	(BL-HS-AR18W)
*ACTA-RELIEF™ 20W BRACKET	\$160.00	(BL-HS-AR20W)

*Bracket only, does not include mounting block.

ADDITIONAL	MSRP	EXTRA PART
1/4" SPACER PLATE	\$39.00/each	(BL-HS-SPCRA)