



F3 Corpus - PP



Dealer Name*:

Account #*:

City, State, Zip:

Phone:

Email:

ATP:

Rep Name:

P.O. Number:

Quote Number:

Client Reference*:

Client Weight:

*Required field

Special considerations that need to be addressed here (e.g., diagnosis):

Chair Model and Wheels

Front-wheel drive F3 power base includes: Corpus seat frame, ComfortRide suspension, Enhanced Steering Performance (ESP), 4-point tie-down hardware, VoltPro charger, and R-net Power Platform 120A electronics with Permobil Connect.

Choose Base

11716	F3 Base Corpus - F3 SPO <i>FWD F3 power base set up for one power seat function, max speed: 6 mph, weight capacity: 300 lb.</i>	K0856	\$ 8,768.00
108993-99-0	F3 Base Corpus - F3 MPO <i>FWD F3 power base set up for multiple power seat functions, max speed: 6 mph, weight capacity: 300 lb.</i>	K0861	\$ 8,768.00
11717	F3 Base Corpus - F3 SPO HD <i>FWD F3 power base set up for one power seat function, max speed: 6 mph, weight capacity: 330 lb.</i>	K0858	\$ 10,343.00
11718	F3 Base Corpus - F3 MPO HD <i>FWD F3 power base set up for multiple power seat functions, max speed: 6 mph, weight capacity: 330 lb.</i>	K0862	\$ 10,343.00

User manual language English Spanish

Support Wheel Options | Choose one

107330-99-0	Support Wheels for F3	No charge
11313	Omit Front Support Wheels <i>User Weight must not exceed 220lbs. Max speed will be limited to 5mph on the F3 Corpus and 6mph on the F5 Corpus. Additional restrictions may apply based on chair configuration.</i>	No charge

Choose Drive Tires

108414-99-0	Flat Free Tire 3"W x 14"T Black - Standard Tread <i>Includes flat free standard tread caster wheels.</i>	No charge
108411-99-0	Pneumatic Tire 3"W x 14"T Black - Standard Tread <i>Includes flat free standard tread caster wheels.</i>	No charge
11333	Flat Free Tire 3"W x 14"T Black - Aggressive Tread <i>Includes flat free aggressive tread caster wheels and mud guards.</i>	\$ 514.00
108412-99-0	Pneumatic/Snow Tire 3"W x 14"T Black - Aggressive Tread <i>Not recommended for year-long use and must be ordered along with standard drive wheels. Does not include caster wheels.</i>	\$ 514.00
11334	Mud Guard Kit F3	\$ 45.00

Base Options

Choose Base Color | Choose both shroud and rim color

Standard Shroud Colors	Aqua	Sky	Violet	Flamingo	Chili	Canary	Spring	No charge
	Arctic	Sunset	Midnight	Gold	Bronze	Platinum		
Premium Shroud Colors	Carbon Fiber		Mossy-Oak		Liquid Fusion		Water Lily	\$ 223.00
	Midnight Patriot		American Liberty					
Standard Rim Colors	Aqua	Sky	Violet	Flamingo	Chili	Canary	Spring	No charge
	Arctic	Sunset	Midnight	Gold	Bronze	Platinum		
Premium Rim Colors	Carbon Fiber		Mossy-Oak		Liquid Fusion		Water Lily	\$ 55.00
	Midnight Patriot		American Liberty					

Batteries | Choose one

108999-99-0	Batteries, Group 24 - Pair	E2363	\$ 1,095.00
109183-99-0	Batteries, Group 34 - Pair	E2359	\$ 1,048.00

Optional Items

11732	LED Light Kit for F3 - Personalized Color Changing Option <i>Color selection requires MyPermobil activation. Not programmed for hazard lights or turn signals.</i>		No charge
109382-99-0	LED Light Kit for F3 - Extended Safety Option <i>Includes hazard lights and turn signals in addition to color changing option. Color selection requires MyPermobil activation.</i>		\$ 1,087.00
109170-99-0	Essentials Box with Reflectors		\$ 535.00
109169-99-0	Essentials Box with Lights <i>Must also select LED Light Kit.</i>		\$ 872.00
10715	Compact Ventilator Tray - Corpus Seat <i>Custom option. 13"T, adjustable width between 8" - 13", and adjustable depth between 2" - 6.5". Additional drive and seating restrictions may apply.</i>	E1029	\$ 1,755.00
10787	Battery Tray for Laptop Ventilator Tray <i>Custom option. 11 3/8" W x 6 1/2" H x 1 5/8" D.</i>		\$ 294.00
10788	O2 Holder for Laptop Ventilator Tray <i>Custom option. Mounts to left or right side of compact laptop ventilator tray. 7"T, adjustable width between 3.5" - 7", and adjustable depth between 3.75" - 5".</i>	E2208	\$ 373.00

Joystick Options

Choose joystick type, handle, and mount.

NOTE: For additional input devices, please refer to the Alternative Drive Systems section or submit a custom request.

Joystick Type

108848-99-0	Color Joystick Module (CJSM) - Power Platform <i>High-resolution 3.5" VGA color display.</i>	No charge
109588-99-0	Permobil Joystick Module (PJSM) - Power Platform <i>Streamlined profile with 2.8" color display.</i>	\$ 1,182.00
107829-99-0	3.5mm Y-Splitter Cable for PJSM/CJSM <i>Allows 2 mechanical switches to control up to 4 separate functions through assignable buttons programming.</i>	\$ 39.00

Joystick Handle

10045	Cone Shaped Handle - Standard	No charge
100960-99-0	Chin Cup Handle	E2324 \$ 139.00
100961-99-0	T Handle	E2323 \$ 196.00
100962-99-0	Large Ball	E2323 \$ 139.00
100963-99-0	Softball	E2323 \$ 183.00
100964-99-0	Mushroom Handle	E2323 \$ 183.00
100965-99-0	Stick Handle	E2323 \$ 183.00
105442-99-0	1.75" Black Foam Ball	E2323 \$ 139.00
105846-99-0	Bodypoint 3" U-Shaped Handle w/ Flex-Shaft	E2323 \$ 183.00
105847-99-0	Bodypoint 4" U-Shaped Handle w/ Flex-Shaft	E2323 \$ 183.00
105848-99-0	Bodypoint Rubber Dome Handle	E2323 \$ 183.00

Joystick Mount

Fixed	Mount on Left	Mount on Right	No charge
Swing-Away	Mount on Left	Mount on Right	E1028/ E1032 \$ 409.00
Slimline Retractable	Mount on Left	Mount on Right	E1028/ E1032 \$ 550.00
11232	Adjustable Height Panel Bracket <i>Provides 1.77" of vertical adjustment and 25° of angle adjustment.</i>		\$ 299.00
105210-99-0	Mount Standard Joystick on Corpus Back Hinge <i>Can be chosen instead of attendant control.</i>		\$ 376.00
	Mount on Left	Mount on Right	
108268-99-0	Bodypoint Tri-Lock Midline Mount Assembly <i>Must select specific joystick attachment separately.</i>		\$ 702.00
	Mount on Left	Mount on Right	
108280-99-0	Bodypoint Tri-Lock Mount for R-net Joystick		\$ 292.00

Other Options

10411	Expandable Controller - R-net <i>Required when ordering power tilt, power recline, and power elevating legrests together OR when ordering an alternative drive control.</i>	E2377	\$ 764.00
10513	Harness for Expandable Controller <i>Must be ordered with any expandable controller.</i>	E2313	\$ 678.00
109532-99-0	mo-vis Scoot Attendant Control <i>Fully proportional handlebar with configurable thumb throttles, access to wheelchair functions and status and a foldable, adjustable mounting bracket. Depending on chair configuration, custom mounting may be required.</i>	E2331	\$ 4,200.00
105408-99-0	Attendant Control R-net	E2331	\$ 1,499.00
110286-99-0	Remote Stop System - F/M Series <i>Consists of a 315MHz radio frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair.</i>		\$ 1,381.00
10582	Output Module (IOM) for R-net Control System <i>Enables user to control an external device (e.g., communication device, door opener, etc.) through the drive input (e.g., JSM, head array, etc.). May require specialty cable or wiring to be ordered separately.</i>		\$ 1,412.00

Power Seating

Please Note: Power seat adjustability may be limited to prevent contact with base and/or backrest-mounted items.

Adjustable Seat Height | *If ordering power adjustable seat height, elevator travel is used to set seat heights greater than 17.5". Please note: Wheelchair configuration can affect seat-to-floor height.*

11738	ActiveHeight Power Adjustable Seat Height - 12" Travel <i>Provides 12" of power seat elevation</i>	E2298	\$ 3,786.00
11117	Seat Tube Fixed		No charge
Seat-to-Floor Height	17.5" 18.5" 19.5"		

Posterior Tilt/Recline Packages | Choose one

11726	Corpus Seat w/ No Tilt & 180° Power Recline	E1004	\$ 9,555.00
11727	Corpus Seat w/ 50° Tilt	E1002	\$ 8,525.00
11728	Corpus Seat w/ 50° Tilt & 180° Power Recline	E1007	\$ 13,492.00

ActiveReach Packages | *Chest support belt is strongly recommended. Power legrest recommended for full functionality. Seating and drive capabilities may be limited in certain configurations.*

11141	5° ActiveReach		No charge
107090-99-0	10° ActiveReach <i>ActiveHeight and power tilt required.</i>		No charge
107091-99-0	20° ActiveReach <i>ActiveHeight, power tilt/recline and Corpus VS elevating legrests required.</i>		\$ 2,967.00
11100	30° ActiveReach w/Transfer <i>ActiveHeight, power tilt/recline and Corpus VS elevating legrests required.</i>		\$ 3,335.00

Seat Function Control Options | ICS switchbox availability may be limited to certain configurations.

10418	Seat Function Control Switch for R-net <i>Does not include operation of power seat functions through drive control.</i>		No charge
10419	Single Seat Function Control Kit for R-net <i>Allows operation of a single power seat function through drive control.</i>	E2310	\$ 2,110.00
10585	Multiple Seat Function Control Kit for R-net <i>Allows operation of a two or more power seat functions through drive control.</i>	E2311	\$ 3,994.00
102798-99-0	Push Buttons on ICS Switchbox		No charge
102799-99-0	Toggle Switches on ICS Switchbox		\$ 357.00

Switchbox Mounting

10626	Mount ICS Switchbox onto R-net Joystick Bracket		No charge
-------	--	--	-----------

104651-99-0	Mount ICS Switchbox onto Armrest Mount on Left Mount on Right	No charge
106256-99-0	Mount ICS Switchbox on Back <i>No charge when ordering Omni2 Control Module.</i> Mount on Left Mount on Right	\$ 172.00
105409-99-0	Do Not Install ICS Switchbox, Send as Part	No charge

Power Seating Options

104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Accommodates up to 8 momentary switches.</i>	\$ 1,396.00
110441-99-0	Display Memory Seating <i>With this selection, three memory seat positions will be included (M home, M1, M2). Individual customization and additional memory positions are available through QuickConfig programming.</i>	No charge
110156-99-0	Enable Independent Repositioning Mode (IRM) <i>With this selection, a sequence for IRM (X3) will be included. IRM allows for sequential movement of power seat functions into a saved posterior position and returns to a saved upright position.</i>	No charge

Backrest Options

Installed backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Accessible seat depth will be dependent upon backrest selection.

Corpus Backrest Kits

105067-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No charge
105066-99-0	Fixed Backrest Kit - No Mechanical Shear	\$ 44.00

Corpus Backrests | *Corpus Ergo Backs include lateral wedges and foam lumbar pads. Solid Backs do not include foam lumbar pads.*

	20" Tall <i>Power Recline will be limited to 150°</i>	23" - 28" Tall <i>Height adjustable in 1" increments</i>		
Corpus Ergo (C3GERGBCA)				
14" W	14 x 20	14 x 23-28	E2620	\$ 1,192.00
16" W	16 x 20	16 x 23-28	E2620	\$ 1,192.00
18" W	18 x 20	18 x 23-28	E2620	\$ 1,192.00
20" W	20 x 20	20 x 23-28	E2620	\$ 1,192.00
Solid Corpus Ergo				
16" W	16 x 20	16 x 23-28		No charge
18" W	18 x 20	18 x 23-28		No charge
20" W	20 x 20	20 x 23-28		No charge

Set Back Height:	20"	23"	24"	25"	26"	27"	28"
-------------------------	-----	-----	-----	-----	-----	-----	-----

Corpus Backrest Options

11335	Corpus Ergo Back - Leatherette	No charge
11336	Corpus Ergo Back - Stretch-Air	No charge
10725	Split Backrest Cushion - Corpus <i>Custom option.</i>	\$ 377.00
10786	Tailored Backrest Cushion - Corpus <i>Custom option. Sides of cushion are trimmed to provide additional clearance for lateral trunk supports.</i>	\$ 141.00
11339	Foam Lumbar Pads for Solid Back - Corpus	\$ 45.00
109029-99-0	BodiLink Premium 3" x 5" Lateral Supports - Pair <i>Must also select hardware.</i>	E0956 \$ 446.00
109030-99-0	BodiLink Premium 4" x 6" Lateral Supports - Pair <i>Must also select hardware.</i>	E0956 \$ 446.00
109031-99-0	Swing-Away Adj BodiLink Lateral Hardware Corpus - Pair	E1028/ E1034 \$ 625.00

ROHO Backrests | Includes backrest shell, cushion, air inserts, UT rail for headrest and Direct Mount adapter.

	16" Tall Installed height between 22" – 25"	18" Tall Installed height between 22" – 26"	20" Tall Installed height between 24" – 27"		
ROHO AGILITY Max Contour					
14" W	14 x 16	14 x 18	14 x 20	E2620	\$ 874.00
16" W	16 x 16	16 x 18	16 x 20	E2620	\$ 874.00
18" W	18 x 16	18 x 18	18 x 20	E2620	\$ 874.00
20" W	20 x 16	20 x 18	20 x 20	E2620	\$ 874.00
Set Back Height:	22" 23" 24" 25" 26" 27"				

Comfort Backrests | Includes backrest shell, cushion and Compass Power Mount.

Note: User weight may be limited based upon the size of the backrest selected.

	14" Tall Installed height between 17" – 23"	16" Tall Installed height between 19" – 25"	18" Tall Installed height between 21" – 27"	20" Tall Installed height between 23" – 29"	
Acta-Back Stretch-Air					
14" W		14 x 16	14 x 18		E2613 \$ 668.00
16" W		16 x 16	16 x 18	16 x 20	E2613 \$ 751.00
18" W		18 x 16	18 x 18	18 x 20	E2613 \$ 751.00
20" W			20 x 18	20 x 20	E2613 \$ 775.00
Acta-Embrace Comfort-Tek					
12" W	12 x 14				E2613 \$ 694.00
14" W		14 x 16	14 x 18		E2613 \$ 694.00

Comfort Backrest Options

11714	BodiLink Premium 3" x 5" Lateral Supports - Pair Must also select hardware.	E0956	\$ 446.00
11715	BodiLink Premium 4" x 6" Lateral Supports - Pair Must also select hardware.	E0956	\$ 446.00
109903-99-0	BodiLink GT2 Swing-Away Lateral Hardware - Pair Must select with BodiLink lateral trunk supports on Comfort backrests.	E1028/ E1034	\$ 625.00

Aftermarket Backrest Options

105782-99-0	Tubular Backrest Frame - Corpus (1" Canes) Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Up to several inches of seat depth may be lost when mounting this system.		\$ 557.00
106072-99-0	Direct Backrest Frame - Corpus Used to mount backrests ranging in width from 14" to 20". Some seat depth may be lost when mounting this system.		\$ 557.00

Seat Sizing

Please select both width and depth options below.

Seat Pan Width | Choose one

Set at:	17" W	19" W	21" W	23" W	No charge
----------------	-------	-------	-------	-------	-----------

Seat Pan Depth | Choose one

Kit 1, Set at:	14"	15"	16"	No charge
Kit 2, Set at:	17"	18"	19"	No charge
Kit 3, Set at:	20"	21"	22"	No charge

Seat Cushion Options

Corpus Ergo Seat Cushions | *Corpus Ergo Air cushion features ROHO Dry Floatation technology.*

	15" D	17" D	19" D	21" D	23" D		
Corpus Ergo <i>Leatherette</i>							
17" W	17 x 15	17 x 17	17 x 19	17 x 21		E2605	\$ 632.00
19" W	19 x 15	19 x 17	19 x 19	19 x 21	19 x 23	E2605	\$ 632.00
21" W		21 x 17	21 x 19	21 x 21	21 x 23	E2605	\$ 632.00
23" W			23 x 19	23 x 21	23 x 23	E2605	\$ 632.00
Corpus Ergo <i>Stretch-Air</i>							
17" W		17 x 17	17 x 19	17 x 21		E2605	\$ 632.00
19" W		19 x 17	19 x 19	19 x 21	19 x 23	E2605	\$ 632.00
21" W		21 x 17	21 x 19	21 x 21	21 x 23	E2605	\$ 632.00
23" W			23 x 19	23 x 21	23 x 23	E2605	\$ 632.00
Corpus Ergo Air <i>C3GRHOSCA</i>							
17" W		17 x 17	17 x 19	17 x 21	17 x 23	E2622	\$ 667.00
19" W		19 x 17	19 x 19	19 x 21		E2622	\$ 667.00
21" W			21 x 19	21 x 21		E2622	\$ 667.00

ROHO Neoprene Seat Cushions

	8 cells 14.75" D	9 cells 16.5" D	10 cells 18.25" D	11 cells 20" D	12 cells 21.75" D	13 cells 23.5" D		
Quadro Select <i>High Profile</i>								
9 cells / 16.5" W	9 x 8	9 x 9	9 x 10	9 x 11	9 x 12		E2624	\$ 532.00
10 cells / 18.25" W		10 x 9	10 x 10	10 x 11	10 x 12	10 x 13	E2624	\$ 532.00
11 cells / 20" W		11 x 9	11 x 10	11 x 11	11 x 12	11 x 13	E2624	\$ 532.00
12 cells / 21.75" W			12 x 10	12 x 11	12 x 12	12 x 13	E2624	\$ 532.00
13 cells / 23.5" W		13 x 9	13 x 10	13 x 11	13 x 12	13 x 13	E2625	\$ 592.00
Single Compartment <i>High Profile, Sensor-Ready (* = non-SR)</i>								
9 cells / 16.5" W	9 x 8	9 x 9	9 x 10	9 x 11	9 x 12*		E2622	\$ 532.00
10 cells / 18.25" W		10 x 9	10 x 10	10 x 11	10 x 12	10 x 13*	E2622	\$ 532.00
11 cells / 20" W		11 x 9	11 x 10	11 x 11	11 x 12	11 x 13*	E2622	\$ 532.00
12 cells / 21.75" W			12 x 10	12 x 11	12 x 12	12 x 13*	E2622	\$ 532.00
13 cells / 23.5" W			13 x 10	13 x 11			E2623	\$ 592.00

ROHO Hybrid Seat Cushions

	14" D	15" D	16" D	17" D	18" D	20" D		
Hybrid Select								
14" W	14 x 14		14 x 16				E2624	\$ 585.00
15" W		15 x 15	15 x 16	15 x 17			E2624	\$ 585.00
16" W		16 x 15	16 x 16	16 x 17	16 x 18	16 x 20	E2624	\$ 585.00
17" W			17 x 16	17 x 17	17 x 18	17 x 20	E2624	\$ 585.00
18" W			18 x 16	18 x 17	18 x 18	18 x 20	E2624	\$ 585.00
20" W			20 x 16	20 x 17	20 x 18	20 x 20	E2624	\$ 585.00

Comfort Seat Cushions | *User weight may be limited based upon the size of the cushion selected.*

	14" D	16" D	18" D	20" D	22" D	
Embrace Anti-Thrust <i>Comfort-Tek</i>						
12" W	12 x 14					E2607 \$ 465.00
14" W	14 x 14	14 x 16				E2607 \$ 465.00
M2 Anti-Thrust <i>Comfort-Tek</i>						
14" W	14 x 14	14 x 16				E2607 \$ 474.00
16" W		16 x 16	16 x 18	16 x 20		E2607 \$ 474.00
18" W		18 x 16	18 x 18	18 x 20	18 x 22	E2607 \$ 474.00
20" W		20 x 16	20 x 18	20 x 20	20 x 22	E2607 \$ 474.00
22" W			22 x 18	22 x 20	22 x 22	E2608 \$ 508.00

Cushion Options

107505-99-0	ROHO Smart Check	\$ 185.00
11113	Incontinence Cover for Ergo Cushion <i>Sent as a part.</i>	\$ 203.00
11173	Incontinence Cover for ROHO Cushion <i>Sent as a part.</i>	\$ 82.00
105118-99-0	Omit Seat Cushion	No charge

Armrest Options

Corpus armrests are infinitely adjustable from 8" - 15" (measured from the seat pan to the top of the armrest pad); however, minimum height may be limited for certain armrest pad options.

Armrest Set-Up | *Width must be wider than backrest to ensure proper fit.*

106543-99-0	Narrow Armbar Kit (13"W - 19"W) <i>For 13" width, armrests will be mounted to inside of bracket and custom will be required for VS chest support.</i> <i>Set width between armrests to:</i> 13" 15" 17" 19"	No charge
106544-99-0	Wide Armbar Kit (19"W - 23"W) <i>Set width between armrests to:</i> 19" 21" 23"	No charge
105069-99-0	Height Adjustable Armrests - NC <i>No charge when ordering power tilt and/or recline. Set armrest height to:</i> 8" 9" 10" 11" 12" 13" 14" 15"	No charge
105081-99-0	Height Adjustable Armrests <i>Set armrest height to:</i> 8" 9" 10" 11" 12" 13" 14" 15"	\$ 391.00
105071-99-0	Medial Swivel <i>Medial swivel cannot be enabled when selecting fixed option.</i> Left Fixed Left Enabled Right Fixed Right Enabled	No charge

Choose Armrest Pads | *Must select left and right. Priced per side.*

CL Pads	Corpus Leatherette						No charge
	13" Left	16" Left	18" Left	13" Right	16" Right	18" Right	
CF Pads	Corpus Black Fabric						No charge
	13" Left	16" Left	18" Left	13" Right	16" Right	18" Right	
CG Pads	Corpus Gel						\$ 202.00
	13" Left	16" Left	18" Left	13" Right	16" Right	18" Right	
NL Pads	2" Narrow (Custom)						No charge
	14" Left	14" Right					
ATA	Arm Trough Adapters						\$ 377.00
	1" Pair	7/8" Pair	Install		Send as Parts		

Optional Items | *Priced per side/support. Not available with narrow armrest pads.*

Armrest Pouch	Left	Right				\$ 267.00
Armrest Supports	Left Outer	Left Inner	Right Outer	Right Inner		\$ 385.00
Elbow Supports (Must also select hardware)	Left Elbow Support	Right Elbow Support				\$ 191.00
Elbow Support Hardware	Left	Right				\$ 280.00

Legrest Options

Legrest Type | *Complete column for legrest type chosen*

	Corpus <i>Includes mechanical articulation Footplate range is 6" - 17"</i>	Corpus VS <i>Only available with > or = 20° ActiveReach; Includes power articulation Footplate range is 9" - 17"</i>		
Type	Power	Power	E1012, K0040	\$ 4,102.00
	Manual			No charge
Footplate Style	2pc Wide - 9"D x 7.5"W	2pc Wide - 9"D x 7.5"W		No charge
	2pc Narrow - 9"D x 6"W	2pc Narrow - 9"D x 6"W		No charge
	1pc Narrow - 9.75"D x 12.5"W	1pc Wide - 9"D x 16.5"W		\$294.00/ \$317.00
Footplate Height <i>Height measured from seat pan with legrest at 90°.</i>	Set footplate height at: (Round to nearest inch)	Set footplate height at: (Round to nearest inch)		
Calf Supports <i>*Not compatible with Swing-Away knee supports ^Obstructs headlights</i>				
B Pads	3.5"W x 4.5"T	3.5"W x 4.5"T		No charge
J Pads	5"W x 6.5"T	5"W x 6.5"T		No charge
A Pads	5.5"W x 6"T	5.5"W x 6"T*		No charge
H Pads (curved)	7"W x 7.5"T^	7"W x 7.5"T*^		\$ 419.00
D Pads	8"W x 8"T^	8"W x 8"T*^		\$ 419.00
	Omit Calf Supports	Omit Calf Supports (Note: Swing-away knee support bases include fabric covers)		
Knee Support/ Abductor <i>Choose one for Corpus VS Legrest.</i>		Single Post	E0957	\$ 292.00
		Swing-Away - Large	E0957	\$ 475.00
		Swing-Away - Small	E0957	\$ 475.00
Knee Support Hardware <i>Required with Knee Support/Abductor</i>		Adj/Removable Knee Support Hardware	E1028	\$ 306.00

Footplate Options

106530-99-0	Stealth Gel Pads - Narrow Footplates - Pair	\$ 470.00
106225-99-0	Stealth Gel Pads - Wide Footplates - Pair	\$ 481.00

Accessories

Belt Options | *Must select one of the following*

104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No charge
104520-99-0	Bodypoint Padded Hip Belt - Medium	\$ 167.00
104514-99-0	Bodypoint Non-Padded Hip Belt - Large	\$ 142.00
104516-99-0	Bodypoint Padded Hip Belt - Large	\$ 223.00
106531-99-0	Bodypoint Padded Hip Belt - Extra Large	\$ 265.00
100519-99-0	Positioning Belt Retractable	\$ 332.00
100600-99-0	Positioning Belt Retractable Long	\$ 345.00
108267-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Medium	\$ 224.00
108266-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Large	\$ 236.00

Head Support Options

107470-99-0	Permobil 7"W Headrest <i>Must also select headrest mounting hardware.</i>	E0955	\$ 434.00
107469-99-0	Permobil 10"W Headrest <i>Must also select headrest mounting hardware.</i>	E0955	\$ 434.00
107468-99-0	Permobil 14"W Headrest <i>Must also select headrest mounting hardware.</i>	E0955	\$ 434.00
107480-99-0	Adj/Removable Headrest Link Hardware - Standard	E1028/ E1033	\$ 299.00
109178-99-0	Extra Link for Link Headrest Hardware <i>Sent as a part.</i>		\$ 55.00
108721-99-0	Adj/Removable BodiLink Head Support Hardware - MAX 8 <i>Required with more than 2 swing-away headrest accessories or chin boom.</i>	E1028/ E1033	\$ 365.00
109625-99-0	Adj/Removable BodiLink Head Support Hardware - MAX 8 Flip-Down <i>Flip-down joint allows head support to swing back 120°.</i>	E1028/ E1033	\$ 501.00
109624-99-0	Extra 4" Link for BodiLink MAX 8 Hardware <i>Sent as a part.</i>		\$ 80.00
11305	Swing-Away Egg Switch Assembly Mount on Upper Left Mount on Upper Right	E1028/ E1032	\$ 640.00
11306	Swing-Away Microlight Assembly Mount on Upper Left Mount on Upper Right	E1028/ E1032	\$ 730.00
11307	Swing-Away Facial Pad Assembly 4"L x 2"T Mount on Upper Left Mount on Upper Right	E1028	\$ 519.00
11308	Swing-Away Facial Pad Assembly 6"L x 2"T Mount on Upper Left Mount on Upper Right	E1028	\$ 519.00
11276	Provide Headrest Receiver for AGILITY Backrest <i>Select for AGILITY backs and third party backs ordered off chair.</i>		No charge
106331-99-0	Permobil 10"W Headrest, Leatherette <i>Only available with Corpus backrests. Must also select hardware.</i>	E0955	\$ 397.00
10552	Adj/Removable Headrest Hardware	E1028/ E1033	\$ 306.00
107499-99-0	Universal Headrest Adapter		\$ 302.00
Extended UniTrack Rail	10" 13" 15"		\$ 285.00

Leg Support Options | *Lateral pelvic thigh supports (LPTS)*

103070-99-0	UT Amp. Support Pad (H), Mesh Cover - 1 Each <i>7"W x 7"T curved pad. Must also select hardware.</i>	E1020	\$ 181.00
103071-99-0	UT Amp. Support Pad (D), Mesh Cover - 1 Each <i>8"W x 8"T flat pad. Must also select hardware.</i>	E1020	\$ 181.00
10553	Adj Removable - Amp Support Hardware - 1 Each Mount on Left Mount on Right		\$ 306.00
109040-99-0	BodiLink Premium LPTS Pad (B), Stretch-Air - Pair <i>4.5"W x 3.5"T pad. Must also select hardware.</i>	E0953/ E0956	\$ 359.00
109043-99-0	BodiLink Premium LPTS Pad (A), Stretch-Air - Pair <i>6"W x 5.5"T pad. Must also select hardware.</i>	E0953/ E0956	\$ 359.00
109041-99-0	BodiLink Premium LPTS Pad (C), Stretch-Air - Pair <i>8"W x 3.5"T pad. Must also select hardware.</i>	E0953/ E0956	\$ 359.00
109042-99-0	BodiLink Premium LPTS Pad (G), Stretch-Air - Pair <i>12"W x 3.5"T pad. Must also select hardware.</i>	E0953/ E0956	\$ 359.00
109044-99-0	Adj/Removable BodiLink LPTS Hardware Med Corpus - Pair <i>Height adjustable.</i>	E1028/ E1034	\$ 689.00
109194-99-0	Adj/Removable Permobil Lateral Thigh/Hip Hardware - Pair	E1028/ E1034	\$ 602.00
106762-99-0	Abductor Kit Corpus - Leatherette <i>Not compatible with Corpus VS Legrest. Must also select hardware.</i>	E0957	\$ 423.00
10555	Adj Removable - Abductor Hardware - 1 Each	E1028	\$ 306.00

Chest Support Options

105840-99-0	Bodypoint Monoflex Chest Support Belt Kit - Medium <i>18.25"L x 4"W Pad Size.</i>	E0960	\$ 204.00
105839-99-0	Bodypoint Monoflex Chest Support Belt Kit - Large <i>21.5"L x 4.5"W Pad Size.</i>	E0960	\$ 204.00
105850-99-0	Bodypoint Monoflex Chest Support Belt Kit - X-Large <i>24.5"L x 5"W Pad Size.</i>	E0960	\$ 204.00

Other Options

105260-99-0	Transfer Handles		\$ 550.00
105292-99-0	Push Handles, Corpus Seat		\$ 489.00
106210-99-0	RAM Self-Leveling Cup Holder		\$ 93.00
106806-99-0	RAM X-Grip Phone Holder		\$ 93.00
107100-99-0	RAM X-Grip 10" Tablet Holder w/ Swing-Away Mount		\$ 467.00
108995-99-0	USB Charger - Power Platform <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets.</i> Mount on Left Armrest Mount on Right Armrest		\$ 231.00
105805-99-0	Medical Necessities Bag Hooks		\$ 111.00
102059-99-0	Medical Necessities Bag		\$ 321.00
108231-99-0	Permobil Adjustable Tray - 11"W <i>Must also select hardware.</i>	E0950	\$ 284.00
108230-99-0	Permobil Adjustable Tray - 13"W <i>Must also select hardware.</i>	E0950	\$ 284.00
108228-99-0	Permobil Adjustable Tray - 15"W <i>Must also select hardware.</i>	E0950	\$ 284.00
108227-99-0	Permobil Adjustable Tray - 17"W <i>Must also select hardware.</i>	E0950	\$ 284.00
108229-99-0	Permobil Adjustable Tray - 19"W <i>Must also select hardware.</i>	E0950	\$ 284.00
Permobil Adjustable Tray Hardware	Left Side Adjustable Right Side Adjustable		\$ 208.00

100578-99-0	Small Angle Adjustable Tray <i>Must also select hardware.</i>	E0950	\$ 339.00
105262-99-0	Joystick Cutout Tray - 16" W <i>Must also select hardware.</i>	E0950	\$ 438.00
105263-99-0	Joystick Cutout Tray - 18" W <i>Must also select hardware.</i>	E0950	\$ 438.00
Small/Cutout Adjustable Tray Hardware	Mount on Left Mount on Right		\$ 306.00
105224-99-0	Oxygen Holder, Corpus	E2208	\$ 301.00
105223-99-0	Crutch Holder, Corpus	E2207	\$ 301.00

Alternative Drive Systems

Omni2 Control Module required for all alternative drive systems.

Alt Proportional Joysticks

109525-99-0	mo-vis Micro Joystick for Omni <i>Operating Force = 8 g (0.28 oz). Neutral to Full Deflection = 0.13" (3.3 mm).</i>	E2312	\$ 3,200.00
109527-99-0	mo-vis Multi Joystick for Omni <i>Operating Force = 50 g (1.8 oz). Neutral to Full Deflection = 0.31" (8 mm).</i>	E2312	\$ 3,200.00
103894-99-0	Compact Joystick Lite for Omni - (9-pin connector) <i>Operating Force = 50 g (1.8 oz). Neutral to Full Deflection = 0.7" (17mm).</i>	E2373	\$ 2,003.00
103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 220 g (7.8 oz). Neutral to Full Deflection = 1.1" (28mm).</i>	E2373	\$ 1,771.00

Alt Proportional Joystick Mounts | *Please see Alt Drive Accessories section for additional mounting options.*

107857-99-0	Swing-Away Chin Boom for Micro/Multi Joystick - Left	E1028/ E1032	\$ 465.00
107947-99-0	Swing-Away Chin Boom for Micro/Multi Joystick - Right	E1028/ E1032	\$ 465.00
107946-99-0	Swing-Away Chin Boom for Compact Lite Joystick - Left	E1028/ E1032	\$ 465.00
107948-99-0	Swing-Away Chin Boom for Compact Lite Joystick - Right	E1028/ E1032	\$ 465.00
104690-99-0	Alt Joystick Swing-Away Mount Kit - Armrest Mount on Left Mount on Right	E1028/ E1032	\$ 310.00
107135-99-0	Alt Joystick / Omni Combo Swing-Away Mount Kit - Armrest Mount on Left Mount on Right	E1028/ E1032	\$ 586.00
106657-99-0	PERMOfix UniTrack Swing-Away Mount Kit		\$ 786.00
108268-99-0	Bodypoint Tri-Lock Midline Mount Assembly <i>Must select specific joystick attachment separately.</i> Mount on Left Mount on Right		\$ 702.00
108279-99-0	Bodypoint Tri-Lock Mount for Compact Joystick		\$ 292.00
10842	PERMOFix Mount for mo-vis/Compact Joystick <i>Permofix mount for attaching Compact and mo-vis joysticks to Tri-Lock midline mount.</i>		\$ 292.00
104143-99-0	Alt Joystick Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	E1028/ E1032	\$ 966.00

Alt Proportional Joystick Accessories

104904-99-0	Dual Switches for Compact Joystick - Fixed Mount Mount Left of JS Mount Right of JS Mount front of JS Mount back of JS		\$ 288.00
104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mount		\$ 518.00

Total Control Head Array

109558-99-0	Permobil Total Control Head Array <i>Features an occipital pad sensor and two lateral arms configured with proximity sensors. Must select swing away hardware and mounting hardware below to complete assembly. An optional lateral arm with egg switch is strongly recommended.</i>	E2330	\$ 4,819.00
11575	Swing-Away Hardware for Side Sensor - Left	E1028	\$ 263.00
11576	Swing-Away Hardware for Side Sensor - Right	E1028	\$ 263.00
11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the notes section.</i>		No charge
109563-99-0	Adj/Removable BodiLink HA Hardware - MAX 8	E1028/ E1033	\$ 365.00
109564-99-0	Adj/Removable BodiLink HA Hardware - MAX 8 Flip-Down <i>Flip-down joint allows head array to swing back 120°.</i>	E1028/ E1033	\$ 501.00
10976	(Optional) Add Left Lateral Arm w Proximity Sensor	E1028	\$ 688.00
10977	(Optional) Add Left Lateral Arm w Egg Switch	E1028	\$ 563.00
10978	(Optional) Add Left Lateral Arm w Pad (No Sensor or Switch)	E1028	\$ 459.00
10979	(Optional) Add Right Lateral Arm w Proximity Sensor	E1028	\$ 688.00
10980	(Optional) Add Right Lateral Arm w Egg Switch	E1028	\$ 563.00
10981	(Optional) Add Right Lateral Arm w Pad (No Sensor or Switch)	E1028	\$ 459.00

ASL FUSION Head Array

109582-99-0	ASL FUSION Proportional and Digital Head Array <i>Includes occipital (ASL 935) and lateral (ASL 933) pads, and pulse mechanical switch (ASL 320BK) for reset/mode changes. Must select hardware listed below to complete assembly.</i>	E2328	\$ 6,715.00
11580	Swing-Away Hardware for Lateral Sensor - Left	E1028	\$ 342.00
11581	Swing-Away Hardware for Lateral Sensor - Right	E1028	\$ 342.00
11582	Swing-Away Hardware for Reset/Mode Switch - Right	E1028	\$ 342.00
11583	Adj/Removable Headrest Mounting Hardware <i>Includes ASL 606 and BodiLink receiver.</i>	E1028/ E1033	\$ 342.00
109583-99-0	ASL FUSION Display, Programmer and Attendant Control w/ Mount <i>Includes ASL 165 and quick release mount.</i>	K0108	\$ 730.00

Sip & Puff Systems

10573	Sip & Puff System for Omni Mount on Left Mount on Right	E2325, E2326	\$ 1,937.00
10574	Sip & Puff System with Fiber Optic Opti-Stop Mount on Left Mount on Right	E2325, E2326	\$ 2,577.00
10762	Adjustable Mounting Hardware - Sip & Puff <i>Required with any sip & puff system.</i>		\$ 310.00

Proportional Foot Control

10767	Proportional Foot Control <i>Not available with ActiveReach packages.</i> Mount on Left Mount on Right	E2328	\$ 3,690.00
10766	Adjustable Mounting Hardware - Foot Control <i>Required with proportional foot control.</i>		\$ 310.00

Alt Drive System Setup

Omni2 Options | Required for all alternative drive systems.

109589-99-0	Omni2 Control Module for R-net - Corpus - Power Platform <i>3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology. Must select mounting hardware below. Specify alternative input device being used:</i>		\$ 2,061.00
103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet Mount on Left Mount on Right		\$ 310.00

PERMOfix Mounting Solutions

106655-99-0	PERMOfix UniTrack Fixed Mount Kit		\$ 312.00
106656-99-0	PERMOfix UniTrack Removable Mount Kit		\$ 549.00
105907-99-0	Proximity Sensor + Swivel Assembly <i>Swivel assembly allows sensor to be mounted with PERMOfix. Since Proximity Sensors require power, CANNOT be used as On/Off switches.</i>	E1028	\$ 397.00
105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix.</i>		\$ 268.00
105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix.</i>		\$ 162.00
106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Microlight, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix.</i>		\$ 142.00

Switches and Cables

101769-99-0	Microlight Switch Qty:		\$ 217.00
101767-99-0	Buddy Button Qty:		\$ 183.00
106293-99-0	Egg Switch Black Qty:		\$ 183.00
101771-99-0	Mini Cup Switch Qty:		\$ 183.00
107500-99-0	Proximity Switch Qty:		\$ 723.00
10772	Dome Switch, Backrest Mount w/ Gooseneck Qty: <i>Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights.</i>		\$ 573.00
	Mount on Left	Mount on Right	
10846	Dome Switch, Backrest Swing-Away Mount w/ Gooseneck Qty: <i>Swing-away, mounts at top of backrest shell. Custom charge may be applied if ordered with short backrest heights.</i>		\$ 847.00
	Mount on Left	Mount on Right	
105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long		\$ 36.00
105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long		\$ 39.00
105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long		\$ 47.00
105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long		\$ 51.00
105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>		\$ 39.00

Additional Order Instructions:

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of these items and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that any parts ordered for service or future alterations may carry different pricing and warranty. For US customers, HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.