

2026 order form

Use dimension chart to select appropriate size.

Scan for sizing and pricing:
<https://www.rifton.com/mobilestander>



The **conversion kit** converts the K110 mini stander to the K130 small. It includes the deck and complete body support of the K130. Request from customer service.

Conversion kit K116 \$850



* **Mobile Stander** frame includes body support with padding, seat pad, straps and fleece.

Specify color:
pink, blue or lime

CHOOSE

User dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Height	25–32	30–40	37–52	48–65
Key user dimension: height Select the appropriate Mobile Stander by the user's overall height. Choose the model that allows for growth. Important: User's weight must not exceed the maximum working load.				
Item dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Overall length and width	25 x 25	25 x 25½	32 x 29	39 x 33½
Overall height	22½–28	23½–30	31–39	37½–49
Width without large wheels	19½	19½	24½	28¾
Wheel diameter	20	20	27	35
Width of body support	7–10	7–10	9–12	10–14
Height of body support	16½–22	22–29	29–38	36–48
Item weight (lb)	28	30	40	65
Max. working load (lb)	50	50	120	175

		K110 mini	K130 small	K150 medium	K170 large
CHOOSE		K101 \$3,664 ☐	K103 \$4,088 ☐	K105 \$4,968 ☐	K107 \$6,053 ☐
	Pink	☐	Pink ☐	Pink ☐	Pink ☐
	Blue	☐	Blue ☐	Blue ☐	Blue ☐
	Lime	☐	Lime ☐	Lime ☐	Lime ☐
Wheels (pair) (include brakes)		K113 \$2,169 ☐	K113 \$2,169 ☐	K115 \$2,435 ☐	K117 \$2,938 ☐
Mounting bracket		K133 \$115 ☐	K133 \$115 ☐	K173 \$129 ☐	K173 \$129 ☐
Arm prompts (pair) (require mounting bracket)		K512 \$964 ☐	K512 \$964 ☐	K512 \$964 ☐	K512 \$964 ☐
CHOOSE	Single	K578 \$177 ☐	K578 \$177 ☐	K578 \$177 ☐	K578 \$177 ☐
	Pair	K579 \$378 ☐	K579 \$378 ☐	K579 \$378 ☐	K579 \$378 ☐
Communication tray (requires mounting bracket)		K518 \$738 ☐	K518 \$738 ☐	K518 \$738 ☐	K518 \$738 ☐
CHOOSE	Small	N/A	N/A	K135 \$260 ☐	K135 \$260 ☐
	Large	N/A	N/A	K155 \$287 ☐	K155 \$287 ☐
CHOOSE	Small	K114 \$201 ☐	K114 \$201 ☐	N/A	N/A
	Medium	N/A	N/A	K154 \$214 ☐	N/A
	Large	N/A	N/A	N/A	K174 \$228 ☐
Total		K110 \$ _____	K130 \$ _____	K150 \$ _____	K170 \$ _____

* Required