

2026 order form

Scan for sizing and pricing:
https://www.rifton.com/



Use dimension chart to select appropriate size.

			K610 mini			K620 small			K630 medium			K640 large			K650 XL				
* Color			CHOOSE	Red	☐	Red	☐	Red	☐	Red	☐	Red	☐	Red	☐				
				Blue	☐	Blue	☐	Blue	☐	Blue	☐	Blue	☐						
				Black	☐	Black	☐	Black	☐	Black	☐	Black	☐						
				Lime	☐	Lime	☐	Lime	☐	Lime	☐	Lime	☐						
				Pink	☐	Pink	☐	Pink	☐	Pink	☐	Pink	☐						
* Upper frame		CHOOSE	Dynamic	N/A		K621	\$2,391	☐	K631	\$3,166	☐	K641	\$3,563	☐	K651	\$4,009	☐		
		Standard	N/A		K622	\$1,273	☐	K632	\$1,665	☐	K642	\$1,962	☐	K652	\$2,865	☐			
* Base	Standard base		CHOOSE	With odometer	K613	\$1,850	☐	K623	\$1,500	☐	K633	\$1,643	☐	K643	\$1,961	☐	K653	\$2,481	☐
			Without odometer	K614	\$1,531	☐	K624	\$800	☐	K634	\$963	☐	K644	\$1,366	☐	K654	\$2,141	☐	
	Utility base (larger wheels)		CHOOSE	With odometer	N/A		N/A		K635	\$2,163	☐	K645	\$2,562	☐	N/A				
			Without odometer	N/A		N/A		K636	\$1,844	☐	K646	\$2,163	☐	N/A					
Arm supports (pair)			CHOOSE	Small arm prompts	K512	\$964	☐	K512	\$964	☐	K512	\$964	☐	K512	\$964	☐	K512	\$964	☐
				Large arm prompts	N/A		N/A		K532	\$1,104	☐	K532	\$1,104	☐	K532	\$1,104	☐		
				Arm platforms	N/A		N/A		K697	\$842	☐	K697	\$842	☐	K697	\$842	☐		
				Arm platforms w/ handgrips	N/A		N/A		K696	\$968	☐	K696	\$968	☐	K696	\$968	☐		
Handbrake for utility base (requires arm supports)				N/A		N/A		K687	\$651	☐	K688	\$651	☐	N/A					
Handlebars (pair)			CHOOSE	Small	K625	\$486	☐	K625	\$486	☐	K625	\$486	☐	K625	\$486	☐	K625	\$486	☐
				Large	N/A		N/A		K675	\$511	☐	K675	\$511	☐	K675	\$511	☐		
Hand loops (pair)				K579	\$378	☐	K579	\$378	☐	K579	\$378	☐	K579	\$378	☐	K579	\$378	☐	
Chest prompt			CHOOSE	Small	K514	\$1,115	☐	K514	\$1,115	☐	K514	\$1,115	☐	N/A		N/A			
				Medium	N/A		K524	\$1,183	☐	K524	\$1,183	☐	K524	\$1,183	☐	K524	\$1,183	☐	
				Large	N/A		N/A		K534	\$1,255	☐	K534	\$1,255	☐	K534	\$1,255	☐		
Chest pad				N/A		K629	\$728	☐	K639	\$782	☐	K649	\$838	☐	N/A				
Multi-position saddle			CHOOSE	Small	N/A		K672	\$2,242	☐	N/A		N/A		N/A					
				Medium	N/A		N/A		K673	\$2,154	☐	N/A		N/A					
				Large	N/A		N/A		K674	\$2,242	☐	K674	\$2,242	☐	K674	\$2,242	☐		
MPS seat cover			CHOOSE	Small/Medium	N/A		K627	\$174	☐	K627	\$174	☐	N/A		N/A				
				Large	N/A		N/A		K628	\$174	☐	K628	\$174	☐	K628	\$174	☐		
Pelvic support (includes handholds)			CHOOSE	Small	K556	\$797	☐	K556	\$797	☐	K556	\$797	☐	N/A		N/A			
				Medium	N/A		K557	\$847	☐	K557	\$847	☐	K557	\$847	☐	K557	\$847	☐	
				Large	N/A		N/A		K558	\$891	☐	K558	\$891	☐	K558	\$891	☐		
Hip positioner (includes handholds)			CHOOSE	Small	K513	\$797	☐	K513	\$797	☐	K513	\$797	☐	K513	\$797	☐	K513	\$797	☐
				Large	N/A		N/A		K533	\$891	☐	K533	\$891	☐	K533	\$891	☐		
Hip positioner pad			CHOOSE	Small	K541	\$206	☐	K541	\$206	☐	K541	\$206	☐	K541	\$206	☐	K541	\$206	☐
				Large	N/A		N/A		K542	\$270	☐	K542	\$270	☐	K542	\$270	☐		
Thigh prompts (pair)			CHOOSE	Small	N/A		K515	\$797	☐	K515	\$797	☐	K515	\$797	☐	N/A			
				Large	N/A		N/A		K535	\$964	☐	K535	\$964	☐	K535	\$964	☐		
Ankle prompts for standard base (pair)					K681	\$363	☐	K682	\$384	☐	K683	\$418	☐	K684	\$483	☐	K684	\$483	☐
Ankle prompts for utility base (pair)					N/A		N/A		K693	\$418	☐	K694	\$483	☐	N/A				
Additional handholds (pair)					K519	\$219	☐	K519	\$219	☐	K519	\$219	☐	K519	\$219	☐	K519	\$219	☐
Communication tray					K518	\$738	☐	K698	\$771	☐	K698	\$771	☐	K698	\$771	☐	K698	\$771	☐
Attendant guide bar					K589	\$270	☐	K679	\$315	☐	K689	\$315	☐	K689	\$315	☐	N/A		
Rifton accessories tote					K522	\$310	☐	K522	\$310	☐	K522	\$310	☐	K522	\$310	☐	K522	\$310	☐
Total					K610 \$		K620 \$		K630 \$		K640 \$		K650 \$						

* Required

Pacer dimensions

Key user dimensions (inches)		K610 mini	K620 small	K630 medium	K640 large	K650 XL
Elbow height		15½–20½	18½–27½	24–35	32–47	34–49
 Key user dimension: elbow height Measure the vertical distance from the bent elbow to the floor while the user is standing upright. Choose the size that allows for growth.						
Item dimensions (inches)		K610 mini	K620 small	K630 medium	K640 large	K650 XL
Floor to top of arm prompt		15½–20½	18½–27½	24–35	32–47	34–49
Floor to top of MPS (not tilted)*		N/A	10–20	9½–23½	17–35	19½–37½
Overall width	Standard base:	20½	23	26	28	31½
	Utility base:	N/A	N/A	30	32	N/A
Overall length	Standard base:	22½	26	30	36	39¾
	Utility base:	N/A	N/A	36	40½	N/A
Frame height	Standard & utility base:	12½	16–21	21½–27	28½–38½	31–41
Base height (without upper frame)	Standard & utility base:	N/A	11	14½	16	18½
Frame weight (lb)	Standard base:	10¼	11½	15	15½	22½
	Utility base:	N/A	N/A	18½	19¼	N/A
	Standard upper:	N/A	6	7½	9½	11
	Dynamic upper:	N/A	7½	11½	13½	15
Dynamic upper movement	Vertical:	N/A	1¼	2	2	2
	Horizontal:	N/A	1½	2	2	2
Maximum working load (lb)		50	75	150	200	250

*This measurement is taken at the lowest frame and MPS height, and the highest frame and MPS height.

Pelvic support

User dimensions (inches)

Indicated by the red arrow below

Key user dimension:


half of thigh circumference

measured at the groin.


K556 small	K557 med.	K558 large
4 – 8	6 – 10	8 – 14

Product dimensions (inches)		K556 small	K557 med.	K558 large
a		6½	8½	12
b		12	16	20
c		7	9	12

Hip positioner

Dimensions (inches) – without pad		
		
	K513 small	K533 large
a Crotch length	2½	3½
b Width	1¼	1¾
c Length	9	11

MPS saddle

Dimensions (inches)		
		
	K672/K673 small/med	K674 large
a Width at rear	3½	5¼
b Width at front	3	4½
c Length	9½	10¾



All five sizes are available in five colors: red, black, pink, lime and blue.