Order form

Comfort Unite 2 back support





Dealer Name*: _____ Account #*: _____

City, State, Zip: _____

Phone: _____ Email: _____

ATP: _____ Rep Name: _____

P.O. Number: _____ Quote Number: _____

Client Reference*: _____ Client Weight: _____

Required field

Special considerations that need to be addressed here (e.g., diagnosis):

_____

Notes | *Refer to reference pages at end of order form for additional information and images.*

Unite 2 back support with foam only insert

Part number: U2-BSXXWYYL.

Length	12”	14”	16”	18”	20”	22”	24”
22”			☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,408.00	☐ $ 1,408.00
20”		☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,408.00	☐ $ 1,408.00
18”	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,408.00	☐ $ 1,408.00
16”	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,408.00	
14”	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00		
12”	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00	☐ $ 1,106.00			

Width

12”W x 16”L, and 12”W x 18”L not available with power mount hardware. 12L, 14L, and 12W not available with direct mount hardware. 24” wide back supports require 4 point, power mount hardware or direct mount hardware.

Cover | *Required, choose one.*

☐ C-COMFORT-TEK	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ C-DIAMOND-SPACER	**Diamond Spacer** *For comfort and heat dissipation. NOT a fluid resistant material.*	NC
☐ C-GLIDEWEAR-DIAMOND	**Diamond Spacer with GlideWear** *For skin protection and shear reduction. NOT a fluid resistant material*	$ 85.00

Unite 2 back support with ROHO Insert

Part number: U2R-BSXXWYYL

Length	14"	16"	18"	20"	22"	24"
22"		☐ $ 1,254.00	☐ $ 1,254.00	☐ $ 1,511.00	☐ $ 1,511.00	☐ $ 1,511.00
20"	☐ $ 1,254.00	☐ $ 1,254.00	☐ $ 1,254.00	☐ $ 1,511.00	☐ $ 1,511.00	☐ $ 1,511.00
18"	☐ $ 1,254.00	☐ $ 1,254.00	☐ $ 1,254.00	☐ $ 1,254.00	☐ $ 1,511.00	☐ $ 1,511.00
16"	☐ $ 1,254.00	☐ $ 1,254.00	☐ $ 1,254.00	☐ $ 1,254.00	☐ $ 1,511.00	
14"	☐ $ 1,254.00	☐ $ 1,254.00	☐ $ 1,254.00	☐ $ 1,254.00		
			Width			

14L not available with direct mount hardware. 24" wide back supports require 4 point, power mount or direct mount hardware.

Cover | *Required, choose one.*

☐ C-COMFORT-TEK	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ C-DIAMOND-SPACER	**Diamond Spacer** *For comfort and heat dissipation. NOT a fluid resistant material.*	NC
☐ C-GLIDEWEAR-DIAMOND	**Diamond Spacer with GlideWear** *For skin protection and shear reduction. NOT a fluid resistant material*	$ 85.00

Unite 2 back support

Mounting Hardware | *Required, choose one.*

Compass 4 Quick Release	☐ Quick Release 2 Point (C-Q-COMPASS) ☐ Quick Release 4 Point (C-Q-COMPASS4)	NC
Compass 4 Fixed	☐ Fixed 2 Point (C-NQ-COMPASS) ☐ Fixed 4 Point (C-NQ-COMPASS4)	NC
Compass Power Mount	☐ Permobil w/ Corpus 3G (C-PMB-PERMOBIL) ☐ Quantum w/ Tru Balance 3 or 4 (C-PMB-QUANTUM) ☐ Quickie w/ 3.7 Power (C-PMB-QUICKIE) ☐ Avid Rehab w/ Power T & R (C-PMB-AVID) ☐ ROVI w/ Motion Concepts (C-PMB-ROVI) *Using Power Mount hardware on ROVI requires mounting to the recessed planar interface plate available when ordering the Motion Concepts seating system. Drilling is required for mounting.*	$ 125.00
Permobil Direct Mount	☐ Permobil Direct Mount (C-DIRECT-MOUNT) *Includes Unitrack Assembly*	$ 125.00

Embroidery | *Optional, not available with GlideWear cover. Includes corresponding color zipper pull*

Thread color	☐ Red (C-EMB-RED)	☐ Purple (C-EMB-PURPLE)	☐ Pink (C-EMB-PINK)	$ 60.00
	☐ Blue (C-EMB-BLUE)	☐ Orange (C-EMB-ORANGE)	☐ Green (C-EMB-GREEN)	
	☐ Silver (C-EMB-SILVER)			

Font

☐ Helvetica (C-EMB-HELVETICA) ☐ Garamond (C-EMB-GARAMOND)

☐ Impress (C-EMB-IMPRESS) ☐ Kaufmann (C-EMB-KAUFMANN)

Font styles

HELVETICA
GARAMOND
Impress
Kaufmann

Word one ☐ Input is case sensitive:

Word two (optional) ☐ Input is case sensitive:

Other Accessories | *Optional*

☐ C-PV-SHIELD	**Privacy Shield** *Connects with hook & loop attachment*	$ 85.00

BodiLink lateral trunk supports & hardware

Optional, if selecting must choose pad shape & size, pad style, pad covers, and mounting hardware.

Basic pad shape & size | *Priced per pad.*

	Left *Must enter the number of pads needed*	**Right** *Must enter the number of pads needed*	
Supportive			
3”L x 4”D	☐ BL-LTSP1S-3L4D-LH	☐ BL-LTSP1S-3L4D-RH	$ 161.00
3”L x 5”D	☐ BL-LTSP1S-3L5D-LH	☐ BL-LTSP1S-3L5D-RH	$ 161.00
4”L x 5”D	☐ BL-LTSP1S-4L5D-LH	☐ BL-LTSP1S-4L5D-RH	$ 161.00
4”L x 6”D	☐ BL-LTSP1S-4L6D-LH	☐ BL-LTSP1S-4L6D-RH	$ 161.00
5”L x 7”D	☐ BL-LTSP1S-5L7D-LH	☐ BL-LTSP1S-5L7D-RH	$ 161.00
6”L x 6”D	☐ BL-LTSP1S-6L6D-LH	☐ BL-LTSP1S-6L6D-RH	$ 161.00
Active			
3”L x 4”D	☐ BL-LTSP1A-3L4D-LH	☐ BL-LTSP1A-3L4D-RH	$ 161.00
3”L x 5”D	☐ BL-LTSP1A-3L5D-LH	☐ BL-LTSP1A-3L5D-RH	$ 161.00
4”L x 5”D	☐ BL-LTSP1A-4L5D-LH	☐ BL-LTSP1A-4L5D-RH	$ 161.00
4”L x 6”D	☐ BL-LTSP1A-4L6D-LH	☐ BL-LTSP1A-4L6D-RH	$ 161.00
5”L x 7”D	☐ BL-LTSP1A-5L7D-LH	☐ BL-LTSP1A-5L7D-RH	$ 161.00
6”L x 6”D	☐ BL-LTSP1A-6L6D-LH	☐ BL-LTSP1a-6L6D-RH	$ 161.00

Basic pad style | *Must match number of pads selected.*

☐ LTS-P1-MOLD — **Moldable** — NC
☐ Left ☐ Right

☐ LTS-P1-RIGID — **Rigid** — NC
☐ Left ☐ Right

Basic pad cover options | *Must match number of pads selected.*

☐ LTS-P1-COMFORT — **COMFORT-TEK®** — NC
For fluid protection & an easily cleaned surface
☐ Left ☐ Right

☐ LTS-P1-STRETCH — **STRETCH-AIR®** — NC
For comfort & heat dissapation
☐ Left ☐ Right

Premium pad shape & size | *Priced per pad.*

	Left *Enter number of pads required*	Right *Enter number of pads required*	
Supportive			
3”L x 4” D	☐ BL-LTSP2S-3L4D-LH	☐ BL-LTSP2S-3L4D-RH	$ 218.00
3” L x 5.5” D	☐ BL-LTSP2S-3L5D-LH	☐ BL-LTSP2S-3L5D-RH	$ 218.00
4” L x 5” D	☐ BL-LTSP2S-4L5D-LH	☐ BL-LTSP2S-4L5D-RH	$ 218.00
4” L x 6.5” D	☐ BL-LTSP2S-4L6D-LH	☐ BL-LTSP2S-4L6D-RH	$ 218.00
5” L x 6.5” D	☐ BL-LTSP2S-5L6D-LH	☐ BL-LTSP2S-5L6D-RH	$ 218.00
6” L x 6.5” D	☐ BL-LTSP2S-6L6D-LH	☐ BL-LTSP2S-6L6D-RH	$ 218.00
Active			
3”L x 4” D	☐ BL-LTSP2A-3L4D-LH	☐ BL-LTSP2A-3L4D-RH	$ 218.00
3” L x 5.5” D	☐ BL-LTSP2A-3L5D-LH	☐ BL-LTSP2A-3L5D-RH	$ 218.00
4” L x 5” D	☐ BL-LTSP2A-4L5D-LH	☐ BL-LTSP2A-4L5D-RH	$ 218.00
4” L x 6.5” D	☐ BL-LTSP2A-4L6D-LH	☐ BL-LTSP2A-4L6D-RH	$ 218.00
5” L x 6.5” D	☐ BL-LTSP2A-5L6D-LH	☐ BL-LTSP2A-5L6D-RH	$ 218.00
6” L x 6.5” D	☐ BL-LTSP2A-6L6D-LH	☐ BL-LTSP2A-6L6D-RH	$ 218.00

Premium pad style | *Must match number of pads selected.*

☐ LTS-P2-MOLD	**Moldable** ☐ Left ☐ Right	NC
☐ LTS-P2-RIGID	**Rigid** ☐ Left ☐ Right	NC

Premium pad cover options | *Must match number of pads selected. Priced per selection.*

☐ LTS-P2-COMFORT	**COMFORT-TEK** *For fluid protection & an easily cleaned surface* ☐ Left ☐ Right	NC
☐ LTS-P2-STRETCH	**STRETCH-AIR** *For comfort & heat dissapation* ☐ Left ☐ Right	NC
☐ LTS-P2-GLIDE	**GlideWear®** *For skin protection & shear reduction. NOT fluid resistant* ☐ Left ☐ Right	$ 75.00
☐ LTS-P2-COMFORT-FM	**COMFORT-TEK with foam** ☐ Left ☐ Right	$ 30.00
☐ LTS-P2-STRETCH-FM	**STRETCH-AIR with foam** ☐ Left ☐ Right	$ 30.00
☐ LTS-P2-GLIDE-FM	**GlideWear with foam** ☐ Left ☐ Right	$ 103.00

GT2 mounting hardware | *Used with basic and premium pads. Both fixed and swing away hardware styles are adjustable.*

	Left *Must match number of pads selected*	Right *Must match number of pads selected*	
Swing-away Slot Mount	☐ BL-LTS-GT2SAS-LH	☐ BL-LTS-GT2SAS-RH	$ 368.00
Fixed Slot Mount	☐ BL-LTS-GT2FXS-LH	☐ BL-LTS-GT2FXS-RH	$ 202.00

Additional hardware options | *Priced per selection.*

Extra GT2 hardware link	☐ Left (BL-LTS-GT2LINK)	☐ Right (BL-LTS-GT2LINK)	$ 58.00

Long mounting bracket	☐ Left (BL-LTS-GT2LSLOT-LH)		$ 58.00
	☐ Right (BL-LTS-GT2LSLOT-RH)		
Hardware cover	☐ Left (BL-LTS-HCOVER)	☐ Right (BL-LTS-HCOVER)	$ 52.00

BodiLink premium head support & hardware

Optional, if selecting must choose pad size, foam selection, cover, and mounting hardware.

Pad size and foam selection

Pad size	Soft	Extra soft	
6”W x 3”L	☐ BL-HSP1-6W3L	☐ BL-HSP2-6W3L	$ 260.00
8”W x 4”L	☐ BL-HSP1-8W4L	☐ BL-HSP2-8W4L	$ 260.00
10”W x 5”L	☐ BL-HSP1-10W5L	☐ BL-HSP2-10W5L	$ 260.00
12”W x 5”L	☐ BL-HSP1-12W5L	☐ BL-HSP2-12W5L	$ 260.00
14”W x 5”L	☐ BL-HSP1-14W5L	☐ BL-HSP2-14W5L	$ 260.00
18”W x 5”L	☐ BL-HSP1-18W5L	☐ BL-HSP2-18W5L	$ 260.00

Cover options

☐ COMFORT-TEK	**COMFORT-TEK** *For fluid protection & an easily cleaned surface*	NC
☐ STRETCH-AIR	**STRETCH-AIR** *For comfort & heat dissapation*	NC
☐ GLIDEWEAR	**GlideWear** *For skin protection & shear reduction. NOT fluid resistant*	$ 82.00

PT hardware, removable

☐ BL-HS-PT1RMMID6	**Mid 6**	$ 373.00
☐ BL-HS-PT1RMLAT6	**Lat 6**	$ 489.00
☐ BL-HS-PT1RMMID8	**Mid 8**	$ 373.00
☐ BL-HS-PT1RMLAT8	**Lat 8**	$ 489.00
☐ BL-HS-PT1RMMAX8	**Max 8**	$ 545.00

Stem length option | *For removable PT hardware only*

☐ HS-STEM-7	**Standard 7” stem**	NC
☐ HS-STEM-3	**3”stem (recommended for use with power mounting where there may be constraints)**	NC

PT hardware, flip down | *Flip down hardware has 7” stem.*

☐ BL-HS-PT1FDMID6	**Mid 6**	$ 419.00
☐ BL-HS-PT1FDLAT6	**Lat 6**	$ 530.00
☐ BL-HS-PT1FDMID8	**Mid 8**	$ 419.00
☐ BL-HS-PT1FDLAT8	**Lat 8**	$ 530.00
☐ BL-HS-PT1FDMAX8	**Max 8**	$ 591.00

Mounting option

☐ HS-MB-POSTA	**Posterior A (For Compass 4 and Compass Power Mount hardware)**	NC
☐ HS-MB-PW2C	**Permobil C (For Permobil Direct Mount only)**	$155.00

Additional hardware options | *Priced per selection.*

Extra link for PT hardware	☐ Extra 2” link (BL-HS-PT1LINK2)	☐ Extra 3” link (BL-HS-PT1LINK3)	$ 114.00
	☐ Extra 4” link (BL-HS-PT1LINK4)		

☐	BL-HS-52WINGA1	**Tool-less wing lock for PT hardware** *Enter desired number.*	$ 47.00
☐	BL-HS-SPCRA	**1/4 Spacer Plate**	$ 87.00

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- The client's name or code
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, , am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of seating and positioning products and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____


Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.


References

Section 1 | *Hardware options*


Compass 4 quick release




Compass 4 fixed




Compass power mount
Power mount hardware has 8” of height adjustment in ½” increments without using tools.




4-point




Permobil Direct Mount
For use on Permobil Corpus power wheelchairs only



UniTrack Assembly
Included with Direct Mount




Section 2 | *Accessories*




Privacy shield
Provides coverage between the bottom of the back support to the seat base of their wheelchair. Attaches directly to the canes of a wheelchair.

Section 3 | *BodiLink lateral trunk support hardware*

GT2 mounting hardware, slot mount, swing-away




GT2 mounting hardware, slot mount, fixed




Section 4 | *BodiLink lateral trunk support pads*


Premium pads, supportive




Premium pads, active



Basic pads, supportive




Basic pads, active




Section 5 | *Bodilink lateral trunk support additional hardware options*

Extra hardware links




Long mounting bracket

Adds 1” of additional adjustment per bracket




Hardware cover



Section 6 | *Bodilink head support pads*

Basic pad sizes




Length(L) refers to the actual dimension of the support from bottom to top edge



Basic head support pads can be formed to any shape using the moldable steel insert inside the pad.

Premium pad sizes



Length(L) refers to the actual dimension of the support from bottom to top edge



Premium head support pad angles can be formed by adjusting the plates in the pad

Section 7 | *BodiLink head support hardware options*

GT hardware, removable

Range of motion


GT model	Total # of links	Total link length	Max. anterior reach	Max. posterior reach	Max. lateral reach
LAT	5	10.5”	10”	2”	7.5”
MID	3	8.5”	8”	2”	N/A



BodiLink head support hardware: PT hardware, removable

Range of motion


PT model	Total # of links	Total link length	Max anterior reach	Max posterior reach	Max lateral reach
MID 6	2	6”	5.5”	1”	N/A
LAT 6	3	6”	5.5”	1”	6”
MID 8	2	8”	7.5”	1”	N/A
LAT 8	3	8”	7.5”	1”	8”
MAX 8	4	8”	7.5”	3”	8”




PT hardware, flip down hardware, removable
Flip down hardware has 7" stem



BodiLink head support hardware: PT hardware, removable stem length option
3" stem is recommended for use with power mounting where there may be special constraints with power chair functions



Mounting options




Mounting block Posterior A
For Comfort Unite, Embrace and Visco backs as well as Comfort Relief with Compass power mount




Mounting block Permobil C
For use with Permobil Direct Mount and Unitrack assembly

Extra hardware

¼" spacer plate



Extra links: 2", 3", 4"



Tool-less wing lock

