Order form

Comfort cushions



Dealer Name*: Account #*:

City, Province, Postal Code*:

Phone: Email:

ATP: Rep Name:

P.O. Number: Quote Number:

Client Reference*: Client Weight*:

***Required field**

Special considerations that need to be addressed here (e.g., diagnosis):

NOTE | *Refer to reference pages at end of order form for additional imformation and images.*

Curve Cushion

Highlighted box indicates quick ship is available in this size

Depth	8”	9”	10”	11”	12”	13”	14”	15”	16”	17”	18”	19”	20”	21”	22”	23”	24”
20”								☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $304.00	☐ $304.00	☐ $304.00	☐ $304.00
19”								☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $304.00	☐ $304.00	☐ $304.00	☐ $304.00
18”							☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $304.00	☐ $304.00	☐ $304.00	☐ $304.00
17”						☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $304.00	☐ $304.00	☐ $304.00	☐ $304.00
16”					☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $304.00	☐ $304.00	☐ $304.00	☐ $304.00
15”			☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00						
14”			☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00								
13”		☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00										
12”	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00										
11”	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00										
10”	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00										
9”	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00	☐ $177.00												
8”	☐ $177.00	☐ $177.00	☐ $177.00														

Width

Cover and liner options

☐	CU-FV-WWDD	**COMFORT-TEK ®** *For fluid protection, an easily cleaned surface*	NC
☐	CU-SV-WWDD	**STRETCH-AIR ®** *For comfort and heat dissipation*	NC
☐	INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover and liner options | *Optional*

One extra outer cover	☐ COMFORT-TEK (CU-FV-WWDD-CV)	☐ STRETCH-AIR (CU-XV-WWDD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (CU-FV-WWDD-CV)	☐ STRETCH-AIR (CU-XV-WWDD-CV)	$264.00

☐ BINC-CU-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-CU-xxWyyD	**Two extra incontinence liners**	$252.00

Ascent Cushion

Depth	16”	17”	18”	19”	20”	21”	22”	23”	24”
20”	☐ $420.00	☐ $420.00	☐ $420.00	☐ $420.00	☐ $420.00	☐ $578.00	☐ $578.00	☐ $578.00	☐ $578.00
19”	☐ $420.00	☐ $420.00	☐ $420.00	☐ $420.00	☐ $420.00	☐ $578.00	☐ $578.00	☐ $578.00	☐ $578.00
18”	☐ $420.00	☐ $420.00	☐ $420.00	☐ $420.00	☐ $420.00	☐ $578.00	☐ $578.00	☐ $578.00	☐ $578.00
17”	☐ $420.00	☐ $420.00	☐ $420.00	☐ $420.00	☐ $420.00	☐ $578.00	☐ $578.00	☐ $578.00	☐ $578.00
16”	☐ $420.00	☐ $420.00	☐ $420.00	☐ $420.00	☐ $420.00	☐ $578.00	☐ $578.00	☐ $578.00	☐ $578.00

Width

Cover options

☐ AC-GF-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ AC-GS-WWDD	**STRETCH-AIR** *For comfort and heat dissipation*	NC
☐ INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover and liner options | *Optional*

One extra outer cover	☐ COMFORT-TEK (AC-xxWyyD-CV)	☐ STRETCH-AIR (AC-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (AC-xxWyyD-CV)	☐ STRETCH-AIR (AC-xxWyyD-CV)	$264.00

☐ BINC-AC-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-AC-xxWyyD	**Two extra incontinence liners**	$252.00

Comfort Embrace™ Cushion Zero-Elevation

Highlighted box indicates quick ship is available in this size

Depth	10”	11”	12”	13”	14”	15”	16”	17”	18”	19”	20”	21”	22”	23”	24”
24”									☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
23”									☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
22”							☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
21”						☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
20”					☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
19”					☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
18”			☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
17”		☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
16”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
15”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00				
14”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00					
13”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00							
12”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00							
11”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00									
10”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00									

Width

Cover and liner options

☐ EM-FZ-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ EM-SZ-WWDD	**STRETCH-AIR** *For comfort and heat dissipation. Not avaible with quick ship option.*	NC
☐ INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (EM-Z-xxWyyD-CV) ☐ STRETCH-AIR (EM-Z-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (EM-Z-xxWyyD-CV) ☐ STRETCH-AIR (EM-Z-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00

☐ BINC-EMZ-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-EMZ-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐ GROWTH-NOTCH — **Growth notches** — NC

If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.

☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)

☐ RAIL-CUT — **Rail cuts** — $58.00

A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.

☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5)

☐ 2" high (RAIL-CUT-2)

☐ COLOR — **Colored piping on cover** — $58.00

Accent colors can be applied to the cushion via colored piping.

☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK)

☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)

Comfort Embrace™ Cushion Zero-Elevation with GlideWear

Depth \ Width	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
24"									☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
23"									☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
22"							☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
21"						☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
20"					☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
19"					☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
18"			☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
17"		☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
16"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
15"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00				
14"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00					
13"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00							
12"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00							
11"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00									
10"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00									

Cover and liner options

☐	EM-FQZ-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐	EM-SQZ-WWDD	**STRETCH-AIR** *For comfort and heat dissipation*	NC
☐	INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (EM-Z-xxWyyD-CV) ☐ STRETCH-AIR (EM-Z-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (EM-Z-xxWyyD-CV) ☐ STRETCH-AIR (EM-Z-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00

☐	BINC-EMZ-xxWyyD	**One extra incontinence liner**	$126.00
☐	BINC-EMZ-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐	GROWTH-NOTCH	**Growth notches** *If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.* ☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)	NC
☐	RAIL-CUT	**Rail cuts** *A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.* ☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5) ☐ 2" high (RAIL-CUT-2)	$58.00
☐	COLOR	**Colored piping on cover** *Accent colors can be applied to the cushion via colored piping.* ☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK) ☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)	$58.00

Comfort Embrace™ Cushion Anti-Thrust

Highlighted box indicates quick ship is available in this size

Depth	10”	11”	12”	13”	14”	15”	16”	17”	18”	19”	20”	21”	22”	23”	24”
24”									☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
23”									☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
22”							☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
21”						☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
20”					☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
19”					☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
18”			☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
17”		☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
16”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $718.00	☐ $718.00	☐ $718.00	☐ $718.00
15”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00				
14”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00					
13”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00							
12”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00							
11”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00									
10”	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00	☐ $508.00									

Width

Cover and liner options

☐ EM-FA-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ EM-SA-WWDD	**STRETCH-AIR** *For comfort and heat dissipation. Not availible with quick ship option.*	NC
☐ INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (EM-A-xxWyyD-CV) ☐ STRETCH-AIR (EM-A-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (EM-A-xxWyyD-CV) ☐ STRETCH-AIR (EM-A-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00
☐ BINC-EMA-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-EMA-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐ GROWTH-NOTCH	**Growth notches** *If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.* ☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)	NC
☐ RAIL-CUT	**Rail cuts** *A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.* ☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5) ☐ 2" high (RAIL-CUT-2)	$58.00
☐ COLOR	**Colored piping on cover** *Accent colors can be applied to the cushion via colored piping.* ☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK) ☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)	$58.00

Comfort Embrace™ Cushion Anti-Thrust with GlideWear

Depth	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
24"									☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
23"									☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
22"							☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
21"						☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
20"					☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
19"					☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
18"			☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
17"		☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00
16"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00
15"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00				
14"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00					
13"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00							
12"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00							
11"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00									
10"	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00	☐ $536.00									

Width

Cover and liner options

☐ EM-FQA-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ EM-SQA-WWDD	**STRETCH-AIR** *For comfort and heat dissipation*	NC
☐ INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (EM-A-xxWyyD-CV) ☐ STRETCH-AIR (EM-A-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (EM-A-xxWyyD-CV) ☐ STRETCH-AIR (EM-A-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00
☐ BINC-EMA-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-EMA-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional*

☐ GROWTH-NOTCH	**Growth notches** *If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.* ☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)	NC
☐ RAIL-CUT	**Rail cuts** *A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.* ☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5) ☐ 2" high (RAIL-CUT-2)	$58.00
☐ COLOR	**Colored piping on cover** *Accent colors can be applied to the cushion via colored piping.* ☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK) ☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)	$58.00

Comfort M2 Cushion Zero-Elevation

Highlighted box indicates quick ship is available in this size

Depth

Depth	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
24"									☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
23"									☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
22"							☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
21"						☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
20"					☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
19"					☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
18"			☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
17"		☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
16"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
15"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00				
14"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00					
13"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00							
12"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00							
11"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00									
10"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00									

Width

Cover and liner options

☐	M2-F-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐	M2-S-WWDD	**STRETCH-AIR** *For comfort and heat dissipation. Not availible with quick ship option.*	NC
☐	INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (M2-Z-xxWyyD-CV) ☐ STRETCH-AIR (M2-Z-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (M2-Z-xxWyyD-CV) ☐ STRETCH-AIR (M2-Z-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00

☐	BINC-M2-xxWyyD	**One extra incontinence liner**	$126.00
☐	BINC-M2-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐	GROWTH-NOTCH	**Growth notches** *If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.* ☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)	NC
☐	RAIL-CUT	**Rail cuts** *A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.* ☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5) ☐ 2" high (RAIL-CUT-2)	$58.00
☐	COLOR	**Colored piping on cover** *Accent colors can be applied to the cushion via colored piping.* ☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK) ☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)	$58.00

Comfort M2 Cushion Zero-Elevation with GlideWear

Depth	10”	11”	12”	13”	14”	15”	16”	17”	18”	19”	20”	21”	22”	23”	24”
24”									☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
23”									☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
22”							☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
21”						☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
20”					☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
19”					☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
18”			☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
17”		☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
16”	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
15”	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00				
14”	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00					
13”	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00							
12”	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00							
11”	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00									
10”	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00									

Width

Cover options

☐ M2-FQ-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ M2-SQ-WWDD	**STRETCH-AIR** *For comfort and heat dissipation*	NC
☐ INCON-LINER	**Incontinence liner** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (M2-Z-xxWyyD-CV) ☐ STRETCH-AIR (M2-Z-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (M2-Z-xxWyyD-CV) ☐ STRETCH-AIR (M2-Z-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00
☐ BINC-M2-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-M2-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐ GROWTH-NOTCH	**Growth notches** *If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.* ☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)	NC
☐ RAIL-CUT	**Rail cuts** *A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.* ☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5) ☐ 2" high (RAIL-CUT-2)	$58.00
☐ COLOR	**Colored piping on cover** *Accent colors can be applied to the cushion via colored piping.* ☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK) ☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green(GREEN)	$58.00

Comfort M2 Cushion Anti-Thrust

Highlighted box indicates quick ship is available in this size

Depth	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
24"									☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
23"									☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
22"							☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
21"						☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
20"					☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
19"					☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
18"			☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
17"		☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
16"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $762.00	☐ $762.00	☐ $762.00	☐ $762.00
15"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00				
14"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00					
13"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00							
12"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00							
11"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00									
10"	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00	☐ $580.00									

Width

Cover and liner options | *Required, choose one.*

☐ M2-FA-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ M2-SA-WWDD	**Stretch Air** *For comfort and heat dissipation. Not availible with quick ship option.*	NC
☐ INCON-LINER	**Incontinence liner** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (M2-A-xxWyyD-CV) ☐ STRETCH-AIR (M2-A-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (M2-A-xxWyyD-CV) ☐ STRETCH-AIR (M2-A-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00

☐ BINC-M2A-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-M2A-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐ GROWTH-NOTCH — **Growth notches** — NC

If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.

☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)

☐ RAIL-CUT — **Rail cuts** — $58.00

A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.

☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5)
☐ 2" high (RAIL-CUT-2)

☐ COLOR — **Colored piping on cover** — $58.00

Accent colors can be applied to the cushion via colored piping.

☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK)
☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)

Comfort M2 Cushion Anti-Thrust with GlideWear

Depth	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
24"									☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
23"									☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
22"							☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
21"						☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
20"					☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
19"					☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
18"			☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
17"		☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
16"	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $789.00	☐ $789.00	☐ $789.00	☐ $789.00
15"	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00				
14"	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00					
13"	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00							
12"	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00							
11"	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00									
10"	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00	☐ $607.00									

Width

Cover and liner options

☐ M2-FQA-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ M2-SQA-WWDD	**Stretch Air** *For comfort and heat dissipation*	NC
☐ INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (M2-A-xxWyyD-CV) ☐ STRETCH-AIR (M2-A-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (M2-A-xxWyyD-CV) ☐ STRETCH-AIR (M2-A-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00
☐ BINC-M2A-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-M2A-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐ GROWTH-NOTCH	**Growth notches** *If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.* ☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)	NC
☐ RAIL-CUT	**Rail cuts** *A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.* ☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5) ☐ 2" high (RAIL-CUT-2)	$58.00
☐ COLOR	**Colored piping on cover** *Accent colors can be applied to the cushion via colored piping.* ☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK) ☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)	$58.00

Saddle 7 Series Cushion without Quadrangel

Part number: 53S-WWDD-B

Includes Solace cover and Kwik Strap

	22"	24"	26"	28"	30"
20"	☐ $818.00	☐ $818.00	☐ $818.00	☐ $818.00	☐ $818.00
18"	☐ $818.00	☐ $818.00	☐ $818.00	☐ $818.00	☐ $818.00

Width

Additional cover options | *Optional*

☐ 53S-WWDD-B-CV	**One extra Solace outer cover**	$132.00
☐ 53S-WWDD-B-CV	**Two extra Solace outer covers**	$264.00

Saddle 7 Series Cushion with Quadragel

Part number: 53SG-WWDD-B

Includes Solace cover and Kwik Strap

	22”	24”	26”	28”	30”
20”	☐ $1,133.00	☐ $1,133.00	☐ $1,133.00	☐ $1,133.00	☐ $1,133.00
18”	☐ $1,133.00	☐ $1,133.00	☐ $1,133.00	☐ $1,133.00	☐ $1,133.00

Width

Additional cover options | *Optional*

☐ 53S-WWDD-B-CV	**One extra Solace outer cover**	$132.00
☐ 53S-WWDD-B-CV	**Two extra Solace outer covers**	$264.00

Comfort M2 Cushion with ATI Zero-Elevation

Depth	10”	11”	12”	13”	14”	15”	16”	17”	18”	19”	20”	21”	22”	23”	24”
24”									☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
23”									☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
22”							☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
21”						☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
20”					☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
19”					☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
18”			☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
17”		☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
16”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
15”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00				
14”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00					
13”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00							
12”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00							
11”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00									
10”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00									

Width

Cover options

☐ M2-FM-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ M2-SM-WWDD	**STRETCH-AIR** *For comfort and heat dissipation*	NC
☐ INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (M2-Z-xxWyyD-CV) ☐ STRETCH-AIR (M2-Z-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (M2-Z-xxWyyD-CV) ☐ STRETCH-AIR (M2-Z-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00
☐ BINC-M2-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-M2-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐ GROWTH-NOTCH — **Growth notches** — NC

If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.

☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)

☐ RAIL-CUT — **Rail cuts** — $58.00

A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.

☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5)

☐ 2" high (RAIL-CUT-2)

☐ COLOR — **Colored piping on cover** — $58.00

Accent colors can be applied to the cushion via colored piping.

☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK)

☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green(GREEN)

Comfort M2 Cushion with ATI Zero-Elevation with GlideWear

Depth \ Width	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
24"									☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
23"									☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
22"							☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
21"						☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
20"					☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
19"					☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
18"			☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
17"		☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
16"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
15"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00				
14"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00					
13"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00							
12"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00							
11"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00									
10"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00									

Cover options

☐ M2-FQM-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ M2-SQM-WWDD	**STRETCH-AIR** *For comfort and heat dissipation*	NC
☐ INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (M2-Z-xxWyyD-CV) ☐ STRETCH-AIR (M2-Z-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (M2-Z-xxWyyD-CV) ☐ STRETCH-AIR (M2-Z-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00
☐ BINC-M2-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-M2-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐ GROWTH-NOTCH	**Growth notches** *If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.* ☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)	NC
☐ RAIL-CUT	**Rail cuts** *A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.* ☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5) ☐ 2" high (RAIL-CUT-2)	$58.00
☐ COLOR	**Colored piping on cover** *Accent colors can be applied to the cushion via colored piping.* ☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK) ☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)	$58.00

Comfort M2 Cushion with ATI Anti-Thrust

Depth	10”	11”	12”	13”	14”	15”	16”	17”	18”	19”	20”	21”	22”	23”	24”
24”									☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
23”									☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
22”							☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
21”						☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
20”					☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
19”					☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
18”			☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
17”		☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
16”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
15”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00				
14”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00					
13”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00							
12”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00							
11”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00									
10”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00									

Width

Cover options

☐ M2-FAM-WWDD48	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ M2-SAM-WWDD49	**STRETCH-AIR** *For comfort and heat dissipation*	NC
☐ INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (M2-A-xxWyyD-CV) ☐ STRETCH-AIR (M2-A-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (M2-A-xxWyyD-CV) ☐ STRETCH-AIR (M2-A-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00
☐ BINC-M2A-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-M2A-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐	GROWTH-NOTCH	**Growth notches** *If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.* ☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)	NC
☐	RAIL-CUT	**Rail cuts** *A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.* ☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5) ☐ 2" high (RAIL-CUT-2)	$58.00
☐	COLOR	**Colored piping on cover** *Accent colors can be applied to the cushion via colored piping.* ☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK) ☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)	$58.00

Comfort M2 Cushion with ATI Anti-Thrust with GlideWear

Depth	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
24"									☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
23"									☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
22"							☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
21"						☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
20"					☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
19"					☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
18"			☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
17"		☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00
16"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $993.00	☐ $993.00	☐ $993.00	☐ $993.00
15"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00				
14"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00					
13"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00							
12"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00							
11"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00									
10"	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00	☐ $801.00									

Width

Cover options

☐	M2-FQAM-WWDD51	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐	M2-SQAM-WWDD52	**STRETCH-AIR** *For comfort and heat dissipation*	NC
☐	INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (M2-A-xxWyyD-CV) ☐ STRETCH-AIR (M2-A-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (M2-A-xxWyyD-CV) ☐ STRETCH-AIR (M2-A-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00

☐ BINC-M2A-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-M2A-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐ GROWTH-NOTCH — **Growth notches** — NC

If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.

☐ 1" Deep(GROWTH-NOTCH-1) ☐ 2" Deep(GROWTH-NOTCH-2)

☐ RAIL-CUT — **Rail cuts** — $58.00

A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.

☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5)
☐ 2" high (RAIL-CUT-2)

☐ COLOR — **Colored piping on cover** — $58.00

Accent colors can be applied to the cushion via colored piping.

☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK)
☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)

Comfort Embrace™ Cushion with ATI Zero-Elevation

Depth	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
24"									☐ $772.00	☐ $772.00	☐ $772.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
23"									☐ $772.00	☐ $772.00	☐ $772.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
22"							☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
21"						☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
20"					☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
19"					☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
18"			☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
17"		☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
16"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
15"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00				
14"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00					
13"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00							
12"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00							
11"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00									
10"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00									

Width

Cover options

☐ EM-FZM-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ EM-SZM-WWDD	**STRETCH-AIR** *For comfort and heat dissipation*	NC
☐ INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (EM-Z-xxWyyD-CV) ☐ STRETCH-AIR (EM-Z-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (EM-Z-xxWyyD-CV) ☐ STRETCH-AIR (EM-Z-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00
☐ BINC-EMA-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-EMA-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐ GROWTH-NOTCH	**Growth notches** *If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.* ☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)	NC
☐ RAIL-CUT	**Rail cuts** *A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.* ☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5) ☐ 2" high (RAIL-CUT-2)	$58.00
☐ COLOR	**Colored piping on cover** *Accent colors can be applied to the cushion via colored piping.* ☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK) ☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)	$58.00

Comfort Embrace™ Cushion with ATI Zero-Elevation with GlideWear

Depth	10”	11”	12”	13”	14”	15”	16”	17”	18”	19”	20”	21”	22”	23”	24”
24”									☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
23”									☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
22”							☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
21”						☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
20”					☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
19”					☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
18”			☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
17”		☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
16”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
15”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00				
14”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00					
13”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00							
12”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00							
11”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00									
10”	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00									

Width

Cover options

☐ EM-FQZM-WWDD57	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ EM-SQZM-WWDD58	**STRETCH-AIR** *For comfort and heat dissipation*	NC
☐ INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (EM-Z-xxWyyD-CV) ☐ STRETCH-AIR (EM-Z-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (EM-Z-xxWyyD-CV) ☐ STRETCH-AIR (EM-Z-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00
☐ BINC-EMZ-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-EMZ-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐ GROWTH-NOTCH	**Growth notches** *If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.* ☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)	NC
☐ RAIL-CUT	**Rail cuts** *A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.* ☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5) ☐ 2" high (RAIL-CUT-2)	$58.00
☐ COLOR	**Colored piping on cover** *Accent colors can be applied to the cushion via colored piping.* ☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK) ☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)	$58.00

Comfort Embrace™ Cushion with ATI Anti-Thrust

Depth	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
24"									☐ $745.00	☐ $745.00	☐ $745.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
23"									☐ $745.00	☐ $745.00	☐ $745.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
22"							☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
21"						☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
20"					☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
19"					☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
18"			☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
17"		☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $938.00	☐ $938.00	☐ $938.00
16"	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $938.00	☐ $938.00	☐ $938.00	☐ $938.00
15"	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00				
14"	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00					
13"	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00							
12"	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00							
11"	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00									
10"	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00	☐ $745.00									

Width

Cover options

☐ EM-FAM-WWDD	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ EM-SAM-WWDD	**STRETCH-AIR** *For comfort and heat dissipation*	NC
☐ INCON-LINER	**Incontinence liner (on chair)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (EM-A-xxWyyD-CV) ☐ STRETCH-AIR (EM-A-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (EM-A-xxWyyD-CV) ☐ STRETCH-AIR (EM-A-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00

☐ BINC-EMA-xxWyyD	**One extra incontinence liner**	$126.00
☐ BINC-EMA-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐ GROWTH-NOTCH — **Growth notches** — NC

If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.

☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)

☐ RAIL-CUT — **Rail cuts** — $58.00

A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.

☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5)
☐ 2" high (RAIL-CUT-2)

☐ COLOR — **Colored piping on cover** — $58.00

Accent colors can be applied to the cushion via colored piping.

☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK)
☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)

Comfort Embrace™ Cushion with ATI Anti-Thrust with GildeWear

Depth	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"
24"									☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
23"									☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
22"							☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
21"						☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
20"					☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
19"					☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
18"			☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
17"		☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
16"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $965.00	☐ $965.00	☐ $965.00	☐ $965.00
15"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00				
14"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00					
13"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00							
12"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00							
11"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00									
10"	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00	☐ $772.00									

Width

Cover options

☐	EM-FQAM-WWDD63	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐	EM-SQAM-WWDD64	**STRETCH-AIR** *For comfort and heat dissipation*	NC
☐	INCON-LINER	**Incontinence liner (on cushion)** *For extra incontinence protection*	$56.00

Additional cover options | *Optional*

One extra outer cover	☐ COMFORT-TEK (EM-A-xxWyyD-CV) ☐ STRETCH-AIR (EM-A-xxWyyD-CV)	$132.00
Two extra outer covers	☐ COMFORT-TEK (EM-A-xxWyyD-CV) ☐ STRETCH-AIR (EM-A-xxWyyD-CV)	$264.00
One extra GlideWear cover	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$158.00
Two extra GlideWear covers	☐ COMFORT-TEK with GlideWear (GLIDEWEAR-C) ☐ STRETCH-AIR with GlideWear (GLIDEWEAR-C)	$316.00

☐	BINC-EMA-xxWyyD	**One extra incontinence liner**	$126.00
☐	BINC-EMA-xxWyyD	**Two extra incontinence liners**	$252.00

Cushion options | *Optional. Refer to reference pages for details and images.*

☐	GROWTH-NOTCH	**Growth notches** *If growth notches are selected, a 1 1/4" wide cut will be taken from each side to allow cushion to fit between chair canes.* ☐ 1" Deep (GROWTH-NOTCH-1) ☐ 2" Deep (GROWTH-NOTCH-2)	NC
☐	RAIL-CUT	**Rail cuts** *A standard 1.5" wide cut will be made and you decide how high the cut will need to be. All cuts run the entire depth of the cushion.* ☐ 1" high (RAIL-CUT-1.0) ☐ 1.5" high (RAIL-CUT-1.5) ☐ 2" high (RAIL-CUT-2)	$58.00
☐	COLOR	**Colored piping on cover** *Accent colors can be applied to the cushion via colored piping.* ☐ Red (RED) ☐ Purple (PURPLE) ☐ Pink (PINK) ☐ Orange (ORANGE) ☐ Blue (BLUE) ☐ Green (GREEN)	$58.00

Cushion accessories

Refer to reference pages for details and images.

☐	RSSP-1-WWDD	**Removable solid seat pan & hardware**	$690.00
Rigid insert	☐ Rigid insert glued (RIDWWDD) ☐ Rigid insert not glued (RIDWWDD)		$96.00

BodiLink® pelvic/thigh support hardware and pads

Refer to reference pages for details and images.

Basic pad shape & size | *Priced per pad.*

	Left *Must enter number of pads*	Right *Must enter number of pads*	
3”L x 4”D	☐ BL-LPTSP1Z-3L4D-LH	☐ BL-LPTSP1Z-3L4D-RH	$139.00
3”L x 5”D	☐ BL-LPTSP1Z-3L5D-LH	☐ BL-LPTSP1Z-3L5D-RH	$139.00
4”L X 4”D	☐ BL-LPTSP1Z-4L4D-LH	☐ BL-LPTSP1Z-4L4D-RH	$139.00
4”L x 6”D	☐ BL-LPTSP1Z-4L6D-LH	☐ BL-LPTSP1Z-4L6D-RH	$139.00
4”L x 8”D	☐ BL-LPTSP1Z-4L8D-LH	☐ BL-LPTSP1Z-4L8D-RH	$139.00
4”L x 10”D	☐ BL-LPTSP1Z-4L10D-LH	☐ BL-LPTSP1Z-4L10D-RH	$139.00
4”L x 12”D	☐ BL-LPTSP1Z-4L12D-LH	☐ BL-LPTSP1Z-4L12D-RH	$139.00
4”L x 14”D	☐ BL-LPTSP1Z-4L14D-LH	☐ BL-LPTSP1Z-4L14D-RH	$139.00
4”L x 16”D	☐ BL-LPTSP1Z-4L16D-LH	☐ BL-LPTSP1Z-4L16D-RH	$139.00
5”L x 7”D	☐ BL-LPTSP1Z-5L7D-LH	☐ BL-LPTSP1Z-5L7D-RH	$139.00

Basic pad cover options | *Required, must enter number of covers needed and match number of pads selected. Includes foam insert.*

☐ LPTS-P1-COMFORT	**COMFORT-TEK®** *For fluid protection & an easily cleaned surface* ☐ Left ☐ Right	NC
☐ LPTS-P1-STRETCH	**STRETCH-AIR®** *For comfort & heat dissipation* ☐ Left ☐ Right	NC
☐ LPTS-P1-GLIDE	**GlideWear®** *For skin protection & shear reduction. NOT fluid resistant.* ☐ Left ☐ Right	$36.00

Premium pad shape & size | *Priced per pad.*

	Left *Must enter the number of pads needed*	Right *Must enter the number of pads needed*	
3.5”L x 4”D	☐ BL-LPTSP2Z-3L4D-LH	☐ BL-LPTSP2Z-3L4D-RH	$150.00
3.5”L x 8”D	☐ BL-LPTSP2Z-3L8D-LH	☐ BL-LPTSP2Z-3L8D-RH	$150.00
3.5”L x 12”D	☐ BL-LPTSP2Z-3L12D-LH	☐ BL-LPTSP2Z-3L12D-RH	$150.00
5.5” x 6”D	☐ BL-LPTSP2Z-5L6D-LH	☐ BL-LPTSP2Z-5L6D-RH	$150.00

Premium pad cover options | *Required, must enter number of covers needed and match number of pads selected. Includes foam insert.*

☐ LPTS-P2-COMFORT-FM	**COMFORT-TEK** *For fluid protection & an easily cleaned surface* ☐ Left ☐ Right	
☐ LPTS-P2-STRETCH-FM	**STRETCH-AIR** *For comfort & heat dissapation* ☐ Left ☐ Right	
☐ LPTS-P2-GLIDE-FM	**GlideWear** *For skin protection & shear reduction. NOT fluid resistant.* ☐ Left ☐ Right	$36.00

GT mounting hardware, slot mount only

Fixed	☐ Left (BL-LPTS-GT2FXS-LH)	☐ Right (BL-LPTS-GT2FXS-RH)		$218.00
Removable	☐ Left (BL-LPTS-GT2RMS-LH)	☐ Right (BL-LPTS-GT2RMS-RH)		$260.00

TT mounting hardware | *Required, must match number of pads selected. If selecting power mount hardware, must choose power chair mounting option.*

	Small *Recommened cushion thickness: 2.0”-3.0”*	**Medium** *Recommened cushion thickness: 3.0”-4.5”*	**Large** *Recommened cushion thickness: 4.0”-6.5”*	
TT hardware, slot mount				
Fixed-Left	☐ BL-LPTS-TT1FXSL1-LH	☐ BL-LPTS-TT1FXSL2-LH	☐ BL-LPTS-TT1FXSL3-LH	$218.00
Fixed-Right	☐ BL-LPTS-TT1FXSL1-RH	☐ BL-LPTS-TT1FXSL2-RH	☐ BL-LPTS-TT1FXSL3-RH	$218.00
Removable-Left	☐ BL-LPTS-TT1RMSL1-LH	☐ BL-LPTS-TT1RMSL2-LH	☐ BL-LPTS-TT1RMSL3-LH	$311.00
Removable-Right	☐ BL-LPTS-TT1RMSL1-RH	☐ BL-LPTS-TT1RMSL2-RH	☐ BL-LPTS-TT1RMSL3-RH	$311.00
TT hardware, power mount				
Fixed-Left	☐ BL-LPTS-TT1FXPWL1-LH	☐ BL-LPTS-TT1FXPWL2-LH	☐ BL-LPTS-TT1FXPWL3-LH	$282.00
Fixed-Right	☐ BL-LPTS-TT1FXPWL1-RH	☐ BL-LPTS-TT1FXPWL2-RH	☐ BL-LPTS-TT1FXPWL3-RH	$282.00
Removable-Left	☐ BL-LPTS-TT1RMPWL1-LH	☐ BL-LPTS-TT1RMPWL2-LH	☐ BL-LPTS-TT1RMPWL3-LH	$368.00
Removable-Right	☐ BL-LPTS-TT1RMPWL1-RH	☐ BL-LPTS-TT1RMPWL2-RH	☐ BL-LPTS-TT1RMPWL3-RH	$368.00

Power mount option

- ☐ Quantum w/ Corpus 3G (LPTS-TT-PW1)
- ☐ Permobil w/ Corpus 3G (LPTS-TT-PW2)
- ☐ Quickie w/ 3.7 Power (LPTS-TT-PW3)
- ☐ Rovi w/ Motion Concepts (LPTS-TT-PW4)
- ☐ Avid Rehab w. Power T & R (LPTS-TT-PW5)

Additional options

☐ Bl-LPTS-TT1LARM4 **TT long extension arm, 4”** $58.00
Price per arm

☐ Left ☐ Right ☐ Extra left ☐ Extra right

Bodilink® medial knee/thigh support hardware and pads

Refer to reference pages for details and images.

Pad shape & size | *Priced per pad.*

	Wedge	**Oval**	
2.5”W x 3.5”D	☐ BL-MKTSP1-3W4D	☐ BL-MKTSP2-3W4D	$203.00
3.5”W x 5”D	☐ BL-MKTSP1-4W5D	☐ BL-MKTSP2-4W5D	$203.00
4”W x 6”D	☐ BL-MKTSP1-4W6D		$203.00
5”W x 6”D		☐ BL-MKTSP2-5W6D	$203.00

Pad cover options | *Required, choose one.*

☐	MKTSP-COMFORT	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	
☐	MKTSP-STRETCH	**STRETCH-AIR** *For comfort and heat dissipation*	
☐	MKTSP-GLIDE	**GlideWear** *For skin protection & shear reduction. NOT fluid resistant*	$36.00

Swing away hardware

☐	BL-MKTS-ST1SASL1	**Small** *Recommended inferior thigh thickness: 1.0”-3.0”*	$280.00
☐	BL-MKTS-ST1SASL2	**Medium** *Recommended inferior thigh thickness: 2.0”-5.0”*	$280.00
☐	BL-MKTS-ST1SASL3	**Large** *Recommended inferior thigh thickness: 3.0”-7.0”*	$280.00

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- The client's name or code and the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____


Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.

References


Section 1 | *Cushion options*

Growth notches
If growth notches are selected, a 1 ¼" wide cut will be taken from each side to allow the cushion to fit between the wheelchair canes




Top view

Rail cuts
A standard 1 ½" wide cut will be made, the height of the cut is chosen by customer. All cuts run the entire depth of the cushion.




Front

Color piping
Accent colors can be applied to the cushion via colored piping.



Color location

Section 2 | *Cushion accessories*

Removable seat pan kit *(includes hardware)*
Includes a slotted aluminum pan and attaching hardware to accommodate 7/8" or 1" tubing. Kit includes two different cross bars to accommodate various wheelchair frame types. This will fit both folding and non-folding seat rails. Lateral and medial thigh support hardware can easily be attached.




Rigid insert
3/16" rigidizing board.




Section 3 | *Bodilink lateral/pelvic thigh supports*


Basic pads
Length (L) refers to the actual size dimension of the support from top to bottom edge. Depth (D) refers to the actual size dimension from anterior to posterior edge




Premium pads
Length (L) refers to the actual size dimension of the support from top to bottom edge. Depth (D) refers to the actual size dimension from anterior to posterior edge.




GT hardware



Bodilink lateral/pelvic thigh support hardware: Slot mount




Hardware size VS. cushion thickness




Hardware size	Max cushion thickness clearance	Recommended cushion thickness
Small	3.0"	2.0"- 3.0"
Medium	4.5"	3.0"- 4.5"
Large	6.5"	4.0"- 6.5"

Bodilink lateral/pelvic thigh support hardware: Power mount




Hardware size	Max cushion thickness clearance	Recommended cushion thickness
Small	3.5”	2.0”- 3.5”
Medium	5.5”	3.0”- 5.5”
Large	7.5”	4.0”- 7.5”

TT long extension arm 4”





Section 4 | *Bodilink medial knee/thigh supports*

Premium pad size/shape



Swing-away hardware





Size	Recommended inferior thigh thickness
Small	1.00”- 3.00”
Medium	2.00”- 5.00”
Large	3.00”- 7.00”