Order form

Custom Comfort Unite 1 back support



Dealer Name*: _____ Account #*: _____

City, Province, Post Code: _____

Phone: _____ Email: _____

ATP: _____ Rep Name: _____

P.O. Number: _____ Quote Number: _____

Client Reference*: _____ Client Weight: _____

Required field

Special considerations that need to be addressed here (e.g., diagnosis):

_____

Notes | *Refer to reference pages at end of order form for additional information and images.*

Seating system measurements

Seating system measurements are for reference only and do not affect selections made throughout the order form.

A - Back support width

B- Back support length

C - Width between lateral trunk supports


D - Width between lateral thigh supports

E - Seat width

F - Seat depth

G - Effective seat depth

H - Lateral trunk support height



Custom Unite 1 back support

Part number: C-U1-BSXXWYYL.

Length	12”	14”	16”	18”	20”	22”	24”	26”	28”
22”		☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00
20”		☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00
18”	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00
16”	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00
14”	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00
12”	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00	☐ $ 1,400.00		

Width

12”W x 16”L, and 12”W x 18”L not available with power mount hardware. 12L, 14L, and 12W not available with direct mount hardware. 24”, 26”, and 28” wide back supports require 4 point or power mount hardware.

Cover material selection | *Required, choose one. For back supports with custom ROHO insert, the four-way stretch cover needs to be selected.*

☐	CCD-COMFORT-TEK	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐	CCD-DIAMOND-SPACER	**Diamond Spacer** *For comfort and heat dissipation, not a fluid resistant material*	NC
☐	CCD-GLIDEWEAR-DIAMOND	**Diamond Spacer with GlideWear** *For skin protection and shear reduction. not a fluid resistant material*	$ 85.00
☐	CCD-FOUR-WAY	**Four-Way Stretch** *For use with custom ROHO air insert only. For optimal immersion and envelopment.*	NC

Custom foam selection


Select only one type of foam per layer. | *1/4” Crosslink foundation layer will come standard on all back supports. If Foam in Place Kit is selected, a middle and top layer needs to be selected. If Custom ROHO insert is being selected, do not select foam layers. For other foam combinations, call Customer Experience at 800-736-0925.*

Standard foam configuration:


12”-14” wide
0.25” Crosslink base layer
0.75” Comfort Medium foam

16”-28” wide
0.25” Crosslink base layer
1.25” Comfort Medium foam


Rear layer
REQUIRED
H.R. MSRP $0.00
Comfort MSRP $160.00
Visco MSRP $225.00



Middle layer
OPTIONAL
H.R. MSRP $0.00
Comfort MSRP $160.00
Visco MSRP $225.00



Topper layer
OPTIONAL
H.R. MSRP $0.00
Comfort MSRP $160.00
Visco MSRP $225.00



Foam type	Rear 1/4”	Rear 1/2”	Rear 3/4”	Rear 1”	Rear 1 1/4”	Middle 1/4”	Middle 1/2”	Middle 3/4”	Middle 1”	Middle 1 1/4”	Topper 1/4”	Topper 1/2”	Topper 3/4”	Topper 1”	Topper 1 1/4”
H.R. Medium	☐ CCD-REAR-HRM-0.25	☐ CCD-REAR-HRM-0.50	☐ CCD-REAR-HRM-0.75	☐ CCD-REAR-HRM-1.00	☐ CCD-REAR-HRM-1.25	☐ CCD-MID-HRM-0.25	☐ CCD-MID-HRM-0.50	☐ CCD-MID-HRM-0.75	☐ CCD-MID-HRM-1.00	☐ CCD-MID-HRM-1.25	☐ CCD-TOP-HRM-0.25	☐ CCD-TOP-HRM-0.50	☐ CCD-TOP-HRM-0.75	☐ CCD-TOP-HRM-1.00	☐ CCD-TOP-HRM-1.25
H.R. Soft	☐ CCD-REAR-HRS-0.25	☐ CCD-REAR-HRS-0.50	☐ CCD-REAR-HRS-0.75	☐ CCD-REAR-HRS-1.00	☐ CCD-REAR-HRS-1.25	☐ CCD-MID-HRS-0.25	☐ CCD-MID-HRS-0.50	☐ CCD-MID-HRS-0.75	☐ CCD-MID-HRS-1.0	☐ CCD-MID-HRS-1.25	☐ CCD-TOP-HRS-0.25	☐ CCD-TOP-HRS-0.50	☐ CCD-TOP-HRS-0.75	☐ CCD-TOP-HRS-1.00	☐ CCD-TOP-HRS-1.25
Comfort Medium	☐ CCD-REAR-COM-0.25	☐ CCD-REAR-COM-0.50	☐ CCD-REAR-COM-0.75	☐ CCD-REAR-COM-1.00	☐ CCD-REAR-COM-1.25	☐ CCD-MID-COM-0.25	☐ CCD-MID-COM-0.50	☐ CCD-MID-COM-0.75	☐ CCD-MID-COM-1.00	☐ CCD-MID-COM-1.25	☐ CCD-TOP-COM-0.25	☐ CCD-TOP-COM-0.50	☐ CCD-TOP-COM-0.75	☐ CCD-TOP-COM-1.00	☐ CCD-TOP-COM-1.25
Comfort Soft	☐ CCD-REAR-COS-0.25	☐ CCD-REAR-COS-0.50	☐ CCD-REAR-COS-0.75	☐ CCD-REAR-COS-1.00	☐ CCD-REAR-COS-1.25	☐ CCD-MID-COS-0.25	☐ CCD-MID-COS-0.50	☐ CCD-MID-COS-0.75	☐ CCD-MID-COS-1.00	☐ CCD-MID-COS-1.25	☐ CCD-TOP-COS-0.25	☐ CCD-TOP-COS-0.50	☐ CCD-TOP-COS-0.75	☐ CCD-TOP-COS-1.00	☐ CCD-TOP-COS-1.25
Visco Medium		☐ CCD-REAR-VIS-M-0.50		☐ CCD-REAR-VIS-M-1.00			☐ CCD-MID-VIS-M-0.50		☐ CCD-MID-VIS-M-1.00			☐ CCD-TOP-VIS-M-0.50		☐ CCD-TOP-VIS-M-1.00	
Visco Soft		☐ CCD-REAR-VIS-S-0.50		☐ CCD-REAR-VIS-S-1.00			☐ CCD-MID-VIS-S-0.50		☐ CCD-MID-VIS-S-1.00			☐ CCD-TOP-VIS-S-0.50		☐ CCD-TOP-VIS-S-1.00	
Visco X-Soft		☐ CCD-REAR-VIS-XS-0.50		☐ CCD-REAR-VIS-XS-1.00			☐ CCD-MID-VIS-XS-0.50		☐ CCD-MID-VIS-XS-1.00			☐ CCD-TOP-VIS-XS-0.50		☐ CCD-TOP-VIS-XS-1.00	
Visco XX-Soft		☐ CCD-REAR-VIS-XXS-0.50		☐ CCD-REAR-VIS-XXS-1.00			☐ CCD-MID-VIS-XXS-0.50		☐ CCD-MID-VIS-XXS-1.00			☐ CCD-TOP-VIS-XXS-0.50		☐ CCD-TOP-VIS-XXS-1.00	

Foam in place kit

Kit includes |


1) Back cushion cover with an optional 3" or 4" deep back, with a stretchy, pleated contact surface to allow for foam growth.
2) Two aluminum stays sewn into the rear of the cover which can be molded to support spinal curvature prior to pouring FIP chemicals.
3) Middle foam layer is cut to fit the shape of the back support and is placed loosely inside the cover. This will give the user some padding prior to pouring the foam.
4) Top foam layer is cut to fit the shape of the back support and is placed inside the front of the cover. This will give the final mold a plush feel for the user and will also smooth out any potential inconsistencies after pouring the foam.
5) Two plastic bags are included to contain the foam chemicals as they expand.

***Note:** Price of Visco layers will be added to the order separately.*

Bag pour size recommendations:

3" border thickness for frame widths 10" - 16", 1 bag
3" border thickness for frame widths 17"-22", 1 1/2 bag
4" border thickness for frame widths 10" - 14", 1 bag
4" border thickness for frame widths 15" - 20", 1 1/2 bag

MSRP:
1 Bag FIP:$650.00
1 1/2 Bag FIP:$800.00



Pour size	Extra soft memory		Soft memory foam		Medium soft foam	
	1 bag	1 1/2 bag	1 bag	1 1/2 bag	1 bag	1 1/2 bag
3" border thickness	☐ CCD-FIP3-1-XS	☐ CCD-FIP3-1.5-XS	☐ CCD-FIP3-1-S	☐ CCD-FIP3-1.5-S	☐ CCD-FIP3-1-M	☐ CCD-FIP3-1.5-M
4" border thickness	☐ CCD-FIP4-1-XS	☐ CCD-FIP4-1.5-XS	☐ CCD-FIP4-1-S	☐ CCD-FIP4-1.5-S	☐ CCD-FIP4-1-M	☐ CCD-FIP4-1.5-M

Custom ROHO Insert

Optional | *Completion of the additional Custom Unite 1 ROHO insert order form is required if selecting the Custom ROHO insert. Custom ROHO inserts include a four-way stretch cover.*

Select size range
Note: Shell size must match this order form.

☐ Custom ROHO insert, SMALL (12"- 14" wide)
(CCD-U1-ROHO1214)
MSRP: $1,500.00

☐ Custom ROHO insert, MEDIUM (16"- 20" wide)
(CCD-U1-ROHO1620)
MSRP: $1,500.00

☐ Custom ROHO insert, LARGE (22" - 26" wide)
(CCD-U1-ROHO2226)
MSRP: $1,500.00



Available custom ROHO insert sizes

Unite 1 Custom ROHO Insert								
LENGTH	22								
	20								
	18								
	16								
	14								
	12								
		12	14	16	18	20	22	24	26
		WIDTH							


Back support options

Moldable aluminum stays | *This back support comes standard with two aluminum stays. A standard alignment of the moldable stays will be chosen based on back support width if you choose not to select a specific width between the stays. For 12W back using custom ROHO insert, Max width in between stays is 1". For 14W back using custom ROHO insert, Max width in between stays is 3".*

Note: The stays are accessible while the user is seated against the back support through the zipper located at the top of the back support.

MSRP: $0.00

- ☐ A = 1" (CCD-STAY-1.0)
- ☐ A = 1 1/2" (CCD-STAY-1.5)
- ☐ A = 2" (standard alignment) (CCD-STAY-2.0)
- ☐ A = 2 1/2" (CCD-STAY-2.5)
- ☐ A = 3" (CCD-STAY-3.0)
- ☐ A = 3 1/2" (CCD-STAY-3.5)
- ☐ A = 4" (CCD-STAY-4.0)
- ☐ A = 4 1/2" (CCD-STAY-4.5)
- ☐ Specify custom spacing (CCD-STAY-OTHER)
- ☐ None (CCD-STAY-NONE)



Mounting Hardware | *Required, choose one.*

Compass 4 Quick Release
☐ Quick Release 2 Point (CCD-Q-COMPASS) ☐ Quick Release 4 Point (CCD-Q-COMPASS4) — NC

Compass 4 Fixed
☐ Fixed 2 Point (CCD-NQ-COMPASS) ☐ Fixed 4 Point (CCD-NQ-COMPASS4) — NC

Compass Power Mount — $ 125.00
- ☐ Permobil w/ Corpus 3G (CCD-PMB-PERMOBIL)
- ☐ Quantum w/ Tru Balance 3 or 4 (CCD-PMB-QUANTUM)
- ☐ Quickie w/ 3.7 Power (CCD-PMB-QUICKIE)
- ☐ Avid Rehab w/ Power T & R (CCD-PMB-AVID)
- ☐ ROVI w/ Motion Concepts (CCD-PMB-ROVI)

Using Power Mount hardware on ROVI requires mounting to the recessed planar interface plate available when ordering the Motion Concepts seating system. Drilling is required for mounting.

Permobil Direct Mount
☐ Permobil Direct Mount (CCD-DIRECT-MOUNT) — $ 125.00

Includes Unitrack assembly

Color Piping | *Optional. If no color is selected, the standard silver reflective piping is used.*

Choose one — NC
☐ Remove (CCD-REMOVE) ☐ Red (CCD-RED) ☐ Purple (CCD-PURPLE)
☐ Pink (CCD-PINK) ☐ Blue (CCD-BLUE) ☐ Green (CCD-GREEN)
☐ Orange (CCD-ORANGE)

Embroidery | *Optional, not available with GlideWear or four-way stretch cover. Includes corresponding color zipper pull.*

Thread color — NC
☐ Red (CCD-EMB-RED) ☐ Purple (CCD-EMB-PURPLE) ☐ Pink (CCD-EMB-PINK)
☐ Blue (CCD-EMB-BLUE) ☐ Orange (CCD-EMB-ORANGE) ☐ Green (CCD-EMB-GREEN)
☐ Silver (CCD-EMB-SILVER)

Font — NC
☐ Helvetica (CCD-EMB-HELVETICA) ☐ Garamond (-EMB-GARAMOND)
☐ Impress (CCD-EMB-IMPRESS) ☐ Kaufmann (CCD-EMB-KAUFMANN)

Font styles
HELVETICA
GARAMOND
Impress
Kaufmann

Word one — NC
☐ Input is case sensitive:

Word two (optional) — NC
☐ Input is case sensitive:

Other Accessories | *Optional*

☐	CCD-PV-SHIELD	**Privacy Shield** *Connects with hook & loop attachment*	$ 85.00

BodiLink lateral trunk supports & hardware

Optional, if selecting must choose pad shape & size, pad style, pad covers, and mounting hardware.

Basic pad shape & size | *Priced per pad.*

	Left *Must enter the number of pads needed*		**Right** *Must enter the number of pads needed*		
Supportive					
3”L x 4”D	☐	BL-LTSP1S-3L4D-LH	☐	BL-LTSP1S-3L4D-RH	$ 161.00
3”L x 5”D	☐	BL-LTSP1S-3L5D-LH	☐	BL-LTSP1S-3L5D-RH	$ 161.00
4”L x 5”D	☐	BL-LTSP1S-4L5D-LH	☐	BL-LTSP1S-4L5D-RH	$ 161.00
4”L x 6”D	☐	BL-LTSP1S-4L6D-LH	☐	BL-LTSP1S-4L6D-RH	$ 161.00
5”L x 7”D	☐	BL-LTSP1S-5L7D-LH	☐	BL-LTSP1S-5L7D-RH	$ 161.00
6”L x 6”D	☐	BL-LTSP1S-6L6D-LH	☐	BL-LTSP1S-6L6D-RH	$ 161.00
Active					
3”L x 4”D	☐	BL-LTSP1A-3L4D-LH	☐	BL-LTSP1A-3L4D-RH	$ 161.00
3”L x 5”D	☐	BL-LTSP1A-3L5D-LH	☐	BL-LTSP1A-3L5D-RH	$ 161.00
4”L x 5”D	☐	BL-LTSP1A-4L5D-LH	☐	BL-LTSP1A-4L5D-RH	$ 161.00
4”L x 6”D	☐	BL-LTSP1A-4L6D-LH	☐	BL-LTSP1A-4L6D-RH	$ 161.00
5”L x 7”D	☐	BL-LTSP1A-5L7D-LH	☐	BL-LTSP1A-5L7D-RH	$ 161.00
6”L x 6”D	☐	BL-LTSP1A-6L6D-LH	☐	BL-LTSP1a-6L6D-RH	$ 161.00

Basic pad style | *Must match number of pads selected.*

☐	LTS-P1-MOLD	**Moldable**	NC
		☐ Left ☐ Right	NC
☐	LTS-P1-RIGID	**Rigid**	NC
		☐ Left ☐ Right	NC

Basic pad cover options | *Must match number of pads selected.*

☐	LTS-P1-COMFORT	**COMFORT-TEK®** *For fluid protection & an easily cleaned surface*	NC
		☐ Left ☐ Right	
☐	LTS-P1-STRETCH	**STRETCH-AIR®** *For comfort & heat dissapation*	NC
		☐ Left ☐ Right	

Premium pad shape & size | *Priced per pad.*

	Left *Enter number of pads required*	Right *Enter number of pads required*	
Supportive			
3”L x 4” D	☐ BL-LTSP2S-3L4D-LH	☐ BL-LTSP2S-3L4D-RH	$ 218.00
3” L x 5.5” D	☐ BL-LTSP2S-3L5D-LH	☐ BL-LTSP2S-3L5D-RH	$ 218.00
4” L x 5” D	☐ BL-LTSP2S-4L5D-LH	☐ BL-LTSP2S-4L5D-RH	$ 218.00
4” L x 6.5” D	☐ BL-LTSP2S-4L6D-LH	☐ BL-LTSP2S-4L6D-RH	$ 218.00
5” L x 6.5” D	☐ BL-LTSP2S-5L6D-LH	☐ BL-LTSP2S-5L6D-RH	$ 218.00
6” L x 6.5” D	☐ BL-LTSP2S-6L6D-LH	☐ BL-LTSP2S-6L6D-RH	$ 218.00
Active			
3”L x 4” D	☐ BL-LTSP2A-3L4D-LH	☐ BL-LTSP2A-3L4D-RH	$ 218.00
3” L x 5.5” D	☐ BL-LTSP2A-3L5D-LH	☐ BL-LTSP2A-3L5D-RH	$ 218.00
4” L x 5” D	☐ BL-LTSP2A-4L5D-LH	☐ BL-LTSP2A-4L5D-RH	$ 218.00
4” L x 6.5” D	☐ BL-LTSP2A-4L6D-LH	☐ BL-LTSP2A-4L6D-RH	$ 218.00
5” L x 6.5” D	☐ BL-LTSP2A-5L6D-LH	☐ BL-LTSP2A-5L6D-RH	$ 218.00
6” L x 6.5” D	☐ BL-LTSP2A-6L6D-LH	☐ BL-LTSP2A-6L6D-RH	$ 218.00

Premium pad style | *Must match number of pads selected.*

☐ LTS-P2-MOLD	**Moldable**	NC
	☐ Left ☐ Right	NC
☐ LTS-P2-RIGID	**Rigid**	NC
	☐ Left ☐ Right	NC

Premium pad cover options | *Must match number of pads selected. Priced per selection.*

☐ LTS-P2-COMFORT	**COMFORT-TEK** *For fluid protection & an easily cleaned surface*	NC
	☐ Left ☐ Right	NC
☐ LTS-P2-STRETCH	**STRETCH-AIR** *For comfort & heat dissapation*	NC
	☐ Left ☐ Right	NC
☐ LTS-P2-GLIDE	**GlideWear®** *For skin protection & shear reduction. NOT fluid resistant*	$ 75.00
	☐ Left ☐ Right	
☐ LTS-P2-COMFORT-FM	**COMFORT-TEK with foam**	$ 30.00
	☐ Left ☐ Right	
☐ LTS-P2-STRETCH-FM	**STRETCH-AIR with foam**	$ 30.00
	☐ Left ☐ Right	
☐ LTS-P2-GLIDE-FM	**GlideWear with foam**	$ 103.00
	☐ Left ☐ Right	

GT2 mounting hardware | *Both fixed and swing away hardware styles are adjustable.*

	Left *Must match number of pads selected*	Right *Must match number of pads selected*	
Swing-away			
Slot Mount	☐ BL-LTS-GT2SAS-LH	☐ BL-LTS-GT2SAS-RH	$ 368.00
Track Mount	☐ BL-LTS-GT2SAT-LH	☐ BL-LTS-GT2SAT-RH	$ 368.00
Fixed			
Slot Mount	☐ BL-LTS-GT2FXS-LH	☐ BL-LTS-GT2FXS-RH	$ 202.00
Track Mount	☐ BL-LTS-GT2FXT-LH	☐ BL-LTS-GT2FXT-RH	$ 202.00

Additional hardware options | *Priced per selection.*

Extra GT2 hardware link	☐ Left (BL-LTS-GT2LINK)	☐ Right (BL-LTS-GT2LINK)	$ 58.00
Long mounting bracket	☐ Left (BL-LTS-GT2SLOT-LH)	☐ Right (BL-LTS-GT2SLOT-RH)	$ 58.00
Hardware cover	☐ Left (BL-LTS-HCOVER)	☐ Right (BL-LTS-HCOVER)	$ 52.00

BodiLink premium head support & hardware

Optional, if selecting must choose pad size, foam selection, cover, and mounting hardware.

Pad size and foam selection

Pad size	Soft	Extra soft	
6”W x 3”L	☐ BL-HSP1-6W3L	☐ BL-HSP2-6W3L	$ 260.00
8”W x 4”L	☐ BL-HSP1-8W4L	☐ BL-HSP2-8W4L	$ 260.00
10”W x 5”L	☐ BL-HSP1-10W5L	☐ BL-HSP2-10W5L	$ 260.00
12”W x 5”L	☐ BL-HSP1-12W5L	☐ BL-HSP2-12W5L	$ 260.00
14”W x 5”L	☐ BL-HSP1-14W5L	☐ BL-HSP2-14W5L	$ 260.00
18”W x 5”L	☐ BL-HSP1-18W5L	☐ BL-HSP2-18W5L	$ 260.00

Cover options

☐ COMFORT-TEK	**COMFORT-TEK** *For fluid protection & an easily cleaned surface*	NC
☐ STRETCH-AIR	**STRETCH-AIR** *For comfort & heat dissapation*	NC
☐ GLIDEWEAR	**GlideWear** *For skin protection & shear reduction. NOT fluid resistant*	$ 87.00

PT hardware, removable

☐ BL-HS-PT1RMMID6	**Mid 6**	$ 373.00
☐ BL-HS-PT1RMLAT6	**Lat 6**	$ 489.00
☐ BL-HS-PT1RMMID8	**Mid 8**	$ 373.00
☐ BL-HS-PT1RMLAT8	**Lat 8**	$ 489.00
☐ BL-HS-PT1RMMAX8	**Max 8**	$ 545.00

Stem length option | *For removable PT hardware only*

☐ HS-STEM-7	**Standard 7” stem**	NC
☐ HS-STEM-3	**3”stem (recommended for use with power mounting where there may be constraints)**	NC

PT hardware, flip down | *Flip down hardware has 7” stem.*

☐ BL-HS-PT1FDMID6	**Mid 6**	$ 419.00
☐ BL-HS-PT1FDLAT6	**Lat 6**	$ 530.00
☐ BL-HS-PT1FDMID8	**Mid 8**	$ 419.00
☐ BL-HS-PT1FDLAT8	**Lat 8**	$ 530.00
☐ BL-HS-PT1FDMAX8	**Max 8**	$ 591.00

Mounting option

☐ HS-MB-POSTA	**Posterior A (For Compass 4 and Compass Power Mount hardware)**	NC
☐ HS-MB-PW2C	**Permobil C (For Permobil Direct Mount only)**	$155.00

Additional hardware options | *Priced per selection.*

Extra link for PT hardware	☐ Extra 2” link (BL-HS-PT1LINK2) ☐ Extra 3” link (BL-HS-PT1LINK3) ☐ Extra 4” link (BL-HS-PT1LINK4)		$ 114.00
☐	BL-HS-52WINGA1	**Tool-less wing lock for PT hardware** *Enter desired number.*	$ 47.00
☐	BL-HS-SPCRA	**1/4 Spacer Plate**	$ 87.00

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- The client's name or code
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of seating and positioning products and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____


Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.


References

Section 1 | *Hardware options*


Compass 4 quick release




Compass 4 fixed




Compass power mount
Power mount hardware has 8” of height adjustment in ½” increments without using tools.




4-point




Permobil Direct Mount
For use on Permobil Corpus power wheelchairs only



UniTrack Assembly
Included with Direct Mount




Section 2 | *Accessories*




Privacy shield
Provides coverage between the bottom of the back support to the seat base of their wheelchair. Attaches directly to the canes of a wheelchair.

Section 3 | *BodiLink lateral trunk support hardware*


GT2 mounting hardware, slot mount, swing-away



GT2 mounting hardware, slot mount, fixed




GT2 mounting hardware, track mount, swing away
Track mount hardware is only available when using Compass 4 hardware. If using compass power mount hardware, select slot mounting hardware.




Section 4 | *BodiLink lateral trunk support pads*


Premium pads, supportive




Premium pads, active



Basic pads, supportive




Basic pads, active




Section 5 | *Bodilink lateral trunk support additional hardware options*

Extra hardware links




Long mounting bracket

Adds 1” of additional adjustment per bracket




Hardware cover



Section 6 | *Bodilink head support pads*

Basic pad sizes




Length(L) refers to the actual dimension of the support from bottom to top edge



Basic head support pads can be formed to any shape using the moldable steel insert inside the pad.

Premium pad sizes



Length(L) refers to the actual dimension of the support from bottom to top edge



Premium head support pad angles can be formed by adjusting the plates in the pad

Section 7 | *BodiLink head support hardware options*

GT hardware, removable

Range of motion


GT model	Total # of links	Total link length	Max. anterior reach	Max. posterior reach	Max. lateral reach
LAT	5	10.5”	10”	2”	7.5”
MID	3	8.5”	8”	2”	N/A



BodiLink head support hardware: PT hardware, removable

Range of motion


PT model	Total # of links	Total link length	Max anterior reach	Max posterior reach	Max lateral reach
MID 6	2	6”	5.5”	1”	N/A
LAT 6	3	6”	5.5”	1”	6”
MID 8	2	8”	7.5”	1”	N/A
LAT 8	3	8”	7.5”	1”	8”
MAX 8	4	8”	7.5”	3”	8”




PT hardware, flip down hardware, removable
Flip down hardware has 7" stem



BodiLink head support hardware: PT hardware, removable stem length option
3" stem is recommended for use with power mounting where there may be special constraints with power chair functions



Mounting options




Mounting block Posterior A
For Comfort Unite, Embrace and Visco backs as well as Comfort Relief with Compass power mount




Mounting block Permobil C
For use with Permobil Direct Mount and Unitrack assembly

Extra hardware

¼" spacer plate



Extra links: 2", 3", 4"



Tool-less wing lock

