Order form

Custom Comfort Embrace back support



Dealer Name*: _____ Account #*: _____

City, Province, Post Code: _____

Phone: _____ Email: _____

ATP: _____ Rep Name: _____

P.O. Number: _____ Quote Number: _____

Client Reference*: _____ Client Weight: _____

Required field

Special considerations that need to be addressed here (e.g., diagnosis):

_____

Notes | *Refer to reference pages at end of order form for additional information and images.*

Seating system measurements

Seating system measurements are for reference only and do not affect selections made throughout the order form.

A - Back support width

B- Back support length

C - Width between lateral trunk supports


D - Width between lateral thigh supports

E - Seat width

F - Seat depth

G - Effective seat depth

H - Lateral trunk support height



Comfort Embrace Back Support

Part number: C-E-BSXXWYYL.

Length	10”	12”	14”	16”	18”	20”	22”	24”	26”	28”
22”		☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00
21”		☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00
20”		☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00
19”		☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00
18”	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00
17”	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00
16”	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00
15”	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00
14”	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	
13”	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	
12”	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00			
11”	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00			
10”	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00	☐ $ 1,500.00			

Width

10” length and/or width not available with power mount hardware. 10” length not available with 4 point hardware if using head support. 24” wide back supports require 4 point or power mount hardware.


Cover material selection | *Required, choose one.*

☐ C-COMFORT-TEK	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	NC
☐ C-DIAMOND-SPACER	**Diamond Spacer** *For comfort and heat dissipation. NOT a fluid resistant material.*	NC

Custom foam selection


Select only one type of foam per layer. | *1/4” Crosslink foundation layer will come standard on all back supports. If Foam in Place Kit is selected, a middle and top layer needs to be selected. For other foam combinations, call Customer Experience at 800-736-0925.*

Rear layer
REQUIRED
H.R. MSRP $0.00
Comfort MSRP $160.00
Visco MSRP $225.00




Foam type	1/4”	1/2”	3/4”	1”	1 1/4”
H.R. Medium	☐ CCD-REAR-HRM-0.25	☐ CCD-REAR-HRM-0.50	☐ CCD-REAR-HRM-0.75	☐ CCD-REAR-HRM-1.00	☐ CCD-REAR-HRM-1.25
H.R. Soft	☐ CCD-REAR-HRS-0.25	☐ CCD-REAR-HRS-0.50	☐ CCD-REAR-HRS-0.75	☐ CCD-REAR-HRS-1.00	☐ CCD-REAR-HRS-1.25
Comfort Medium	☐ CCD-REAR-COM-0.25	☐ CCD-REAR-COM-0.50	☐ CCD-REAR-COM-0.75	☐ CCD-REAR-COM-1.00	☐ CCD-REAR-COM-1.25
Comfort Soft	☐ CCD-REAR-COS-0.25	☐ CCD-REAR-COS-0.50	☐ CCD-REAR-COS-0.75	☐ CCD-REAR-COS-1.00	☐ CCD-REAR-COS-1.25
Visco Medium		☐ CCD-REAR-VIS-M-0.50		☐ CCD-REAR-VIS-M-1.00	
Visco Soft		☐ CCD-REAR-VIS-S-0.50		☐ CCD-REAR-VIS-S-1.00	
Visco X-Soft		☐ CCD-REAR-VIS-XS-0.50		☐ CCD-REAR-VIS-XS-1.00	
Visco XX-Soft		☐ CCD-REAR-VIS-XXS-0.50		☐ CCD-REAR-VIS-XXS-1.00	

Middle layer
OPTIONAL
H.R. MSRP $0.00
Comfort MSRP $160.00
Visco MSRP $225.00



Foam type	1/4”	1/2”	3/4”	1”	1 1/4”
H.R. Medium	☐ CCD-MID-HRM-0.25	☐ CCD-MID-HRM-0.50	☐ CCD-MID-HRM-0.75	☐ CCD-MID-HRM-1.00	☐ CCD-MID-HRM-1.25
H.R. Soft	☐ CCD-MID-HRS-0.25	☐ CCD-MID-HRS-0.50	☐ CCD-MID-HRS-0.75	☐ CCD-MID-HRS-1.0	☐ CCD-MID-HRS-1.25
Comfort Medium	☐ CCD-MID-COM-0.25	☐ CCD-MID-COM-0.50	☐ CCD-MID-COM-0.75	☐ CCD-MID-COM-1.00	☐ CCD-MID-COM-1.25
Comfort Soft	☐ CCD-MID-COS-0.25	☐ CCD-MID-COS-0.50	☐ CCD-MID-COS-0.75	☐ CCD-MID-COS-1.00	☐ CCD-MID-COS-1.25
Visco Medium		☐ CCD-MID-VIS-M-0.50		☐ CCD-MID-VIS-M-1.00	
Visco Soft		☐ CCD-MID-VIS-S-0.50		☐ CCD-MID-VIS-S-1.00	
Visco X-Soft		☐ CCD-MID-VIS-XS-0.50		☐ CCD-MID-VIS-XS-1.00	
Visco XX-Soft		☐ CCD-MID-VIS-XXS-0.50		☐ CCD-MID-VIS-XXS-1.00	

Topper layer
OPTIONAL
H.R. MSRP $0.00
Comfort MSRP $160.00
Visco MSRP $225.00



Foam type	1/4”	1/2”	3/4”	1”	1 1/4”
H.R. Medium	☐ CCD-TOP-HRM-0.25	☐ CCD-TOP-HRM-0.50	☐ CCD-TOP-HRM-0.75	☐ CCD-TOP-HRM-1.00	☐ CCD-TOP-HRM-1.25
H.R. Soft	☐ CCD-TOP-HRS-0.25	☐ CCD-TOP-HRS-0.50	☐ CCD-TOP-HRS-0.75	☐ CCD-TOP-HRS-1.00	☐ CCD-TOP-HRS-1.25
Comfort Medium	☐ CCD-TOP-COM-0.25	☐ CCD-TOP-COM-0.50	☐ CCD-TOP-COM-0.75	☐ CCD-TOP-COM-1.00	☐ CCD-TOP-COM-1.25
Comfort Soft	☐ CCD-TOP-COS-0.25	☐ CCD-TOP-COS-0.50	☐ CCD-TOP-COS-0.75	☐ CCD-TOP-COS-1.00	☐ CCD-TOP-COS-1.25
Visco Medium		☐ CCD-TOP-VIS-M-0.50		☐ CCD-TOP-VIS-M-1.00	
Visco Soft		☐ CCD-TOP-VIS-S-0.50		☐ CCD-TOP-VIS-S-1.00	
Visco X-Soft		☐ CCD-TOP-VIS-XS-0.50		☐ CCD-TOP-VIS-XS-1.00	
Visco XX-Soft		☐ CCD-TOP-VIS-XXS-0.50		☐ CCD-TOP-VIS-XXS-1.00	

Foam in place kit

Kit includes |


1) Back cushion cover with an optional 3" or 4" deep back, with a stretchy, pleated contact surface to allow for foam growth.
2) Two aluminum stays sewn into the rear of the cover which can be molded to support spinal curvature prior to pouring FIP chemicals.
3) Middle foam layer is cut to fit the shape of the back support and is placed loosely inside the cover. This will give the user some padding prior to pouring the foam.
4) Top foam layer is cut to fit the shape of the back support and is placed inside the front of the cover. This will give the final mold a plush feel for the user and will also smooth out any potential inconsistencies after pouring the foam.
5) Two plastic bags are included to contain the foam chemicals as they expand.

***Note:** Price of Visco layers will be added to the order separately.*

Bag pour size recommendations:

3" border thickness for frame widths 10" - 16", 1 bag
3" border thickness for frame widths 17"-22", 1 1/2 bag
4" border thickness for frame widths 10" - 14", 1 bag
4" border thickness for frame widths 15" - 20", 1 1/2 bag

MSRP:
1 Bag FIP: $650.00
1 1/2 Bag FIP: $800.00



Pour size	Extra soft memory		Soft memory foam		Medium soft foam	
	1 bag	1 1/2 bag	1 bag	1 1/2 bag	1 bag	1 1/2 bag
3" border thickness	☐ CCD-FIP3-1-XS	☐ CCD-FIP3-1.5-XS	☐ CCD-FIP3-1-S	☐ CCD-FIP3-1.5-S	☐ CCD-FIP3-1-M	☐ CCD-FIP3-1.5-M
4" border thickness	☐ CCD-FIP4-1-XS	☐ CCD-FIP4-1.5-XS	☐ CCD-FIP4-1-S	☐ CCD-FIP4-1.5-S	☐ CCD-FIP4-1-M	☐ CCD-FIP4-1.5-M


Back support options

Moldable aluminum stays | *This back support comes standard with two patented aluminum stays that can be molded to support spinal curvature. A standard alignment of the moldable stays will be chosen based on back support width if you choose not to select a specific width between the stays.*

Note: The stays are accessible while the user is seated against the back support through the zipper located at the top of the back support.

MSRP: $0.00

- ☐ A = 1" (CCD-STAY-1.0)
- ☐ A = 1 1/2" (CCD-STAY-1.5)
- ☐ A = 2" (standard alignment) (CCD-STAY-2.0)
- ☐ A = 2 1/2" (CCD-STAY-2.5)
- ☐ A = 3" (CCD-STAY-3.0)
- ☐ A = 3 1/2" (CCD-STAY-3.5)
- ☐ A = 4" (CCD-STAY-4.0)
- ☐ A = 4 1/2" (CCD-STAY-4.5)
- ☐ Specify custom spacing (CCD-STAY-OTHER)
- ☐ None (CCD-STAY-NONE)




Back support wing style |

MSRP $0.00 / MSRP $110.00 Extra wing

Wing size recommendation chart

Back Length → Wing Length
10"L - 12"L Back → 5"L Wing
13"L - 16"L Back → 6"L Wing
17"L - 22"L Back → 8"L Wing



Note:
- Contoured wings create a height and width adjustable back support
- Multiple wings may be ordered if more than one lateral is to be used per side
- Wings are need to mount lateral trunk support hardware

Wing Length	Left	Extra Left	Right	Extra Right
5" L Contoured Wing	☐ CD-LEFTWING-5-CON	☐ CD-ADD-L-WING-5-CON	☐ CD-RIGHTWING-5-CON	☐ CD-ADD-R-WING-5-CON
6" L Contoured Wing	☐ CD-LEFTWING-6-CON	☐ CD-ADD-L-WING-6-CON	☐ CD-RIGHTWING-6-CON	☐ CD-ADD-R-WING-6-CON
8" L Contoured Wing	☐ CD-LEFTWING-8-CON	☐ CD-ADD-L-WING-8-CON	☐ CD-RIGHTWING-8-CON	☐ CD-ADD-R-WING-8-CON

Mounting Hardware | *Required, choose one.*

Compass 4 Quick Release
☐ Quick Release 2 Point (CCD-Q-COMPASS) ☐ Quick Release 4 Point (CCD-Q-COMPASS4) — NC

Compass 4 Fixed
☐ Fixed 2 Point (CCD-NQ-COMPASS) ☐ Fixed 4 Point (CCD-NQ-COMPASS4) — NC

Compass Power Mount — $ 125.00
- ☐ Permobil w/ Corpus 3G (CCD-PMB-PERMOBIL)
- ☐ Quantum w/ Tru Balance 3 or 4 (CCD-PMB-QUANTUM)
- ☐ Quickie w/ 3.7 Power (CCD-PMB-QUICKIE)
- ☐ Avid Rehab w/ Power T & R (CCD-PMB-AVID)
- ☐ ROVI w/ Motion Concepts (CCD-PMB-ROVI)

Using Power Mount hardware on ROVI requires mounting to the recessed planar interface plate available when ordering the Motion Concepts seating system. Drilling is required for mounting.

Color Piping | *Optional. If no color is selected, the standard silver reflective piping is used.*

Choose one — $ 125.00
☐ Remove (CCD-REMOVE) ☐ Red (CCD-RED) ☐ Purple (CCD-PURPLE)
☐ Pink (CCD-PINK) ☐ Blue (CCD-BLUE) ☐ Green (CCD-GREEN)
☐ Orange (CCD-ORANGE)

Embroidery | *Optional. Includes corresponding color zipper pull.*

Thread color — NC
☐ Red (CCD-EMB-RED) ☐ Purple (CCD-EMB-PURPLE) ☐ Pink (CCD-EMB-PINK)
☐ Blue (CCD-EMB-BLUE) ☐ Orange (CCD-EMB-ORANGE) ☐ Green (CCD-EMB-GREEN)
☐ Silver (CCD-EMB-SILVER)

Font — NC
☐ Helvetica (CCD-EMB-HELVETICA) ☐ Garamond (CCD-EMB-GARAMOND)
☐ Impress (CCD-EMB-IMPRESS) ☐ Kaufmann (CCD-EMB-KAUFMANN)

Font styles
HELVETICA
GARAMOND
Impress
Kaufmann


Word one — NC
☐ Input is case sensitive:

Word two (optional) — NC
☐ Input is case sensitive:

Growth Kit | *Optional. Includes additional foam, cover and hardware.*

Size — $ 550.00
☐ 2" Growth plate (CCD-2"GK) ☐ 4" Growth plate (CCD-4"GK)

KIT INCLUDES:
- A 2" or 4" growth plate that bolts onto the bottom of existing shell to grow the back support in length.
- An additional 2" or 4" larger back support cushion and cover will be supplied based on selected growth plate.
- Mounting hardware will be included in kit.



Strap guide extension bracket

☐ CHG2-EB1	**Extra small** *Up to 1 3/4" with webbing*	$ 135.00
☐ CHG2-EB2	**Small** *Up to 2 1/4" with webbing*	$ 135.00

Other Accessories | *Optional*

☐ C-PV-SHIELD	**Privacy Shield** *Connects with hook & loop attachment*	$ 85.00

BodiLink lateral trunk supports & hardware

Optional, if selecting must choose pad shape & size, pad style, pad covers, and mounting hardware.

Basic pad shape & size | *Priced per pad.*

	Left *Must enter the number of pads needed*	**Right** *Must enter the number of pads needed*	
Supportive			
3"L x 4"D	☐ BL-LTSP1S-3L4D-LH	☐ BL-LTSP1S-3L4D-RH	$ 161.00
3"L x 5"D	☐ BL-LTSP1S-3L5D-LH	☐ BL-LTSP1S-3L5D-RH	$ 161.00
4"L x 5"D	☐ BL-LTSP1S-4L5D-LH	☐ BL-LTSP1S-4L5D-RH	$ 161.00
4"L x 6"D	☐ BL-LTSP1S-4L6D-LH	☐ BL-LTSP1S-4L6D-RH	$ 161.00
5"L x 7"D	☐ BL-LTSP1S-5L7D-LH	☐ BL-LTSP1S-5L7D-RH	$ 161.00
6"L x 6"D	☐ BL-LTSP1S-6L6D-LH	☐ BL-LTSP1S-6L6D-RH	$ 161.00
Active			
3"L x 4"D	☐ BL-LTSP1A-3L4D-LH	☐ BL-LTSP1A-3L4D-RH	$ 161.00
3"L x 5"D	☐ BL-LTSP1A-3L5D-LH	☐ BL-LTSP1A-3L5D-RH	$ 161.00
4"L x 5"D	☐ BL-LTSP1A-4L5D-LH	☐ BL-LTSP1A-4L5D-RH	$ 161.00
4"L x 6"D	☐ BL-LTSP1A-4L6D-LH	☐ BL-LTSP1A-4L6D-RH	$ 161.00
5"L x 7"D	☐ BL-LTSP1A-5L7D-LH	☐ BL-LTSP1A-5L7D-RH	$ 161.00
6"L x 6"D	☐ BL-LTSP1A-6L6D-LH	☐ BL-LTSP1a-6L6D-RH	$ 161.00

Basic pad style | *Must match number of pads selected.*

☐ LTS-P1-MOLD	**Moldable**	NC
	☐ Left ☐ Right	NC
☐ LTS-P1-RIGID	**Rigid**	NC
	☐ Left ☐ Right	NC

Basic pad cover options | *Must match number of pads selected.*

☐ LTS-P1-COMFORT	**COMFORT-TEK®** *For fluid protection & an easily cleaned surface*	NC
	☐ Left ☐ Right	
☐ LTS-P1-STRETCH	**STRETCH-AIR®** *For comfort & heat dissapation*	NC
	☐ Left ☐ Right	

Premium pad shape & size | *Priced per pad.*

	Left *Enter number of pads required*	Right *Enter number of pads required*	
Supportive			
3”L x 4” D	☐ BL-LTSP2S-3L4D-LH	☐ BL-LTSP2S-3L4D-RH	$ 218.00
3” L x 5.5” D	☐ BL-LTSP2S-3L5D-LH	☐ BL-LTSP2S-3L5D-RH	$ 218.00
4” L x 5” D	☐ BL-LTSP2S-4L5D-LH	☐ BL-LTSP2S-4L5D-RH	$ 218.00
4” L x 6.5” D	☐ BL-LTSP2S-4L6D-LH	☐ BL-LTSP2S-4L6D-RH	$ 218.00
5” L x 6.5” D	☐ BL-LTSP2S-5L6D-LH	☐ BL-LTSP2S-5L6D-RH	$ 218.00
6” L x 6.5” D	☐ BL-LTSP2S-6L6D-LH	☐ BL-LTSP2S-6L6D-RH	$ 218.00
Active			
3”L x 4” D	☐ BL-LTSP2A-3L4D-LH	☐ BL-LTSP2A-3L4D-RH	$ 218.00
3” L x 5.5” D	☐ BL-LTSP2A-3L5D-LH	☐ BL-LTSP2A-3L5D-RH	$ 218.00
4” L x 5” D	☐ BL-LTSP2A-4L5D-LH	☐ BL-LTSP2A-4L5D-RH	$ 218.00
4” L x 6.5” D	☐ BL-LTSP2A-4L6D-LH	☐ BL-LTSP2A-4L6D-RH	$ 218.00
5” L x 6.5” D	☐ BL-LTSP2A-5L6D-LH	☐ BL-LTSP2A-5L6D-RH	$ 218.00
6” L x 6.5” D	☐ BL-LTSP2A-6L6D-LH	☐ BL-LTSP2A-6L6D-RH	$ 218.00

Premium pad style | *Must match number of pads selected.*

☐ LTS-P2-MOLD	**Moldable**	NC
	☐ Left ☐ Right	NC
☐ LTS-P2-RIGID	**Rigid**	NC
	☐ Left ☐ Right	NC

Premium pad cover options | *Must match number of pads selected. Priced per selection.*

☐ LTS-P2-COMFORT	**COMFORT-TEK** *For fluid protection & an easily cleaned surface* ☐ Left ☐ Right	NC
☐ LTS-P2-STRETCH	**STRETCH-AIR** *For comfort & heat dissapation* ☐ Left ☐ Right	NC
☐ LTS-P2-GLIDE	**GlideWear®** *For skin protection & shear reduction. NOT fluid resistant* ☐ Left ☐ Right	$ 75.00
☐ LTS-P2-COMFORT-FM	**COMFORT-TEK with foam** ☐ Left ☐ Right	$ 30.00
☐ LTS-P2-STRETCH-FM	**STRETCH-AIR with foam** ☐ Left ☐ Right	$ 30.00
☐ LTS-P2-GLIDE-FM	**GlideWear with foam** ☐ Left ☐ Right	$ 103.00

GT2 mounting hardware | *Used with basic and premium pads. Both fixed and swing away hardware styles are adjustable.*

	Left *Must match number of pads selected*	Right *Must match number of pads selected*	
Swing-away			
Slot Mount	☐ BL-LTS-GT2SAS-LH	☐ BL-LTS-GT2SAS-RH	$ 368.00
Slot Mount w/ Winter/Summer	☐ BL-LTS-GT2SAWSE-LH	☐ BL-LTS-GT2SAWSE-RH	$ 626.00
Fixed			
Slot Mount	☐ BL-LTS-GT2FXS-LH	☐ BL-LTS-GT2FXS-RH	$ 202.00
Slot Mount w/ Winter/Summer	☐ BL-LTS-GT2FXWSE-LH	☐ BL-LTS-GT2FXWSE-RH	$ 426.00

Additional hardware options | *Priced per selection.*

Extra GT2 hardware link	☐ Left (BL-LTS-GT2LINK)	☐ Right (BL-LTS-GT2LINK)	$ 58.00
Long mounting bracket	☐ Left (BL-LTS-GT2SLOT-LH)	☐ Right (BL-LTS-GT2SLOT-RH)	$ 58.00
Hardware cover	☐ Left (BL-LTS-HCOVER)	☐ Right (BL-LTS-HCOVER)	$ 52.00
Winter/summer hardware	☐ Left (BL-LTS-HWSAE-LH)	☐ Right (BL-LTS-HWSAE-RH)	$ 225.00

BodiLink premium head support & hardware

Optional, if selecting must choose pad size, foam selection, cover, and mounting hardware.

Pad size and foam selection

Pad size	Soft	Extra soft	
6”W x 3”L	☐ BL-HSP1-6W3L	☐ BL-HSP2-6W3L	$ 260.00
8”W x 4”L	☐ BL-HSP1-8W4L	☐ BL-HSP2-8W4L	$ 260.00
10”W x 5”L	☐ BL-HSP1-10W5L	☐ BL-HSP2-10W5L	$ 260.00
12”W x 5”L	☐ BL-HSP1-12W5L	☐ BL-HSP2-12W5L	$ 260.00
14”W x 5”L	☐ BL-HSP1-14W5L	☐ BL-HSP2-14W5L	$ 260.00
18”W x 5”L	☐ BL-HSP1-18W5L	☐ BL-HSP2-18W5L	$ 260.00

Cover options

☐ COMFORT-TEK	**COMFORT-TEK** *For fluid protection & an easily cleaned surface*	NC
☐ STRETCH-AIR	**STRETCH-AIR** *For comfort & heat dissapation*	NC
☐ GLIDEWEAR	**GlideWear** *For skin protection & shear reduction. NOT fluid resistant*	$ 87.00

PT hardware, removable

☐ BL-HS-PT1RMMID6	**Mid 6**	$ 373.00
☐ BL-HS-PT1RMLAT6	**Lat 6**	$ 489.00
☐ BL-HS-PT1RMMID8	**Mid 8**	$ 373.00
☐ BL-HS-PT1RMLAT8	**Lat 8**	$ 489.00
☐ BL-HS-PT1RMMAX8	**Max 8**	$ 545.00

Stem length option | *For removable PT hardware only.*

☐ HS-STEM-7	**Standard 7” stem**	NC
☐ HS-STEM-3	**3”stem (recommended for use with power mounting where there may be constraints)**	NC

PT hardware, flip down | *Flip down hardware has 7” stem.*

☐	BL-HS-PT1FDMID6	**Mid 6**	$ 419.00
☐	BL-HS-PT1FDLAT6	**Lat 6**	$ 530.00
☐	BL-HS-PT1FDMID8	**Mid 8**	$ 419.00
☐	BL-HS-PT1FDLAT8	**Lat 8**	$ 530.00
☐	BL-HS-PT1FDMAX8	**Max 8**	$ 591.00

Mounting option

☐	HS-MB-POSTA	**Posterior A**	NC

Additional hardware options | *Priced per selection.*

Extra link for PT hardware	☐	Extra 2” link (BL-HS-PT1LINK2) ☐ Extra 3” link (BL-HS-PT1LINK3) ☐ Extra 4” link (BL-HS-PT1LINK4)	$ 114.00
☐	BL-HS-52WINGA1	**Tool-less wing lock for PT hardware** *Enter desired number.*	$ 47.00
☐	BL-HS-SPCRA	**1/4 Spacer Plate**	$ 87.00

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- The client's name or code
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, , am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of seating and positioning products and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____


Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.


References

Section 1 | *Hardware options*

Compass 4 quick release




Compass 4 fixed




Compass power mount


Power mount hardware has 8” of height adjustment in ½” increments without using tools.



4-point




Section 2 | *Accessories*



Privacy shield


Provides coverage between the bottom of the back support to the seat base of their wheelchair. Attaches directly to the canes of a wheelchair.




Strap guide extension bracket

Section 3 | *BodiLink lateral trunk support hardware*

GT2 mounting hardware, slot mount, swing-away




GT2 mounting hardware, slot mount, fixed




Section 4 | *BodiLink lateral trunk support pads*


Premium pads, supportive




Premium pads, active



Basic pads, supportive




Basic pads, active




Section 5 | *Bodilink lateral trunk support additional hardware options*

Extra hardware links




Long mounting bracket

Adds 1” of additional adjustment per bracket



Hardware cover




Winter/Summer hardware



Section 6 | *Bodilink head support pads*

Basic pad sizes




Length(L) refers to the actual dimension of the support from bottom to top edge



Basic head support pads can be formed to any shape using the moldable steel insert inside the pad.

Premium pad sizes



Length(L) refers to the actual dimension of the support from bottom to top edge



Premium head support pad angles can be formed by adjusting the plates in the pad

Section 7 | *BodiLink head support hardware options*

GT hardware, removable

Range of motion


GT model	Total # of links	Total link length	Max. anterior reach	Max. posterior reach	Max. lateral reach
LAT	5	10.5”	10”	2”	7.5”
MID	3	8.5”	8”	2”	N/A



BodiLink head support hardware: PT hardware, removable

Range of motion


PT model	Total # of links	Total link length	Max anterior reach	Max posterior reach	Max lateral reach
MID 6	2	6”	5.5”	1”	N/A
LAT 6	3	6”	5.5”	1”	6”
MID 8	2	8”	7.5”	1”	N/A
LAT 8	3	8”	7.5”	1”	8”
MAX 8	4	8”	7.5”	3”	8”




PT hardware, flip down hardware, removable
Flip down hardware has 7" stem



BodiLink head support hardware: PT hardware, removable stem length option
3" stem is recommended for use with power mounting where there may be special constraints with power chair functions




Mounting options




Mounting block Posterior A
For Comfort Unite, Embrace and Visco backs as well as Comfort Relief with Compass power mount

Extra hardware

¼" spacer plate



Extra links: 2", 3", 4"



Tool-less wing lock

