Order form

Bodilink® head support



Dealer Name*: _____ Account #*: _____

City, State, Zip: _____

Phone: _____ Email: _____

ATP: _____ Rep Name: _____

P.O. Number: _____ Quote Number: _____

Client Reference*: _____ Client Weight: _____

***Required field**

Special considerations that need to be addressed here (e.g., diagnosis):

_____

NOTE | *Refer to reference page at end of order form for additional information and images.*

Basic head support pad options

Must select pad size and cover option.

Pad size

☐	BL-HSP3-6W5L	**6”W x 5”L**	$ 260.00
☐	BL-HSP3-8W5L	**8”W x 5”L**	$ 260.00
☐	BL-HSP3-10W6L	**10”W x 6”L**	$ 260.00
☐	BL-HSP3-12W6L	**12”W x 6”L**	$ 260.00
☐	BL-HSP3-14W6L	**14”W x 6”L**	$ 260.00
☐	BL-HSP3-18W6L	**18”W x 6”L**	$ 260.00

Cover options | *Required, choose one.*

☐	COMFORT-TEK®	**COMFORT-TEK®** *For fluid protection & an easily cleaned surface*	NC
☐	STRETCH-AIR®	**STRETCH-AIR®** *For comfort & heat dissapation*	NC

Premium head support pad options

Must select pad size, foam and cover option.

Pad size and foam selection

Pad size	Soft	Extra soft	
6”W x 3”L	☐ BL-HSP1-6W3L	☐ BL-HSP2-6W3L	$ 260.00
8”W x 4”L	☐ BL-HSP1-8W4L	☐ BL-HSP2-8W4L	$ 260.00
10”W x 5”L	☐ BL-HSP1-10W5L	☐ BL-HSP2-10W5L	$ 260.00
12”W x 5”L	☐ BL-HSP1-12W5L	☐ BL-HSP2-12W5L	$ 260.00
14”W x 5”L	☐ BL-HSP1-14W5L	☐ BL-HSP2-14W5L	$ 260.00
18”W x 5”L	☐ BL-HSP1-18W5L	☐ BL-HSP2-18W5L	$ 260.00

Cover options | *Required, choose one.*

☐ COMFORT-TEK	**COMFORT-TEK** *For fluid protection & an easily cleaned surface*	NC
☐ STRETCH-AIR	**STRETCH-AIR** *For comfort & heat dissapation*	NC
☐ GLIDEWEAR	**GlideWear®** *For skin protection & shear reduction. NOT fluid resistant*	$ 82.00

Hardware options

Required, choose one.

GT hardware, removable

☐ BL-HS-GT3RMLAT	**Lateral**	$ 350.00
☐ BL-HS-GT3RMMID	**Mid-line**	$ 311.00

PT hardware, removable

☐ BL-HS-PT1RMMID6	**Mid 6**	$ 373.00
☐ BL-HS-PT1RMLAT6	**Lat 6**	$ 489.00
☐ BL-HS-PT1RMMID8	**Mid 8**	$ 373.00
☐ BL-HS-PT1RMLAT8	**Lat 8**	$ 489.00
☐ BL-HS-PT1RMMAX8	**Max 8**	$ 545.00

Stem length option | *Must select when choosing removable PT hardware.*

☐ HS-STEM-7	**Standard 7”**	NC
☐ HS-STEM-3	**3” stem**	NC

PT hardware, flip down | *Includes 7” stem.*

☐ BL-HS-PT1FDMID6	**Mid 6**	$ 419.00
☐ BL-HS-PT1FDLAT6	**Lat 6**	$ 530.00
☐ BL-HS-PT1FDMID8	**Mid 8**	$ 419.00
☐ BL-HS-PT1FDLAT8	**Lat 8**	$ 530.00
☐ BL-HS-PT1FDMAX8	**Max 8**	$ 591.00

Mounting hardware options

Required, choose one.

Mounting option

☐ HS-MB-POSTA	**Posterior A**	NC
☐ HS-MB-POSTB	**Posterior B**	NC
☐ HS-MB-ANTA	**Anterior A**	NC
☐ HS-MB-PW2A	**Permobil A**	$ 155.00
☐ HS-MB-PW2B	**Permobil B**	$ 155.00
☐ HS-MB-PW2C	**Permobil C**	$ 155.00
☐ HS-MB-PW1A	**Quantum A**	$ 155.00
☐ HS-MB-AR14W	**Comfort Relief 14W**	$ 161.00
☐ HS-MB-AR16W	**Comfort Relief 16W**	$ 161.00
☐ HS-MB-AR18W	**Comfort Relief 18W**	$ 161.00
☐ HS-MB-AR20W	**Comfort Relief 20W**	$ 161.00

Additional hardware options

Optional, priced per selection.

Extra link for PT hardware	☐ Extra 2” link (BL-HS-PT1LINK2)	☐ Extra 3” link (BL-HS-PT1LINK3)	$ 114.00
	☐ Extra 4” link (BL-HS-PT1LINK4)		
☐ BL-HS-52WINGA1	**Tool-less wing lock for PT hardware** *Entered desired number*		$ 47.00
☐ BL-HS-SPCRA	**1/4 spacer plate**		$ 87.00

Extra mounting hardware options

Optional

☐ BL-HS-MBPOSTA	**Posterior A**	$ 119.00
☐ BL-HS-MBPOSTB	**Posterior B**	$ 119.00
☐ BL-HS-MBANTA	**Anterior A**	$ 119.00
☐ BL-HS-MBPW2A	**Permobil A**	$ 119.00
☐ BL-HS-MBPW2B	**Permobil B**	$ 119.00
☐ BL-HS-MBPW2C	**Permobil C**	$ 119.00
☐ BL-HS-MBPW1A	**Quantum A**	$ 119.00
☐ BL-HS-AR14W	**Comfort Relief 14W**	$ 161.00
☐ BL-HS-AR16W	**Comfort Relief 16W**	$ 161.00
☐ BL-HS-AR18W	**Comfort Relief 18W**	$ 161.00
☐ BL-HS-AR20W	**Comfort Relief 20W**	$ 161.00

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- The client's name or code.
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, , am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of seating and positioning products and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____


Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.

References

Section 1 | *Basic pad option*

Pad size




Length(L) refers to the actual dimension of the support from bottom to top edge



Basic head support pads can be formed to any shape using the moldable steel insert inside the pad.

Section 2 | *Premium pad options*

Pad size



Length(L) refers to the actual dimension of the support from bottom to top edge




Premium head support pad angles can be formed by adjusting the plates in the pad

Section 3 | *Hardware options*

GT hardware, removable

Range of motion


GT model	Total # of links	Total link length	Max. anterior reach	Max. posterior reach	Max. lateral reach
LAT	5	10.5”	10”	2”	7.5”
MID	3	8.5”	8”	2”	N/A



PT hardware, removable

Range of motion


PT model	Total # of links	Total link length	Max anterior reach	Max posterior reach	Max lateral reach
MID 6	2	6”	5.5”	1”	N/A
LAT 6	3	6”	5.5”	1”	6”
MID 8	2	8”	7.5”	1”	N/A
LAT 8	3	8”	7.5”	1”	8”
MAX 8	4	8”	7.5”	3”	8”



PT hardware, flip down hardware, removable
Flip down hardware has 7" stem



PT hardware, removable stem length option
3" stem is recommended for use with power mounting where there may be special constraints with power chair functions




Note:
Mid and Lat 8 hardware is recommended for Comfort Relief back supports


Max 8 hardware provides the greatest amount or articulation and range of motion

Section 4 | *Head support mounting options*


Mounting options




Mounting block Posterior A or B
A: For Comfort Unite, Embrace and Visco backs as well as Comfort Relief with Compass power mount.
B: For OEM wood/t-nutback supports




Mounting block Anterior A
For OEM back supports



Mounting block Permobil A
For Permobil Corpus 3G Back Supports (23"-28"H)




Mounting block Quantum A
For Quantum Tru-Comfort back supports




Mounting block Permobil B
For Permobil Corpus 3G Back Supports (20"H)



Mounting block Permobil C
For use with Permobil Direct Mount and Unitrack rail




Comfort Relief Bracket with mounting block
18W Shown Above For Non-Power Mount Only.



¼" spacer plate

Section 5 | *Hardware*

Extra links: 2", 3", 4"



Tool-less wing lock

