Order form

Comfort Saddle 7 Series cushion



Dealer Name*: Account #*:

City, Province, Postal Code*:

Phone: Email:

ATP: Rep Name:

P.O. Number: Quote Number:

Client Reference*: Client Weight*:

***Required field**

Special considerations that need to be addressed here (e.g., diagnosis):

Saddle 7 Series Cushion without Quadragel

Part number: 53S-WWDD-B

Includes Solace cover and Kwik Strap

	22”	24”	26”	28”	30”
20”	☐ $818.00	☐ $818.00	☐ $818.00	☐ $818.00	☐ $818.00
18”	☐ $818.00	☐ $818.00	☐ $818.00	☐ $818.00	☐ $818.00
	Width				

Additional cover options | *Optional*

☐	53S-WWDD-B-CV	**One extra Solace outer cover**	$132.00
☐	53S-WWDD-B-CV	**Two extra Solace outer covers**	$264.00

Saddle 7 Series Cushion with Quadragel

Part number: 53SG-WWDD-B

Includes Solace cover and Kwik Strap

	22”	24”	26”	28”	30”
20”	☐ $1,133.00	☐ $1,133.00	☐ $1,133.00	☐ $1,133.00	☐ $1,133.00
18”	☐ $1,133.00	☐ $1,133.00	☐ $1,133.00	☐ $1,133.00	☐ $1,133.00
	Width				

Additional cover options | *Optional*

☐	53S-WWDD-B-CV	**One extra Solace outer cover**	$132.00
☐	53S-WWDD-B-CV	**Two extra Solace outer covers**	$264.00

Cushion accessories

Refer to reference pages for details and images.

☐	RSSP-1-WWDD	**Removable solid seat pan & hardware**	$690.00
Rigid insert	☐ Rigid insert glued (RIDWWDD)	☐ Rigid insert not glued (RIDWWDD)	$96.00

BodiLink® pelvic/thigh support hardware and pads

Refer to reference pages for details and images.

Basic pad shape & size | *Priced per pad.*

	Left *Must enter number of pads*		**Right** *Must enter number of pads*		
3”L x 4”D	☐	BL-LPTSP1Z-3L4D-LH	☐	BL-LPTSP1Z-3L4D-RH	$139.00
3”L x 5”D	☐	BL-LPTSP1Z-3L5D-LH	☐	BL-LPTSP1Z-3L5D-RH	$139.00
4”L X 4”D	☐	BL-LPTSP1Z-4L4D-LH	☐	BL-LPTSP1Z-4L4D-RH	$139.00
4”L x 6”D	☐	BL-LPTSP1Z-4L6D-LH	☐	BL-LPTSP1Z-4L6D-RH	$139.00
4”L x 8”D	☐	BL-LPTSP1Z-4L8D-LH	☐	BL-LPTSP1Z-4L8D-RH	$139.00
4”L x 10”D	☐	BL-LPTSP1Z-4L10D-LH	☐	BL-LPTSP1Z-4L10D-RH	$139.00
4”L x 12”D	☐	BL-LPTSP1Z-4L12D-LH	☐	BL-LPTSP1Z-4L12D-RH	$139.00
4”L x 14”D	☐	BL-LPTSP1Z-4L14D-LH	☐	BL-LPTSP1Z-4L14D-RH	$139.00
4”L x 16”D	☐	BL-LPTSP1Z-4L16D-LH	☐	BL-LPTSP1Z-4L16D-RH	$139.00
5”L x 7”D	☐	BL-LPTSP1Z-5L7D-LH	☐	BL-LPTSP1Z-5L7D-RH	$139.00

Basic pad cover options | *Required, must enter number of covers needed and match number of pads selected. Includes foam insert.*

☐	LPTS-P1-COMFORT	**COMFORT-TEK®** *For fluid protection & an easily cleaned surface* ☐ Left ☐ Right	NC
☐	LPTS-P1-STRETCH	**STRETCH-AIR®** *For comfort & heat dissipation* ☐ Left ☐ Right	NC
☐	LPTS-P1-GLIDE	**GlideWear®** *For skin protection & shear reduction. NOT fluid resistant.* ☐ Left ☐ Right	$36.00

Premium pad shape & size | *Priced per pad.*

	Left *Must enter the number of pads needed*		**Right** *Must enter the number of pads needed*		
3.5”L x 4”D	☐	BL-LPTSP2Z-3L4D-LH	☐	BL-LPTSP2Z-3L4D-RH	$150.00
3.5”L x 8”D	☐	BL-LPTSP2Z-3L8D-LH	☐	BL-LPTSP2Z-3L8D-RH	$150.00
3.5”L x 12”D	☐	BL-LPTSP2Z-3L12D-LH	☐	BL-LPTSP2Z-3L12D-RH	$150.00
5.5” x 6”D	☐	BL-LPTSP2Z-5L6D-LH	☐	BL-LPTSP2Z-5L6D-RH	$150.00

Premium pad cover options | *Required, must enter number of covers needed and match number of pads selected. Includes foam insert.*

☐	LPTS-P2-COMFORT-FM	**COMFORT-TEK** *For fluid protection & an easily cleaned surface* ☐ Left ☐ Right	
☐	LPTS-P2-STRETCH-FM	**STRETCH-AIR** *For comfort & heat dissapation* ☐ Left ☐ Right	
☐	LPTS-P2-GLIDE-FM	**GlideWear** *For skin protection & shear reduction. NOT fluid resistant.* ☐ Left ☐ Right	$36.00

GT mounting hardware, slot mount only

Fixed	☐	Left (BL-LPTS-GT2FXS-LH)	☐	Right (BL-LPTS-GT2FXS-RH)	$218.00
Removable	☐	Left (BL-LPTS-GT2RMS-LH)	☐	Right (BL-LPTS-GT2RMS-RH)	$260.00

TT mounting hardware | *Required, must match number of pads selected. If selecting power mount hardware, must choose power chair mounting option.*

	Small *Recommened cushion thickness: 2.0"-3.0"*	**Medium** *Recommened cushion thickness: 3.0"-4.5"*	**Large** *Recommened cushion thickness: 4.0"-6.5"*	
TT hardware, slot mount				
Fixed-Left	☐ BL-LPTS-TT1FXSL1-LH	☐ BL-LPTS-TT1FXSL2-LH	☐ BL-LPTS-TT1FXSL3-LH	$218.00
Fixed-Right	☐ BL-LPTS-TT1FXSL1-RH	☐ BL-LPTS-TT1FXSL2-RH	☐ BL-LPTS-TT1FXSL3-RH	$218.00
Removable-Left	☐ BL-LPTS-TT1RMSL1-LH	☐ BL-LPTS-TT1RMSL2-LH	☐ BL-LPTS-TT1RMSL3-LH	$311.00
Removable-Right	☐ BL-LPTS-TT1RMSL1-RH	☐ BL-LPTS-TT1RMSL2-RH	☐ BL-LPTS-TT1RMSL3-RH	$311.00
TT hardware, power mount				
Fixed-Left	☐ BL-LPTS-TT1FXPWL1-LH	☐ BL-LPTS-TT1FXPWL2-LH	☐ BL-LPTS-TT1FXPWL3-LH	$282.00
Fixed-Right	☐ BL-LPTS-TT1FXPWL1-RH	☐ BL-LPTS-TT1FXPWL2-RH	☐ BL-LPTS-TT1FXPWL3-RH	$282.00
Removable-Left	☐ BL-LPTS-TT1RMPWL1-LH	☐ BL-LPTS-TT1RMPWL2-LH	☐ BL-LPTS-TT1RMPWL3-LH	$368.00
Removable-Right	☐ BL-LPTS-TT1RMPWL1-RH	☐ BL-LPTS-TT1RMPWL2-RH	☐ BL-LPTS-TT1RMPWL3-RH	$368.00

Power mount option

- ☐ Quantum w/ Corpus 3G (LPTS-TT-PW1)
- ☐ Permobil w/ Corpus 3G (LPTS-TT-PW2)
- ☐ Quickie w/ 3.7 Power (LPTS-TT-PW3)
- ☐ Rovi w/ Motion Concepts (LPTS-TT-PW4)
- ☐ Avid Rehab w. Power T & R (LPTS-TT-PW5)

Additional options

☐	Bl-LPTS-TT1LARM4	**TT long extension arm, 4"** *Price per arm* ☐ Left ☐ Right ☐ Extra left ☐ Extra right	$58.00

Bodilink® medial knee/thigh support hardware and pads

Refer to reference pages for details and images.

Pad shape & size | *Priced per pad.*

	Wedge	Oval	
2.5”W x 3.5”D	☐ BL-MKTSP1-3W4D	☐ BL-MKTSP2-3W4D	$203.00
3.5”W x 5”D	☐ BL-MKTSP1-4W5D	☐ BL-MKTSP2-4W5D	$203.00
4”W x 6”D	☐ BL-MKTSP1-4W6D		$203.00
5”W x 6”D		☐ BL-MKTSP2-5W6D	$203.00

Pad cover options | *Required, choose one.*

☐	MKTSP-COMFORT	**COMFORT-TEK** *For fluid protection, an easily cleaned surface*	
☐	MKTSP-STRETCH	**STRETCH-AIR** *For comfort and heat dissipation*	
☐	MKTSP-GLIDE	**GlideWear** *For skin protection & shear reduction. NOT fluid resistant*	$36.00

Swing away hardware

☐	BL-MKTS-ST1SASL1	**Small** *Recommended inferior thigh thickness: 1.0”-3.0”*	$280.00
☐	BL-MKTS-ST1SASL2	**Medium** *Recommended inferior thigh thickness: 2.0”-5.0”*	$280.00
☐	BL-MKTS-ST1SASL3	**Large** *Recommended inferior thigh thickness: 3.0”-7.0”*	$280.00

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- The client's name or code and the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, , am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____


Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.


References

Section 1 | *Cushion accessories*

Removable seat pan kit *(includes hardware)*
Includes a slotted aluminum pan and attaching hardware to accommodate 7/8" or 1" tubing. Kit includes two different cross bars to accommodate various wheelchair frame types. This will fit both folding and non-folding seat rails. Lateral and medial thigh support hardware can easily be attached.




Rigid insert
3/16" rigidizing board.




Section 2 | *Bodilink lateral/pelvic thigh supports*

Premium pad size/shape
Length(L) refers to the actual size dimension of the support from top to bottom edge. Depth(D) refers to the actual size dimension from anterior to posterior edge




Bodilink lateral/pelvic thigh support hardware: Slot mount




Hardware size	Max cushion thickness clearance	Recommended cushion thickness
Small	3.0"	2.0"- 3.0"
Medium	4.5"	3.0"- 4.5"
Large	6.5"	4.0"- 6.5"

Bodilink lateral/pelvic thigh support hardware: Power mount




Hardware size	Max cushion thickness clearance	Recommended cushion thickness
Small	3.5"	2.0"- 3.5"
Medium	5.5"	3.0"- 5.5"
Large	7.5"	4.0"- 7.5"

TT long extension arm 4"





Section 3 | *Bodilink medial knee/thigh supports*

Premium pad size/shape



Swing-away hardware





Size	Recommended inferior thigh thickness
Small	1.00”- 3.00”
Medium	2.00”- 5.00”
Large	3.00”- 7.00”